




## 4R57=WR/4R57=WR-ST

|           |  |     |
|-----------|--|-----|
| <b>DE</b> | Gebrauchsanweisung (Fachpersonal) .....                      | 3   |
| <b>EN</b> | Instructions for use (qualified personnel) .....             | 18  |
| <b>FR</b> | Instructions d'utilisation (Personnel spécialisé) .....      | 32  |
| <b>IT</b> | Istruzioni per l'uso (Personale tecnico specializzato) ..... | 47  |
| <b>ES</b> | Instrucciones de uso (Personal técnico especializado) .....  | 62  |
| <b>PT</b> | Manual de utilização (Pessoal técnico) .....                 | 77  |
| <b>NL</b> | Gebruiksaanwijzing (Vakmensen) .....                         | 92  |
| <b>SV</b> | Bruksanvisning (Fackpersonal) .....                          | 108 |
| <b>DA</b> | Brugsanvisning (Faguddannet personale) .....                 | 122 |
| <b>NO</b> | Bruksanvisning (Fagpersonell) .....                          | 136 |
| <b>FI</b> | Käyttöohje (Ammattihenkilöstö) .....                         | 150 |
| <b>PL</b> | Instrukcja użytkowania (Personel fachowy) .....              | 165 |
| <b>HU</b> | Használati utasítás (szakszemélyzet) .....                   | 180 |
| <b>CS</b> | Návod k použití (Odborný personál) .....                     | 194 |
| <b>RO</b> | Instrucțiuni de utilizare (Personal de specialitate) .....   | 209 |
| <b>HR</b> | Upute za uporabu (Stručno osoblje) .....                     | 224 |
| <b>SL</b> | Navodila za uporabo (Strokovno osebje) .....                 | 238 |
| <b>SK</b> | Návod na používanie (Odborný personál) .....                 | 251 |
| <b>BG</b> | Инструкция за употреба (Специалисти) .....                   | 265 |
| <b>TR</b> | Kullanma talimatı (Uzman personel) .....                     | 281 |
| <b>EL</b> | Οδηγίες χρήσης (Τεχνικό προσωπικό) .....                     | 295 |
| <b>RU</b> | Руководство по применению (Квалифицированный персонал) ...   | 311 |
| <b>JA</b> | 取扱説明書(有資格担当者) .....  | 326 |
| <b>ZH</b> | 使用说明书(专业人员) .....  | 340 |
| <b>KO</b> | 사용 설명서 (전문기술요원) .....  | 352 |

## INFORMATION

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### INFORMATION

Date of last update: 2020-10-23

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- ▶ Please keep this document for your records.

The product "4R57=WR/4R57=WR-ST rotation adapter" is referred to as the product/adapter/rotation adapter below.

These instructions for use provide you with important information on the use, adaptation and handling of the product.

Only put the product into use in accordance with the information contained in the accompanying documents supplied.

## 2 Product description

### 2.1 Function

The product is used as a component for modular lower limb prostheses.

Through incorporation of the rotation adapter above the knee joint, the lower leg may be rotated medially or laterally relative to the socket with the knee flexed.

This makes everyday activities such as putting on shoes and changing socks much easier and allows the user to sit comfortably. The sitting position can be varied up to sitting cross-legged.

The rotation mechanism is activated by pressing the release button (see fig. 1, (1)). The position of the button is marked by a deeper profile of the protective cap. Upon reaching the initial position, the rotation movement is locked automatically.

## 2.2 Combination possibilities

This prosthetic component is compatible with Ottobock's system of modular connectors. Functionality with components of other manufacturers that have compatible modular connectors has not been tested.

### Prohibited combination possibilities

The rotation adapter **cannot be combined** with the following Ottobock components:

- 2R49, 2R50, 4R22, 4R37, 4R51, 4R55, 4R69, 4R72=32, 4R95, 4R98

## 3 Intended use

### 3.1 Indications for use

The product is to be used **solely** for lower limb exoprosthetic fittings above the prosthetic knee joint.

### 3.2 Conditions of use

The product was developed for everyday use and must not be used for unusual activities. These unusual activities include, for example, extreme sports (free climbing, parachuting, paragliding, etc.).

Permissible ambient conditions are described in the technical data (see page 30).

Our components perform optimally when paired with appropriate components based upon weight and mobility grades identifiable by our MOBIS classification information and which have appropriate modular connectors.

The product is recommended for mobility grade 1 (indoor walker), mobility grade 2 (restricted outdoor walker), mobility grade 3 (unrestricted outdoor walker) and mobility grade 4 (unrestricted outdoor walker with particularly high demands). Approved for a body weight of **max. 166 kg**.

### 3.3 Indications

- For patients with a transfemoral amputation, hip disarticulation
- For unilateral or bilateral amputation

### 3.4 Contraindications

- All conditions which contradict or go beyond the specifications listed in the section on "Safety" and "Indications for use".

### 3.5 Qualification

The fitting of a patient with the product may only be carried out by an O&P professional.

## 4 Safety

### 4.1 Explanation of warning symbols

 **WARNING** Warning regarding possible serious risks of accident or injury.

 **CAUTION** Warning regarding possible risks of accident or injury.

 **NOTICE** Warning regarding possible technical damage.

### 4.2 Structure of the safety instructions

#### **CAUTION**

##### **The heading describes the source and/or the type of hazard**

The introduction describes the consequences in case of failure to observe the safety instructions. Consequences are presented as follows if more than one consequence is possible:

- > E.g.: Consequence 1 in case of failure to observe the hazard
- > E.g.: Consequence 2 in case of failure to observe the hazard
- ▶ This symbol identifies activities/actions that must be observed/carried out in order to avert the hazard.

### 4.3 General safety instructions

#### **WARNING**

##### **Non-observance of safety notices**

Personal injury/damage to the product due to using the product in certain situations.

- ▶ Observe the safety notices and the stated precautions in this accompanying document.

#### **WARNING**

##### **Fitting in the presence of children**

Injury due to swallowing small parts.

- ▶ Never leave children unattended during the fitting.

#### **CAUTION**

##### **Remaining in areas outside the allowable temperature range**

Falling due to breakage of load-bearing components or malfunction of the product.

- ▶ Avoid remaining in areas with temperatures outside of the permissible range (see page 30).

**⚠ CAUTION**

**Independent manipulation of the product**

Falling due to breakage of load-bearing components or malfunction of the product.

- ▶ Manipulations to the product other than the tasks described in these instructions for use are not permitted.
- ▶ Only Ottobock Service Centres may open the product or replace damaged components.

**NOTICE**

**Use with damaged or missing protective cap of the release button**

Limited functionality of the product

- ▶ Replace a damaged or missing protective cap for the release button.
- ▶ Install a new protective cap as described in the section "Maintenance" (see page 28).

#### 4.4 Information on Alignment/Adjustment

**⚠ CAUTION**

**Incorrect alignment or assembly**

Falling due to damage to the prosthesis components.

- ▶ Observe the alignment and assembly instructions.

**⚠ CAUTION**

**Collision of the rotation adapter with the components (hydraulic unit, frame) of the knee joint**

- > Falling due to unexpected behaviour of the knee joint as the result of malfunction.
- > Falling due to breakage of load-bearing components.
- ▶ In case of contact between the rotation adapter and the hydraulic unit and/or frame of the knee joint at maximum flexion, check whether the contact was caused by extreme tilting of the rotation adapter. In this case install the rotation adapter without tilting.

- ▶ If contact between the rotation adapter and the hydraulic unit and/or frame of the knee joint occurs anyway, the knee joint has to be equipped with a flexion stop or the existing flexion stop has to be replaced (see the instructions for use of the knee joint).

**⚠ CAUTION**

**Failure to observe the torque values of the screw connections**

Falling due to breakage of load-bearing components caused by screw connections coming loose.

- ▶ Observing the specified tightening torques is essential for safe operation.
- ▶ The specified tightening torques **must not** be exceeded.

**⚠ CAUTION**

**Incorrectly secured screws**

Falling due to breakage of load-bearing components caused by screw connections coming loose.

- ▶ After completing all settings, the set screws in the rotation adapter must be secured with adhesive before they are tightened to the specified torque.

**⚠ CAUTION**

**Use of unsuitable prosthetic components**

Falling due to unexpected behaviour of the product or breakage of load-bearing components.

- ▶ Use the product only in combination with components approved by Ottobock (see instructions for use and catalogues).
- ▶ Note the prohibited combination possibilities in the section "Combination possibilities".
- ▶ If the product is to be used in water, verify that each prosthetic component is waterproof.

**⚠ CAUTION**

**Improper assembly of the screw connections**

Falling due to breakage or loosening of the screw connections.

- ▶ Clean the threads before every installation.
- ▶ Apply the specified tightening torque values for installation (see the section "Technical data" see page 30).

- ▶ Observe the instructions for securing the screw connections and the use of the correct length.

#### NOTICE

##### Improper product care

Damage to the product due to the use of incorrect cleaning agents and/or disinfectants.

- ▶ Only clean the product with a damp cloth and mild soap (e.g. 453H10=1 Ottobock DermaClean).
- ▶ Disinfect the product with 634A58 isopropyl alcohol when needed.
- ▶ Do not use disinfectants based on hydrogen peroxide and silver salt to disinfect the product.
- ▶ Only use fresh water with a temperature of no more than 65°C for cleaning.
- ▶ If dirt cannot be removed, the product has to be sent to an authorised Ottobock Service Centre.

## 4.5 Patient information

#### CAUTION

##### Penetration of dirt and humidity into the product

- > Falling due to unexpected product behaviour as the result of a malfunction.
- > Falling due to breakage of load-bearing components.
- ▶ Ensure that no solid particles or foreign objects can penetrate into the product.
- ▶ The product is waterproof and resistant to corrosion, and can therefore be operated in fresh water as well as salt water and chlorinated water. Do not use the product under extreme conditions like diving or jumping into water. The product is designed for underwater use.
- ▶ After using the product in salt water or chlorinated water, rinse the product with fresh water, dry it with a lint-free cloth and allow it to air dry fully.
- ▶ If the product comes into contact with **solutions other than fresh water, salt water or chlorinated water, promptly rinse the product with fresh water** and allow it to dry.
- ▶ In case of a malfunction after drying, the product must be inspected by an authorised Ottobock Service Centre. The O&P professional is your contact person.
- ▶ The product is not resistant to penetration from water jets or steam.

**⚠ CAUTION**

**Overloading due to unusual activities**

Falling due to breakage of load-bearing components.

- ▶ The product was developed for everyday use and must not be used for unusual activities. These unusual activities include, for example, extreme sports (free climbing, paragliding, etc.).
- ▶ Careful handling of the product and its components not only increases their service life but, above all, ensures your personal safety!
- ▶ If the product and its components have been subjected to extreme loads (e.g. due to a fall, etc.), then the product must be inspected for damage immediately by an O&P professional. If necessary, he or she will forward the product to an authorised Ottobock Service Centre.

**⚠ CAUTION**

**Remaining in areas outside the allowable temperature range**

Falling due to breakage of load-bearing components or malfunction of the product.

- ▶ Avoid remaining in areas with temperatures outside of the permissible range (see page 30).

**⚠ CAUTION**

**Independent manipulation of the product**

Falling due to breakage of load-bearing components or malfunction of the product.

- ▶ Manipulations to the product other than the tasks described in these instructions for use are not permitted.
- ▶ Replacing a damaged protective cap may only be carried out by an O&P professional.
- ▶ Only Ottobock Service Centres may open the product or replace damaged components (except for the protective cap).

**⚠ CAUTION**

**Accidental activation of the release button**

Falling due to twisting of the prosthesis under load.

- ▶ Only activate the release button of the rotation adapter deliberately and with due consideration of the respective situation.

**⚠ CAUTION****Use with the rotation adapter unlocked**

Falling due to twisting of the prosthesis under load.

- ▶ Before loading the prosthesis, check whether the rotation adapter is correctly locked in the initial position.
- ▶ If correct locking is no longer possible, continued use of the rotation adapter is prohibited. The rotation adapter must be inspected by an authorised Ottobock Service Centre. The O&P professional is your contact person.

**NOTICE****Use with damaged or missing protective cap of the release button**

Limited functionality of the product.

- ▶ Have a damaged or missing protective cap replaced by the O&P professional.
- ▶ Do not perform the replacement independently.

**NOTICE****Improper product care**

Damage to the product due to the use of incorrect cleaning agents and/or disinfectants.

- ▶ Only clean the product with a damp cloth and mild soap (e.g. 453H10=1 Ottobock DermaClean).
- ▶ Disinfect the product with 634A58 isopropyl alcohol when needed.
- ▶ Do not use disinfectants based on hydrogen peroxide and silver salt to disinfect the product.
- ▶ Only use fresh water with a temperature of no more than 65°C for cleaning.
- ▶ If dirt cannot be removed, the product has to be sent to an authorised Ottobock Service Centre. The O&P professional is your contact person.

## 5 Scope of Delivery and Accessories

### 5.1 Scope of delivery

- 1 pc. 4R57=WR rotation adapter or  
1 pc. 4R57=WR-ST rotation adapter
- 1 pc. Instructions for use (qualified personnel)

## 5.2 Accessories

The following components are not included in the scope of delivery and must be ordered separately:

- 4X46=ST lamination dummy (for 4R57\*=ST)

## 5.3 Spare parts

- 4Y492 protective cap

# 6 Preparing the product for use

## 6.1 Required tools

- 710D4 torque wrench (up to 30 Nm)
- 636K13 adhesive (Loctite 241)

## 6.2 Mounting the adapter

### 6.2.1 Rotation adapter with pyramid adapter

The pyramid adapter of the rotation adapter is fixed to the pyramid receiver with the set screws of the proximal prosthetic component. All four set screws of the rotation adapter have to be unscrewed for replacement or disassembly. The adjustment possibilities for static corrections in this area are minor.

#### INFORMATION

Deviating from the specification in the corresponding instructions for use of the proximal prosthetic component, the **set screws** of the prosthetic component must be tightened to a **maximum of 10 Nm**.

The following points have to be observed during assembly:

- The alignment of the release button must be medial.
- After activating the release button, it must be possible to rotate the housing of the rotation adapter 360° without touching the socket and/or components of the knee joint (such as the frame or hydraulic unit).
- A collision with the frame or hydraulic unit at maximum flexion of the knee joint must be prevented. If necessary, flexion stops have to be installed or a tilt of the rotation adapter has to be eliminated.

#### Fitting:

- ▶ Use the torque wrench to tighten the set screws of the pyramid receiver to **10 Nm**.

#### Definitive mounting:

- 1) Secure the set screws of the pyramid receiver with adhesive.
- 2) Use the torque wrench to tighten the set screws to **10 Nm**.

## 6.2.2 Rotation adapter with threaded connector

### CAUTION

#### **Insufficient screw-in depth in the lamination anchor**

Falling due to breakage of load-bearing components.

- ▶ After fully screwing the rotation adapter into the lamination anchor, unscrew the rotation adapter **out of the thread by no more than 3/4 turn**.

### INFORMATION

#### **4X46=ST lamination dummy for lamination process**

In order to properly screw the rotation adapter with threaded connector into the lamination anchor, the 4X46=ST lamination dummy must be used for laminating. It must be ordered separately.

- 1) Screw the rotation adapter completely into the thread of the lamination anchor.
- 2) Unscrew the rotation adapter **out of the thread by no more than 3/4 turn**.
  - The alignment of the release button must be medial
  - After activating the release button, it must be possible to rotate the housing of the rotation adapter by 360° without touching the lamination anchor.
- 3) Secure the cap screw of the threaded clamping of the lamination anchor using adhesive and tighten using the torque wrench (**10 Nm**).

## 6.3 Installing prosthetic components

The pyramid adapter of the prosthetic component is fixed with the set screws of the pyramid receiver of the rotation adapter.

### **Fitting:**

- ▶ Use the torque wrench to tighten the set screws of the pyramid receiver to **10 Nm**.

### **Definitive mounting:**

- 1) Secure the set screws with the adhesive.
- 2) Pre-tighten the set screws with the torque wrench to **10 Nm** and then tighten to **12 Nm**.

## Alignment

The set screws in the pyramid receiver can be used to make static adjustments during alignment, trial fittings and after the prosthesis is finished.

## Replacement and disassembly

The set position of the prosthetic component can be maintained during replacement or disassembly. In order to do this, unscrew the two set screws that are screwed in the furthest but are not opposite each other.

## Replacing set screws

Replace any set screws that are protruding or recessed too much with suitable ones. For approved set screws, see the section "Technical data" (see page 30).

## 7 Handling

- 1) To release, press the release button on the side (see fig. 1 (1)) and rotate the prosthesis. The position of the button is marked by a deeper profile of the protective cap.
  - Subsequently the prosthesis can be rotated from the initial position without holding the button.
- 2) Upon reaching the initial position, the rotation movement is locked automatically.

## 8 Maintenance

Annual safety inspections are recommended in the interest of the patient's safety and in order to maintain operational safety.

### 8.1 Replacing the protective cover for the release button

#### CAUTION

#### **Non-compliance with safety data sheets**

Injury due to skin, eye and respiratory tract irritation.

- ▶ Read the appropriate safety data sheets before working with 634A58 isopropyl alcohol, 636W145 adhesive and 636P4 primer.
- ▶ Wear suitable protective gloves during processing.

#### **INFORMATION**

#### **Inspection for contamination or damage**

Check the following points after removing the protective cap:

- ▶ Condition of the groove (a) in the housing (Gummikappe abziehen).
- ▶ Penetration of foreign bodies in the area around the release button.

The protective cap is not a seal to prevent the penetration of liquids, but protects against coarse contaminants (mud, sand etc.). If the protective cap comes loose during the service life or if the cap is damaged, it has to be replaced with a new protective cap.

> **Recommended tools:**

4Y492 protective cap, tweezers with rounded tips, Loctite 454 adhesive (636W145), Loctite 770 primer (636P4), isopropyl alcohol (634A58), protective gloves.

- 1) Remove the prosthesis.
- 2) Use tweezers with rounded tips to remove the glued protective cap from the housing of the rotation adapter.
- 3) Remove the protective cap.
- 4) Clean the area around the release button and the groove in the housing and check for damage.
- 5) Use sandpaper (220 grit or finer) to sand adhesive residue off the rotation adapter (see fig. 4).
- 6) Clean the inside of the new protective cap (see fig. 2) and the exterior gluing surface of the rotation adapter (see fig. 3) with isopropyl alcohol.
- 7) After the marked areas are dry, apply one coat of primer to the new protective cap (see fig. 2) and the rotation adapter (see fig. 3).
- 8) Allow the primer to dry for at least 60 seconds and no more than 60 minutes.
- 9) Apply a few drops of adhesive within the marked areas (towards the groove in the housing) of the rotation adapter (see fig. 4).
- 10) Slide on the new protective cap with proper alignment (see fig. 5) and lightly press the protective cap into place for about 15 seconds at the positions of the adhesive dots.

## 9 Cleaning

- 1) Clean the product with a damp cloth and soap.
- 2) Rinse the product with clean, fresh water.
- 3) Dry the product with a soft cloth.
- 4) **INFORMATION: Do not store the product while moist.**  
Allow to air dry in order to remove residual moisture.

## 10 Disposal

This product may not be disposed of with regular domestic waste in all jurisdictions. Disposal that is not in accordance with the regulations of the country where the product is used may have a detrimental impact on health and the environment. Please observe the information provided by the responsible

authorities in the country of use regarding return, collection and disposal procedures.

## 11 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

### 11.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

### 11.2 Trademarks

All product names mentioned in this document are subject without restriction to the respective applicable trademark laws and are the property of the respective owners.

All brands, trade names or company names may be registered trademarks and are the property of the respective owners.

Should trademarks used in this document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

### 11.3 CE conformity

Otto Bock Healthcare Products GmbH hereby declares that the product is in compliance with European requirements for medical devices.

The full text of the regulations and requirements is available at the following Internet address: <http://www.ottobock.com/conformity>

## 12 Technical data

| Environmental conditions                    |   |
|---|---|
| Storage and transport in original packaging | -25 °C/-13 °F to +70 °C/+158 °F<br>Max. 93% relative humidity, non-condensing |
| Storage and transport without packaging     | -25 °C/-13 °F to +70 °C/+158 °F<br>Max. 93% relative humidity, non-condensing |
| Operation                                   | -10 °C/+14 °F to +60 °C/+140 °F<br>Max. 93% relative humidity, non-condensing |

| <b>Environmental conditions</b> |  |
|---------------------------------|--|
| Waterproof                      | Fresh water, salt water or chlorinated water |
| Maximum water depth             | 3 m  |

| <b>Product</b>  | <b>4R57=WR</b>   | <b>4R57=WR-ST</b> |
|---|------------------|-------------------|
| Weight [g]  | 214              | 253               |
| Proximal connection   | Pyramid adapter  | Thread            |
| Torque moment for the set screw/clamp screw of the proximal connection [Nm] | 10               |                   |
| Distal connection   | Pyramid receiver |                   |
| Max. body weight [kg]   | 166              |                   |
| System height [mm]  | 25               | 47                |
| Material  | Coated steel     |                   |
| Rotation range [°]  | 360              |                   |
| Service life [years]  | 5                |                   |

| <b>Approved set screws</b> | <b>506G3= M8x12-ZNNI</b> | <b>506G3= M8x16 ZN</b> |
|----------------------------|--------------------------|------------------------|
| Length [mm]                | 12                       | 16                     |
| Maximum torque [Nm]        | 12                       |                        |

### **Torque values of the screw connections**

Using a torque wrench, tighten the corresponding screws alternately in several cycles until the specified tightening torque is reached.

| <b>Screw connection</b>                                | <b>Tightening torque</b> |
|--|--------------------------|
| Proximal prosthesis components with threaded connector | 10 Nm / 89 lbf. In.      |
| Proximal prosthesis components with pyramid adapter    | 10 Nm / 89 lbf. In.      |
| Pyramid receiver of the rotation adapter               | 12 Nm / 133 lbf. In.     |

## 13 Appendices

### 13.1 Symbols Used

Declaration of conformity according to the applicable European directives

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