## 8E33=\*, 8E34=\*

<b>DE</b> Gebrauchsanweisung (Benutzer) .....	3
<b>EN</b> Instructions for use (user) .....	16
<b>FR</b> Instructions d'utilisation (Utilisateur) .....	29
<b>IT</b> Istruzioni per l'uso (Utilizzatore) .....	43
<b>ES</b> Instrucciones de uso (Usuario) .....	56
<b>PT</b> Manual de utilização (Usuário) .....	70
<b>NL</b> Gebruiksaanwijzing (Gebruiker) .....	84
<b>SV</b> Bruksanvisning (Användare) .....	98
<b>JA</b> 取扱説明書(ユーザー用) .....	111

## INFORMATION

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### INFORMATION

Date of last update: 2020-10-19

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Obtain instruction from the qualified personnel in the safe use of the product.
- ▶ Please contact the qualified personnel if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

The 8E33=9-1, 8E34=9-1 System Electric Greifers are referred to as the product/prosthesis/terminal device in the following.

These instructions for use provide you with important information on the use, adaptation and handling of the product.


Only put the product into use in accordance with the information contained in the accompanying documents supplied.

## 2 Product description

### 2.1 Design

The product consists of the following components:

- 1. Fingertips**  
The position can be rotated by  $\pm 60^\circ$  to adapt to various tasks.  
For other types of fingertips, see the section "Accessories" (see page 22)
- 2. Safety lever**  
To unlock the prehensile arm drive (see page 26)
- 3. Dial**  
For manual opening/closing
- 4. On/off switch**  
To turn the electronics and flexion joint on/off  
Position of the rocker switch towards the prosthetic socket: electronics switched on and flexion brake released  
Position of the rocker switch towards the terminal device: electronics switched off and flexion joint blocked
- 5. Quick-disconnect wrist unit** (only with 8E33\*)  
Mechanical connection to the lamination ring

### 2.2 Function

The product is a myoelectrically controlled terminal device of the MyoBock prosthesis system. Various control programmes make optimal adaptation to the individual needs and abilities possible.

The terminal device is based on the Ottobock DMC system (DMC=Dynamic Mode Control): This system, which has been developed by Ottobock, uses two independent measurement and control systems in order to optimally control the gripping speed and gripping force in accordance with the muscle signal.

## 3 Intended use

### 3.1 Indications for use

The product is intended **exclusively** for upper limb exoprosthetic fittings.

### 3.2 Conditions of use

The product was developed for everyday use and must not be used for unusual activities. These unusual activities include, for example, sports which involve excessive strain and/or shocks to the wrist joint (pushups, downhill racing, mountain biking, ...) or extreme sports (free climbing, paragliding, etc.). Furthermore, the product should not be used to operate motor vehicles, heavy equipment (e.g. construction machines), industrial machines or motor-driven equipment.

Permissible ambient conditions are described in the technical data (see page 27).

The product is intended **exclusively** for use on **one** user. Use of the product by another person is not approved by the manufacturer.

### 3.3 Indications

- Amputation level below-elbow, above-elbow and shoulder disarticulation
- For unilateral or bilateral amputation
- Dysmelia of the forearm or upper arm

- The user must be able to understand usage instructions and safety notices and put them into practice.
- The user must fulfil the physical and mental requirements for perceiving optical/acoustic signals and/or mechanical vibrations

### 3.4 Contraindications

- All conditions which contradict or go beyond the specifications listed in the section on "Safety" and "Indications for use".

### 3.5 Qualification

The fitting of a user with the product may only be carried out by O&P professionals who have been authorised through the corresponding Ottobock training.

## 4 Safety

### 4.1 Explanation of warning symbols

 <b>WARNING</b>	Warning regarding possible serious risks of accident or injury.
 <b>CAUTION</b>	Warning regarding possible risks of accident or injury.
 <b>NOTICE</b>	Warning regarding possible technical damage.

### 4.2 Structure of the safety instructions

 <b>WARNING</b> <b>The heading describes the source and/or the type of hazard</b> The introduction describes the consequences in case of failure to observe the safety instructions. Consequences are presented as follows if more than one consequence is possible: > E.g.: Consequence 1 in the event of failure to observe the hazard > E.g.: Consequence 2 in the event of failure to observe the hazard ► This symbol identifies activities/actions that must be observed/carried out in order to avert the hazard.
---

### 4.3 General safety instructions

 <b>WARNING</b> <b>Non-observance of safety notices</b> Personal injury/damage to the product due to using the product in certain situations. ► Observe the safety notices and the stated precautions in this accompanying document.
---

 <b>WARNING</b> <b>Operating a vehicle with the prosthesis</b> Accidents due to unexpected actions of the prosthesis. ► The prosthesis should not be used for the operation of motor vehicles and heavy equipment (e.g. construction machines).
--

 <b>WARNING</b> <b>Using the prosthesis while operating machines</b> Injury due to unexpected actions of the prosthesis. ► The prosthesis should not be used to operate industrial machines or motor-driven equipment.
---

**⚠ WARNING**

**Operating the prosthesis near active implanted systems**

Interference with active implantable systems (e.g., pacemaker, defibrillator, etc.) due to electromagnetic emissions of the prosthesis.

- ▶ When operating the prosthesis in the immediate vicinity of active implantable systems, ensure that the minimum distances stipulated by the manufacturer of the implant are observed.
- ▶ Make sure to observe any operating conditions and safety instructions stipulated by the manufacturer of the implant.

**⚠ CAUTION**

**Signs of wear on the product**

Injury due to faulty control or malfunction of the product

- ▶ In the interest of personal safety and in order to maintain operating reliability, the prosthesis has to be inspected by a O&P professional in case of noticeable restrictions of functionality.
- ▶ Note that the functionality of the terminal device may be limited if the battery charge level gets too low.

**⚠ WARNING**

**Use of damaged power supply unit, adapter plug or battery charger**

Risk of electric shock due to contact with exposed, live components.

- ▶ Do not open the power supply unit, adapter plug or battery charger.
- ▶ Do not expose the power supply unit, adapter plug or battery charger to extreme loading conditions.
- ▶ Immediately replace damaged power supply units, adapter plugs or battery chargers.

**⚠ CAUTION**

**Independent manipulation of the product**

Injury due to product damage or malfunction.

- ▶ Manipulations to the product other than the tasks described in these instructions for use are not permitted.
- ▶ The product and any damaged components may only be opened and repaired by authorised, qualified Ottobock personnel.

#### 4.4 Information on Proximity to Certain Areas

**⚠ CAUTION**

**Proximity to sources of strong magnetic or electrical interference (e.g. theft prevention systems, metal detectors)**

Injury due to unexpected behaviour of the prosthesis system caused by interference with internal data communication.

- ▶ Avoid remaining in the vicinity of visible or concealed theft prevention systems at the entrance/exit of stores, metal detectors/body scanners for persons (e.g. in airports) or other sources of strong magnetic and electrical interference (e.g. high-voltage lines, transmitters, transformer stations, computer tomographs, magnetic resonance tomographs, etc.).
- ▶ When walking through theft prevention systems, body scanners or metal detectors, watch for unexpected behaviour of the prosthetic system.

**⚠ CAUTION**

**Distance to HF communication devices is too small (e.g. mobile phones, Bluetooth devices, WiFi devices)**

Injury due to unexpected behaviour of the prosthesis system caused by interference with internal data communication.

- ▶ Therefore, keeping the following minimum distances to these HF communication devices is recommended:
  - Mobile phone GSM 850 / GSM 900: 0.99 m
  - Mobile phone GSM 1800 / GSM 1900 / UMTS: 0.7 m
  - DECT cordless phones incl. base station: 0.35 m
  - WiFi (routers, access points ...): 0.22 m
  - Bluetooth devices (third-party products not approved by Ottobock): 0.22 m

**⚠ CAUTION**

**Wearing in extreme ambient temperatures**

Injury due to faulty control or malfunction of the Axon-Bus prosthetic system.

- ▶ Avoid wearing the prosthesis in areas with temperatures outside of the permissible range (see page 27).

#### 4.5 Information on Use

**⚠ WARNING**

**Operating the prosthesis near active implanted systems**

Interference with active implantable systems (e.g., pacemaker, defibrillator, etc.) due to electromagnetic emissions of the prosthesis.

- ▶ When operating the prosthesis in the immediate vicinity of active implantable systems, ensure that the minimum distances stipulated by the manufacturer of the implant are observed.
- ▶ Make sure to observe any operating conditions and safety instructions stipulated by the manufacturer of the implant.

**⚠ CAUTION**

**Using the prosthesis with pointed or sharp objects (e.g. knives in the kitchen)**

Injury due to unintentional movements.

- ▶ Use extra caution when using the prosthesis for handling pointed or sharp objects.

**⚠ CAUTION**

**Risk of pinching between the prehensile arms/fingertips**

Injury due to pinching of body parts.

- ▶ Ensure that no body parts are between the prehensile arms/fingertips when using the product.

**⚠ CAUTION**

**Overloading due to unusual activities**

Injury due to unexpected product behaviour as a result of a malfunction.

- ▶ The product was developed for everyday use and must not be used for unusual activities. These unusual activities include, for example, sports with excessive strain and/or shocks to the wrist joint (pushups, downhill racing, mountain biking, ...) or extreme sports (free climbing, paragliding, etc.).

- ▶ Careful handling of the product and its components not only increases their service life but, above all, ensures your personal safety!
- ▶ If the product and its components have been subjected to extreme loads (e.g. due to a fall, etc.), then the product must be inspected for damage immediately by an O&P professional. If necessary, he or she will forward the product to an authorised Ottobock Service Centre.

**⚠ CAUTION**

**Mechanical stress on the product**

Injury due to faulty control or malfunction of the product components.

- ▶ Do not subject the product to mechanical vibrations or impacts.
- ▶ Check the product for visible damage before each use.

**⚠ CAUTION**

**Penetration of dirt and humidity into the prosthetic components**

Injury due to faulty control or malfunction of the prosthetic components.

- ▶ Ensure that neither solid particles nor liquids can penetrate into the product or prosthetic components (e.g. the terminal device).

**NOTICE**

**Improper care of the prosthesis system**

Damage to the prosthesis system due to the use of incorrect cleaning agents.

- ▶ Clean the prosthesis system with a damp cloth and mild soap only (e.g. 453H10=1 Ottobock DermaClean).

**⚠ CAUTION**

**Changing prosthetic components when switched on**

Injury due to faulty control or malfunction of the prosthesis system.

- ▶ Before changing prosthetic components (such as the terminal device), remove the rechargeable battery from the battery receptacle or switch off the prosthesis system.

#### 4.6 Information on the Power Supply/Battery Charging

**⚠ CAUTION**

**Charging the product with soiled or damaged contacts**

Injury due to unexpected behaviour of the product caused by insufficient charging.

- ▶ Make sure that the contacts are always clean and free of grease.
- ▶ Clean the electrical contacts of the charging plug and charging receptacle regularly using cotton swabs and a mild soap solution.
- ▶ Take care to avoid damaging the contacts with pointed or sharp objects.

**NOTICE**

**Use of incorrect power supply unit/battery charger**

Damage to product due to incorrect voltage, current or polarity.

- ▶ Use only power supply units/battery chargers approved for this product by Ottobock (see instructions for use and catalogues).

## 4.7 Notices regarding battery charger with magnetic charging plug

### WARNING

#### **Charging the prosthesis without taking it off**

Risk of electric shock due to defects in the power supply unit or in the battery charger.

- ▶ For safety reasons, remove the prosthesis prior to charging it.

### WARNING

#### **Storing/transporting the battery charger near active implanted systems**

Interference with active implantable systems (e.g. pacemaker, defibrillator, etc.) due to the magnetic field of the charging plug.

- ▶ When storing/transporting the battery charger in the immediate vicinity of active implantable systems, ensure that the minimum distances stipulated by the manufacturer of the implant are observed.
- ▶ Make sure to observe any operating conditions and safety notices stipulated by the manufacturer of the implant.

### NOTICE

#### **Contact of the charging plug with magnetic data carriers**

Wiping of the data carrier.

- ▶ Do not place the charging plug on credit cards, diskettes, audio or video cassettes.

## 4.8 Notice for using the torch

### WARNING

#### **Use of the incorrect battery type and/or failure to note the polarity**

Explosion due to escaping gases and their self-ignition caused by short circuit and excessive temperature.

- ▶ Only use type CR2032 batteries.
- ▶ When inserting the two batteries, note the polarity information on the cover of the battery compartment.

### CAUTION

#### **Looking directly into the beam of light from the torch**

Eye damage due to looking directly into the beam of light.

- ▶ Do not use the torch for examinations, for instance of the eyes (pupil test).
- ▶ Do not aim the beam of light directly at the eyes of people or animals.
- ▶ Do not look directly into the beam of light.
- ▶ If the beam of light is aimed at the eyes, be sure to close your eyes and immediately move the head out of the beam of light.

## 5 Scope of Delivery and Accessories

### 5.1 Scope of delivery

- 2 pc. CR2032 lithium battery
- 1 pc. 8E33=9-1 System Electric Greifer DMC VariPlus or  
1 pc. 8E34=9-1 System Electric Greifer DMC VariPlus
- 1 pc. Instructions for use (user)

## 5.2 Accessories

### Fingertips

- 9S138 tip set
- 9S145-1 Rubber gripping pad set
- 9S146=PAA-1 Tip padding (pair)
- 9S149 fingertip blank set
- 9S234=PAA tip pair
- 9S278=PAA tip pair

## 6 Charging the battery

The following points must be observed when charging the battery:

- The capacity of a fully charged battery is sufficient for one full day.
- We recommend charging the product every day when used on a daily basis.
- The battery should be charged for at least 3 hours prior to initial use.

### 6.1 757L20 battery charger and 757B2\* EnergyPack

#### NOTICE

#### Deep discharge of the battery

Irreparable damage to the product and resulting loss of functionality in case of irregular and incomplete charging of the battery.

- ▶ In case of long-term storage, fully charge the battery every 4 to 6 months.

- 1) Plug the charging station with power supply unit into the wall socket.  
→ The green LED on the power supply unit and the green LED in the middle of the LED row light up.
- 2) Insert one or two rechargeable batteries into the charging station.
- 3) The green LED for the respective charging socket lights up ● and the charging process starts.
- 4) When the rechargeable battery is fully charged, the green LED for the respective charging socket flashes ✨.
- 5) Remove the rechargeable battery after the charging process is complete.

#### 6.1.1 Display of the current charge level

When inserting/engaging the rechargeable battery in the prosthesis, a charge level indicator is activated on the battery for a few seconds.

LED indicator	Event
●	Charge level above 50% (illuminated in green)
● and ●	Charge level below 50% (flashing alternating green and orange)
●	Charge level below 5% (illuminated in orange)

## 6.1.2 Troubleshooting

Charging socket	LED for the charging socket	Event
Empty	● or ● or ●	<b>The battery charger is defective</b> The battery charger and power supply must be inspected by an authorised Ottobock Service Centre.
Empty or rechargeable battery inserted		
Rechargeable battery inserted		<b>The rechargeable battery is charging</b>
		<b>The rechargeable battery is fully charged</b>
		<b>The temperature of the rechargeable battery is too high</b> Remove the rechargeable battery and let it cool down
		<b>The inserted rechargeable battery is defective</b> The rechargeable battery must be inspected by an authorised Ottobock Service Centre.

## 6.2 757L35 battery charger and 757B35=\* rechargeable battery

### 6.2.1 Connect battery charger to the product

- 1) Connect the charging plug to the charging receptacle on the prosthetic socket.  
→ A correct connection between the battery charger and the prosthesis is indicated by feedback (Status signals).
- 2) The charging process starts.  
→ The prosthesis is switched off automatically.
- 3) Disconnect the prosthesis when the charging process is complete.

#### INFORMATION

The product cannot be used as long as the charging plug is connected to the charging receptacle. The product is deactivated for the duration of charging.

### 6.2.2 Safety Shutoff

The purpose of the prosthesis safety shutoff is to protect the battery; it is triggered in case of:

- Excessively high or low temperature
- Overvoltage or undervoltage
- Short circuit

After a short circuit, the charging plug must be plugged into the charging receptacle and then unplugged in order for the electronics to be activated.

### 6.2.3 Display of the current charge level

The battery charge level can be queried at any time.

- 1) With the prosthesis switched on, press the charging receptacle button and hold for less than one second.

2) The LED display on the charging receptacle provides information on the current battery charge level.


Charging receptacle	Event
	Battery fully charged (illuminated in green)
	Battery 50% charged (illuminated in yellow)
	Battery drained (illuminated in orange)

### 6.2.4 Beep signals

Beep signal	Additional display	Event
1 x long	-	<ul style="list-style-type: none"> <li>Turn off prosthesis on charging receptacle</li> <li>Charging begins (charging plug connected to charging receptacle)</li> <li>Charging ends (charging plug disconnected from charging receptacle)</li> </ul>
2 x short	LED on charging receptacle lights up briefly	Switch on prosthesis
3 x short	-	Battery voltage too low, prosthesis shuts off automatically

## 7 Use

### 7.1 Switching the terminal device on and off

#### Switching off

Move the rocker switch to position 1 until it engages.

The electronics of the terminal device are switched off and the flexion joint is blocked.


#### Switching on

Move the rocker switch to position 2 until it engages.

The electronics of the terminal device are switched on and the flexion brake is released. When the flexion brake is released, the terminal device can be moved by 45 ° in both directions.

### 7.2 Torch

(for terminal devices with built-in torch)

1. **2x powerful LEDs**  
to light the working environment

2. **On/off switch**  
to turn the LEDs on/off

**Switching on:** Press the button to turn the LEDs on. If this button is pressed immediately after switching on, the brightness can be changed in 3 levels by pressing it again.

**Switching off:** About 5 seconds after switching on, press the button again to turn the LEDs off.

After about 15 minutes, the torch is switched off automatically.

3. **Battery compartment**  
See the next section for battery replacement


## Inserting the batteries (2x CR2032)

### INFORMATION

With a bilateral amputation, torch battery replacement depends on the prosthesis, the anatomical conditions and the patient's experience with the prosthesis system. In principle the support of another person can also be obtained for battery replacement.

To insert the batteries, complete the following steps:

- 1) Turn the battery compartment cover counter-clockwise (e.g. with a coin) until it is released from the battery compartment.
- 2) Remove the battery compartment cover.
- 3) Insert 2 CR2032 button cells with the + pole facing the cover.
- 4) Install the cover. In doing so, make sure the rubber seal is not damaged.
- 5) Close the cover by turning it clockwise.

## 7.3 Opening the terminal device in an emergency

This safety function allows the prehensile arms to be opened regardless of the control signals.

- Swing the unlock lever away from the dial (arrow) and open the prehensile arms manually.

## 7.4 Changing the 8E33=9\* terminal device

### Separating the terminal device from the socket

- 1) Turn the prosthesis off with the button in the charging receptacle (prosthetic socket) or the switch on the elbow joint.
- 2) Rotate the terminal device once around its own axis until a slight resistance is felt (about 360°).
- 3) Overcome this resistance and pull the terminal device off the socket.

### Attaching the terminal device to the socket

- 1) Insert the quick-disconnect wrist unit into the lamination ring and press it in firmly.
- 2) Turn the terminal device slightly to the left or right.
- 3) Verify proper attachment by pulling on the terminal device.

## 8 Cleaning

### NOTICE

#### Improper product care

Damage to the product due to the use of incorrect cleaning agents.

- Only clean the product with a damp cloth and mild soap (e.g. 453H10=1-N Ottobock DermaClean).

- 1) Clean the product with a damp cloth and mild soap (e.g. 453H10=1-N Ottobock DermaClean) when needed.  
Make sure that no liquids get into the product and product components.

2) Dry the product with a lint-free cloth and allow it to air dry fully.

## 9 Maintenance

Performing regular maintenance (service inspections) every 24 months is recommended to prevent injuries and maintain the quality of the product.

In general, all products are subject to compliance with the maintenance intervals during the warranty period. This is the only way to maintain full warranty cover.

Additional services such as repairs may be provided in the course of maintenance. These additional services may be provided free of charge or can be billable according to an advance cost estimate, depending on the extent and validity of the warranty.

The following components must be taken to the O&P professional on the agreed maintenance dates:

The prosthesis, battery charger and power supply unit.

## 10 Legal information

### 10.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

### 10.2 Trademarks

All product names mentioned in this document are subject without restriction to the respective applicable trademark laws and are the property of the respective owners.

All brands, trade names or company names may be registered trademarks and are the property of the respective owners.

Should trademarks used in this document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

### 10.3 CE conformity

Otto Bock Healthcare Products GmbH hereby declares that the product is in compliance with European requirements for medical devices.

The product meets the requirements of the RoHS Directive 2011/65/EU on the restriction of the use of certain hazardous substances in electrical and electronic devices.

The full text of the regulations and requirements is available at the following Internet address: <http://www.ottobock.com/conformity>

## 11 Technical data

Ambient conditions	
Storage (with and without packaging)	+5 °C/+41 °F to +40 °C/+104 °F Max. 85% relative humidity, non-condensing
Transport (with and without packaging)	-20 °C/-4 °F to +60 °C/+140 °F Max. 90% relative humidity, non-condensing
Operation	-5 °C/+23 °F to +45 °C/+113 °F Max. 95% relative humidity, non-condensing
Charging the battery	+5 °C/+41 °F to +40 °C/+104 °F Max. 85% relative humidity, non-condensing

Reference number	8E33=9-1	8E34=9-1
Operating voltage [V]	7.4	

Reference number	8E33=9-1	8E34=9-1
Opening width [mm]	95	
Gripping force [N]	0–160	
Weight [g]	540	520
Expected lifetime given compliance with the recommended maintenance intervals [years]	5	
Fingertip tightening torque [Ncm]	40–50	
Proportional speed [mm/sec]	8–200	
Flexion angle of the wrist joint with the flexion brake released [°]	±45	

Torch	
Illuminants [number, type]	2, LED
Total maximum luminosity of the two LEDs [mcd]	16400–24000
Battery [number, type]	2x CR2032
Maximum operating time with new batteries at highest brightness and at room temperature [h]	8
Brightness levels [number]	3
Automatic switch-off	After 15 minutes

Rechargeable battery reference number	757B20	757B21	757B35 =1	757B35 =3	757B35 =5
Temperature range for charging [°C]	+5 to +40				
Capacity [mAh]	900 <sup>1</sup> / 950 <sup>1</sup>	680 <sup>1</sup> / 800 <sup>1</sup>	600	1150	3450
Rechargeable battery lifetime [years]	2				
Product behaviour during the charging process	The product is non-functional				
Operating time of the product with fully charged battery [grip cycles]	Approx. 2500– 3000	Approx. 2000– 2500 (depend- ing on the bat- tery capacity)	Approx. 1000– 2000	Approx. 4000	Approx. 10000
Charging times (rechargeable battery fully dis- charged) [hours]	Approx. 3.5	Approx. 3	Approx. 2.5		Approx. 3
Nominal voltage [V]	Approx. 7.4				
Battery technology	Li-Ion		Li-Po	Li-Ion	
Approved battery chargers	757L20		757L35		

<sup>1</sup> See labeling on battery

Power supply	
Reference number	757L16*
Storage and transport in original packaging	-40 °C/-40 °F to +70 °C/+158 °F
Storage and transport without packaging	-40 °C/-40 °F to +70 °C/+158 °F 10 % to 93 % relative humidity, non-condens- ing

<b>Power supply</b>	
Operation	0 °C/+32 °F to +40 °C/+104 °F max. 90 % relative humidity, non-condensing
Input voltage	100 V~ to 240 V~
Mains frequency	50 Hz to 60 Hz
Output voltage	12 V ===

## 12 Symbols Used

In some jurisdictions it is not permissible to dispose of these products with unsorted household waste. Disposal that is not in accordance with the regulations of your country may have a detrimental impact on health and the environment. Please observe the instructions of your national authority pertaining to return and collection.

Declaration of conformity according to the applicable European directives

Serial number (YYYY WW NNN)  
 YYYY – year of manufacture  
 WW – week of manufacture  
 NNN – sequential number

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