
## AxonRotation 9S503

<b>DE</b> Gebrauchsanweisung (Fachpersonal) .....	3
<b>EN</b> Instructions for use (qualified personnel) .....	7
<b>FR</b> Instructions d'utilisation (Personnel spécialisé) .....	12
<b>IT</b> Istruzioni per l'uso (Personale tecnico specializzato) .....	17
<b>ES</b> Instrucciones de uso (Personal técnico especializado) .....	22
<b>PT</b> Manual de utilização (Pessoal técnico) .....	27
<b>NL</b> Gebruiksaanwijzing (Vakmensen) .....	32
<b>SV</b> Bruksanvisning (Fackpersonal) .....	37
<b>DA</b> Brugsanvisning (Faguddannet personale) .....	41
<b>NO</b> Bruksanvisning (Fagpersonell) .....	46
<b>FI</b> Käyttöohje (Ammattihenkilöstö) .....	50
<b>CS</b> Návod k použití (Odborný personál) .....	55
<b>RU</b> Руководство по применению (Квалифицированный персонал) .....	60
<b>JA</b> 取扱説明書 (有資格担当者) .....	65

1

## INFORMATION

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### INFORMATION

Date of last update: 2020-12-02

► Please read this document carefully before using the product and observe the safety notices.

- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

The "9S503 AxonRotation" Axon-Bus system component will be referred to simply as the product below.

These instructions for use provide you with important information on the use, adaptation and handling of the product.

## **2 Product description**

### **2.1 Function**

In combination with an Axon-Bus gripping component, e.g. the 8E500 Michelangelo Hand, the product facilitates fitting with components in the modular Ottobock Axon-Bus prosthesis system. Pronation/supination of the Axon-Bus gripping component are possible from the neutral position and are carried out actively.

These functions support the patient during his or her daily activities and provide the highest rehabilitation value.

The product promotes a physiologically correct body posture and reduces the need for unnatural, compensating body movements.

The product can only be used in combination with the Axon-Bus system component "13E500 AxonMaster" version 1.4 or higher.

Power is supplied to the product by the battery which comes integrated in the socket.

Control settings for the product can only be carried out when all Axon-Bus components are connected. For further information, see the instructions for use provided with the Axon-Bus components.

To prevent mechanical damage, the coupling (connection between the motor and gripping component) slips when rotary movement is obstructed.

### **2.2 Combination possibilities**

This product is to be combined exclusively with components in the Axon-Bus prosthetic system. Components of the Ottobock MyoBock system or components from other manufacturers cannot be used with this product.

## **3 Intended use**

### **3.1 Indications for use**

The product is to be used exclusively for exoprosthetic fittings of the upper limbs and enables active pronation and supination of the Axon-Bus gripping component.

### **3.2 Conditions of use**

In combination with myoelectrically controlled Axon-Bus gripping components in the modular Ottobock Axon-Bus prosthesis system, the product can be used from the transradial or transhumeral amputation level or, in case of dysmelia, for forearm or upper arm fittings.

The product was developed for everyday use and must not be used for unusual activities such as extreme sports (free climbing, paragliding, etc.).

Permissible ambient conditions are described in the technical data (see page 11).

The product is intended **exclusively** for use on **one** patient. Use of the product by another person is not approved by the manufacturer.

### 3.3 Contraindications

- All conditions which contradict or go beyond the specifications listed in the section on "Safety" and "Indications for use".

### 3.4 Qualification

The product may be fitted only by qualified personnel authorised by Ottobock after completing the corresponding training.

## 4 Safety

### 4.1 Explanation of warning symbols

Warning regarding possible risks of accident or injury.

Warning regarding possible technical damage.

### 4.2 Structure of the safety instructions

#### The heading describes the source and/or the type of hazard

The introduction describes the consequences in case of failure to observe the safety instructions. Consequences are presented as follows if more than one consequence is possible:

- > E.g.: Consequence 1 in case of failure to observe the hazard
- > E.g.: Consequence 2 in case of failure to observe the hazard
- ▶ This symbol identifies activities/actions that must be observed/carried out in order to avert the hazard.

### 4.3 General safety instructions

#### Non-observance of safety instructions

Personal injury/damage to the product due to using the product in certain situations.

- ▶ Observe the safety instructions and the stated precautions in this accompanying document.

#### Independent manipulation of the product

Injury due to malfunction and resulting unexpected prosthesis actions.

- ▶ Manipulations to the product other than the tasks described in these instructions for use are not permitted.
- ▶ The product and any damaged components may only be opened and repaired by authorised, qualified Ottobock personnel.

#### Insufficient distance to HF communication devices (e.g. mobile phones, Bluetooth devices, WiFi devices)

Falling due to unexpected behaviour of the product caused by interference with internal data communication.

- ▶ Therefore, keeping a minimum distance of 30 cm to HF communication devices is recommended.

## CAUTION

### **Proximity to sources of strong magnetic or electrical interference (e.g. theft prevention systems, metal detectors)**

Injury due to unexpected behaviour of the product caused by interference with internal data communication.

- ▶ Avoid remaining in the vicinity of visible or concealed theft prevention systems at the entrance/exit of stores, metal detectors/body scanners for persons (e.g. in airports) or other sources of strong magnetic and electrical interference (e.g. high-voltage lines, transmitters, transformer stations, computer tomographs, magnetic resonance tomographs, etc.).
- ▶ When walking through theft prevention systems, body scanners or metal detectors, watch for unexpected behaviour of the product.

## NOTICE

### **Improper product care**

Damage to the product due to the use of incorrect cleaning agents.

- ▶ Clean the product and the contacts with a damp cloth and mild soap only.

## **5 Scope of delivery**

- 9S503 AxonRotation (active rotation unit)
- 4 pcs. M3 x 8 Torx screws
- 1 pc. Torx offset screwdriver
- 1 pc. Instructions for use (qualified personnel)

## **6 Preparing the product for use**

- 1) Use silicone grease to seal the plug connection.
- 2) Connect the product with the connecting cable.
- 3) Position the product so that the arrow mark points in the direction of the ulna (outside of the arm). (see fig. 1)
- 4) Screw the product to the lamination ring. (see fig. 1)

**INFORMATION: Tighten the screws hand-tight and do not use a thread lock such as Loctite.**

## **7 Maintenance**

Performing regular maintenance (service inspections) every 24 months is recommended to prevent injuries and maintain the quality of the product.

In general, all products are subject to compliance with the maintenance intervals during the warranty period. This is the only way to maintain full warranty cover.

Additional services such as repairs may be provided in the course of maintenance. These additional services may be provided free of charge or can be billable according to an advance cost estimate, depending on the extent and validity of the warranty.

The following components must always be sent in for maintenance and repairs:

The product, battery charger and power supply. The shipping container for the loaner unit you receive must be reused for sending back the components requiring inspection.

## **8 Legal information**

### **8.1 Liability**

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for

damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

## 8.2 Trademarks

All product names mentioned in this document are subject without restriction to the respective applicable trademark laws and are the property of the respective owners.

All brands, trade names or company names may be registered trademarks and are the property of the respective owners.

Should trademarks used in this document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

## 8.3 CE conformity

Otto Bock Healthcare Products GmbH hereby declares that the product is in compliance with applicable European ergänzen!European requirements for medical devices.

The product meets the requirements of the RoHS Directive 2011/65/EU on the restriction of the use of certain hazardous substances in electrical and electronic devices.

The full text of the regulations and requirements is available at the following Internet address: <http://www.ottobock.com/conformity>

## 9 Technical data

<b>Ambient conditions</b>	
Storage (with and without packaging)	+5 °C/+41 °F to +40 °C/+104 °F Max. 85% relative humidity, non-condensing
Transport (with and without packaging)	-20 °C/-4 °F to +60 °C/+140 °F Max. 90% relative humidity, non-condensing
Operation	-5 °C/+23 °F to +45 °C/+113 °F Max. 95% relative humidity, non-condensing

<b>General</b>	
Reference number	9S503
Product service life	5 years
Active pronation from the neutral position	160°
Active supination from the neutral position	160°
Maximum torque	1.5 Nm
Idle speed	25 RPM
Rotation angle	150°/sec.
Weight	approx. 140 g/4.94 oz

<b>Load limits</b>	
Maximum static vertical load	700 N
Maximum dynamic vertical load	200 N

## 10 Appendix

### 10.1 Symbols Used

Declaration of conformity according to the applicable European directives

In some jurisdictions it is not permissible to dispose of these products with unsorted household waste. Disposal that is not in accordance with the regulations of your country may have a detrimental impact on health and the environment. Please observe the instructions of your national authority pertaining to return and collection.

MD

Medical device

SN

Serial number (YYYY WW NNN)

YYYY – year of manufacture

WW – week of manufacture

NNN – sequential number

Manufacturer

IP33

Protection against penetration of solid foreign objects with a diameter > 2.5 mm, protection against water spray

## 10.2 Axon

The term "Axon" stands for **A**daptive **eX**change **O**f **N**europlacement data. The Axon-Bus is an Ottobock innovation for the field of exoprosthetics: a data transmission system, derived from safety-related bus systems in the aviation and automobile industries. For the user this means enhanced safety and reliability because of a considerably reduced sensitivity to electromagnetic interference in comparison with conventional systems.

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