


3R90, 3R92

DE Gebrauchsanweisung	4
EN Instructions for use	10
FR Instructions d'utilisation	16
IT Istruzioni per l'uso	22
ES Instrucciones de uso	28
PT Manual de utilização	34
SV Bruksanvisning	41
NL Gebruiksaanwijzing	46
DA Brugsanvisning	53
NO Bruksanvisning	58
FI Käyttöohje	64
PL Instrukcja użytkowania	70
HU Használati utasítás	77
CS Návod k použití	83
RO Instrucțiuni de utilizare	88
HR Upute za uporabu	95
EL Οδηγίες χρήσης	101
RU Руководство по применению	108
TR Kullanma talimatı	114
ZH 使用说明书	120
JA 取扱説明書	125

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Last update: 2021-01-18


- Please read this document carefully before using the product and observe the safety notices.
- Instruct the user in the safe use of the product.
- Please contact the manufacturer if you have questions about the product or in case of problems.
- Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- Please keep this document in a safe place.

1 Description

1.1 Indications for use


The 3R90 and 3R92 Modular Knee Joints with Friction Brake are to be used **exclusively** for the prosthetic fitting of the lower extremity.

1.2 Field of application according to the Ottobock Mobility System (MOBIS):

3R90: Recommended for amputees with **mobility grade 1 and 2** (indoor walkers and restricted outdoor walkers).

Approved for a patient weight of up to 125 kg / 275 lbs.

3R92: Recommended for amputees with **mobility grade 2 and 3** (restricted outdoor walkers and non-restricted outdoor walkers).

Approved for a patient weight of up to 125 kg / 275 lbs.

The 3R90 and 3R92 Modular Knee Joints with Friction Brake are not suitable for bilateral trans-femoral amputees.

1.3 Safety instructions

Attention!

These knee joints with friction brake are contraindicated for patients who are chronically dependent on mobility devices (e.g. anterior walkers, crutches) as well as for insecure patients whose gait cycle does not permit them to apply the brake mechanism as required, i.e. to lock the knee joint at heel contact and unlock it when the forefoot is bearing weight.

The Ottobock MOBIS is not intended to be used as a guide to obtaining reimbursement for prosthetic components in the USA.

Please avoid exposing prosthetic components to surroundings that corrode metal parts, for example, freshwater, saltwater, acids and other liquids. Using this medical product in such environmental conditions will render all claims against Otto Bock HealthCare null and void.

Please inform your patients.

Reuse on another patient

Fall due to loss of functionality as well as damage to the product

- Only use the product for a single patient.
- Inform the patient.

1.4 Construction and function

These single-axis knee joints are equipped with a novel, weight-dependent braking mechanism that offers the patient a high degree of security. At heel contact, the brake responds, stabilizes the prosthesis and ensures a high degree of stance phase security. During loading of the forefoot, the brake deactivates to allow the patient to initiate the swing phase easily and harmoniously without raising his or her hip. When used correctly, this makes it possible for patients to walk much more naturally and exert less energy in comparison to conventional knee joints with friction brake.

Instructions:

Use parallel bars for first trial walk with the knee joints. Instruct the patient that lifting the hip is not necessary for disengaging the brake and starting the swing phase. Unlike conventional knee joints with friction brake such as the 3R49 Modular Knee Joint, the brake disengages as early as when the forefoot is loaded. Thus the securing effect of the brake is missing at the point of forefoot contact e.g. when the prosthesis is in the heel off phase. **The brake also disengages/switches off when the forefoot is loaded unintentionally, which may under circumstances result in a fall, if this is not conveyed to the patient.**

The brake also offers the desired security when standing in a slightly knee flexed position. Due to the construction, the braking system is durable and self-adapting, that is, adjusting the effectiveness of the brake is generally not necessary.

For joint 3R90, a mechanical extension assist with a progressively acting spring combination serves to control the swing phase. The joint 3R92 is equipped with a progressively acting dual-chamber pneumatic system.

A flexible plastic cover encases the mechanical parts of the lower joint section.

1.5 Lifetime

This product was tested by the manufacturer with 3 million load cycles. Depending on the user's activity level, this corresponds to a maximum lifetime of 5 years.

1.6 Combination possibilities

CAUTION

Failure to observe the manufacturer requirements regarding combination possibilities

Injuries, malfunctions or product damage due to **unallowable combination** of prosthesis components

- Based on the instructions for use of all prosthetic components used, verify that they may be combined with each other and are approved for the patient's field of application.
- Contact the manufacturer with questions.

INFORMATION

When fabricating a prosthesis for the lower limbs, all of the prosthetic components used have to meet the requirements in regard to the patient's body weight and activity level.

Unallowable combinations	
Prosthetic hip joints	7E4, 7E5, 7E7, 7E9, 7E10

1.7 Prosthetic feet

Connection to the prosthetic foot is provided by employing the included tube adapter, which is to be shortened according to the leg length.

	<p>Attention!</p> <p>Do not clamp the tube into a vise! Cut the tube using 719R2 Tube Cutter placed at a right angle to the tube. Afterwards carefully debur the tube inside and outside using the 718R1 Tube Inside and Outside Reamer. Insert the tube until it is fully seated inside the integrated tube clamp. Tighten the cap screw using the 710D1 Torque Wrench.</p> <p>Torque: 10 Nm (7.4 lbf. ft/88.5 lbf. in) (Fig. 3)</p>
--	---

	<p>Attention!</p> <p>For trial fitting: Tighten the 506G3 Set Screws of the tube adapter which are side-by-side to a torque of 10 Nm using the 710D4 Torque Wrench. To finish the prosthesis: Apply 636K13 Loctite to the set screws, pre-tighten the side-by-side set screws to a torque of 10 Nm, then tighten all four set screws to a torque of 15 Nm (11.1 lbf. ft/133 lbf. in).</p>
--	--

2 Technical data

Article number	3R90	3R92
Proximal connection	Pyramid adapter	
Distal connection	Tube clamp	
Knee flexion angle	135°	
Weight	700 g (24.7 oz)	840 g (28.2 oz)
System height	745 g	895 g
Proximal system height to alignment reference point	8 mm (0.31 in)	
Distal system height to alignment reference point	89 mm (3.5 in)	146 mm (5.75 in)
Max. user weight	125 kg (275 lbs)	
Mobility Grade	1, 2	2, 3

3 Alignment instructions

The advantages of the 3R90/3R92 can only be made optimal use of in case of a correct alignment. The optimal residual limb position must be anticipated when positioning the socket connector.

Plumb lines in the frontal and sagittal planes drawn from the proximal center of the socket and marked during plaster cast taking and trial fitting of the test socket will facilitate correct positioning of the lamination anchor/socket adapter.

To align the prosthesis please proceed in two steps:

First make the bench alignment using an alignment tool such as 743L200 L.A.S.A.R. Assembly. For **subsequent** optimization of the alignment, please use 743L100 L.A.S.A.R. Posture.

3.1 Bench alignment with alignment tool (Fig. 7)

- 1 **Position** the middle of the foot approx. **30 mm anterior** to the alignment reference line. **This applies to all feet recommended for use with 3R90 / 3R92, independent of the values in the instruction manuals of these feet!**
- 2 Add 5 mm to the required heel height. Set correct outward rotation of the foot.
- 3 Clamp knee joint using Assembly parts. Place the knee center **10 mm** posterior to the alignment reference line. Pay attention to the knee ground distance and outward rotation of the knee (approx. 5°). The ring disc – divided into 5 mm sections – makes it easier to read (Fig. 1).
- 4 Connect knee and foot. Either use the adjustable 2R45=34 Tube Adapter or shorten tube adapter to the correct length.
- 5 Mark the center of the socket proximal and distal on the lateral side. Connect both marks with a line.
- 6 Now position the socket such that the upper mark passes through the alignment reference line. Socket flexion normally will be somewhere between 3° and 5°; however, the individual situation (e.g. hip joint contractures) must be taken into account and, if necessary, more flexion should be provided. Pay attention to the ischial tuberosity to ground distance!

Attention!

If the residual limb flexion is not taken into account, the joint will be positioned too far anterior. This will lead to malfunction and premature wear.

- 7 Connect socket and knee using a corresponding adapter (e.g. 4R111, 4R41, 4R55, 4R51 Socket Adapter). Use the 4R112 Sliding Adapter for alignment correction (Fig. 10). If adapters with 4-hole connection are used like 5R1=*, 5R6=*, etc., you can use the 4R118 Sliding Adapter for a permanent extra posterior placement of 10-25 mm (Fig. 9).

3.1.1 Static alignment optimization **ONLY** using 743L100 L.A.S.A.R. Posture (Fig. 8)

The static alignment can be substantially improved using L.A.S.A.R. Posture. In order to ensure appropriate stability combined with easy swing phase initiation, please proceed as follows:

- 1 To make the load line visible, the amputee stands on the L.A.S.A.R. Posture with the prosthetic side on the force plate and with the contralateral side leg on the height compensation panel. The prosthetic side should at least be loaded with 35% of the body weight, please check on the display.
- 2 Now adapt the alignment by only adjusting the plantar flexion of the foot. The load line / laser line should be **40 mm** anterior to the knee joint center (See fig. 8).
- 3 After step 2, dynamic optimization can take place between parallel bars.

3.2 Adjustment and final assembly

3.2.1 Fine adjustment of the brake

The sensitivity of the brake has a factory setting, so that in general no adjustments to the patient's requirements are necessary. If need be, however, the sensitivity of the brake can be adjusted by slightly turning the setting screw (Fig. 4).

Turn to the right = brake becomes less sensitive

Turn to the left = brake becomes more sensitive

Attention!

The sensitivity of the brake may only be reduced to a point where a sufficient braking effect is still ensured.

3.2.2 Setting the swing phase control

The knee joints are delivered as follows:

- 3R90 in the highest spring catch – the least effect
- 3R92 with the least pneumatic resistance

3.2.3 Mechanical extension assist of the 3R90

The effect of the mechanical extension assist can be regulated at 5 levels by the tension of the spring extension assist. If necessary, the spring effect, as delivered, must be increased.

To increase the spring effect: (Fig. 5)

Slightly press the setting unit and turn to the right.

To reduce the spring effect:

Slightly press the setting unit and turn to the left.

3.2.4 Pneumatic of the 3R92

The pneumatic swing phase control makes the gait of the patient more harmonious. At the same time, movement resistances prevent the prosthetic lower leg from swinging too far when flexing and ensure dampened extension. Along with the patient's habits, the foot weight and the length of the lower leg, which have a pendular weight, also have an influence on his or her gait.

Before changing the factory settings, be sure to pay attention to the following instructions:

During trial walking, first adjust the flexion resistance:

Increase the flexion resistance, as required, by turning the left valve screw (F) to the right (Fig. 6).

Afterwards, set the extension resistance by turning the right valve screw (E) to the right.

Carefully fine-tune flexion and extension.

Status	Measure	Procedure
Excessive heel rise	Increase resistance	Turn the left screw (F) to the right.
Insufficient heel rise	Decrease resistance	Turn the left screw (F) to the left.
Extension impact too hard	Increase resistance	Turn the right screw (E) to the right.
Insufficient extension movement	Decrease resistance	Turn the right screw (E) to the left.

Attention!

The resistance for extension may only be adjusted to the point that the complete extension of the leg is still always possible.

3.3 Cosmetic foam cover

Use the 3S107 Foam Cover for the 3R90/3R92 Joints.

Attention!

To eliminate noise in the cosmetic foam cover, do not use talcum powder! Talcum powder reduces the lubrication of the mechanical parts, which may lead to a malfunction and thus increase the risk of failure. Using this medical product after application of talcum powder will render all claims against Otto Bock HealthCare null and void.

Note

An alternative for eliminating noises is a silicone spray (519L5) that can be sprayed directly on the friction surface of the cosmetic foam cover.

3.4 Maintenance instructions

Ottobock recommends the readjustment of the knee joint's settings after the patient has spent a period of time getting used to the prosthesis. This period of time varies depending on individual patient characteristics.

Please check the wear and functionality of the knee joint at least once a year and make readjustments if need be. Pay special attention to the resistance to movement and to the development of unusual noise. Complete flexion and extension must be ensured.

Do not disassemble the joint. If there are potential malfunctions, please send the entire joint to Ottobock Service.

Attention – Please inform your patient!

Due to different environmental conditions and conditions of application, the function of the knee joint can be impaired. To avoid the risk of accident and possible injury to the patient, the knee joint must no longer be used after a noticeable change of function has occurred. This noticeable change of function can be for example stiffness, non-attainment of the extension stop, decreasing swing phase control or stance phase stability, abnormal noises, etc.

Measure

Have the knee joint examined by your prosthetist.

4 Disposal

In some jurisdictions it is not permissible to dispose of the product with unsorted household waste. Improper disposal can be harmful to health and the environment. Observe the information provided by the responsible authorities in your country regarding return, collection and disposal procedures.

5 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

5.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

5.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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