Myo Plus

DE Gebrauchsanweisung (Benutzer)	3
EN Instructions for use (user)	19
FR Instructions d'utilisation (Utilisateur)	36
IT Istruzioni per l'uso (Utilizzatore)	54
ES Instrucciones de uso (Usuario)	69
PT Manual de utilização (Usuário)	85
NL Gebruiksaanwijzing (Gebruiker)	101
SV Bruksanvisning (Användare)	117
DA Brugsanvisning (Bruger)	133
NO Bruksanvisning (Bruker)	148
FI Käyttöohje (Käyttäjä)	164
TR Kullanma talimatı (Kullanıcı)	179
RU Руководство по применению (Пользователь)	195
JA 取扱説明書(ユーザー用)	212

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INFORMATION

Date of last update: 2021-01-20

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Obtain instruction from the qualified personnel in the safe use of the product.
- ▶ Please contact the qualified personnel if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

The product "Myo Plus" is referred to as the product below.


These instructions for use provide you with important information on the use, adaptation and handling of the product.

Only put the product into use in accordance with the information contained in the accompanying documents supplied.

2 Product description

2.1 Design

The product consists of the following components:

1. Myo Plus TR
2. Remote electrode with electrode dome

2.2 Function

The product is used to control a myoelectric prosthesis.

The product measures the control signals and assigns them to the prosthesis movements.

By calibration via the Myo Plus app, the control unit learns to assign the sampled muscle signals to the various movement types. This calibration can be carried out directly by the user and repeated at regular intervals.

2.3 Combination possibilities

This product can be combined with the following Ottobock components:

Terminal devices

- System Electric Hand DMC plus: 8E38=6, 8E39=6
- Sensor Hand Speed: 8E38=8, 8E39=8
- MyoHand VariPlus Speed: 8E38=9, 8E39=9, 8E41=9
- Transcarpal Hand DMC plus: 8E44=6
- System Electric Greifer DMC VariPlus: 8E33=9, 8E34=9
- DMC VariPlus System Electric Greifer: 8E33=9-1, 8E34=9-1
- bebionic hand EQD: 8E70=*
- bebionic hand Short Wrist: 8E71=*
- bebionic hand Flex: 8E72=*

3 Intended use

3.1 Indications for use

Myo Plus is to be used exclusively for exoprosthetic fittings of the upper limbs.

3.2 Conditions of use

The product was developed for everyday use and must not be used for unusual activities such as extreme sports (free climbing, paragliding, etc.). Furthermore, the product must not be used for the operation of motor vehicles, heavy equipment (e.g. construction machines), industrial machines, firearms or motor-driven equipment.

The product is intended **exclusively** for use on **one** user. Use of the product by another person is not approved by the manufacturer.

3.3 Indications

- For users with unilateral or bilateral amputation.
- For users with hand and forearm amputation.

3.4 Contraindications

- For users with upper arm amputation.
- For users with partial hand amputation.
- For users with shoulder amputation.
- For users who are not able to generate an adequately separable pattern.

3.5 Qualification

Fitting a patient with the product may only be carried out by O&P professionals, therapists or nursing staff who have been authorised by Ottobock upon completion of a corresponding training course. Furthermore, the O&P professional must have the technical qualifications required for the alignment of a prosthesis together with all necessary settings and adjustments.

The user must be instructed about the use of the product by authorised, qualified personnel.

4 Safety

4.1 Explanation of warning symbols

 WARNING	Warning regarding possible serious risks of accident or injury.
 CAUTION	Warning regarding possible risks of accident or injury.
 NOTICE	Warning regarding possible technical damage.

4.2 Structure of the safety instructions

 WARNING
The heading describes the source and/or the type of hazard
The introduction describes the consequences in case of failure to observe the safety instructions. Consequences are presented as follows if more than one consequence is possible:
> E.g.: Consequence 1 in the event of failure to observe the hazard
> E.g.: Consequence 2 in the event of failure to observe the hazard
▶ This symbol identifies activities/actions that must be observed/carried out in order to avert the hazard.

4.3 General safety instructions

 WARNING
Non-observance of safety notices
Personal injury/damage to the product due to using the product in certain situations.
▶ Observe the safety notices and the stated precautions in this accompanying document.

 WARNING
Use of the product while operating a vehicle
Accident due to unexpected behaviour of the product.
▶ The product must not be used for the operation of motor vehicles and heavy equipment (e.g. construction machines).

 WARNING
Using the prosthesis while operating machines
Injury due to unexpected actions of the prosthesis.
▶ The prosthesis should not be used to operate industrial machines or motor-driven equipment.

 WARNING
Use of the product while handling a firearm.
Injury or fatal injuries due to unexpected product behaviour.
▶ Do not handle firearms while using the product.

⚠ WARNING

Operating the product near medical devices critical for safety

Interference with life-sustaining medical devices (e.g., pacemaker, defibrillator, heart-lung machine, etc.) due to electromagnetic interference of the product.

- ▶ When operating the product in the immediate vicinity of life-sustaining medical devices, ensure that the minimum distances stipulated by the manufacturer are observed.
- ▶ Make sure to observe the operating conditions stipulated by the manufacturer and the safety notices.

⚠ WARNING

Use of inappropriate components

Injury due to malfunction of the product.

- ▶ Operate the product only with accessories recommended by the manufacturer. You can find the list of components in the section "Scope of delivery and accessories".

⚠ WARNING

Changes or modifications to the product made independently

Injury due to faulty operation or malfunction of the product.

- ▶ Have any changes or modifications to the product carried out only by authorised, qualified Ottobock personnel.
- ▶ The battery may only be handled by authorised, qualified Ottobock personnel (replacement by the user is not permitted).
- ▶ The product and any damaged components may only be opened and repaired by authorised, qualified Ottobock personnel.

⚠ CAUTION

Insufficient distance to HF communication devices (e.g. mobile phones, Bluetooth devices, WiFi devices)

Injury due to unexpected behaviour of the product caused by interference with internal data communication.

- ▶ Therefore, keeping a minimum distance of 30 cm from HF communication devices is recommended.

⚠ CAUTION

Operating the product in very close proximity to other electronic devices

Injury due to unexpected behaviour of the product caused by interference with internal data communication.

- ▶ Do not operate the product in the immediate vicinity of other electronic devices.
- ▶ Do not stack the product with other electronic devices during operation.
- ▶ If simultaneous operation cannot be avoided, monitor the product and verify proper use in the existing setup.

⚠ CAUTION

Penetration of dirt and humidity into the product

Injury due to faulty control or malfunction of the product.

- ▶ Ensure that neither solid particles nor liquids can penetrate into the product.

⚠ CAUTION

Mechanical loads

Injury due to faulty control or malfunction of the product.

- ▶ Watch for loosened electrode domes.
- ▶ Do not subject the product to mechanical vibrations or impacts.
- ▶ Check the product and its housing for visible damage before each use (e.g., cracking or breakage). Should the product be damaged or not function properly, please contact your O&P professional immediately.

⚠ CAUTION

Changing terminal devices when switched on

Injury due to faulty control or malfunction of the product.

- ▶ Power down the product by pushing the button in the charging receptacle before changing terminal devices.

⚠ CAUTION

Proximity to sources of strong magnetic or electrical interference (e.g. theft prevention systems, metal detectors)

Injury due to unexpected behaviour of the product caused by interference with internal data communication.

- ▶ Avoid remaining in the vicinity of visible or concealed theft prevention systems at the entrance/exit of stores, metal detectors/body scanners for persons (e.g. in airports) or other sources of strong magnetic and electrical interference (e.g. high-voltage lines, transmitters, transformer stations, computer tomographs, magnetic resonance tomographs, etc.).
- ▶ When walking through theft prevention systems, body scanners or metal detectors, watch for unexpected behaviour of the product.

⚠ CAUTION

Insufficient skin contact of the electrodes

Injury due to faulty control or malfunction of the product.

- ▶ The electrodes are to be placed on intact skin only and with as much electrode-skin contact as possible.
- ▶ Make sure that skin contact of the electrodes continues to be given when carrying heavy loads.
- ▶ Make sure that skin contact of the electrodes continues to be given when the residual limb volume changes.

⚠ CAUTION

Use of damaged electrodes

Injury due to unexpected product behaviour as a result of a malfunction.

- ▶ Switch off the system and contact your O&P professional.

⚠ CAUTION

Remaining in areas outside the allowable temperature range

Injury due to faulty control or malfunction of the product.

- ▶ Avoid remaining in areas outside the allowable temperature range (see page 31).

⚠ CAUTION

Incorrect electrode settings due to muscle fatigue

Injury due to faulty control or malfunction of the product.

- ▶ Make sure you take breaks during calibration.

⚠ CAUTION

Skin irritation due to inadequate cleaning of the product

Skin irritation due to contact with soiled electrode domes.

- ▶ Clean the product only as described in the section "Cleaning and care" (Cleaning and care).

NOTICE

Failure to observe the system requirements for the installation of the Myo Plus app

Malfunction of the device.

- ▶ Install the Myo Plus app only on the operating systems listed in the section "System requirements" (see page 27).

NOTICE

Mechanical damage to the product

Change in or loss of functionality due to damage.

- ▶ Use caution when working with the product.
- ▶ If the product is damaged, check it for proper function and readiness for use.
- ▶ Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

5 Scope of Delivery and Accessories

Items in the scope of delivery or accessories marked with ■ are application parts according to the IEC 60601-1:2005/A1:2012 standard.

- 1 pc. 13E520 Myo Plus TR
- 3 pc - 8 pc 13E400=*/13E401=* Remote electrodes

The remote electrodes can be combined with any of the following electrode domes:

- 13Z161 electrode dome, flat ■ (6 pc. per package)
- or
- 13Z162 electrode dome, medium ■ (6 pc. per package)
- or
- 13Z163 electrode dome, high ■ (6 pc. per package)
 - 1 pc. 646C107 Bluetooth PIN card
 - Myo Plus app

Adjustment app available for download from the Apple App Store or Google Play Store. Enter the following search terms: Ottobock, Myo Plus.

- 1 pc. Instructions for use (user)

6 Charging the battery

The following points must be observed when charging the battery:

- The capacity of a fully charged battery is sufficient for one full day.
- We recommend charging the product every day when used on a daily basis.
- The battery should be charged for at least 3 hours prior to initial use.

6.1 757L20 battery charger and 757B2* EnergyPack

NOTICE

Deep discharge of the battery

Irreparable damage to the product and resulting loss of functionality in case of irregular and incomplete charging of the battery.

- ▶ In case of long-term storage, fully charge the battery every 4 to 6 months.

- 1) Plug the charging station with power supply unit into the wall socket.
→ The green LED on the power supply unit and the green LED in the middle of the LED row light up.
- 2) Insert one or two rechargeable batteries into the charging station.
- 3) The green LED for the respective charging socket lights up ● and the charging process starts.
- 4) When the rechargeable battery is fully charged, the green LED for the respective charging socket flashes ✨.
- 5) Remove the rechargeable battery after the charging process is complete.

6.1.1 Display of the current charge level

When inserting/engaging the rechargeable battery in the prosthesis, a charge level indicator is activated on the battery for a few seconds.

LED indicator	Event
●	Charge level above 50% (illuminated in green)
● and ●	Charge level below 50% (flashing alternating green and orange)
●	Charge level below 5% (illuminated in orange)

6.1.2 Troubleshooting

Charging socket	LED for the charging socket	Event
Empty	● or ● or ●	The battery charger is defective The battery charger and power supply must be inspected by an authorised Ottobock Service Centre.
Empty or rechargeable battery inserted	☀	
Rechargeable battery inserted	●	The rechargeable battery is charging
	☀	The rechargeable battery is fully charged
	●	The temperature of the rechargeable battery is too high Remove the rechargeable battery and let it cool down
	●	The inserted rechargeable battery is defective The rechargeable battery must be inspected by an authorised Ottobock Service Centre.

6.2 757L35 battery charger and 757B35=* rechargeable battery

6.2.1 Connect battery charger to the product

- 1) Connect the charging plug to the charging receptacle on the prosthetic socket.
→ A correct connection between the battery charger and the prosthesis is indicated by feedback (see page 33).
- 2) The charging process starts.
→ The prosthesis is switched off automatically.
- 3) Disconnect the prosthesis when the charging process is complete.

INFORMATION

The product cannot be used as long as the charging plug is connected to the charging receptacle. The product is deactivated for the duration of charging.

6.2.2 Safety Shutoff

The purpose of the prosthesis safety shutoff is to protect the battery; it is triggered in case of:

- Excessively high or low temperature
- Overvoltage or undervoltage
- Short circuit

After a short circuit, the charging plug must be plugged into the charging receptacle and then unplugged in order for the electronics to be activated.

6.2.3 Display of the current charge level

The battery charge level can be queried at any time.

- 1) With the prosthesis switched on, press the charging receptacle button and hold for less than one second.
- 2) The LED display on the charging receptacle provides information on the current battery charge level.

Charging receptacle	Event
	Battery fully charged (illuminated in green)
	Battery 50% charged (illuminated in yellow)
	Battery drained (illuminated in orange)

6.2.4 Beep signals

Beep signal	Additional display	Event
1 x long	-	<ul style="list-style-type: none">• Turn off prosthesis on charging receptacle• Charging begins (charging plug connected to charging receptacle)• Charging ends (charging plug disconnected from charging receptacle)
2 x short	LED on charging receptacle lights up briefly	Switch on prosthesis
3 x short	-	Battery voltage too low, prosthesis shuts off automatically

7 Use

INFORMATION

Prior to use, make sure the terminal device is connected to the prosthetic socket before switching on. If the terminal device is connected to the prosthetic socket after switching on, rotation and the closing and opening function are active, but the various grip patterns cannot be recognised.

Perform the following steps to put the product into operation:

- 1) Put on the product.
 - 2) Switch on the prosthesis.
- Recalibrate the prosthesis if existing movements are to be adapted or new movements added.

Calibration procedure

The following steps have to be completed to start calibration:

- 1) Tap the navigation menu.
- 2) Tap the "Settings" menu option.
- 3) Set the pause time between the movements.

INFORMATION: The pause time should be adapted according to previous experience with the product (inexperienced user – maximum, experienced user or difficult arm position – minimum).

- 4) Bring the arm to the selected arm position.
- 5) Tap "→" to start recording the calibration data.
- 6) Perform all 6 sets.

INFORMATION: Perform all recordings in the illustrated arm position.

8 Myo Plus app

With the Myo Plus app, the user can change the behaviour of the product to a limited extent and is able to use the frequently used functions. The user can adapt the product to conditions on the current day and adjust the movements through calibration. A menu item in the app allows the user to train their muscle signal.

8.1 System Requirements

See the information in the Apple App Store or Google Play Store regarding compatibility with mobile devices and versions.

8.2 Initial connection between Myo Plus app and component

The following points should be observed before the initial connection:

- Bluetooth on the mobile device must be switched on.
- The Bluetooth ID and Bluetooth PIN of the component being connected must be known. They are found on the enclosed Bluetooth PIN card. The Bluetooth ID starts with the letters "BT ID".

8.2.1 Starting the Myo Plus app for the first time

INFORMATION

Keep mobile app up to date.


- Please contact the manufacturer if you suspect cybersecurity problems.

INFORMATION

Bluetooth on the prosthesis must be turned on in order to use the MyoPlus app. After turning on Bluetooth, it remains active for approx. 5 minutes. The app must be started and the connection established during this time.


- 1) Tap the icon for the Myo Plus app .
→ The end user licence agreement (EULA) is displayed.
 - 2) If rights are requested the first time the app is started, these have to be granted. Otherwise, the app cannot be started.
 - 3) Accept the end user licence agreement (EULA). If the end user license agreement (EULA) is not accepted, the Myo Plus app cannot be used.
 - 4) The "Bluetooth connection" menu item is opened.
 - 5) Tap  to establish a connection.
 - 6) Select the desired component from the list.
 - 7) Enter the Bluetooth PIN code and tap "Connect".
→ The  icon is displayed when the connection has been established.
- Once the connection has been established, the data are read from the component. Then the radar chart is displayed.

8.3 Control elements of the Myo Plus app

1.  Access online help
2. Display and name of the currently selected motion sequence
3.  Connection to component has been established
4. Display of the selected screen
5. Display of the currently selected menu item
6.  Access the navigation menu

8.3.1 Navigation menu

Tap the ☰ symbol in the menus to display the navigation menu. Additional settings for the connected component can be configured in this menu.

The navigation menu is divided into 3 main menu items:

1. My Myo Plus
2. O&P professional mode
3. General information

See the online help for further information about the respective submenus.

9 Cleaning and Care

Cleaning the electrode domes

- 1) Clean the electrode domes with a cleaning cloth and 634A58 isopropyl alcohol after each application.
- 2) Dry the electrode domes with a cloth.

10 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

10.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

10.2 Trademarks

All product names mentioned in this document are subject without restriction to the respective applicable trademark laws and are the property of the respective owners.

All brands, trade names or company names may be registered trademarks and are the property of the respective owners.

Should trademarks used in this document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

10.3 CE conformity

Otto Bock Healthcare Products GmbH hereby declares that the product is in compliance with applicable European ergänzen!European requirements for medical devices.

The product meets the requirements of the RoHS Directive 2011/65/EU on the restriction of the use of certain hazardous substances in electrical and electronic devices.

This product meets the requirements of the 2014/53/EU directive.

The full text of the regulations and requirements is available at the following Internet address:
<http://www.ottobock.com/conformity>

10.4 Local Legal Information

Legal information that applies **exclusively** to specific countries is written in the official language of the respective country of use in this chapter.

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:

- 1) This device may not cause harmful interference, and
- 2) This device must accept any interference received, including interference that may cause undesired operation.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/ TV technician for help.

Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

Caution: Exposure to Radio Frequency Radiation.

This device must not be co-located or operating in conjunction with any other antenna or transmitter.

Responsible party:

Otto Bock Health Care, LP
3820 West Great Lakes Drive
Salt Lake City, Utah 84120-7205 USA
Phone + 1-801-956-2400
Fax + 1-801-956-2401

This device complies with RSS 210 of Industry Canada.

Operation is subject to the following two conditions:

- (1) this device may not cause interference, and
- (2) this device must accept any interference, including interference that may cause undesired operation of this device.

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Caution: Exposure to Radio Frequency Radiation.

The installer of this radio equipment must ensure that the antenna is located or pointed such that it does not emit RF field in excess of Health Canada limits for the general population; consult Safety Code 6, obtainable from Health Canada's website

<http://www.hc-sc.gc.ca/rpb>.

Responsible party:

Otto Bock Healthcare Canada Ltd.

5470 Harvester Road

L7L 5N5 Burlington, Ontario

Canada

Phone + 1-800-665-3327

Caution: Federal law (USA) restricts this device to sale by or on the order of a practitioner licensed by law of the State in which he/she practices to use or order the use of the device.

11 Technical data

Ambient conditions	
Storage in original packaging	+5 °C/41 °F to +40 °C/104 °F Max. 85% relative humidity, non-condensing
Transport in original packaging	-25 °C/-13 °F to +35 °C/95 °F 15% to 90% relative humidity, non-condensing +35 °C/95 °F to +70 °C/158 °F Water vapour pressure up to 50 hPa
Storage between subsequent applications	-25 °C/-13 °F to +35 °C/95 °F 15% to 90% relative humidity, non-condensing +35 °C/95 °F to +70 °C/158 °F Water vapour pressure up to 50 hPa
Operation	+5 °C/41 °F to +40 °C/104 °F 15% to 90% relative humidity, non-condensing; air pressure 533 hPa to 1060 hPa

Myo Plus TR	
Reference number	13E520
Dimensions	67 x 27 x 9.2 mm/2.64 x 1.1 x 0.36 inch
Weight	15 g/0.53 oz
Operating voltage	6 V–11.1 V DC
Current draw	Max. 25 mA
Power supply	757B35=3, =5 MyoEnergy Integral or 757B20, 757B21 Energy Pack
Lifetime	5 years
Operating time	MyoBock hand: approx. 14 hours bebionic hand: approx. 8 hours

Data transfer	
Wireless technology	Bluetooth Smart Ready
Range	min. 3 m/9.84 ft
Frequency range	2402 MHz to 2480 MHz
Modulation	GFSK, $\pi/4$ DQPSK, 8DPSK

Data transfer	
Data rate (over the air)	2178 kbps (asymmetrical)
Maximum output power (EIRP):	+8.5 dBm
Remote electrode	
Reference number	13E400=*, 13E401=*
Product service life	5 years
Frequency bandwidth	80 - 500 Hz
Sensitivity range	1800x
Myo Plus app	
Reference number	560X18-ANDR=V* / 560X18-IOS-V*
Supported operating system	See the information in the Apple App Store or Google Play Store regarding compatibility with mobile devices and versions.

12 Appendices

12.1 Symbols Used

Type BF applied part

Declaration of conformity according to the applicable European directives

In some jurisdictions it is not permissible to dispose of these products with unsorted household waste. Disposal that is not in accordance with the regulations of your country may have a detrimental impact on health and the environment. Please observe the instructions of your national authority pertaining to return and collection.

Manufacturer

Compliance with the requirements under "FCC Part 15" (USA)

Non-ionising radiation

Compliance with the requirements under the "Radiocommunications Act" (AUS)

Lot number (PPPP YYYY WW)
PPPP – plant

YYYY – year of manufacture

WW – week of manufacture

REF

Article number

BT ID: YYWWNNN

Bluetooth ID number

12.2 Status signals

Charge level of the 757B35=* MyoEnergy Integral battery


The battery charge level can be queried at any time.

- 1) With the prosthesis switched on, press the charging receptacle button and hold for less than one second.
- 2) The LED display on the charging receptacle provides information on the current battery charge level (Status signals).

Charging receptacle	Event
	Battery fully charged (illuminated in green)
	Battery 50% charged (illuminated in yellow)
	Battery drained (illuminated in orange)

Charge level of the 757B20/757B21 EnergyPack battery

The charge level is activated for a few seconds when the battery is inserted into the prosthetic socket.

LED indicator	Event
	Charge level above 50% (illuminated in green)
	Charge level below 50% (flashing alternating green and orange)
	Charge level below 5% (illuminated in orange)

12.3 Operating States

12.3.1 Beep signals

The following beep signals are set by default:

Beep signal	Additional display	Event
1 x long	-	<ul style="list-style-type: none">• Turn off product on charging receptacle• Charging begins (charging plug connected to charging receptacle)• Charging ends (charging plug disconnected from charging receptacle)
2 x short	LED on charging receptacle lights up briefly	Switching on the product
3 x short	-	Battery voltage too low, product shuts off automatically

12.4 Troubleshooting

Event	Cause	Required action
Unexpected prosthesis behaviour in various arm positions, e.g. during overhead use	Use in this arm position was not taken into account during calibration, or was incorrectly configured in the app.	<ul style="list-style-type: none"> Note arm position during calibration Take the weight into account during calibration Adjust the values in the app under "Settings > Advanced Settings"
Unexpected prosthesis behaviour in case of different weight settings, e.g. when loading the prosthesis with excess weight		
Hand cannot be opened any more	Faulty control	<ul style="list-style-type: none"> Use mechanical overload protection: Pull fingers open with great force Rotate hand with great force
Hand cannot be rotated any more		
-	General fault	<ul style="list-style-type: none"> Switch the prosthesis off and back on Consult O&P professional

12.5 Directives and manufacturer's declaration

12.5.1 Electromagnetic environment

This product is designed for operation in the following electromagnetic environments:

- Operation in a professional healthcare facility (e.g. hospital, etc.)
- Operation in areas of home healthcare (e.g. use at home, use outdoors)

Observe the safety notices in the section "Information on proximity to certain areas" (Information on Proximity to Certain Areas).

Electromagnetic emissions

Interference measurements	Compliance	Electromagnetic environment directive
HF emissions according to CISPR 11	Group 1/class B	The product uses HF energy exclusively for its internal functioning. Its HF emissions are therefore very low, and interference with neighbouring electronic devices is unlikely.
Harmonics according to IEC 61000-3-2	Not applicable – power below 75 W	–
Voltage fluctuations/flicker according to IEC 61000-3-3	Product meets the requirements of the standard.	–

Electromagnetic interference immunity

Phenomenon	EMC basic standard or test procedure	Interference immunity test level
Electrostatic discharge	IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air,
High-frequency electromagnetic fields	IEC 61000-4-3	10 V/m 80 MHz to 2.7 GHz 80% AM at 1 kHz

Phenomenon	EMC basic standard or test procedure	Interference immunity test level
Magnetic fields with rated power frequencies	IEC 61000-4-8	30 A/m 50 Hz or 60 Hz
Electrical fast transients/bursts	IEC 61000-4-4	± 2 kV 100 kHz repetition rate
Surges Line against line	IEC 61000-4-5	± 0.5 kV, ± 1 kV
Conducted interference induced by high-frequency fields	IEC 61000-4-6	3 V 0.15 MHz to 80 MHz 6 V in ISM and amateur frequency bands between 0.15 MHz and 80 MHz 80% AM at 1 kHz
Voltage drops	IEC 61000-4-11	0% U_T ; 1/2 period At 0, 45, 90, 135, 180, 225, 270 and 315 degrees
		0% U_T ; 1 period and 70% U_T ; 25/30 periods Single phase: at 0 degrees
Voltage interruptions	IEC 61000-4-11	0% U_T ; 250/300 periods

Interference resistance against wireless communication devices

Test frequency [MHz]	Frequency band [MHz]	Radio service	Modulation	Maximum power [W]	Distance [m]	Interference immunity test level [V/m]
385	380 to 390	TETRA 400	Pulse modulation 18 Hz	1.8	0.3	27
450	430 to 470	GMRS 460, FRS 460	FM ± 5 kHz deviation 1 kHz sine	1.8	0.3	28
710	704 to 787	LTE band 13, 17	Pulse modulation 217 Hz	0.2	0.3	9
745						
780						
810	800 to 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, GSM 800/900, LTE band 5	Pulse modulation 18 Hz	2	0.3	28
870						
930						

Test frequency [MHz]	Frequency band [MHz]	Radio service	Modulation	Maximum power [W]	Distance [m]	Interference immunity test level [V/m]
1,720	1,700 to 1,990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE band 1, 3, 4, 25; UMTS	Pulse modulation 217 Hz	2	0.3	28
1,845						
1,970						
2,450	2,400 to 2,570	Bluetooth WLAN 802.1-1 b/g/n, RFID 2450 LTE band 7	Pulse modulation 217 Hz	2	0.3	28
5,240	5,100 to 5,800	WLAN 802.1-1 a/n	Pulse modulation 217 Hz	0.2	0.3	9
5,500						
5,785						

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