


## 3WR95

<b>DE</b> Gebrauchsanweisung .....	5
<b>EN</b> Instructions for use .....	14
<b>FR</b> Instructions d'utilisation .....	23
<b>IT</b> Istruzioni per l'uso .....	32
<b>ES</b> Instrucciones de uso .....	42
<b>PT</b> Manual de utilização .....	52
<b>NL</b> Gebruiksaanwijzing .....	61
<b>SV</b> Bruksanvisning .....	71
<b>DA</b> Brugsanvisning .....	79
<b>NO</b> Bruksanvisning .....	88
<b>FI</b> Käyttöohje .....	96
<b>PL</b> Instrukcja użytkowania .....	105
<b>HU</b> Használati utasítás .....	114
<b>CS</b> Návod k použití .....	123
<b>RO</b> Instrucțiuni de utilizare .....	131
<b>HR</b> Upute za uporabu .....	141
<b>TR</b> Kullanma talimatı .....	150
<b>EL</b> Οδηγίες χρήσης .....	159
<b>RU</b> Руководство по применению .....	169
<b>JA</b> 取扱説明書 .....	179
<b>ZH</b> 使用说明书 .....	187

**INFORMATION**

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### INFORMATION

Date of last update: 2021-01-25

- Please read this document carefully before using the product and observe the safety notices.
- Instruct the user in the safe use of the product.
- Please contact the manufacturer if you have questions about the product or in case of problems.
- Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- Please keep this document for your records.

## 1 Product description

### 1.1 Function

The Aqua knee is a knee joint component intended for use in a waterproof modular prosthesis and can be used to walk and stand in wet areas. The upper and lower parts of the joint are made of aluminium alloy and connected to one another by a specially designed axis. A small hydraulic unit is integrated into the joint to control the swing phase. It generates dynamic resistance to movement which optimises the gait pattern by preventing excessive swing through during flexion and an excessively hard strike at the end of the extension phase.

The movement resistance (damping behaviour) can be controlled by changing the flow cross-section. The flexion and extension movements can be adjusted independently.

- The control geometry of the hydraulics is designed to create a high flexion resistance even with a high degree of flexion.
- A special feature of the hydraulics reduces the flexion angle at high walking speeds.
- The extension resistance can be set to a low value, since hydraulic end position damping ensures a soft transition to the elastic stop.

Unlike mechanical gait controllers, the hydraulics automatically adjust to different walking speeds which results in harmonious motion. In the 3WR95, the flexion angle of approx. 135° is limited by the stop of the modular adapter against the lower part of the joint (Fig. 1).

## 1.2 Combination possibilities

### INFORMATION

In a prosthesis, all prosthetic components have to meet the patient's requirements regarding the amputation level, body weight, activity level, environmental conditions and field of application.

## 2 Application

### 2.1 Indications for use

#### ⚠ CAUTION

##### Reuse on another patient

Fall due to loss of functionality as well as damage to the product

- Only use the product for a single patient.

The 3WR95 Aqua Knee with integrated miniature hydraulics is intended **exclusively** for the exo-prosthetic fitting of the lower extremities.

### 2.2 Area of application

The 3WR95 Knee Joint can be used for transfemoral amputees and patients with disarticulated knees.

**Approved for a body weight up to 150 kg (330 lbs).**

#### ⚠ CAUTION

**Risk of falls due to improper operation.** The amputee must have the necessary muscular and cognitive abilities to control a freely moving knee joint without stance phase control. If he or she is not able to do so (this usually applies to patients with mobility grades 1 and 2), then the knee joint must always be locked in order to walk.

### 2.3 Conditions of use

The Aqua Knee is a knee joint component intended for use in waterproof walking devices and can be used to walk and stand in wet areas. The waterproof walking device is not intended to replace the initial prosthesis and should only be used in wet areas. It may also not be used for extreme activities such as extreme sports (water skiing, kite surfing, etc.).

#### ⚠ CAUTION

**Overheating of the hydraulic unit.** A continuous high level of activity can cause the hydraulic unit to overheat. Please note that burns may result from contact with overheated components.

## 2.4 Environmental conditions

### CAUTION

#### Use under unallowable environmental conditions

Fall due to damaged product

- Do not expose the product to unallowable environmental conditions (see the table „Unallowable environmental conditions“ in this section).
- If the product was exposed to unallowable environmental conditions, take suitable steps (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop etc.).

#### Allowable environmental conditions

Temperature range for use: -10 °C to +60 °C

Storage and transport conditions: -10 °C to +60 °C, 250 mbar to 1013 mbar, condensing

Acid resistance: salt water, chlorine water, soapsuds, perspiration, urine

Moisture: fresh water, salt water, chlorine water, soapsuds, steam

Allowable immersion depth: 3 m

#### Unallowable environmental conditions

Mechanical vibrations or impacts

Dust, sand, highly hygroscopic particles (e.g. talcum)

Sauna

## 2.5 Reuse and lifetime

### CAUTION

#### Reuse on another patient

Fall due to loss of functionality as well as damage to the product

- Only use the product for a single patient.

### CAUTION

#### Exceeding the lifetime

Fall due to change in or loss of functionality and damage to the product

- Ensure that the maximum lifetime defined in this section is not exceeded.

This product was tested by the manufacturer with 3 million load cycles. Depending on the user's activity level, this corresponds to a maximum lifetime of 5 years.

## 2.6 Required qualification

### 2.6.1 Prosthetist qualifications

Fitting of the 3WR95 Knee Joint may only be performed by prosthetists.


### 2.6.2 Qualification of the therapist or nursing staff

The therapist or the nursing staff must be trained in the use of the 3WR95 Knee Joint (and accessories) by the authorized prosthetist.

### 3 Safety

#### 3.1 Explanation of the ANSI symbols

 Warning regarding possible risks of accident or injury.

 Warning regarding possible technical damage.

 Additional information on the fitting/use.

#### 3.2 General safety instructions

##### CAUTION

**Non-observance of safety instructions.** Non-observance of the following safety instructions can lead to malfunction of the 3WR95 Knee Joint, resulting in a risk of injury to the patient and in destruction of the 3WR95 Knee Joint.

##### CAUTION

**Failure to comply with torque values.** Loss of functionality due to structural failure may result if the tightening torques required for the screw connections are not observed. This can cause the patient to fall.

Observing the specified tightening torques is essential for safe operation.

##### CAUTION

**Incorrect use of the knee joint.** Any sort of excessive strain or overloading or incorrect handling can lead to defects and resulting malfunction of the 3WR95 Knee Joint, including loss of function due to structural failure. This can cause the patient to fall.

Instruct the patient in the proper handling of the 3WR95 Knee Joint and explain the following patient information.

#### 3.3 Safety instructions for patients / users

##### CAUTION

**Risk of burns in case of heat exposure.** The 3WR95 Knee Joint is not suitable for use in saunas or for activities during which the prosthesis is exposed to high heat.

##### CAUTION

**Mechanical overloading.** External mechanical loads, such as impacts and vibration, can cause defects and result in malfunction of the 3WR95 Knee Joint, including loss of function due to structural failure. This can cause the patient to fall.

Do not subject prosthetic components to mechanical vibrations or impacts.

##### CAUTION

**Risk of loss of function due to incorrect care.** The 3WR95 Knee Joint must be cleaned with clean fresh water following contact with saltwater, water containing chlorine/soap, or contact with abrasive substances (e.g. sand).

It must also be dried carefully immediately after every use.

Contact with abrasive substances (e.g. sand) can lead to premature wear.

### CAUTION

**Manipulation of prosthetic components.** Any changes or modifications to prosthetic components which are made on the user's own initiative may result in defects and malfunction of the 3WR95 knee joint, including loss of function due to structural failure. This can cause the patient to fall.

- Any changes or modifications to the 3WR95 may result in a loss of functionality.
- Opening and repairing the 3WR95 Knee Joint and/or the repair of damaged components may only be performed by a certified Ottobock Service provider.

### CAUTION

**Environment and usage conditions.** The environment and usage conditions may compromise the functionality of the 3WR95 Knee Joint. To reduce the risk for the patient, the 3WR95 Knee Joint must not be used after noticeable functional changes occur. Noticeable functional changes may include poor response, incomplete extension, diminished swing phase control or stance phase stability, abnormal noises, etc.

Please contact your practitioner to have the prosthesis inspected if you discover damage or noticeable functional changes.

### CAUTION

**Incorrect use of the knee joint.** Any sort of excessive strain or overloading or incorrect handling can lead to defects and resulting malfunction of the 3WR95 Knee Joint, including loss of function due to structural failure.

This can cause the patient to fall.

The 3WR95 Knee Joint may also not be used for extreme activities such as extreme sports (water skiing, kite surfing, etc.).

Careful handling of the prosthesis and its components not only increases their service life, but is primarily in the interest of your own safety. If the prosthesis has been subjected to extreme loads (e.g. by falling or the like), then it must be inspected by a prosthetist immediately for possible damage. The contact person for this is the responsible prosthetist, who will send the prosthesis to Ottobock Service if necessary.

### CAUTION

**Risk of pinching where the joint bends.** Ensure that fingers and other body parts are not in this area when bending the joint.

## 4 Delivery

### 4.1 Scope of Delivery

1 pc. 3WR95 Aqua Knee

1 pc. 647G818=all\_INT Instructions for Use

## 5 Preparation for use

### 5.1 Alignment/assembly

The three-dimensional alignment of the prosthesis socket and the modular components affects the static and dynamic function of the prosthesis. The position of the axes affects the function of the joint.

The functions of the 3WR95 can only be put to optimal use if alignment was done correctly. **The**

**residual limb position must be taken into account for positioning of the socket connection.**

Plumb lines in the frontal and sagittal planes, which are marked during plaster casting and trial fitting from the hip pivot point, simplify proper positioning of the lamination anchor or socket adapter.

### **Perform alignment in 2 steps:**

1. First basic alignment is done using the alignment tool (e.g. 743L200 L.A.S.A.R. Assembly).
2. This is followed by static alignment optimisation using the 743L100 L.A.S.A.R. Posture.

#### **5.1.1 Basic alignment using the alignment tool (Fig. 6)**

1. Move the mid-foot forward approx. **30 mm relative to the alignment reference line.**
2. Set the effective heel height of the foot and **add 5 mm.** Set the outer position of the foot.
3. Clamp the 3WR95 Knee Joint in place. For basic alignment, **set the knee joint axis of rotation (alignment reference point) 15 mm back from the alignment reference line.** The joint should be positioned horizontally at this time. Pay attention to the knee-ground distance and the knee outer rotation (adapter inserts provide for a rotation of approx. 5°)! Recommended position of the alignment reference point: 20 mm above the medial tibial plateau.
4. Use the tube adapter to connect the foot to the modular knee joint.
5. Mark the lateral center of the socket with a centered, proximal and a centered, distal dot. Mark a line through both points from the edge to the end of the socket.
6. Position the socket so that the proximal center point of the socket coincides with the alignment reference line. Adjust the socket flexion to 3° – 5°, but take the individual situation (e.g. hip joint contractures) and the ischial tuberosity-to-ground distance into account!

#### **⚠ CAUTION**

**Errors during prosthesis alignment.** Errors during prosthesis alignment can lead to malfunction of the 3WR95 Knee Joint and to early wear. This can cause the patient to fall.

The joint will be positioned too far to the anterior if socket flexion is not taken into account. Move the lamination anchor to the correct position by foaming up, if necessary.

Position the socket so that the proximal center of the socket coincides with the alignment reference line.

Position the 4WR95=1 or 4WR95=2 Lamination Anchor along the alignment reference line.

Move the 3WR95 Knee Joint axis of rotation (alignment reference point) back 15 mm relative to the alignment reference line.

#### **5.1.2 Static alignment optimisation using the 743L100 L.A.S.A.R. Posture (Fig. 7)**

Significant optimisation of the bench alignment is possible using the L.A.S.A.R. Posture. In order to achieve sufficient safety while providing easy swing phase initiation, please proceed as follows:

1. To measure the load line, the transfemoral amputee steps on the L.A.S.A.R. Posture force measurement plate with the prosthetic fitting side and places the other leg on the height compensation plate. Sufficient weight should be placed on the prosthesis side (> 35 % body weight).
2. Now the alignment should only be adjusted by **changing the plantar flexion** so that the **load line** (laser line) runs approx. **45 mm in front of the knee axis** (alignment reference point).
3. After this step, perform dynamic optimisation during trial walking.

## 5.2 Adjustment

### 5.2.1 Adjusting the swing phase control

To avoid damage, the swing phase control adjustments must be made by the prosthetist using the 4X50 Adjustment Wrench inserted into the joint! **The adjustment wrench should be kept by the prosthetist after the trial fitting.** The flexion and extension resistance is set separately by positioning the pins on the left and right sides of the piston rod. The adjusting pins are **only accessible in the housing slots when the knee joint is flexed.**

### 5.2.2 Adjusting the flexion (Fig. 2)

Bend the knee joint all the way to the stop

Use the adjusting wrench to set the adjusting pin on the side labelled „**Flexion**“, according to the chart.

Turning the key **back** = **low** flexion resistance= **0**

Turning the key **forward** = **highest** flexion resistance= **8**

#### NOTICE

#### Adjusting procedure and removing the adjusting wrench.

The knee joint must be flexed during the adjusting procedure and before the adjusting wrench is removed from the housing.

### 5.2.3 Adjusting the extension (Fig. 3)

Flex the knee joint all the way to the stop

Use the adjusting wrench to set the adjusting pin on the side labelled „**Extension**“, according to the chart.

Turning the key **forward** = **low** extension resistance= **0**

Turning the key **backward** = **highest** extension resistance = **8**

#### ⚠ CAUTION

**Errors during the adjustment procedure.** If extension damping is set too high, the forward swing of the the 3WR95 Knee Joint may be restricted. This can cause the patient to fall.

Only dampen extension to the point where full extension is always attained.

## 5.3 Final assembly

Use the 710D4 Torque Wrench to tighten the set screws. **For the tightening torques, refer to the instructions for use of the adapters you are using .**

## 6 Use

### 6.1 Locking function

The locking function makes it possible to walk and stand in wet areas. To activate the locking function:

- The prosthesis leg must be extended fully while weight is applied to the forefoot (hyperextension).
- To lock the joint, move the locking lever in the direction indicated by the arrow (Fig. 4)

#### ⚠ CAUTION

**Risk of falling due to deactivated locking function.** The behaviour of the prosthesis and the swing phase changes in wet areas. For safety reasons, the amputee must activate the locking function before entering wet areas. The locking function will not be activated if the prosthesis leg is not fully extended when the locking lever is engaged.

- The knee joint is locked

### INFORMATION

Check the position of the locking lever before using the locking function.

- Unlocking the knee joint

The knee joint can be unlocked when necessary (e.g. to sit down). The knee joint is locked.

- The prosthesis leg must be fully extended while weight is applied to the forefoot (hyperextension).
- To unlock the joint, move the locking lever in the direction indicated by the arrow (Fig. 5)
- The knee joint is unlocked

The knee joint remains unlocked until it is placed in a hyperextension position again and locked using the locking lever.

### ⚠ CAUTION

**Risk of falls due to improper operation.** The amputee must have the necessary muscular and cognitive abilities to control a freely moving knee joint without stance phase control. If he or she is not able to do so (this usually applies to patients with mobility grades 1 and 2), then the knee joint must always be locked in order to walk.

## 6.2 Care

### ⚠ CAUTION

**Risk of loss of function due to improper care.** The 3WR95 Knee Joint must be cleaned with clean fresh water following contact with saltwater, water containing chlorine/soap, or contact with abrasive substances (e.g. sand). Do not use lipophilic cleaning agents, high pressure washers, or boiling hot water. The knee joint may not be cleaned in a dishwasher or washing machine. It must also be dried carefully immediately after every use. Contact with abrasive substances (e.g. sand) can lead to premature wear.

## 7 Maintenance/repair

### 7.1 Maintenance

### INFORMATION

We always recommend regular annual safety inspections.

The safety inspections must be performed by a prosthetist.

Do not disassemble the joint!

Ottobock recommends adjusting the knee joint to the patient's requirements again after the patient has taken some time to get used to the prosthesis. Special attention should be paid to movement resistance and abnormal noises.

Full flexion and extension must be ensured.

Replace the stop every time the prosthesis is serviced.

If you have problems, please send the joint in.

## 7.2 Repair

### CAUTION

**Disassembly on your own initiative.** If the 3WR95 is disassembled on your own initiative, this may lead to defects and resulting malfunction of the 3WR95 Knee Joint, including loss of function due to structural failure. This can cause the patient to fall.

Do not disassemble the 3WR95 Knee Joint on your own initiative! Please send the joint to Ottobock Service in case of problems.

## 8 Disposal

### 8.1 Disposal information

#### INFORMATION

This product may not be disposed of with regular domestic waste in all jurisdictions. Disposal that is not in accordance with the regulations of your country may have a detrimental impact on health and the environment. Please observe the information provided by the responsible authorities in your country regarding return and collection processes.

## 9 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

### 9.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

### 9.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

## 10 Appendices

### 10.1 Technical data

Proximal connection	Pyramid adapter
Distal connection	Pyramid adapter
Operating and storage temperature	-10 °C to +60 °C (15 °F to 140 °F)
Weight of the joint	approx. 400 g
System height	62 mm
Proximal system height up to the alignment reference point	6 mm
Distal system height up to the alignment reference point	56 mm
Maximum flexion angle	135°
Maximum patient weight	150 kg / 330 lbs

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