





HR5 1040000-027, HR5 1050000-027

DE Gebrauchsanweisung	4
EN Instructions for use	12
FR Instructions d'utilisation	20
IT Istruzioni per l'uso	28
ES Instrucciones de uso	36
PT Manual de utilização	44
NL Gebruiksaanwijzing	52
SV Bruksanvisning	60
DA Brugsanvisning	68
NO Bruksanvisning	76
FI Käyttöohje	84
PL Instrukcja użytkowania	92
HU Használati utasítás	100
SK Návod na používanie	108
RU Руководство по применению	116
JA 取扱説明書	125

1 Vorwort

INFORMATION

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INFORMATION

Date of last update: 2020-12-16

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Obtain instruction from the qualified personnel in the safe use of the product.
- ▶ Please contact the qualified personnel if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

INFORMATION

- ▶ New information regarding product safety and product recalls as well as the declaration of conformity can be obtained at ccc@ottobock.com or from the manufacturer's service department (see inside or outside of back cover for addresses).
- ▶ You can request this document as a PDF file at ccc@ottobock.com or from the manufacturer's service department (see inside or outside of back cover for addresses). The PDF file can also be displayed in a larger size.

These instructions for use provide you with important information regarding the assembly and proper use of the HR51040000-027 and HR51050000-027 Robby bathing supports.

All users and/or their attendants must be trained by qualified personnel in the use of the product with the aid of these instructions for use. In particular, users and/or attendants must be informed of the residual risks with the aid of the safety notices.

The manufacturer reserves the right to make technical changes to the model described in these instructions for use.

2 Product description

The bathing support is a bathing and showering aid that can be adjusted in various ways to meet the child's needs. The child can be positioned lying down or sitting in the bathing support.

Hook-and-loop straps are used to adjust the tension of the fabric cover to the child's needs, creating a recess to suit the muscle tension.

The head section can swing away to the rear, or the fabric cover on this section can be pulled off to keep the head clear, e.g. while washing hair.

The back support length and seat depth can be adjusted to the child's height by telescoping the frame tubes.

The frame cover consists of two pieces of fabric that overlap in the pelvis area of the back support and are connected with a hook-and-loop closure. The hook-and-loop closure can be opened and the required length of fabric from the overlap can be used, e.g. when adjusting the length of the bathing support.

3 Intended use

The safe use of the product can only be ensured in case of intended use in accordance with the information contained in these instructions for use. The user is ultimately responsible for accident-free operation.

3.1 Indications for use

The product is intended to support and optimise the seating position of children while bathing.

The product is suitable for users with intact skin whose anatomy (such as body dimensions and weight) permits the intended use of the product.

Depending on the product that is used, additional accessories may be needed for the intended use. The prerequisites for safe installation, adjustment and use are described in the instructions for use of the respective product and must be met in all cases.

The products must be replaced exclusively with original parts from the manufacturer. All information in the instructions for use of the products that are used must be observed for combinations with other products. Ottobock assumes no liability for combinations with medical devices and/or accessories from other manufacturers not included in Ottobock's modular system.

Combinations based on a combination agreement that have been evaluated for effectiveness and safety are an exception to this.

3.2 Area of application

The bathing support can be used indoors and outdoors.

It is used mainly for children and adolescents with a lack of torso control and/or hyperextension tendency (e.g. in case of cerebral palsy) starting around the second year of life.

The product is suitable for users with a body weight of up to **30 kg [66 lbs]** (size 1) and up to **60 kg [132 lbs]** (size 2). We assume no liability if the maximum permissible user weight is exceeded.

3.3 Indications

- Minor, pronounced or complete restrictions of movement

3.4 Contraindications

3.4.1 Absolute Contraindications

- None known

3.4.2 Relative Contraindications

- Failure to meet physical or mental requirements

3.5 Further Usage Restrictions

The product is designed for use by **one person** only. Re-use of the product is not permitted due to hygiene considerations.

3.6 Qualification

This product may only be selected and adapted by qualified personnel. Compliance with all manufacturer specifications and all applicable legal provisions is required. Please contact the manufacturer's service department for further information (see inside or outside of rear cover for addresses).

4 Safety

4.1 Explanation of warning symbols

 WARNING	Warning regarding possible serious risks of accident or injury.
 CAUTION	Warning regarding possible risks of accident or injury.
 NOTICE	Warning regarding possible technical damage.

4.2 General safety instructions

WARNING

Improper handling of packaging materials

Risk of suffocation due to neglect of the duty to supervise

- ▶ Packaging materials must be kept out of the reach of children.

WARNING

Independent modification of settings

Serious injuries to the user due to improper changes to the product

- ▶ Do not modify the settings established by the qualified personnel. Only the settings described in the section "Use" in these instructions for use may be adjusted independently.
- ▶ In case of problems with the settings, please contact the qualified personnel who adjusted your product.

WARNING

Exceeding the service life

Serious injuries due to failure to observe the manufacturer's requirements

- ▶ Using the product beyond the specified expected service life leads to increased residual risk.
- ▶ Observe the specified service life.

⚠ WARNING

Inflammable product

Burns due to user error

- ▶ The product is flame retardant. Nevertheless, the possibility that it may catch fire cannot be entirely excluded in case of exposure to an ignition source. Therefore, the utmost caution must be exercised in the vicinity of an open flame.
- ▶ Keep away from all ignition sources, especially lit cigarettes.

⚠ CAUTION

Extreme temperatures

Hypothermia or burns due to contact with components, failure of components

- ▶ Do not expose the product to any extreme temperatures (e.g. direct sunlight, sauna, extreme cold).
- ▶ Do not leave the product in the immediate vicinity of heaters.

⚠ CAUTION

Skin damage

Skin damage or pressure points due to overloading

- ▶ Check your skin for intactness before and during use of the product.
- ▶ Pay attention to diligent skin care and pressure redistribution during interruptions in using the product.
- ▶ If skin damage or other problems occur during use, stop using the product. Consult the qualified personnel.

⚠ CAUTION

Improper use of the product

Falling of the user due to neglect of the duty to supervise

- ▶ The product is a therapy aid and may only be used under supervision. Never leave your child unattended.
- ▶ Pay particular attention to potentially hazardous situations, e.g. while setting up and positioning the product.
- ▶ In the interest of the user's safety, the correct function of the therapy aid must be checked before every use.

⚠ CAUTION

Damage to the product

Health problems for the user due to improper handling

- ▶ Stop using the product immediately in the event of defects or other hazards that can result in personal injury.
- ▶ Contact the qualified personnel.
- ▶ Report damage to the product that may cause health problems for the user to Ottobock.

⚠ CAUTION

Use of the product during diagnostic examinations and therapeutic treatment

Impairment of the examination results or the effectiveness of treatment due to interactions of the product with devices that are used

- ▶ Make sure that examinations and treatments are carried out exclusively under the prescribed conditions.

⚠ CAUTION

Improper use

Product tipping over due to incorrect handling

- ▶ Only use the product on level, solid surfaces.
- ▶ Do not place the product on steep inclines.

⚠ CAUTION

Risk of pinching and crushing on moving parts and mechanical adjustments

Pinching, crushing of limbs due to improper handling

- ▶ Pay attention to your fingers while setting up and folding the product (risk of pinching and crushing).

4.3 Side effects

There are no known side effects given proper use.
Contact a doctor or therapist in case of problems.

4.4 Nameplate

Label	Meaning
	A Manufacturer's product name
	B CE marking
	C Maximum user weight (see section "Technical data")
	D Manufacturer information/address
	E Serial number ¹⁾
	F Manufacturing date ²⁾
	G Symbol for medical device
	H WARNING! Read the instructions for use before using the product. Observe important safety-related information (e.g. warnings, precautions).
	I Manufacturer's reference number for the product variant
	J Serial number (PI) ^{3),1)}
K Global Trade Item Number (DI) ⁴⁾	

1) YYYY = year of manufacture; WW = week of manufacture; PP = production site; XXXX = sequential production number

2) YYYY = year of manufacture; MM = month of manufacture; DD = day of manufacture

3) UDI-PI to GS1 standard; UDI = Unique Device Identifier, PI = Product Identifier

4) UDI-DI to GS1 standard; UDI = Unique Device Identifier, DI = Device Identifier

Additional numbers may be applied to the product for the identification of certain manufacturer data.

5 Delivery

5.1 Scope of delivery

Instructions for use are included with the product.

The scope of delivery (see fig. 1) includes:

Item	Article designation	Quantity
	Bathing support, consisting of:	1
1	Lower leg support	1
2	Seat	1
3	Back support	1
4	Head section	1

5.2 Accessories

Head support

The head support (see fig. 2) provides lateral support for the head. It consists of two foam pads covered with fabric (see fig. 3).

Torso strap

The torso strap (see fig. 4) provides added support for the child when stable positioning is difficult due to involuntary, sudden movements.

Like the head support, the torso strap can be adjusted to the desired height using the lengthwise straps on the cover fabric. Decreasing the tension of the cover fabric creates a recess to provide added support for the child.

Another way to give the child's torso more support and security is to use the head support to guide the torso. Its special shape and size makes it suitable for this area (see fig. 5).

Leg straps

The leg straps (see fig. 6) are used to secure the legs to counteract extension movements that interfere with positioning and the care of the child.

Like the head support and torso strap, the height of the straps is adjustable (see fig. 7). These straps can also be used to hold and guide the thighs when necessary.

5.3 Storage

Store the product in a dry, enclosed room with sufficient air circulation and protection from external influences. Specific information about storage conditions: see page 18.

6 Preparing the product for use

6.1 Installation

INFORMATION

The bathing support is easy and straightforward to adjust when dry.

INFORMATION

The bathing support can be set to a sitting (see fig. 8) or reclining (see fig. 9) position. While the position (see fig. 10) is adjustable, the bathing support may rock if the child makes sudden movements (e.g. due to abrupt spasticity, seizures, athetosis). This can make the child feel insecure, and he or she can slide out if not secured.

- 1) Stabilise the bathing support by placing the head section against your abdomen or chest, and face the foot section away from your body (or the other way around) (see fig. 11). You can use your upper body to make the adjustment and handling easier.
INFORMATION: The adjustment position (see fig. 12) requires too much strength and arm length, and is not recommended.
- 2) Adjust the reclining position using the ratchet joints (see fig. 13). Push the press button to move the joint (see fig. 14). The joint engages in the desired position after the press button is released.
- 3) Position the individual sections of the bathing support relative to each other with the ratchet joints A, B, C and D (see fig. 15). To make an adjustment, press and hold both ratchet joints across from each other at the same time during the adjustment (see fig. 16).
INFORMATION: We recommend using the balls of your thumbs to press the buttons. At the same time, use your fingers (particularly the index and middle fingers) to hold the respective section of the bathing support and move it to the desired position. Once the desired position is reached, release the press buttons on the joints. The joints then engage automatically.

6.2 Settings

Adjusting the head section

- 1) Use the ratchet joints D (see fig. 15) to fold the head section back so you can move the head to the desired position (see fig. 17).
- 2) **Option:** The fabric cover can also be pulled off the head section. The child's head is completely unsupported in this case and should be held by the attendant in case of poor head control (see fig. 18).
INFORMATION: In case of severe hyperextension of the head, the head section can also be adjusted forward to reduce the excessive muscle tension in conjunction with a more pronounced hip flexion setting.
- 3) Adjust the lateral head supports to the correct width: Lift the outer end of the pad you want to adjust and slide the pad to the desired position (see fig. 19). Then fasten the pad to the guide straps again by pressing and securing the hook-and-loop closures. Proceed the same way with the second pad.
- 4) After setting the head support to the correct height, secure it with the lengthwise straps on the cover fabric (see fig. 20).

Adjusting the length of the back section

- 1) Use an Allen key to loosen the two screws in the adjustment holes on the side of the back section (see fig. 21).
- 2) Loosen the hook-and-loop closure of the two overlapping fabric sections in the pelvis area of the back support (see fig. 22).
- 3) Pull up on the head section while holding the seat section with the other hand until the desired position is reached.
- 4) Reinstall the two screws in the same hole on the respective side and tighten them (see fig. 23).
- 5) Pull the two fabric sections back over each other so they are tight and connect them with the hook-and-loop closure (see fig. 24).

Adjusting the length of the seat

- 1) Use an Allen key to loosen the two screws in the adjustment holes on the side of the seat section (see fig. 25).
- 2) Loosen the hook-and-loop closure of the two overlapping fabric sections in the pelvis area of the back support (see fig. 26).

INFORMATION: Because the required fabric comes from the back, the bottom tensioning strap of the back section may have to be loosened depending on how much the seat depth is being adjusted.

- 3) Pull on the foot section to set the desired length.
- 4) Retighten both screws when the desired position is reached.
- 5) Pull the two fabric sections back over each other so they are tight and connect them with the hook-and-loop closure.
- 6) Refasten the bottom tensioning strap of the back section if applicable.

Adjusting the fabric cover

To adjust the tension of the cover fabric to the child's needs, loosen the hook-and-loop connections and adjust the tension of the fabric as needed (see fig. 27).

7 Use

7.1 Using the product

Before using the product, verify that it is securely mounted and in an ergonomically correct position. Inform the qualified personnel if corrections are required.

- 1) Set the bathing support to the ergonomically correct position (see page 16).
- 2) Place the bathing support on a firm, level surface.

INFORMATION: Make sure that the back support height adjustment frame, ratchet joint B (see fig. 15) and the end of the lower leg support are always in contact with the supporting surface (see fig. 28).

- 3) Wipe the bathing support dry after every use and allow to air dry. Set the bathing support with the head section facing down for drying (see fig. 29). The rubber tube protectors prevent slipping and ensure that moisture drains away through the appropriate holes.
- 4) The bathing support can be folded to save space after use (see fig. 30).

7.2 Care

CAUTION

Lack of or improper cleaning

Health hazard due to infections, damage to the product due to user error

- ▶ Clean the product at regular intervals.
- ▶ Do not clean the product with a pressure washer. Doing so can damage the product.

7.2.1 Cleaning

Clean the product regularly depending on the degree of soiling and frequency of use, **at least 1x per month:**

- Use warm soapy water for cleaning. Add some disinfectant if needed. Carefully dab the product with a damp cloth or sponge and wipe with a dry, clean cloth.
- The mesh material can be taken off for care and washed at max. **+60 °C [+140 °F]**. Before washing, fasten the hook-and-loop fasteners of the straps and place the cover fabric into a pillowcase, which makes unfolding easier after washing and prevents tangling with other laundry.
- Do not use any aggressive cleaners, solvents or hard brushes etc.
- Clean the hook-and-loop closure with a dry brush.
- Allow the product to fully air dry after cleaning.
- Do not expose the product to direct heat (e.g. sunshine, stove or radiator heat, hair dryer).

7.2.2 Disinfection

- Thoroughly clean the product before disinfecting.
- Wipe all parts of the product with a disinfectant.
- Only use colourless water-based disinfectants. Observe the instructions for use provided by the manufacturer.

8 Maintenance and repair

⚠ CAUTION

Wear and tear of the padding on the product

Risk of injury due to unallowable continued use, loss of functionality

► Frequent use can result in wearing of the padding surfaces. Have the product replaced immediately in case of damage.

- The function of the product should be checked **before each use**.
- Check the product **monthly** for wear and damage.
- Check screw connections **monthly** for tightness.
- Especially during the initial period of use or after making adjustments, screw connections must be checked for tightness more often. If a screw connection loosens repeatedly, contact the qualified personnel promptly.
- The user is not permitted to make repairs to the product.
- Have damaged or worn parts replaced by qualified personnel as soon as possible.
- Servicing and repairs may only be carried out by qualified personnel or the manufacturer. This will ensure that only Ottobock spare parts are used for repairs.

9 Disposal

9.1 Disposal information

Return the product to the qualified personnel for disposal.

All components of the product must be disposed of properly in accordance with the respective national environmental regulations.

9.2 Information on re-use

The product is not intended for re-use.

10 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

10.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

10.2 Warranty

Further information regarding the warranty terms and conditions may be obtained from the qualified personnel who fitted this product, or from the manufacturer's customer service department (see inside back cover of the instructions for use for addresses).

10.3 Lifetime

Expected lifetime: **4 years**

The design, manufacturing and requirements for the intended use of the product are based on the expected lifetime. These also include the requirements for maintenance, ensuring effectiveness and the safety of the product.

11 Technical data

Article number	HR51040000-027 (size 1)	HR51050000-027 (size 2)
Overall length of the support surface	1150–1330 mm	1470–1730 mm
Seat depth	200–240 mm	310–430 mm
Back support length	630–720 mm	890–940 mm
Width	380 mm	380 mm
Weight	3.8 kg	4.3 kg
Max. allowable user weight	30 kg	60 kg

Ambient conditions

Temperatures and relative humidity	
Temperature during use [°C (°F)]	-10 to +40 (14 to 104)
Transport and storage temperature [°C (°F)]	-10 to +40 (14 to 104)
Relative humidity [%]	10 to 85; non-condensing

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A series of horizontal lines spanning the width of the page, providing a template for writing or editing text.

Lined writing area consisting of 30 horizontal lines.

A series of horizontal lines spanning the width of the page, providing a template for writing or editing text.

Kundenservice/Customer Service

Europe

Otto Bock HealthCare Deutschland GmbH
Max-Näder-Str. 15 · 37115 Duderstadt · Germany
T +49 5527 848-3433 · F +49 5527 848-1460
healthcare@ottobock.de · www.ottobock.de

Otto Bock Healthcare Products GmbH
Brehmstraße 16 · 1110 Wien · Austria
F +43 1 5267985
service-admin.vienna@ottobock.com · www.ottobock.at

Otto Bock Adria d.o.o. Sarajevo
Ramiza Salčina 85
71000 Sarajevo · Bosnia-Herzegovina
T +387 33 255-405 · F +387 33 255-401
obadria@bih.net.ba · www.ottobockadria.com.ba

Otto Bock Bulgaria Ltd.
41 Tzar Boris III' Blvd. · 1612 Sofia · Bulgaria
T +359 2 80 57 980 · F +359 2 80 57 982
info@ottobock.bg · www.ottobock.bg

Otto Bock Suisse AG
Luzerner Kantonsspital 10 · 6000 Luzern 16 · Suisse
T +41 41 455 61 71 · F +41 41 455 61 70
suisse@ottobock.com · www.ottobock.ch

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Otto Bock South Africa (Pty) Ltd
Building 3 Thornhill Office Park · 94 Bekker Road
Midrand · Johannesburg · South Africa
T +27 11 564 9360
info-southafrica@ottobock.co.za
www.ottobock.co.za

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Otto Bock HealthCare Canada
5470 Harvester Road
Burlington, Ontario, L7L 5N5, Canada
T +1 800 665 3327 · F +1 800 463 3659
CACustomerService@ottobock.com
www.ottobock.ca

Oficina Ottobock Habana
Calle 3ra entre 78 y 80.
Edificio Jerusalem · Oficina 112 · Calle 3ra.
Playa, La Habana. Cuba
T +53 720 430 69 · +53 720 430 81
hector.corcho@ottobock.com.br
www.ottobock.com.br

Otto Bock HealthCare Andina Ltda.
Calle 138 No 53-38 · Bogotá · Colombia
T +57 1 8619988 · F +57 1 8619977
info@ottobock.com.co · www.ottobock.com.co

Otto Bock de Mexico S.A. de C.V.
Prolongación Calle 18 No. 178-A
Col. San Pedro de los Pinos
C.P. 01180 México, D.F. · Mexico
T +52 55 5575 0290 · F +52 55 5575 0234
info@ottobock.com.mx · www.ottobock.com.mx

Otto Bock HealthCare LP
11501 Alterra Parkway Suite 600
Austin, TX 78758 · USA
T +1 800 328 4058 · F +1 800 962 2549
USCustomerService@ottobock.com
www.ottobockus.com

Asia/Pacific

Otto Bock Australia Pty. Ltd.
Suite 1.01, Century Corporate Centre
62 Norwest Boulevard
Baulkham Hills NSW 2153 · Australia
T +61 2 8818 2800 · F +61 2 8814 4500
healthcare@ottobock.com.au · www.ottobock.com.au

Beijing Otto Bock Orthopaedic Industries Co., Ltd.
B12E, Universal Business Park
10 Jiuxianqiao Road, Chao Yang District
Beijing, 100015, P.R. China
T +8610 8598 6880 · F +8610 8598 0040
news-service@ottobock.com.cn
www.ottobock.com.cn

Otto Bock Asia Pacific Ltd.
Unit 1004, 10/F, Greenfield Tower, Concordia Plaza
1 Science Museum Road, Tsim Sha Tsui
Kowloon, Hong Kong · China
T +852 2598 9772 · F +852 2598 7886
info@ottobock.com.hk · www.ottobock.com

Otto Bock HealthCare India Pvt. Ltd.
20th Floor, Express Towers
Nariman Point, Mumbai 400 021 · India
T +91 22 2274 5500 / 5501 / 5502
information@indiaottobock.com · www.ottobock.in

Otto Bock Japan K. K.
Yokogawa Building 8F, 4-4-44 Shibaura
Minato-ku, Tokyo, 108-0023 · Japan
T +81 3 3798-2111 · F +81 3 3798-2112
ottobock@ottobock.co.jp · www.ottobock.co.jp

Otto Bock Korea HealthCare Inc.
4F Agaworld Building · 1357-74, Seocho-dong
Seocho-ku, 137-070 Seoul · Korea
T +82 2 577-3831 · F +82 2 577-3828
info@ottobockkorea.com · www.ottobockkorea.com

Otto Bock South East Asia Co., Ltd.
1741 Paholyothin Road
Kwaeng Chatuchark · Khet Chatuchark
Bangkok 10900 · Thailand
T +66 2 930 3030 · F +66 2 930 3311
obsea@ottobock.co.th · www.ottobock.co.th

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Otto Bock Mobility Solutions GmbH