
3R106, 3R106=KD, 3R106=ST, 3R106=HD

DE Gebrauchsanweisung	5
EN Instructions for use	13
FR Instructions d'utilisation	21
IT Istruzioni per l'uso	30
ES Instrucciones de uso	39
PT Manual de utilização	47
NL Gebruiksaanwijzing	56
SV Installationsanvisning	64
DA Installationsvejledning	73
NO Installasjonsanvisning	81
FI Asennusohje	89
PL Instrukcja użytkowania	97
HU Használati utasítás	106
CS Návod k použití	114
TR Kullanma talimatı	122
EL Οδηγίες χρήσης	130
RU Руководство по применению	139
JA 取扱説明書	148
ZH 使用说明书	156




2Z11=KIT




- Einzelteile-Pack 4D3
- Single-Component Pack 4D3
- Pack de pièces détachées 4D3
- Kit componenti 4D3
- Kit de componentes 4D3
- Conjunto de componentes 4D3
- Onderdelenpakket 4D3
- Komponentpakke 4D3
- Zestaw naprawczy 4D3
- 4D3 alkatrészcsomag
- Σετ μεμονωμένων εξαρτημάτων 4D3
- Комплект деталей 4D3
- 交換パーツセット 4D3
- 维修组件 4D3

- ▲ Mindestmenge
- Minimum order quantity
- Quantité minimum
- Quantità minima per l'ordine
- Cantidad mínima
- Quantidade mínima
- Minimikvantitet
- Mindste mængde
- Ilość minimalna
- Minimális mennyiség
- Ελάχιστη ποσότητα
- Минимальное количество
- 最低注文個数あり
- 最低起订量

- Einzelteile
- Single components
- Pièces détachées
- Componentes
- Singoli componenti
- Componentes
- Componentes
- Onderdelen
- Komponenter
- Podzespolny alkatrészek
- Μεμονωμένα εξαρτήματα
- Отдельные детали
- 個別販売
- 单个部件

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English

Last update: 2021-02-19

- Please read this document carefully before using the product and observe the safety notices.
- Instruct the user in the safe use of the product.
- Please contact the manufacturer if you have questions about the product or in case of problems.
- Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- Please keep this document in a safe place.

1 Scope of delivery (Fig. 1)

1.1 4D3 Single Components Set for 3R106 ●

Consisting of: (6) 1 extension stop, (3) 4 truss head screws and (2) 4 two-hole nuts.

1.2 Minimum order quantity ▲

- | | |
|----------------------------------|------------------------------|
| (2) 502R1=M5×16 Two-hole Nut | (4) 506G3=M8×10 Set Screw |
| (3) 501T1=M5×16 Truss Head Screw | (4a) 506G3=M8×12-V Set Screw |

1.3 Single components ■

- | | |
|----------------------------------|--|
| (1) 4G70 Lamination Anchor | (7) 513D83=1.4×12.6×66 Compression Spring (strong) |
| (5) 2R49 Tube Adapter | (7a) 513D83=1.1×12.9×66 Compression Spring (weak) |
| (8) 710H10=2×3 Adjustment Wrench | |


2 Description

2.1 Indications for use

The 3R106 Polycentric Knee Joint is to be **exclusively** used for the prosthetic fitting of lower limb amputees.

2.2 Field of application

Field of application according to the Ottobock MOBIS Mobility System:

Recommended for amputees with **mobility grade 2 and 3** (restricted and unrestricted outdoor walkers).

Approved for a patient weight of up to 100 kg / 220 lbs.

CAUTION

Use under unallowable environmental conditions

Fall due to damaged product

- Do not expose the product to unallowable environmental conditions (see the table „Unallowable environmental conditions“ in this section).
- If the product was exposed to unallowable environmental conditions, take suitable steps (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop etc.).

2.3 Environmental conditions

Allowable environmental conditions

Temperature range for use: -10 °C to +45 °C

Allowable relative humidity 0 % to 96 %, non-condensing

Unallowable environmental conditions

Mechanical vibrations or impacts

Perspiration, urine, fresh water, salt water, acids

Dust, sand, highly hygroscopic particles (e. g. talcum)

2.4 Reuse and lifetime

CAUTION

Reuse on another patient

Fall due to loss of functionality as well as damage to the product

- Only use the product for a single patient.

CAUTION

Exceeding the lifetime

Fall due to change in or loss of functionality and damage to the product

- Ensure that the maximum lifetime defined in this section is not exceeded.

This product was tested by the manufacturer with 3 million load cycles. Depending on the user's activity level, this corresponds to a maximum lifetime of 5 years.

2.5 Construction and function (Fig. 2 + 3)

During the initial design of these knee joints the highest priority was given to the large flexion angle and the new kind of efficient swing phase control. The upper joint section (1) and lower joint section (2) are connected by the anterior links (3) and the posterior link (4) to form a four-bar linkage. In extension, the instantaneous center of rotation is located above the joint and behind

the weight bearing line, thus providing knee stability in the stance phase. To control swing phase, the pneumatic unit is linked to the upper joint section with a connecting rod. The resistances for flexion and extension are independently and individually adjustable. Reversing action of the extension assist occurs at approx. 110°. The distal portion of the pneumatic cylinder is fabricated as a tube clamp and accommodates the provided tube adapter.

Another advantage is the very large knee flexion angle, which goes far beyond the usual range. The joint does not feature a flexion stop bumper. A flexion angle of more than 170° is technically possible, and is limited only by the corresponding socket connectors, socket shapes or by the cosmetic foam cover (Fig. 3).

2.6 Different connection systems

The 3R106 Knee Joint is available in 4 different versions. They differ from each other solely by the type of proximal connection used (see front page):

Order number	Version
3R106	Knee Joint with pyramid connector
3R106=KD	Knee Joint – Version with lamination anchor for knee disarticulation amputees
3R106=ST	Knee Joint – Version with threaded connector for long residual limbs
3R106=HD	Knee Joint – Version with pyramid connector anteriorly angled by 10° for hip disarticulation amputees

Attention!

It is not allowed to loosen, shift, or replace the connection adapter!

For the prosthetic fitting of hip disarticulations, it is of utmost importance that the specially adapted knee joint version 3R106=HD is used.

3 Technical data

Order number	3R106	3R106=KD	3R106=ST	3R106=HD
Proximal connection	Pyramid adapter	Lamination anchor	Thread (M26x1.5)	Pyramid adapter
Distal connection	Tube clamp Ø 30 mm			
Knee flexion angle	170°			
Weight	760 g	755 g	765 g	790 g
System height	162 mm	184 mm	180 mm	164 mm
Proximal system height up to alignment reference point	-6 mm	16 mm	12 mm	-4 mm
Distal system height up to alignment reference point	168 mm			
Max. user weight	100 kg / 220 lbs			
Mobility Grade	2 + 3			

4 Handling

Use 2Z11=KIT

INFORMATION: The connection area of the prosthetic knee joint can be protected against scratches with the protective film of the 2Z11=KIT during alignment in the workshop and testing in the trial fitting area.

- Use the protective film as illustrated in the accompanying document for the 2Z11=KIT.
- Remove the protective film before the patient leaves the fitting area.

4.1 Alignment

The three-dimensional arrangement of the prosthetic socket and the modular components affects the static and dynamic functions of the prosthesis.

The advantages of the 3R160 Knee Joint can only be made optimal use of in case of a correct alignment.

The optimal residual limb position must be anticipated when positioning the socket connector. Plumb lines in the frontal and sagittal planes (drawn from the hip joint's center of rotation and marked during plaster cast taking and trial fitting of the test socket) will facilitate correct positioning of the lamination anchor or socket adapter.

To align the prosthesis please proceed in two steps:

1. First make the bench alignment using an alignment tool (such as 743L200 L.A.S.A.R. Assembly or PRO.S.A. Assembly 743A200).
2. Then use 743L100 L.A.S.A.R. Posture for static alignment optimization.

4.1.1 Bench alignment with alignment tool (the following steps refer to Fig. 7)

- ① Position the middle of the foot **30 mm anterior** to the alignment reference line.
- ② Add 5 mm to the required heel height of the foot. Set correct outward rotation of the foot.
- ③ Clamp the knee joint using the appropriate adapter inserts. For bench alignment, **the alignment reference line should run through the upper anterior axis** (alignment reference point). At that point the pyramid base should be horizontal. Pay attention to the knee ground distance and outward rotation of the knee (adapter inserts provide for a rotation of approx. 5°). Recommended positioning of the alignment reference point: 20 mm above the medial tibial plateau.
- ④ Connect the foot to the modular knee joint using the enclosed tube adapter.

Attention!

1. **Do not clamp the tube into a vise!** Cut the tube using 719R2 Tube Cutter placed at a right angle to the tube.
2. Afterwards carefully debur the inside and outside of the tube using the 718R1 Interior and Exterior Pipe Reamer.
3. It is **indispensable** to insert the tube until it is fully seated inside the integrated tube clamp.
4. Tighten the cap screw using the 710D1 Torque Wrench:
Torque: 10 Nm.
5. **For trial fitting:** Tighten the adjacent 506G3 Set Screws of the tube adapter to a torque of **10 Nm** using the 710D4 Torque Wrench. To finish the prosthesis: Apply 636K13 Loctite to the set screws, pre-tighten the adjacent set screws to a torque of **10 Nm**, and then tighten all four set screws to a torque of **15 Nm**.

For reasons of structural strength, the residual pieces of cut tube adapters must not be used as components in prostheses.

- 5 Mark the center of the socket proximally and distally on the lateral side. Draw a line through both marks from socket brim to the distal end of the socket.
- 6 Now position the socket such that the alignment reference line passes through the proximal center mark of the socket. Set the socket flexion to somewhere between 3° and 5°; however, the individual situation (e.g. hip joint contractures) must be taken into account and, if necessary, more flexion should be provided. Also pay attention to the ischial tuberosity to ground distance.

Attention!

If the residual limb flexion is not taken into account, the joint will be positioned too far to the front. This will lead to malfunction and premature wear. If required, use the 4R118 Adapter Plate for optimal positioning of the prosthetic knee joint.

- 7 Connect the socket and modular knee joint using a corresponding adapter (e.g. 4R111, 4R41, 4R55, 4R51 Socket Adapter).

4.1.2 Static alignment optimization using 743L100 L.A.S.A.R. Posture

(the following steps refer to Fig. 8)

The bench alignment can be substantially improved using L.A.S.A.R. Posture. In order to ensure appropriate stability combined with easy swing phase initiation, please proceed as follows:

- 1 To make the load line visible, the transfemoral amputee stands on the L.A.S.A.R. Posture with the prosthetic side on the force plate and with the other leg on the height compensation panel. The prosthesis side should be sufficiently loaded (> 35 % of the body weight).
- 2 Now adapt the alignment by only **adjusting the plantar flexion of the foot**. The **load line (laser line) should be approx. 35 mm (1.38 in) anterior to the upper front knee axis (see Fig. 8)**.
- 3 After step 2, dynamic optimization can take place between parallel bars.

4.2 Combination possibilities

INFORMATION

In a prosthesis, all prosthetic components have to meet the patient's requirements regarding the amputation level, body weight, activity level, environmental conditions and field of application.

4.3 Adjustments and final assembly

4.3.1 Socket fabrication 3R106=KD

4.3.1.1 Lamination prior to trial fitting

Take a piece of 623T3 Perlon Stockinette twice the length of the plaster model and pull half of it over the isolated socket. Twist the other half of the stockinette on the distal end and reflect it over the socket.

In order to absorb strong forces in the area of the modular knee joint, the socket must be gradually reinforced with 616G13 Fiberglass Stockinette. To do this, pull over the first layer to $\frac{2}{3}$ of the socket length, tie it off and then reflect it up to half of the socket length.

Put 2 layers of 616G12 Carbon-Fiber Cloth on the distal area so as to provide an underlying base with a projecting brim of 3 cm around the 4G70 Adapter, which will be mounted later in the procedure.

Pull on two layers of 623T3 Perlon Stockinette.

The socket lamination is performed in two stages. In the first stage the socket is laminated up to $\frac{2}{3}$ of its length with 617H19 Orthocryl Lamination Resin. In the second stage the proximal part of the socket is laminated with 617H17 Orthocryl Flexible Resin.

Once the first casting has hardened, pull on 2 more layers of 623T3 Perlon Stockinette before laminating the proximal part of the socket with 617H17 Orthocryl Flexible Resin.

4.3.1.2 Applying the lamination anchor

Prior to gluing, prepare the correctly adapted lamination anchor with the residual limb bedding. To do so, fill the coupling opening with 636K8 Plastaband (Fig. 3). Clamp residual limb bedding and joint into the alignment apparatus.

Attention!

Remove soft, flexible or porous materials from beneath the support surface of the lamination anchor. **Mix the putty exclusively from 617H21 Orthocryl Sealing Resin and 639A1 Talcum Powder. Then glue on the lamination anchor.**

Secure with 627B2 Adhesive Tape for trial fitting. Test the stop function. If necessary, form the required contact surface with putty. If need be, glue an extension stop made of Pedilin onto the outer laminate.

4.3.1.3 Socket completion

Once the trial fitting has been completed, fasten the lamination anchor to the socket with truss head screws (Fig. 1, (3)) and two-hole nuts (Fig. 1, (2)). Then laminate. Proceed with further reinforcement as follows: Two layers of 616G12 Carbon Fiber Cloth must be sandwiched between two layers of 623T3 Perlon Stockinette. To achieve this, pull one layer of 623T3 Perlon Stockinette over the whole socket and tie it circularly at the top. Apply two layers of 616G12 Carbon Fiber Cloth over the arms of the 4G70 Adapter. Pull the remaining half of the 623T3 Perlon Stockinette back over the carbon fiber layers. Gradually reinforce the socket again with 616G13 Fiberglass Stockinette (as described in point 4.3.1.1). As a final step, pull over 2 layers of Perlon Stockinette 623T3. Proceed with lamination as in the first casting.

Attention

Deviating from the processing instructions and the recommended materials for the lamination of the anchor can lead to a loosening or breaking of the adapter.

Attention!

After assembly, tighten the set screws (Fig. 1, (3) + (3a)) with the 710D1 Torque Wrench. For trial fitting: Tighten the 506G3 Set Screws to **10 Nm**. For completion: Secure the 506G3 Set Screws with 636K13 Loctite and pre-tighten to **10 Nm**. Then complete assembly by tightening to **15 Nm**.

4.3.2 Adjusting the knee joint during trial walking

As with other polycentric joints, swing phase flexion occurs when the weight is on the ball of the foot.

Alignment, foot design, condition of the residual limb, and the amputee's activity level all influence knee function. The pendulum effect from the weight of the foot and shin length are additional variables. Initial dynamic alignment should begin from the base setting and the function of the knee should be explained to the patient.

The knee is factory adjusted to a basic position for initial dynamic alignment.

Before changing the factory settings, the following guidelines must be observed:

4.3.2.1 Stance phase stability by positioning the joint

Polycentric knee joints are stable at heel strike.

Initiation of swing flexion depends on the position of the joint and thus of the location of the instantaneous center of rotation.

If stance phase stability is **too great**, and swing flexion can only be achieved with a great deal of effort, the instantaneous center of rotation is placed too far posterior and the joint must be tipped forward using the adjustment screws.

If stance phase stability is **insufficient**, please reverse the above procedure.

With the 3R106=KD and 3R106=ST, stance phase stability and initiation of flexion depend on the position of the lamination anchor (See alignment). Later adjustments on the upper (proximal) joint connector are not possible.

4.3.2.2 Adjusting the swing phase

The pneumatic swing phase control ensures a smooth gait. Movement resistances prevent the prosthetic lower leg from swinging too far when flexing and guarantee dampened extension. Along with the patient's habits, the foot weight and the length of the lower leg, which have a pendular weight also have an influence on his or her gait.

Before changing the knee joint's factory settings, read the following instructions:

Upon delivery, the valve screw (Figure 4, Item E) for extension damping is turned out as far as possible (set to the lowest resistance) and the valve screw (Figure 4, Item F) for flexion damping is closed by 2 turns.

Adjusting flexion (F) on the left adjustment screw (from the posterior direction of view) (Fig. 4)

Valve screw of the pneumatic unit **to the right (+)** = greater resistance = less heel rise

Valve screw of the pneumatic unit **to the left (-)** = less resistance = more heel rise

Adjusting extension (E) on the right adjustment screw (from the posterior direction of view)

Valve screw of the pneumatic unit **to the right (+)** = greater resistance = more extension damping

Valve screw of the pneumatic unit **to the left (-)** = less resistance = less extension damping

Attention!

When making the final extension damping adjustment, make sure that the leg still completely extends with each step.

The integrated spring extension assist has no knee stabilizing effect. When fine-tuning the swing phase and the gait pattern, it must always be ensured that the extension stop is fully reached with each step.

Attention – Please inform your patients!

The knee joint function may be affected by both environmental and service conditions. To reduce risk for the patient, the knee joint should not be used when noticeable functional changes occur. Noticeable functional changes may include poor response, incomplete extension, diminished swing phase control or stance phase stability, abnormal noises, etc.

Action to be taken

Have the prosthesis examined at a specialist workshop.

When using the knee joint do not reach into the joint mechanism since there is the risk of injury.

4.4. Replacing the extension assist spring

The knee joint is equipped with an extension assist spring that promotes the extension of the prosthetic lower leg in the swing phase.

If need be, the scope of this assistance can be considerably **reduced** by replacing the extension assist spring.

To replace the spring, undo the end cap at the lower end of the joint (Fig. 5) and remove the heavy-duty compression spring (Fig. 1, (7)). Then insert the light-duty compression spring (Fig. 1, (7a)) (included in delivery). Finally, fully retighten the end cap.

For assembly sequence and arrangement of the components please refer to figure 6.

Attention!

After switching to the light-duty compression spring, ensure that the leg still extends fully at each step. The pneumatic swing phase control will have to be re-adjusted during the gait test (see section 5.2).

4.5 Cosmetic foam cover

The following foam covers are to be used:

Knee Joint	Foam Cover
3R106	3S107
3R106=HD	3S27=L / R 44
3R106=ST	3S107
3R106=KD	6R6 (3S107)

Attention!

The cosmetic foam cover can affect knee joint functions. For this reason, re-adjustments might be necessary.

Attention!

Do not use talcum powder to eliminate noise in the cosmetic foam cover. Talcum powder reduces the lubrication of the mechanical parts, which may lead to a malfunction and thus increase the risk of failure. Using this medical product after application of talcum powder will render all claims null and void.

Note:

To eliminate noise, please apply Silicone Spray (519L5) directly on the contact surfaces of the cosmetic foam cover.

4.6 Maintenance instructions

Attention!

The knee joint must not be lubricated or greased, since this would result in a risk of bearing damage and loss of functionality.

Do not use any aggressive cleaning agents since they could cause damage to bearings, seals and plastic parts.

Do not disassemble the knee joint! If you have a problem, please send the knee joint in for service.

Ottobock recommends readjusting the knee joint's settings once the patient has gotten used to the prosthesis. This period of time varies depending on individual patient characteristics.

Please check the knee joint for wear and proper function at least once a year and make adjustments if necessary. Special attention should be paid to the flexion and extension resistances, bearings, and abnormal noises. Full flexion and extension must be guaranteed.

5 Disposal

In some jurisdictions it is not permissible to dispose of the product with unsorted household waste. Improper disposal can be harmful to health and the environment. Observe the information provided by the responsible authorities in your country regarding return, collection and disposal procedures.

6 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

6.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

6.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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