8E51

DE Gebrauchsanweisung (Fachpersonal)	3
EN Instructions for use (qualified personnel)	12
FR Instructions d'utilisation (Personnel spécialisé)	22
IT Istruzioni per l'uso (Personale tecnico specializzato)	32
ES Instrucciones de uso (Personal técnico especializado)	42
PT Manual de utilização (Pessoal técnico)	52
NL Gebruiksaanwijzing (Vakmensen)	62
SV Bruksanvisning (Fackpersonal)	72
FI Käyttöohje (Ammattihenkilöstö)	81
SK Návod na používanie (Odborný personál)	91
EL Οδηγίες χρήσης (Τεχνικό προσωπικό)	100
RU Руководство по применению (Квалифицированный персонал)	111
JA 取扱説明書 (有資格担当者)	122

INFORMATION

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INFORMATION

Date of last update: 2021-02-22

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.

- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

The product "Electric Hand 2000" is referred to simply as the product/prosthesis/terminal device below.

These instructions for use provide you with important information on the use, adaptation and handling of the product.

Only put the product into use in accordance with the information contained in the accompanying documents supplied.

2 Product description

2.1 Function

The product is a myoelectrically controlled terminal device of the MyoBock prosthesis system.

The thumb and fingers rotate around the same axis on the terminal device and carry out a spatial movement in the process. This makes it possible to grasp objects with natural compensating movements of the forearm and upper arm.

Size 5 is driven by a single motor.

Sizes 5½, 6 and 6½ are driven by two motors with a multi-stage overriding gear mechanism. The first motor opens and closes the hand at high speed and low force. The second motor is activated in order to build up the gripping force required to grasp objects.

2.2 Combination possibilities

This product can be combined with the following Ottobock components:

Use of the 9E420* 7 in 1 Controller

- 9E420=L or 9E420=R 7 in 1 Controller
- 757B35=0 or 757B35=1 or 757B35=3 MyoEnergy Integral
- 757L35 MyoCharge Integral

Wrist joints

- 10V51=1 or 10V51=2 MyolinoWrist 2000

Elbow components

- 12K19=40 elbow component
- 12K12 MovolinoArm Friction

Lamination rings

- 10S16=34 or 10S16=38 lamination ring
- 10S18=40 lamination ring
- 11D12=34 lamination ring

3 Intended use

3.1 Indications for use

The product is to be used exclusively for exoprosthesis fittings of the upper limbs.

3.2 Conditions of use

The product was developed for everyday use and must not be used for unusual activities. These unusual activities include, for example, sports which involve excessive strain and/or shocks to the wrist joint (pushups, downhill racing, mountain biking, ...) or extreme sports (free climbing, paragliding, etc.). Furthermore, the product should not be used to operate motor vehicles, heavy equipment (e.g. construction machines), industrial machines or motor-driven equipment.

The product is intended **exclusively** for use on **one** patient. Use of the product by another person is not approved by the manufacturer.

Permissible ambient conditions are described in the technical data (see page 20).

3.3 Indications

The product can be used for unilateral or bilateral amputations, and in case of dysmelia.

3.4 Contraindications

- All conditions which contradict or go beyond the specifications listed in the section on "Safety" and "Indications for use".

3.5 Qualification

The fitting of a patient with the product may only be carried out by O&P professionals who have been authorised with the corresponding Ottobock training.

4 Safety

4.1 Explanation of warning symbols

 WARNING	Warning regarding possible serious risks of accident or injury.
--	---

 CAUTION	Warning regarding possible risks of accident or injury.
--	---

 NOTICE	Warning regarding possible technical damage.
---	--

4.2 Structure of the safety instructions

 WARNING
--

The heading describes the source and/or the type of hazard

The introduction describes the consequences in case of failure to observe the safety instructions. Consequences are presented as follows if more than one consequence is possible:

- > E.g.: Consequence 1 in the event of failure to observe the hazard
- > E.g.: Consequence 2 in the event of failure to observe the hazard
- ▶ This symbol identifies activities/actions that must be observed/carried out in order to avert the hazard.

4.3 General safety instructions

 WARNING
--

Non-observance of safety notices

Personal injury/damage to the product due to using the product in certain situations.

- ▶ Observe the safety notices and the stated precautions in this accompanying document.

 WARNING
--

Operating a vehicle with the prosthesis

Accidents due to unexpected actions of the prosthesis.

- ▶ The prosthesis should not be used for the operation of motor vehicles and heavy equipment (e.g. construction machines).

 WARNING
--

Using the prosthesis while operating machines

Injury due to unexpected actions of the prosthesis.

- ▶ The prosthesis must not be used to operate industrial machines or motor-driven equipment.

⚠ WARNING

Operating the product near active implanted systems

Interference with active implantable systems (e.g., pacemaker, defibrillator, etc.) due to electromagnetic interference of the product.

- ▶ When operating the product in the immediate vicinity of active implantable systems, ensure that the minimum distances stipulated by the manufacturer of the implant are observed.
- ▶ Make sure to observe any operating conditions and safety instructions stipulated by the manufacturer of the implant.

⚠ WARNING

Use of damaged power supply unit, adapter plug or battery charger

Risk of electric shock due to contact with exposed, live components.

- ▶ Do not open the power supply unit, adapter plug or battery charger.
- ▶ Do not expose the power supply unit, adapter plug or battery charger to extreme loading conditions.
- ▶ Immediately replace damaged power supply units, adapter plugs or battery chargers.

⚠ WARNING

Charging the prosthesis without taking it off

Risk of electric shock due to defects in the power supply unit or in the battery charger.

- ▶ For safety reasons, remove the prosthesis prior to charging it.

⚠ CAUTION

Independent manipulation of the product

Injury due to product damage or malfunction.

- ▶ Manipulations to the product other than the tasks described in these instructions for use are not permitted.
- ▶ The product and any damaged components may only be opened and repaired by authorised, qualified Ottobock personnel.

⚠ CAUTION

Signs of wear to system components

Injury due to faulty control or malfunction of the product.

- ▶ In the interest of the patient's safety and in order to maintain operating reliability, the product should be serviced at regular intervals.

4.4 Information on Alignment/Adjustment

⚠ CAUTION

Operator errors during the adjustment process with the adjustment software

Injury due to unexpected product behaviour.

- ▶ Participation in an Ottobock product training course is mandatory prior to using the product. During this product training course you will receive a password giving you access to the adjustment software. Additional product training courses may become necessary to qualify for software updates.
- ▶ Do not share your unlock PIN.
- ▶ Use the online help which is integrated into the software.

⚠ CAUTION

Insufficient skin contact of the electrodes

Injury due to faulty control or malfunction of the product.

- ▶ The electrodes are to be placed on intact skin only and with as much electrode-skin contact as possible.
- ▶ In the case of strong interference from electronic devices, the position of the electrodes should be checked and changed if necessary.
- ▶ If the interference cannot be eliminated or if you do not achieve the expected results by adjustment or selection of the appropriate control programme, please contact the Ottobock branch responsible for your country.

⚠ CAUTION

Incorrect electrode settings

Injury due to unexpected product behaviour.

- ▶ Set the electrode sensitivity as low as possible in order to reduce interference from powerful electromagnetic radiation (e.g. visible or concealed theft prevention systems at the entrance/exit of stores), metal detectors/body scanners for persons (e.g. in airports) or other sources of strong electromagnetic interference (e.g. high-voltage lines, transmitters, transformer stations, computer tomographs, magnetic resonance tomographs, etc.).

⚠ CAUTION

Incorrect electrode settings due to muscle fatigue

Injury due to faulty control or malfunction of the product.

- ▶ Allow the patient to rest during the adjustment of the electrodes.

⚠ CAUTION

Establishing a connection between the skin and metal components of the prosthesis with the use of carbon fibre

Injury due to unexpected product behaviour as a result of interference conducted by carbon fibres.

- ▶ During fabrication, make sure that no connection is formed between the skin and metal prosthesis components by carbon fibres.

INFORMATION

Do not use silicone spray when putting on the prosthetic glove. This impairs the secure fit of the glove and can limit the functionality of the product (hand).

Observe the instructions for use provided for the prosthetic glove and the process for putting it on and taking it off described in the same.

4.5 Information on Proximity to Certain Areas

⚠ CAUTION

Proximity to sources of strong magnetic or electrical interference (e.g. theft prevention systems, metal detectors)

Injury due to unexpected behaviour of the product caused by interference with internal data communication.

- ▶ Avoid remaining in the vicinity of visible or concealed theft prevention systems at the entrance/exit of stores, metal detectors/body scanners for persons (e.g. in airports) or other sources of strong magnetic and electrical interference (e.g. high-voltage lines, transmitters, transformer stations, computer tomographs, magnetic resonance tomographs, etc.).
- ▶ When walking through theft prevention systems, body scanners or metal detectors, watch for unexpected behaviour of the product.

⚠ CAUTION

Insufficient distance to HF communication devices (e.g. mobile phones, Bluetooth devices, WiFi devices)

Injury due to unexpected behaviour of the product caused by interference with internal data communication.

- ▶ Therefore, keeping a minimum distance of 30 cm from HF communication devices is recommended.

⚠ CAUTION

Remaining in areas outside the allowable temperature range

Injury due to faulty control or malfunction of the product.

- ▶ Avoid remaining in areas outside the allowable temperature range (see page 20).

4.6 Information on Use

⚠ CAUTION

Using the prosthesis with pointed or sharp objects (e.g. knives in the kitchen)

Injury due to unintentional movements.

- ▶ Use extra caution when using the prosthesis for handling pointed or sharp objects.

⚠ CAUTION

Overloading due to unusual activities

Injury due to unexpected product behaviour as a result of a malfunction.

- ▶ The product was developed for everyday use and must not be used for unusual activities. These unusual activities include, for example, sports with excessive strain and/or shocks to the wrist joint (pushups, downhill mountain biking, ...) or extreme sports (free climbing, paragliding, etc.).
- ▶ Careful handling of the product and its components not only increases their service life but, above all, ensures the patient's personal safety!
- ▶ If the product and its components have been subjected to extreme loads (e.g. due to a fall, etc.), then the product must be inspected for damage immediately. If necessary, forward the product to an authorised Ottobock Service Centre.

⚠ CAUTION

Mechanical stress on the product

Injury due to faulty control or malfunction of the product components.

- ▶ Do not subject the product to mechanical vibrations or impacts.
- ▶ Check the product for visible damage before each use.

⚠ CAUTION

Penetration of dirt and humidity into the prosthetic components

Injury due to faulty control or malfunction of the prosthetic components.

- ▶ Ensure that neither solid particles nor liquids can penetrate into the product or prosthetic components (e.g. the terminal device).

NOTICE

Improper product care

Damage to the product due to the use of incorrect cleaning agents.

- ▶ Only clean the product with a damp cloth and mild soap (e.g. 453H10=1-N Ottobock DermaClean).

4.7 Information on the Power Supply/Battery Charging

⚠ CAUTION

Charging the product with soiled or damaged contacts

Injury due to unexpected behaviour of the product caused by insufficient charging.

- ▶ Make sure that the contacts are always clean and free of grease.
- ▶ Clean the electrical contacts of the charging plug and charging receptacle regularly using cotton swabs and a mild soap solution.
- ▶ Take care to avoid damaging the contacts with pointed or sharp objects.

NOTICE

Use of incorrect power supply unit/battery charger

Damage to product due to incorrect voltage, current or polarity.

- ▶ Use only power supply units/battery chargers approved for this product by Ottobock (see instructions for use and catalogues).

NOTICE

Contact of the charging plug with magnetic data carriers

Wiping of the data carrier.

- ▶ Do not place the charging plug on credit cards, diskettes, audio or video cassettes.

INFORMATION

The product cannot be used as long as the charging plug is connected to the charging receptacle. The product is deactivated for the duration of charging.

4.8 Notes on using an elbow joint (e.g. 12K19=40, 12K12)

⚠ CAUTION

Risk of pinching in elbow joint flexion region

Injury due to pinching of body parts.

- ▶ Ensure that fingers and other body parts are not in this area when bending the elbow joint.

⚠ CAUTION

Manual unlocking of elbow lock under load

Injury by release of elbow lock under load.

- ▶ Particular caution should be exercised when unlocking the elbow lock while lifting heavy loads.
- ▶ Be careful when unlocking the lock under such conditions due to the possibility of injury.

⚠ CAUTION

Penetration of dirt and moisture into the prosthetic arm

Injury due to unexpected behaviour of the prosthetic components or due to malfunction.

- ▶ Ensure that neither solid particles nor liquids can penetrate into the prosthetic arm.
- ▶ Do not expose the prosthetic arm, especially the elbow joint, to dripping or splashing water.
- ▶ As a minimum, wear the prosthetic arm and especially the elbow joint under suitably resistant clothing in the rain.

NOTICE

Coating, gluing or painting the prosthesis

Damage or fracture due to chemical processes.

- ▶ The prosthesis must not be coated, glued or painted.

5 Scope of Delivery and Accessories

Scope of Delivery

- 1 pc. Electric Hand 2000
- 1 pc. Instructions for use (qualified personnel)

Required accessories for operation

- 8S20=* Prosthetic glove

6 Preparing the product for use

6.1 Settings

The control settings depend on the controller being used. For more information, see the instructions for use of the controller in question.

6.2 Assembly

6.2.1 Installing O-rings

- 1) Install the outer 623F3=26x1.0 O-ring (1).
- 2) Insert the number of inner O-rings (2) according to the desired friction level (see fig. 1).

6.2.2 Installing controllers

For information on installing the controller, see the instructions for use of the controller being used.

6.2.3 Installing the hand

- 1) Put on the hand (see fig. 2).
- 2) Firmly tighten the wrist nut against the O-ring.

INFORMATION

For daily use, the terminal device has to be worn with a prosthetic glove. The glove protects the mechanism against environmental influences such as moisture, dirt and dust.

For use (putting on and taking off) and care of the prosthetic glove, see the instructions for use included with the prosthetic glove.

7 Cleaning

- 1) Switch off the product before cleaning.
- 2) Clean the product with a damp cloth and mild soap when needed.
Make sure that no liquids get into the product and product components.
- 3) Dry the product with a lint-free cloth and allow it to air dry fully.

8 Maintenance

Performing regular maintenance (service inspections) every 24 months is recommended to prevent injuries and maintain the quality of the product.

In general, all products are subject to compliance with the maintenance intervals during the warranty period. This is the only way to maintain full warranty cover.

Additional services such as repairs may be provided in the course of maintenance. These additional services may be provided free of charge or can be billable according to an advance cost estimate, depending on the extent and validity of the warranty.

The following components must always be sent in for maintenance and repairs:

The product, battery charger and power supply. The shipping container for the loaner unit you receive must be reused for sending back the components requiring inspection.

9 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

9.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

9.2 Trademarks

All product names mentioned in this document are subject without restriction to the respective applicable trademark laws and are the property of the respective owners.

All brands, trade names or company names may be registered trademarks and are the property of the respective owners.

Should trademarks used in this document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

9.3 CE conformity

Otto Bock Healthcare Products GmbH hereby declares that the product is in compliance with applicable European requirements for medical devices.

The product meets the requirements of the RoHS Directive 2011/65/EU on the restriction of the use of certain hazardous substances in electrical and electronic devices.

The full text of the regulations and requirements is available at the following Internet address: <http://www.ottobock.com/conformity>

10 Technical data

Ambient conditions	
Storage (with and without packaging)	+5 °C/+41 °F to +40 °C/+104 °F Max. 85% relative humidity, non-condensing
Transport (with and without packaging)	-20 °C/-4 °F to +60 °C/+140 °F Max. 90% relative humidity, non-condensing
Operation	-5 °C/+23 °F to +45 °C/+113 °F Max. 95% relative humidity, non-condensing

Reference number	8E51			
Size	5	5½	6	6½
Operating voltage [V]	7.4			
Average current consumption [mA]	Approx. 200			

Reference number	8E51			
Electric shut-off "closed" at [mA]	Approx. 400			
Electric shut-off "open" at [mA]	Approx. 400	Approx. 250		
Opening width [mm]	28	37	52	58
Minimum possible gripping force [N]	8	25	35	
Weight without prosthetic glove [g]	86	115	125	130
Expected lifetime given compliance with the recommended maintenance intervals [years]	5			

11 Appendices

11.1 Symbols Used

In some jurisdictions it is not permissible to dispose of these products with unsorted household waste. Disposal that is not in accordance with the regulations of your country may have a detrimental impact on health and the environment. Please observe the instructions of your national authority pertaining to return and collection.

Serial number (YYYY WW NNN)
 YYYY – year of manufacture
 WW – week of manufacture
 NNN – sequential number

Medical device

Manufacturer

Declaration of conformity according to the applicable European directives

11.2 Troubleshooting

Event	Cause/required action
Lack of gripping force	Insufficient contact between hand and switch unit. <ul style="list-style-type: none"> Firmly screw the wrist nut onto the outer O-ring.
Excessive friction in pronation and supination direction	Excessive frictional forces on the O-rings. <ul style="list-style-type: none"> Take out O-ring(s).
Insufficient friction in pronation and supination direction	Insufficient frictional forces on the O-rings, or no O-rings installed. <ul style="list-style-type: none"> Insert O-ring(s).
Hand operates too slowly in open direction	Second motor in operation because the prosthetic glove is too stiff.

Event	Cause/required action
	<ul style="list-style-type: none"> • Warm up the prosthetic glove (max. 60°C/140°F), then open the hand wide. Allow to cool in open position. Battery charge level low <ul style="list-style-type: none"> • Charge battery
Hand operates too slowly in close direction	Second motor in operation because the prosthetic glove is too stiff. <ul style="list-style-type: none"> • Close hand and warm prosthetic glove from the outside (max. 60°C/140°F). Battery charge level low <ul style="list-style-type: none"> • Charge battery
Finger cover damaged	Cut open and pull off finger cover. Install new finger cover.

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