

AxonEnergy Integral 757B501

DE Gebrauchsanweisung (Fachpersonal)	3
EN Instructions for use (qualified personnel)	11
FR Instructions d'utilisation (Personnel spécialisé)	18
IT Istruzioni per l'uso (Personale tecnico specializzato)	27
ES Instrucciones de uso (Personal técnico especializado)	35
PT Manual de utilização (Pessoal técnico)	43
NL Gebruiksaanwijzing (Vakmensen)	51
SV Bruksanvisning (Fackpersonal)	59
DA Brugsanvisning (Faguddannet personale)	67
NO Bruksanvisning (Fagpersonell)	75
FI Käyttöohje (Ammattihenkilöstö)	83
CS Návod k použití (Odborný personál)	90
RU Руководство по применению (Квалифицированный персонал)	98
JA 取扱説明書 (有資格担当者)	107

INFORMATION

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INFORMATION

Date of last update: 2020-12-14

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

These instructions for use provide you with important information on the use, adaptation and handling of the product.

2 Product description

2.1 Function

The product serves to provide power to a prosthesis.

The product consists of several components. These components are inseparable and ready for mounting.

2.2 Components

The product consists of the following components:

Rechargeable battery

The battery (see fig. 1, item 1) consists of two cells with integrated electronics. The electronics protect the battery against short circuits, overvoltage, undervoltage and charging outside the allowable temperature range.

Charging receptacle

The charging receptacle (see fig. 1, item 2) has the following functions:

- Contact for charging the battery
- LED to indicate the current charge level
- LED for providing feedback on operating states
- Button for turning the prosthesis on and off, checking the charge level and emergency opening of the prosthesis
- Beeper for providing feedback on operating states
- Button to activate the Bluetooth function

Communication cable

The communication cable (see fig. 1, item 3) is used to exchange data and connects the respective prosthetic components to the product.

Connection cable

The connection cable (see fig. 1, item 4) is used to connect the 13E500 AxonMaster. When the battery is used in an elbow joint, the forearm cable of the flexion aid (AFB) is connected to this connection cable.

3 Intended use

3.1 Indications for use

The product is to be used exclusively for exoprosthetic fittings of the upper limbs.

3.2 Conditions of use

The product has been developed for use in everyday activities. It must not be used for activities where it would be exposed to unallowable environmental conditions or strains (e.g. extreme sports).

Permissible ambient conditions are described in the technical data (see page 16).

3.3 Contraindications

- All conditions which contradict or go beyond the specifications listed in the section on "Safety" and "Indications for use".

3.4 Qualification

The fitting of a patient with the product may only be carried out by O&P professionals who have been authorised with the corresponding Ottobock training.

4 Safety

4.1 Explanation of warning symbols

 WARNING	Warning regarding possible serious risks of accident or injury.
 CAUTION	Warning regarding possible risks of accident or injury.
 NOTICE	Warning regarding possible technical damage.

4.2 Structure of the safety instructions

 WARNING
The heading describes the source and/or the type of hazard
The introduction describes the consequences in case of failure to observe the safety instructions. Consequences are presented as follows if more than one consequence is possible:
> E.g.: Consequence 1 in the event of failure to observe the hazard
> E.g.: Consequence 2 in the event of failure to observe the hazard
▶ This symbol identifies activities/actions that must be observed/carried out in order to avert the hazard.

4.3 General safety instructions

 CAUTION
Non-observance of safety instructions
Personal injury/damage to the product due to using the product in certain situations.
▶ Observe the safety instructions and the stated precautions in this accompanying document.
 CAUTION
Penetration of dirt and humidity into the product
Injury due to unexpected prosthesis behaviour.
▶ Ensure that neither solid particles nor liquids can penetrate into the product.
▶ Do not use the product for bathing prostheses.

⚠ CAUTION

Changing prosthetic components when switched on

Injury due to faulty control or malfunction of the prosthetic components.

- ▶ Before changing prosthetic components (e.g. terminal device) switch the product off.

⚠ CAUTION

Overloading due to unusual activities

Injury due to unexpected product behaviour as a result of a malfunction.

- ▶ The product was developed for everyday use and must not be used for unusual activities. These unusual activities include, for example, sports with excessive strain and/or shocks to the wrist joint (pushups, downhill mountain biking, ...) or extreme sports (free climbing, paragliding, etc.).
- ▶ Careful handling of the product and its components not only increases their service life but, above all, ensures the patient's personal safety!
- ▶ If the product and its components have been subjected to extreme loads (e.g. due to a fall, etc.), then the product must be inspected for damage immediately. If necessary, forward the product to an authorised Ottobock Service Centre.

⚠ CAUTION

Penetration of perspiration into the product

Injury and chemical burns due to leakage of battery fluids.

- ▶ Wear protective clothing when removing or replacing the battery (disposable gloves and protective goggles at a minimum).
- ▶ Avoid direct skin contact!

⚠ CAUTION

Using the product with pointed or sharp objects

Destruction or damage to the product during installation.

- ▶ Do not use any pointed or sharp objects when installing the product.

4.4 Patient information

⚠ WARNING

Operating the product near active implanted systems

Interference with active implantable systems (e.g., pacemaker, defibrillator, etc.) due to electromagnetic interference of the product.

- ▶ When operating the product in the immediate vicinity of active implantable systems, ensure that the minimum distances stipulated by the manufacturer of the implant are observed.
- ▶ Make sure to observe any operating conditions and safety instructions stipulated by the manufacturer of the implant.

⚠ CAUTION

Mechanical stress on the product

Injury due to faulty control or malfunction of the product components.

- ▶ Do not subject the product to mechanical vibrations or impacts.
- ▶ Check the product for visible damage before each use.

⚠ CAUTION

Remaining, operating and charging in areas outside the allowable temperature range

Injury due to faulty control or malfunction of the product.

- ▶ Avoid remaining, operating and charging in areas outside the allowable temperature range (see page 16).

⚠ CAUTION

Proximity to sources of strong magnetic or electrical interference (e.g. theft prevention systems, metal detectors)

Injury due to unexpected behaviour of the product caused by interference with internal data communication.

- ▶ Avoid remaining in the vicinity of visible or concealed theft prevention systems at the entrance/exit of stores, metal detectors/body scanners for persons (e.g. in airports) or other sources of strong magnetic and electrical interference (e.g. high-voltage lines, transmitters, transformer stations, computer tomographs, magnetic resonance tomographs, etc.).
- ▶ When walking through theft prevention systems, body scanners or metal detectors, watch for unexpected behaviour of the product.

⚠ CAUTION

Independent manipulation of the product

Injury due to malfunction and resulting unexpected prosthesis actions.

- ▶ Manipulations to the product other than the tasks described in these instructions for use are not permitted.
- ▶ The product and any damaged components may only be opened and repaired by authorised, qualified Ottobock personnel.

⚠ CAUTION

Charging the product with soiled or damaged contacts

Injury due to unexpected behaviour of the product caused by insufficient charging.

- ▶ Make sure that the contacts are always clean and free of grease.
- ▶ Clean the electrical contacts of the charging plug and charging receptacle regularly using cotton swabs and a mild soap solution.
- ▶ Take care to avoid damaging the contacts with pointed or sharp objects.

⚠ CAUTION

Insufficient distance to HF communication devices (e.g. mobile phones, Bluetooth devices, WiFi devices)

Injury due to unexpected behaviour of the product caused by interference with internal data communication.

- ▶ Therefore, keeping a minimum distance of 30 cm from HF communication devices is recommended.

5 Scope of delivery

- 1 pc. 757B501 AxonEnergy Integral
- 2 pc. lamination dummy
- 1 pc. lamination dummy – charging receptacle
- 1 pc. drilling template for charging receptacle
- 1 pc. Component holder set
- 1 pc. Instructions for use (qualified personnel)

6 Preparing the product for use

6.1 Assembly

The components are ready for assembly on delivery and cannot be separated.

The product in combination with the component holder can only be assembled in the forearm socket. The charging receptacle must be glued into the outer socket. The battery is attached in the outer socket or between the outer and inner socket, depending on the available space. When little space is available, the battery is laminated into the outer socket.

7 Handling

7.1 Switching the terminal device on and off

NOTICE

Signs of wear and tear on the product components

Injury due to faulty control or malfunction of the product.

- ▶ Excessive pressure on the charging receptacle can lead to increased wear and tear on the button.

- ▶ 1) Push and hold the charging receptacle button until a confirmation sound is heard (at least 1 second).
→ The prosthesis and prosthetic components are switched on.
Repeat this process to switch the prosthesis and prosthetic components off.

7.2 Charging the battery

We recommend charging on a daily basis. If the device is not used for an extended period, the battery charge level should be checked regularly and the battery recharged.

- 1) Connect the charging plug to the charging receptacle of the product.
 - A correct connection between the battery charger and the product is indicated by feedback (see page 17).
 - The prosthesis is switched off automatically.
 - The charging process starts.
- 2) Disconnect the product after the charging process is complete.

7.3 Display of the current charge level

The battery charge level can be queried at any time.

- 1) With the prosthesis switched on, press the charging receptacle button and hold for less than one second.
- 2) The LED display on the charging receptacle provides information on the current battery charge level (Status signals).

7.4 Safety shutoff

The purpose of the safety shutoff is to protect the product; it is triggered in case of:

- Excessively high or low temperature in the battery's environment during the charging process
- Short circuit
- Overvoltage or undervoltage

7.5 Opening the terminal device in an emergency

This safety function allows the Axon-Bus gripping component to be opened regardless of the control signals present.

- 1) With the Axon-Bus prosthetic system switched on, press the charging receptacle button and hold for approximately three seconds until the Axon-Bus gripping component begins to open.
→ A pulsating beep sounds as it opens.
- 2) Releasing the button immediately stops the process of opening the Axon-Bus gripping component and turns the entire Axon-Bus prosthetic system off.

8 Legal information

8.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

8.2 Trademarks

All product names mentioned in this document are subject without restriction to the respective applicable trademark laws and are the property of the respective owners.

All brands, trade names or company names may be registered trademarks and are the property of the respective owners.

Should trademarks used in this document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

8.3 CE conformity

Otto Bock Healthcare Products GmbH hereby declares that the product is in compliance with applicable European requirements for medical devices.

The product meets the requirements of the RoHS Directive 2011/65/EU on the restriction of the use of certain hazardous substances in electrical and electronic devices.

This product meets the requirements of the 2014/53/EU directive.

The full text of the regulations and requirements is available at the following Internet address: <http://www.ottobock.com/conformity>

9 Technical data

Ambient conditions	
Storage (with and without packaging)	+5 °C/+41 °F to +40 °C/+104 °F Max. 85% relative humidity, non-condensing For reasons of quality (service life), it is recommended that the product be stored at low temperatures.
Transport (with and without packaging)	-20 °C/-4 °F to +60 °C/+140 °F Max. 90% relative humidity, non-condensing
Operation	-5 °C/+23 °F to +45 °C/+113 °F Max. 95% relative humidity, non-condensing
Charging the battery	+5 °C/+41 °F to +40 °C/+104 °F Max. 85% relative humidity, non-condensing

Prosthesis battery	757B501
Battery type	Li-Ion
Battery lifetime	2 years
Charging time until battery is fully charged	Approx. 3.5 h
Charging voltage	Approx. 12.3 V
Capacity	1150 mAh
Output voltage	Approx. 11.1 V
Weight	90 g
Dimensions	55 x 35 x 23 mm
Product behaviour during the charging process	The product is non-functional

10 Appendices

10.1 Symbols Used

In some jurisdictions it is not permissible to dispose of these products with unsorted household waste. Disposal that is not in accordance with the regulations of your country may have a detrimental impact on health and the environment. Please observe the instructions of your national authority pertaining to return and collection.

Declaration of conformity according to the applicable European directives

Lot number (PPPP YYYY WW)

PPPP – plant

YYYY – year of manufacture

WW – week of manufacture

Medical device

Manufacturer

10.2 Operating states/error signals

10.2.1 Status signals

Charging receptacle	Event
LED shows green light	Battery fully charged
LED shows yellow light	Battery 50% charged
LED shows orange light	Battery empty

10.2.2 Beep signals

The following beep signals are set by default:

Beep signal	Additional display	Event
1 x long	-	<ul style="list-style-type: none"> Turn off prosthesis on charging receptacle Charging begins (charging plug connected to charging receptacle) Charging ends (charging plug disconnected from charging receptacle)
2 x short	LED on charging receptacle lights up briefly	Switch on prosthesis
3 x short	-	Battery voltage too low, prosthesis shuts off automatically

10.2.3 Warnings/error signals

Charging receptacle	Beeper	Explanation/rectification
LED is flashing red	3x short	<p>Error in the Axon-Bus prosthesis solution</p> <ul style="list-style-type: none"> Turn the Axon-Bus prosthesis solution on/off Check the plug connection in the socket Connect the charger. <p>If the service icon (yellow wrench) or system error icon (red icon) is displayed on the battery charger, send the Axon-Bus prosthesis solution to an Ottobock Service Center.</p>
LED flashes 3x in the colour of the current charge level	-	<p>The connection from the product to the Axon-Master is interrupted</p> <ul style="list-style-type: none"> Check the connection (plug connection, flex cable, Easy Plus, etc.)
LED alternately flashes yellow/red	3x short	<p>Overheating of the Axon-Bus prosthesis solution</p> <ul style="list-style-type: none"> Allow the Axon-Bus prosthesis solution to cool down for a few minutes Do not use the Axon-Bus prosthesis solution at high ambient temperatures (Technical data)

10.3 Axon

The term "Axon" stands for **A**daptive **eX**change **O**f **N**europlacement data. The Axon-Bus is an Ottobock innovation for the field of exoprosthetics: a data transmission system, derived from safety-related bus systems in the aviation and automobile industries. For the user this means enhanced safety and reliability because of a considerably reduced sensitivity to electromagnetic interference in comparison with conventional systems.

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