

## AxonCharge Integral 757L500

<b>DE</b> Gebrauchsanweisung (Fachpersonal) .....	3
<b>EN</b> Instructions for use (qualified personnel) .....	10
<b>FR</b> Instructions d'utilisation (Personnel spécialisé) .....	16
<b>IT</b> Istruzioni per l'uso (Personale tecnico specializzato) .....	23
<b>ES</b> Instrucciones de uso (Personal técnico especializado) .....	30
<b>PT</b> Manual de utilização (Pessoal técnico) .....	37
<b>NL</b> Gebruiksaanwijzing (Vakmensen) .....	44
<b>SV</b> Bruksanvisning (Fackpersonal) .....	51
<b>DA</b> Brugsanvisning (Faguddannet personale) .....	58
<b>NO</b> Bruksanvisning (Fagpersonell) .....	65
<b>FI</b> Käyttöohje (Ammattihenkilöstö) .....	71
<b>CS</b> Návod k použití (Odborný personál) .....	78
<b>RU</b> Руководство по применению (Квалифицированный персонал) .....	85
<b>JA</b> 取扱説明書 (有資格担当者) .....	92

## INFORMATION

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## INFORMATION

Date of last update: 2020-12-02

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- ▶ Please keep this document for your records.

These instructions for use provide you with important information on the use, adaptation and handling of the product.

Only put the product into use in accordance with the information contained in the accompanying documents supplied.

## 2 Product description

### 2.1 Function

The AxonCharge Integral 757L500 is designed to charge the AxonEnergy Integral or 757B501 in the Ottobock Axon-Bus prosthetic system.

Charging is performed automatically after the charging plug has been connected to the charging receptacle of the Axon-Bus prosthetic system. The charging plug is secured to the charging receptacle by the integrated magnet. The special contours of the charging receptacle and charging plug ensure quick, reliable positioning of the two components. LEDs indicate the readiness of the battery charger and the current battery charge level.

## 3 Intended use

### 3.1 Indications for use

The product is to be used **exclusively** for charging batteries.

### 3.2 Conditions of use

The product was developed for everyday activities and can be used for charging unilateral or bilateral upper limb prostheses.

### 3.3 Contraindications

- All conditions which contradict or go beyond the specifications listed in the section on "Safety" and "Indications for use".

### 3.4 Qualification

The fitting of a patient with the product may only be carried out by an O&P professional.

## 4 Safety

### 4.1 Explanation of warning symbols

 <b>WARNING</b>	Warning regarding possible serious risks of accident or injury.
 <b>CAUTION</b>	Warning regarding possible risks of accident or injury.
 <b>NOTICE</b>	Warning regarding possible technical damage.

## 4.2 Structure of the safety instructions

### WARNING

#### **The heading describes the source and/or the type of hazard**

The introduction describes the consequences in case of failure to observe the safety instructions. Consequences are presented as follows if more than one consequence is possible:

- > E.g.: Consequence 1 in the event of failure to observe the hazard
- > E.g.: Consequence 2 in the event of failure to observe the hazard
- ▶ This symbol identifies activities/actions that must be observed/carried out in order to avert the hazard.

## 4.3 General safety instructions

### WARNING

#### **Non-observance of safety notices**

Personal injury/damage to the product due to using the product in certain situations.

- ▶ Observe the safety notices and the stated precautions in this accompanying document.

### WARNING

#### **Use of damaged power supply unit, adapter plug or battery charger**

Risk of electric shock due to contact with exposed, live components.

- ▶ Do not open the power supply unit, adapter plug or battery charger.
- ▶ Do not expose the power supply unit, adapter plug or battery charger to extreme loading conditions.
- ▶ Immediately replace damaged power supply units, adapter plugs or battery chargers.

### WARNING

#### **Operating the product near active implanted systems**

Interference with active implantable systems (e.g., pacemaker, defibrillator, etc.) due to electromagnetic interference of the product.

- ▶ When operating the product in the immediate vicinity of active implantable systems, ensure that the minimum distances stipulated by the manufacturer of the implant are observed.
- ▶ Make sure to observe any operating conditions and safety instructions stipulated by the manufacturer of the implant.

### NOTICE

#### **Penetration of dirt and humidity into the product**

Lack of proper charging functionality due to malfunction.

- ▶ Ensure that neither solid particles nor liquids can penetrate into the product.

### NOTICE

#### **Mechanical stress on the power supply/battery charger**

Lack of proper charging functionality due to malfunction.

- ▶ Do not subject the power supply/battery charger to mechanical vibrations or impacts.
- ▶ Check the power supply/battery charger for visible damage before each use.

**NOTICE****Operating near sources of strong magnetic or electrical interference (e.g. theft prevention systems, metal detectors)**

Lack of proper charging functionality due to malfunction.

- ▶ Avoid operating in the vicinity of visible or concealed theft prevention systems at the entrance/exit of shops, metal detectors/body scanners for persons (e.g. in airports) or other sources of strong magnetic and electrical interference (e.g. high-voltage lines, transmitters, transformer stations, computer tomographs, magnetic resonance tomographs, etc.).

**NOTICE****Improper product care**

Damage to the product due to the use of incorrect cleaning agents.

- ▶ Clean the product with a damp cloth and mild soap only.

**NOTICE****Contact of the charging plug with magnetic data carriers**

Wiping of the data carrier.

- ▶ Do not place the charging plug on credit cards, diskettes, audio or video cassettes.

**NOTICE****Independent changes or modifications carried out to the battery charger**

Lack of proper charging functionality due to malfunction.

- ▶ Have any changes or modifications carried out only by Ottobock authorised, qualified personnel.

**NOTICE****Using the product with soiled or damaged contacts**

Lack of proper charging functionality.

- ▶ Make sure that the contacts are always clean and free of grease.
- ▶ Clean the electrical contacts regularly using cotton swabs and a mild soap solution.
- ▶ Take care to avoid damaging the contacts with pointed or sharp objects.

**NOTICE****Operating the product outside of the permissible temperature range**

Lack of proper charging functionality due to malfunction.

- ▶ Avoid operating in areas outside the permissible temperature range. The section "Technical Data" contains information on the permissible temperature range.

**NOTICE****Distance to HF communication devices is too short (e.g. mobile phones, Bluetooth devices, WiFi devices)**

Lack of proper charging functionality due to malfunction.

- Therefore, keeping the following minimum distances to these HF communication devices is recommended:
- Mobile phone GSM 850 / GSM 900: 0.99 m
  - Mobile phone GSM 1800 / GSM 1900 / UMTS: 0.7 m
  - DECT cordless phones incl. base station: 0.35 m
  - WiFi (routers, access points ...): 0.22 m
  - Bluetooth devices (third-party products not approved by Ottobock): 0.22 m

## 5 Scope of delivery

- 1 pc. 757L500 AxonCharge Integral
- 1 pc. 757L500-1 AxonCharge AC adapter
- 1 pc. 757S500=EU AxonCharge adapter
- 1 pc. 757S500=US AxonCharge adapter
- 1 pc. Instructions for use (qualified personnel)

## 6 Preparing the product for use

- 1) Slide the country-specific plug adapter onto the power supply until it locks into place (see fig. 1).
- 2) Connect the round plug of the power supply to the receptacle on the battery charger so that the plug locks into place. (see fig. 2)

**INFORMATION: Ensure correct polarity (guide lug). Do not use force when connecting the cable plug to the battery charger (see fig. 3).**

- 3) Plug the power supply into the outlet.

## 7 Handling

### 7.1 Starting the Charging Process

- 1) Plug the power supply into the outlet.
- 2) LED 6 (Fig. 5, item 6) lights up green.  
**INFORMATION: If LED 6 (Fig. 5, item 6) flashes red, the charger or power supply is defective. Contact Ottobock Service.**
- 3) Connect the charging plug to the charging receptacle and allow it to snap into place (Fig. 4).
- 4) LED 6 (Fig. 5, item 6) starts flashing green and the beeper sounds.  
→ The charging process is prepared.
- 5) LEDs 2–5 light up (Fig. 5, items 2–5).  
→ The charging process starts.

### 7.2 Display for Charge Condition

The current battery charge level is indicated by the LEDs lighting up in sequence (Fig. 5, items 2–5).

LED display: ● LED illuminated ⚡ LED flashing ○ LED not illuminated

Battery charger	Battery charge level
	Battery empty
	Battery 25% charged
	Battery 50% charged
	Battery 75% charged

Battery charger	Battery charge level
	Battery 100% charged

Once charging is complete the charging plug can be disconnected. If the charging plug remains connected the battery charger will switch to trickle charging. This will not damage the battery or the battery charger.


### INFORMATION

The charging plug may be disconnected at any time. Partially charging the battery does not shorten its lifespan (no memory effect).

## 8 Troubleshooting



### 8.1 Failure

If the charging plug has been connected but charging is not taking place, investigate the following causes:

Battery charger	Cause	Correction
	The charging plug is not properly connected to the charging receptacle	Disconnect and then reconnect the charging plug
	The charging plug or charging receptacle contacts are dirty	Use a cotton swab and mild soap solution to clean the contacts of both components
	The battery or charger's temperature is outside the permissible range	Disconnect the charging plug and allow the battery or the charger to warm up or cool down

### 8.2 Diagnosis and Service Indicators

LED 1 (Fig. 5, item 1) and LED 6 (Fig. 5, item 6) serve as diagnostic and service indicators. Information is exchanged with the AXON-Bus prosthetic system when the battery charger is connected to the charging receptacle. The indicators show when service is required or there is an error in the Axon-Bus prosthetic system.

Battery charger	Probable cause	Remedy
	There are no errors	
	There is an error in the Axon-Bus prosthetic system (battery, Axon-Bus gripping component etc.).	Check the Axon-Bus prosthetic system for error messages with AxonSoft 560X500. If no error messages appear, contact Ottobock Myo-Service.

Battery charger	Probable cause	Remedy
  	Service required	The Axon-Bus components should be serviced by Ottobock Myo-Service. The service intervals depend on the period of use for the AXON-Bus prosthetic system. However, regular servicing is required during the guarantee period. Omitting a scheduled service inspection voids the guarantee.
  	The battery charger or the power supply unit is defective.	Send the product to Ottobock Myo-Service

## 9 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

### 9.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

### 9.2 Trademarks

All product names mentioned in this document are subject without restriction to the respective applicable trademark laws and are the property of the respective owners.

All brands, trade names or company names may be registered trademarks and are the property of the respective owners.

Should trademarks used in this document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

### 9.3 CE conformity

Otto Bock Healthcare Products GmbH hereby declares that the product is in compliance with European requirements for medical devices.

The product meets the requirements of the RoHS Directive 2011/65/EU on the restriction of the use of certain hazardous substances in electrical and electronic devices.

The full text of the regulations and requirements is available at the following Internet address: <http://www.ottobock.com/conformity>

## 10 Technical data

Ambient conditions	
Storage (with and without packaging)	+5 °C/+41 °F to +40 °C/+104 °F Max. 85% relative humidity, non-condensing
Transport (with and without packaging)	-20 °C/-4 °F to +60 °C/+140 °F Max. 90% relative humidity, non-condensing
Operation	-5 °C/+23 °F to +45 °C/+113 °F Max. 95% relative humidity, non-condensing
Charging the battery	+5 °C/+41 °F to +40 °C/+104 °F Max. 85% relative humidity, non-condensing

General information	
Reference number	757L500
Product lifetime	5 years
Max. charging current	1.5 A
Max. output voltage	14 V
Battery charger operating voltage	12 V
Power supply operating voltage	100 V–240 V
Power supply operating frequency	50 Hz–60 Hz

## 11 Appendices

### 11.1 Symbols Used

Manufacturer

Lot number (PPPP YYYY WW)

PPPP – plant

YYYY – year of manufacture

WW – week of manufacture

Declaration of conformity according to the applicable European directives

In some jurisdictions it is not permissible to dispose of these products with unsorted household waste. Disposal that is not in accordance with the regulations of your country may have a detrimental impact on health and the environment. Please observe the instructions of your national authority pertaining to return and collection.

Medical device

### 11.2 Glossary

The term "Axon" stands for **A**daptive **e**xchange of **n**europlacement data. The Axon-Bus is an Ottobock innovation for the field of exoprosthetics: a data transmission system, derived from safety-related bus systems in the aviation and automobile industries. For the user this means enhanced safety and reliability because of a considerably reduced sensitivity to electromagnetic interference in comparison with conventional systems.

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