
21A35=1, 21A35=2, 21A36=1, 21A36=2


DE Gebrauchsanweisung	8
EN Instructions for use	12
FR Instructions d'utilisation	16
IT Istruzioni per l'uso	21
ES Instrucciones de uso	25
PT Manual de utilização	30
NL Gebruiksaanwijzing	34
SV Bruksanvisning	39
DA Brugsanvisning	43
PL Instrukcja użytkowania	47
HU Használati utasítás	52
CS Návod k použití	56
TR Kullanma talimatı	60
EL Οδηγίες χρήσης	65

1 **i** 4 Lieferumfang/4 Scope of delivery

- 21A35=1**
- 10Y31=1
 - 10Y31=2
 - 10Y32=1

- 21A35=2**
- 10Y31=1
 - 10Y31=2
 - 10Y32=2
 - 10Y32=2

- 21A36=1**
- 10Y31=1
 - 10Y31=2
 - 10Y32=1

- 21A36=2**
- 10Y31=1
 - 10Y31=2
 - 10Y32=2
 - 10Y32=2

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INFORMATION

Date of last update: 2021-02-23

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

This document is valid for:

- Triple-control above-elbow harness (21A35=1, 21A35=2)
- Below-elbow harness (21A36=1, 21A36=2)

The product and this document are intended for qualified personnel in an orthopaedic workshop, with expert knowledge of upper limb prosthetic fittings.

INFORMATION

The prosthetic components and prosthesis shown in the illustrations are used as examples to illustrate the general process. The instructions for use of the prosthetic components selected for the patient contain detailed information and have to be used for fabrication of the prosthesis.

1.1 Explanation of warning symbols

 CAUTION	Warning regarding possible risks of accident or injury.
---	---

 NOTICE	Warning regarding possible technical damage.
--	--

1.2 Meanings of pictograms in the illustration

	Numbering for the illustrations
	Reference to the number of the corresponding section
	Numbering for a defined sequence
	Numbering for the parts of an illustration
	Observe the safety notice in the section

2 Product description

The triple-control above-elbow harness and below-elbow harness provide supplementary fixation of the prosthetic socket as well as control of body-powered and hybrid prostheses.

The below-elbow harness only has a prehensile cable, the triple-control above-elbow harness also has a lock cable and flexion cable.

The *=1 versions are supplied with a perlon cable as the pull cable, the *=2 versions with a plastic-coated steel cable.

3 Intended use

3.1 Indications for use

The product is to be used solely for upper limb prosthetic fittings.

3.2 Area of application

Suitable for the following amputation levels:

- 21A35=1, 21A35=2: Transhumeral
- 21A36=1, 21A36=2: Transradial

Suitable for the following amputation sides:

- Unilateral amputation (left or right)

Conditionally suitable for the following amputation sides:

- Bilateral amputation (The product may only be fitted by especially trained, qualified personnel.)

3.3 Lifetime

Expected lifetime: **5 years**

The design, manufacturing and requirements for the intended use of the product are based on the expected lifetime.

4 General safety instructions

 CAUTION

Reuse on other persons and improper cleaning

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ The product may be used by one person only.
- ▶ Clean the product regularly.

⚠ CAUTION

Contact with heat, embers or fire

Risk of injury (such as burns) and risk of product damage.

- ▶ Keep the product away from open flames, embers and other sources of heat.

⚠ CAUTION

Incorrect or excessively tight application

Risk of local pressure and constriction of blood vessels and nerves due to improper or excessively tight application

- ▶ Ensure that the product is applied properly and fits correctly.

NOTICE

Contact with oils, salves, lotions or other products that contain oils or acids

Insufficient stabilisation due to loss of material functionality

- ▶ Do not expose the product to oils, salves, lotions or other products that contain oils or acids.

5 Scope of delivery

The scope of delivery is shown in Figure **1** on page 2.

The article numbers of the product and product versions are shown in bold font, the product components in normal font.

Product components shown in a frame can be ordered together using the article number in the header of the frame.

The supplied quantity is indicated in parenthesis behind an article number.

6 Preparing the product for use

2 (Page 3): 1 Adapting the axilla loop

⚠ CAUTION

Use of the product in case of known allergies to chrome or nickel

Allergic reaction

- ▶ Avoid skin contact with the product in case of known allergies to chrome or nickel.
- ▶ Adapt the axilla loop as shown in the figure, noting the following points in particular:
- ▶ **1 CAUTION! Nickel** – patients sensitive to nickel should avoid wearing the buckle **2** 21Y195=25) in direct contact with their skin.
- ▶ **4** For the best possible force transmission, position the ring on the patient as shown (intersection of the line between the humeral heads and the spine).
- ▶ If required, attach the flexion cable **3** to the ring as shown in the figure (routing of the flexion cable: under the armpit).
- ▶ **5 + 6 INFORMATION:** Complete these points when finishing the definitive prosthesis. These points also apply for the flexion cable.
- ▶ **CAUTION!** Weld the ends of the shortened straps with a lighter flame.

3 (Page 4): 1 Mounting the strap buckle

- ▶ Mount the strap buckle as shown in the figure, noting the following points in particular:
- INFORMATION:** The spacer sleeve **2** allows the strap buckle to adapt to any movement.
- ▶ **1** During positioning, ensure that the straps are routed as shown in the figures (lock cable routing: over the shoulder) and the radii of the pull cables are not too tight.
- ▶ **3** Alternative component: instead of the 29C5 setting nut, the 507S15 lamination disc can also be used; a thread (M4) has to be cut into it.

- ▶ ② + ③ **NOTICE!** Make sure the pull cable is not damaged by the scissors and lighter flame.

4 (Page 5): ① Optional: removing the prehensile cable

- ▶ For a hybrid prosthesis, the pull cable of the prehensile cable can be removed as shown in the figure.

4 (Page 5): ② Mounting the Bowden cable

- ▶ Mount the Bowden cable as shown in the figure, noting the following points in particular:
 - ▶ ① Pull the Bowden cable over the cable and fasten it in the strap buckle by slightly rotating it to the right. Secure with superglue if necessary.
 - ▶ ② Instead of the 29C5 setting nut, the 507S15 lamination disc can also be used; a thread (M4) has to be cut into it.
 - ▶ ③ **NOTICE!** Do not kink the 21A17=3.8* and 21A6.

5 (Page 6): ① Mounting the pull cables

- ▶ Shorten the pull cables to a length suitable for the patient. Then mount the pull cables as shown in the figure, noting the following points in particular:

- ▶ ② Mounting the prehensile cable

Mount the components on the pull cables of the body harness and prosthetic hand (suitable tools: Perlon cable – 711F2; steel cable – 736Y6).

- ▶ ③ Mounting the lock cable

Connect the lock cable to the lock cable of the prosthetic elbow joint.

INFORMATION: The 10Y19=2 screw coupling can be used for a lock cable made of Perlon cable, the 10Y19=1 screw coupling for a lock cable made of steel cable. If necessary, a loop can also be formed on the pull cables using the 21A8 loop sleeve and a centre punch.

- ▶ ④ Mounting the flexion cable

Thread the pull cable through three bore holes and secure it with the 10Y3 clamp with screw.

Alternative mounting for the forearm of the ErgoArm: Use the 21A207 clamping plug set for mounting (not included in the scope of delivery).

INFORMATION: The flexion cable should be taut when the arm hangs straight down.

7 Use

INFORMATION

- ▶ The initial fitting and application of the product must be carried out by qualified personnel.
- ▶ Instruct the patient in the handling and care of the product.
- ▶ Advise the patient that the single-sided muscular work required for control of the body harness may result in muscular strain. Suitable compensating measures are required (such as exercise and therapy). These must be discussed with the attending doctor.

6 (Page 7): ① Practicing the motion sequence

The body harness is designed so that the function of the prosthesis remains unaffected by normal movements.

- ▶ ② Using the prehensile cable

- ① Tighten: push one or both shoulders forward.

- ② Release: push both shoulders back.

- ▶ ③ Using the lock cable

- ④ Shoulder thrust: lower the shoulder girdle.

- ⑤ Alternative movement: bring the fitted upper arm back while bringing the shoulder forward at the same time.

- ▶ ⑥ Using the flexion cable

- ① Tighten: bring the fitted arm forward.

- ② Release: bring the fitted arm back.

Optimising the course of the harness

- ▶ For the best fit of the body harness, readjust the course of the harness using the 21Y195=25 buckles.
- ▶ If necessary, secure the buckle on the strap with a seam.
- ▶ To change the cable travel for grasping, adjust the spiral nut on the Bowden cable (4) on page 3 – 3).

8 Cleaning

NOTICE

Use of improper cleaning agents

Damage to the product due to use of improper cleaning agents

- ▶ Only clean the product with the approved cleaning agents.
- ▶ To prevent skin irritation, wash the axilla pad daily (hand wash with mild detergent/do not iron/do not bleach/do not tumble dry/do not dry clean).
- ▶ Replace the axilla pad in case of noticeable wear.

9 Maintenance

- ▶ Following an individual period for the patient to get accustomed to the prosthesis, check the settings of the prosthesis and adapt them to the patient's requirements again as needed.
- ▶ Inspect the product for proper function and optimum fit at least once a year.

10 Disposal

In some jurisdictions it is not permissible to dispose of the product with unsorted household waste. Improper disposal can be harmful to health and the environment. Observe the information provided by the responsible authorities in your country regarding return, collection and disposal procedures.

11 Legal information

11.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

11.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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