



50K4 Agilium Freestep Classic

DE	Gebrauchsanweisung	4
EN	Instructions for use	10
FR	Instructions d'utilisation	17
IT	Istruzioni per l'uso	24
ES	Instrucciones de uso	30
PT	Manual de utilização	37
NL	Gebruiksaanwijzing	44
SV	Bruksanvisning	51
DA	Brugsanvisning	57
NO	Bruksanvisning	63
FI	Käyttöohje	69
PL	Instrukcja użytkowania	75
HU	Használati utasítás	82
CS	Návod k použití	88
ZH	使用说明书	94

Artikelnummer / Article Number	Seite / Side	Schuhgröße / Shoe Size	Nur bei lateraler Gonarthrose / Only with lateral knee osteoarthritis
50K4	=L	36–39	-LG
50K4	=R	36–39	-LG
50K4	=L	39–42	-LG
50K4	=R	39–42	-LG
50K4	=L	42–45	-LG
50K4	=R	42–45	-LG
50K4	=L	45–48	-LG
50K4	=R	45–48	-LG

	Außengelenkhöhe / Lateral Joint Height (A)	Schuhgröße / Shoe Size (B)
	68 mm / 2.7 inch	36–39
	78 mm / 3.0 inch	39–42
	88 mm / 3.5 inch	42–45
	98 mm / 3.9 inch	45–48

Material	Stahl/Steel, Kunststoff/Plastic, TPE, Textilkaschierter Schaum/Fabric-lined foam
-----------------	---

INFORMATION

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INFORMATION

Date of last update: 2021-02-24

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

These instructions for use contain important information for fitting and setting up the 50K4 Agilium Freestep osteoarthritis brace.

2 Product description

2.1 Components (see Fig. 1)

1	Strap guide loops	5	Allen head screws
2	Strap with padding	6	Brace ankle joint
3	Padding	7	Foot stirrup
4	Connecting tubes	8	Sole
Components not illustrated			
Height adjuster		Allen wrench	

2.2 Design

The brace consists of a foot component and a lower leg component. These components are connected via a freely movable brace ankle joint (sagittal plane) (6).

The foot component consists of a sole (8) that can be cut to size and a foot stirrup (7).

The lower leg component consists of the brace ankle joint, two adjustable connecting tubes (4) and a strap (2) with strap guide loops (1). 2 Allen head screws (5) are used to adjust the lower leg component.

The brace ankle joint, the connecting tubes and the strap guide loops are covered with washable padding (3). The strap is also padded and washable.

3 Intended use

3.1 Indications for use

The brace is intended **exclusively** for orthotic fittings of the lower limbs and **exclusively** for contact with intact skin.

The brace must be used in accordance with the indications.

3.2 Indications

- Unicompartmental knee osteoarthritis
- Conditions/injuries requiring unicompartmental relief (e.g. post-operative treatment following meniscus reconstruction or ligament injuries which require unilateral relief)

The indication must be determined by the physician.

3.3 Contraindications

3.3.1 Absolute Contraindications

None known.

3.3.2 Relative Contraindications

The following indications require consultation with a physician: skin diseases/injuries; inflammation; prominent, swollen scars; reddening and hyperthermia of the treated limb; pronounced varicose veins, especially with impaired return flow; lymphatic flow disorders, including unclear soft tissue swelling distal to the body area where the device will be applied; sensory and circulatory disorders in the legs, e. g. associated with diabetic neuropathy.

- Bicompartamental osteoarthritis of the knee

3.4 Mechanism of Action

The brace reduces or increases the varus moment on the knee joint depending on the chosen brace.

Two mechanisms influence the extent of the varus moment on the knee joint. Initially the rigid bridging of the ankle joint in the frontal plane reduces eversion in the lower ankle joint. At the same time, forces that can be adjusted at the lateral contact point of the brace are transferred in the area of the lower leg.

In the treatment of medial osteoarthritis of the knee (brace bar lateral), this effect shifts the vector of the ground reaction force laterally. Immediate consequences are the reduction of varus moment and relief for the medial compartment. In the treatment of lateral osteoarthritis of the knee (brace bar medial), the force transmission point is shifted in the medial direction. This increases the varus moment, relieving the lateral compartment. Strain can be relieved on either the lateral or medial compartment depending on the brace selection.

4 Safety

4.1 Explanation of warning symbols

 CAUTION	Warning regarding possible risks of accident or injury.
---	---

 NOTICE	Warning regarding possible technical damage.
--	--

4.2 General safety instructions

 CAUTION

Reuse on other persons and improper cleaning

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ The product may be used by one person only.

- ▶ Clean the product regularly.

⚠ CAUTION

Contact with heat, embers or fire

Risk of injury (such as burns) and risk of product damage

- ▶ Keep the product away from open flames, embers and other sources of heat.

NOTICE

Contact with oils, salves, lotions or other products that contain oils or acids

Insufficient stabilization due to loss of material functionality

- ▶ Do not expose the product to oils, salves, lotions or other products that contain oils or acids.

5 Handling

INFORMATION

- ▶ The daily duration of use and period of application are generally determined by the physician.
- ▶ The initial fitting and application of the product must be carried out by qualified personnel.
- ▶ Instruct the patient in the handling and care of the product.
- ▶ Instruct the patient to see a physician immediately if any exceptional changes are noted (e.g. worsening of the complaint).

5.1 Selecting the Size

- 1) Check whether the brace ankle joint height matches the height of the patient's lateral malleolus. If the joint heights do not match, select the brace size based on the height of the patient's lateral malleolus (see sizing table).
- 2) Determine the shoe size and select the appropriate brace size (see sizing table).
- 3) Verify that the lateral contact is at least **3 cm** below the fibular head. If the lateral contact extends closer than **3 cm** to the fibular head, select the next smaller brace size.

5.2 Adaptation

CAUTION

Incorrect alignment, assembly or adjustment

Injury and damage to product components

- ▶ Observe the alignment, assembly and adjustment instructions.
- ▶ Work on the product may only be carried out by qualified personnel.

NOTICE

Improper molding or application

Damage to the product due to overloading of the material and improper fit of the product due to breakage of load-bearing components

- ▶ The product may only be molded by qualified personnel.
- ▶ Do not make any improper changes to the product.
- ▶ Always apply the product according to the information in the instructions.

Adjusting the distance from the brace joint to the ankle

- 1) Position the foot on the sole (8) and make an outline drawing (see fig. 2). Use the removable insole as a template. The pivot point of the brace joint has to line up with the pivot point of the ankle joint.
- 2) Make sure that the ankle joint and brace ankle joint do not touch.
- 3) Cut the sole (8) to the required size (see fig. 3).
- 4) Round off the edges of the sole. Use a silicone sanding cone.
- 5) **Optional:** Take out the removable insole.
- 6) Insert the brace into the shoe (see fig. 5).
- 7) Put the brace on together with the shoe.

Optional: Gluing on the insole

- 1) Roughen the top of the foot component.
- 2) Apply adhesive to the roughened foot component and glue on the insole. In doing so, ensure that the distance between the brace ankle joint and the ankle joint is not reduced.

Optional: Adjusting the correction strength

- > The brace is in the shoe.
- 1) On the contralateral side, insert the height adjuster to compensate for the sole thickness up to **6 mm**.
 - 2) Put the brace on together with the shoe.
 - 3) Loosen the 2 Allen head screws (5) (see fig. 4).

- 4) Set the correction strength. In order to do so, push the rods firmly into the joint up to the marking line and tighten the 2 Allen screws to a torque of **6 Nm** (see fig. 6).

Optional: Optimizing straps

- ▶ Remove the hook-and-loop closure of the strap (2) and shorten the strap.

Conducting the dynamic trial fitting

- ▶ Observe the position of the ankle joint and reduce the setting of the lower leg component in case of a non-physiological heel strike (excessive supination of the foot) (see fig. 7).

5.3 Application

⚠ CAUTION

Using the brace without padding

Risk of local pressure points, superficial injuries or skin irritation due to skin contact with rigid or sharp brace components

- ▶ Only use the brace with undamaged padding.

⚠ CAUTION

Brace is put on too tightly

Risk of local pressure points and restriction of adjacent blood vessels or nerves

- ▶ Check that the brace is correctly positioned and fitted.

NOTICE

Use of a worn or damaged product

Limited effectiveness

- ▶ Before each use, check the product for functional reliability and for possible wear or damage.
- ▶ Do not continue using a product that is no longer functional, or that is worn or damaged.

- 1) Unfasten the hook-and-loop closure of the strap.
- 2) Insert the brace into the shoe.
- 3) Put the brace on together with the shoe.
- 4) Fasten the hook-and-loop closure of the strap.

5.4 Removal

- 1) Unfasten the hook-and-loop closure of the strap.
- 2) Remove the shoe together with the brace.

- 3) **Optional:** Remove the brace from the shoe and re-insert the inner sole of the shoe.

5.5 Cleaning

NOTICE

Use of improper cleaning agents

Damage to the product due to use of improper cleaning agents

- ▶ Only clean the product with the approved cleaning agents.

Clean the brace regularly:

- 1) Remove the strap and the padding.
- 2) Wash the strap and the padding in **30 °C** warm water using a standard mild detergent. Do not use fabric softener. Rinse thoroughly.
- 3) Allow to air dry. Do not expose to direct heat (e.g. from the sun, an oven or radiator).

6 Disposal

Dispose of the product in accordance with national regulations.

7 Legal Information

All legal conditions are subject to the respective national law of the country of use and may vary accordingly.

7.1 Local Legal Information

Legal information that applies **exclusively** to specific countries is written in the official language of the respective country of use after this chapter.

7.2 Liability

The manufacturer shall be liable in the event that the product is used in accordance with the descriptions and instructions in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorized modification of the product.

7.3 CE Conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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