
## 50A8, 50A9 Omo Immobil Sling

<b>DE</b> Gebrauchsanweisung .....	4
<b>EN</b> Instructions for use .....	8
<b>FR</b> Instructions d'utilisation .....	12
<b>IT</b> Istruzioni per l'uso .....	16
<b>ES</b> Instrucciones de uso .....	20
<b>PT</b> Manual de utilização .....	24
<b>NL</b> Gebruiksaanwijzing .....	27
<b>SV</b> Bruksanvisning .....	31
<b>DA</b> Brugsanvisning .....	35
<b>NO</b> Bruksanvisning .....	39
<b>FI</b> Käyttöohje .....	43
<b>PL</b> Instrukcja użytkowania .....	47
<b>HU</b> Használati utasítás .....	51
<b>CS</b> Návod k použití .....	55
<b>HR</b> Upute za uporabu .....	58
<b>SL</b> Navodila za uporabo .....	62
<b>SK</b> Návod na používanie .....	66
<b>BG</b> Инструкция за употреба .....	70
<b>TR</b> Kullanma talimatı .....	74
<b>EL</b> Οδηγίες χρήσης .....	78
<b>RU</b> Руководство по применению .....	82
<b>ZH</b> 使用说明书 .....	86

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### INFORMATION

Date of last update: 2021-02-12

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

These instructions for use provide you with important information on the fitting and application of the 50A8 Omo Immobil Sling and 50A9 Omo Immobil Sling abduction shoulder orthosis.

## 2 Intended use

### 2.1 Indications for use

The orthosis is intended **exclusively** for orthotic fittings of the shoulder and **exclusively** for contact with intact skin.

The orthosis must be used in accordance with the indications.

### 2.2 Indications

Indications must be determined by the physician.

#### 50A8 Omo Immobil Sling

- Acute or chronic shoulder pain
- Irritation in the shoulder region
- Frozen shoulder
- Acromioclavicular separation
- Osteoarthritis of the shoulder
- Subcapital humerus fracture
- Anterior shoulder dislocation
- Bankart lesion

#### 50A9 Omo Immobil Sling, abduction

- Acute or chronic shoulder pain
- Irritation in the shoulder region
- Biceps tendon and muscle injuries
- Posterior shoulder dislocation

## 2.3 Contraindications

### 2.3.1 Absolute Contraindications

Not known.

### 2.3.2 Relative Contraindications

The following complaints require consultation with a physician: skin diseases/injuries, inflammation, prominent scars that are swollen, reddening and hyperthermia in the area treated; lymphatic

flow disorders, including unclear soft tissue swelling distal to the device; sensory and circulatory disorders in the area of the arm and the hand.

## 2.4 Effects

### 50A8 Omo Immobil Sling

The orthosis immobilises and relieves the shoulder joint with interior rotation by holding the arm in place on the body. It has an analgesic effect.

### 50A9 Omo Immobil Sling, abduction

The orthosis immobilises and relieves the shoulder joint in **15°** abduction by holding the arm in place. It has an analgesic effect.

## 3 Safety

### 3.1 Explanation of warning symbols

**CAUTION** Warning regarding possible risks of accident or injury.

**NOTICE** Warning regarding possible technical damage.

### 3.2 General safety instructions

**CAUTION**

#### Reuse on other persons and improper cleaning

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ The product may be used by one person only.
- ▶ Clean the product regularly.

**CAUTION**

#### Contact with heat, embers or fire

Risk of injury (such as burns) and risk of product damage.

- ▶ Keep the product away from open flames, embers and other sources of heat.

**NOTICE**

#### Contact with oils, salves, lotions or other products that contain oils or acids

Insufficient stabilisation due to loss of material functionality

- ▶ Do not expose the product to oils, salves, lotions or other products that contain oils or acids.

**NOTICE**

#### Improper use and changes

Change in or loss of functionality as well as damage to the product

- ▶ Use this product with care and only for its intended purpose.
- ▶ Do not make any improper changes to the product.

## 4 Handling

**INFORMATION**

- ▶ The daily duration of use and period of application are generally determined by the physician.
- ▶ The initial fitting and application of the product must be carried out by qualified personnel.
- ▶ Instruct the patient in the handling and care of the product.

- ▶ Instruct the patient to see a physician immediately if any exceptional changes are noted (e.g. worsening of the complaint).

#### 4.1 Size Selection

- 1) Measure the distance from the elbow to the metacarpophalangeal joints.
- 2) Determine the orthosis size (see size chart).

#### 4.2 Fitting and Application

##### CAUTION

##### **Incorrect or excessively tight application**

Risk of local pressure and constriction of blood vessels and nerves due to improper or excessively tight application

- ▶ Ensure that the product is applied properly and fits correctly.

##### CAUTION

##### **Direct skin contact with the product**

Skin irritation due to chafing or sweat formation

- ▶ Do not wear this product directly on the skin.

##### NOTICE

##### **Use of a worn or damaged product**

Limited effectiveness

- ▶ Before each use, check the product for functional reliability and for possible wear or damage.
- ▶ Do not continue using a product that is no longer functional, or that is worn or damaged.

The shoulder orthosis is delivered closed and ready for trial fitting. It can be used on both sides and is preset for use on the right arm.

##### **Adapting the 50A8 for the left side of the body**

Remove the abdominal closure attached to the arm sling and attach it in the same position on the opposite side, as a mirror image.

##### **Adjusting and applying the 50A8**

- > All hook-and-loop closures on the orthosis are unfastened.
  - > The patient is standing.
- 1) Place the arm in the sling and close the hook-and-loop closures on the forearm and thumb (see fig. 1).
  - 2) Place the shoulder harness diagonally over the shoulder and adjust it (see fig. 2).
  - 3) Close the waist belt (see fig. 3).
  - 4) **Optional:** Shorten the belt.
  - 5) Bring the shoulder pad in the desired position.

##### **Adapting the 50A9 for the left side of the body**

- 1) Open the waist belt.
- 2) Loosen the pommel from the arm sling and fasten it as a mirror image on the opposite side with the double-sided hook-and-loop strap (see fig. 4).
- 3) Open the shoulder harness and switch the closure direction.

##### **Adjusting and applying the 50A9**

- > All hook-and-loop closures on the orthosis are unfastened.
  - > The patient is standing.
- 1) Place the arm into the arm sling, with the elbow lying in the closed corner of the arm sling.

- 2) Position the waist cushion, with the anatomically shaped end forward, on the hip for the affected arm (see fig. 5).
- 3) Adjust and close the hook-and-loop closures on the forearm and thumb (see fig. 5).
- 4) Place the shoulder harness diagonally over the shoulder and adjust it.
- 5) Adjust and close the waist belt (see fig. 6).
- 6) **Optional:** Shorten the belt.
- 7) Bring the shoulder pad in the desired position.

#### INFORMATION

To set the desired internal and external rotation position of the shoulder orthosis, slide the abduction cushion along the hip.

## 5 Cleaning

#### NOTICE

##### Use of improper cleaning agents

Damage to the product due to use of improper cleaning agents

- ▶ Only clean the product with the approved cleaning agents.

#### INFORMATION

To clean the 50A9, remove the foam cushion from the cover and only wash the cover.

- 1) Fasten all hook-and-loop closures.
- 2) Recommendation: Use a laundry bag or net.
- 3) Wash the orthosis in warm water at **30°C** with standard mild detergent. Do not use fabric softener. Rinse thoroughly.
- 4) Allow to air dry. Do not expose to direct heat sources (e.g. sunlight, stove or radiator).

## 6 Disposal

Dispose of the product in accordance with national regulations.

## 7 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

### 7.1 Local Legal Information

Legal information that applies **exclusively** to specific countries is written in the official language of the respective country of use after this chapter.

### 7.2 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

### 7.3 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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## 6 UTILIZATION

UTILIZATION of the device is performed in accordance with the prescriptions, valid in the country of use of the device.

## 7 LEGAL INSTRUCTIONS

All legal instructions apply to the law of the country in which the device is used, therefore these instructions may vary.

### 7.1 Local legal instructions

Legal instructions, which find their application **exclusively** in individual countries, are listed under this chapter in the official language of the corresponding country, in which the device is used.

### 7.2 Liability

The manufacturer bears responsibility in the case where the device is used in accordance with the descriptions and instructions, given in this document. The manufacturer is not responsible for damage, arising as a result of negligence in accordance with the provisions of the document, in particular in the case of improper use or unauthorized modification of the device.

### 7.3 Compliance with CE standards

This device complies with the requirements of the Regulation (EC) 2017/745 on medical devices. Declaration of conformity CE can be downloaded on the manufacturer's website.

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