




**8E38=6, 8E39=6, 8E41=6**

<b>DE</b> Gebrauchsanweisung.....	4
<b>EN</b> Instructions for use.....	10
<b>FR</b> Instructions d'utilisation .....	15
<b>IT</b> Istruzioni per l'uso.....	21
<b>ES</b> Instrucciones de uso .....	27
<b>PT</b> Manual de utilização .....	33
<b>NL</b> Gebruiksaanwijzing .....	39
<b>SV</b> Bruksanvisning .....	46
<b>DA</b> Brugsanvisning .....	51
<b>HU</b> Használati utasítás.....	57
<b>TR</b> Kullanma talimatı .....	63
<b>EL</b> Οδηγίες χρήσης.....	68
<b>JA</b> 取扱説明書 .....	75
<b>ZH</b> 使用说明书 .....	80

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Date of the last update: 2021-04-06

- Please read this document carefully before using the product and observe the safety notices.
- Instruct the user in the safe use of the product.
- Please contact the manufacturer if you have questions about the product or in case of problems.
- Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- Please keep this document for your records.

## 1 Application

Only put the product into use in accordance with the information contained in the accompanying documents supplied. Instruct the patient in the proper use and care of the product. The product may not be transferred to the patient without prior instruction.

The 8E38=6, 8E39=6, 8E41=6 System Electric Hands are to be used **solely** for the exo-prosthetic fitting of upper limb prostheses.

The System Electric Hand DMC plus is designed **exclusively** for use on one patient. Use of the product by other persons is not approved by the manufacturer.

### 1.1 Qualification

Patients may be fitted with the product only by qualified personnel authorised by Ottobock after completing the corresponding training.

## 2 Safety information

Failure to follow the safety instructions given below can lead to damage to or malfunction of the product. Follow the safety instructions and the precautions specified in this accompanying document.

- Do not use silicone spray when putting on the prosthetic glove; silicone spray may compromise secure fit of the glove. As a mounting aid the 633S2 Ottobock Procomfort Gel is recommended.
- The System Electric Hand DMC plus and the connected components might malfunction when near high-tension power lines, transmitters, transformers or other sources of strong electromagnetic radiation (such as security systems for goods in department stores). Setting the electrode gains as low as possible will minimise this risk.
- Before interrupting or establishing electrical connections (such as when removing the Hand from the prosthesis), the system must be disconnected from the power supply. For this purpose, either remove the battery from the battery receptacle or turn the prosthesis off by pressing the button in the charging receptacle.
- If you do not achieve the expected results with the adjustments or with the chosen program, please contact an Ottobock Myo-Service Facility.
- Instruct the patient in the appropriate handling of the System Electric Hand DMC plus as described in Chapter 3 "Patient Information".

## 3 Patient information

- Do not allow foreign particles or liquids to get into the Hand. The Hand should not be subjected to intense smoke, dust, mechanical vibrations, shocks or high temperatures.

- Avoid staying near high-tension power lines, transmitters, transformers or other sources of strong electromagnetic radiation (such as security systems for goods in department stores), as this can lead to malfunction of the System Electric Hand DMC plus.
- The System Electric Hand DMC plus has been developed for everyday use and must not be used for unusual activities such as extreme sports (free climbing, paragliding, etc.). Careful handling of the prosthesis and its components not only increases their service life but, above all, ensures the patient's personal safety! Should the prosthesis be subjected to unusual stresses (such as a fall), immediately contact your prosthetist and have the prosthesis inspected for any damage. If necessary, the responsible prosthetist will pass the prosthesis on to the Ottobock Myo-Service.
- If the Hand is not in use, it should be left in the open position to protect the mechanism and prosthetic glove.
- The System Electric Hand and damaged components may only be opened or repaired by certified Ottobock Myo-Service technicians.
- If the distance to HF communication devices (e.g. mobile phones, Bluetooth devices, WLAN devices) is too small, interference with internal data communications can result in malfunctions of the product. Therefore, keeping the following minimum distances to these HF communication devices is recommended:
  - Mobile phone GSM 850 / GSM 900: 0.99 m
  - Mobile phone GSM 1800 / GSM 1900 / UMTS: 0.7 m
  - DECT cordless phones incl. base station: 0.35 m
  - WLAN (routers, access points,...): 0.22 m
  - Bluetooth devices (third-party products not approved by Ottobock): 0.22 m
- When operating the product in the immediate vicinity of active, implantable systems (e.g. pacemakers, defibrillators, etc.), ensure that the minimum distances stipulated by the manufacturer of the implant are observed. Interference with active, implantable systems may occur due to the electromagnetic radiation emitted by the product. Furthermore, make sure to observe any operating conditions and safety instructions stipulated by the manufacturer of the implant.

### 3.1 Driving vehicles

An upper limb amputee's ability to drive a vehicle is determined on a case-by-case basis. Factors include the type of fitting (amputation level, unilateral or bilateral, residual limb conditions, design of the prosthesis) and the amputee's abilities. Drivers must observe driving laws applicable in their country. For insurance purposes, we strongly recommend having your driving ability examined and confirmed by an authorised agency. For maximum safety and convenience, Ottobock recommends that the vehicle should be adapted in a specialised garage (such as by installing a steering fork). It is imperative that the driver is able to operate the vehicle without any risk with the System Electric Hand DMC plus turned off. Driving with the System Electric Hand DMC plus turned ON may present a risk if the Hand inadvertently opens due to unintentional muscle contraction or other causes.

#### **Caution:**

Before a System Electric Hand DMC plus with a quick-disconnect wrist can be used, the terminal device must be positioned in such a way that a slight turn of it cannot disconnect the Hand from the prosthesis.

These products may not be disposed of with household waste in some jurisdictions. Disposal that is not in accordance with the regulations of your country may have a detrimental impact on health and the environment. Please observe the information provided by the responsible authorities in your country regarding return and collection processes.

### 3.2 Maintenance

Performing regular maintenance (service inspections) every 24 months is recommended to prevent injuries and maintain the quality of the product. The grace period is no more than one month before or three months after maintenance is due. In general, all products are subject to compliance with the maintenance intervals during the warranty period. This is the only way to maintain full warranty cover. Additional services such as repairs may be provided in the course of maintenance. These additional services may be provided free of charge or can be billable according to an advance cost estimate, depending on the extent and validity of the warranty. The following components must always be sent in for maintenance and repairs: The product, battery charger and power supply. The shipping container for the loaner unit you receive must be reused for sending back the components requiring inspection.

### 4 Scope of delivery

- 1 pc. 8E38=6/8E39=6/8E41=6 System Electric Hand DMC plus
- 1 pc. Instructions for Use

### 5 Description and function

The 8E38=6, 8E39=6 und 8E41=6 System Electric Hands are equipped with DMC plus proportional control.

The DMC plus consists of two independent systems which precisely control the speed and grip force of the Hand. Gripping speed and grip force are proportional to the strength of the muscle EMG signal. The gripping process is much more natural, and it is easier to gently grip small or fragile objects.

The DMC plus control is located within a **black housing** with a **black dial**. The Hand's coaxial bushing or hand cable is marked with a **brown ring**.

The System Electric Hands equipped with DMC plus are controlled by the electrodes 13E125, 13E200 or 13E202.

**The prosthetist may choose from two different control modes by inserting or removing a Function Plug from the System Electric Hand DMC plus. The DMC plus Hand will be delivered with the Function Plug inserted.**

### 6 Switching the System Electric Hand DMC plus ON and OFF

#### 6.1 Size 7

The push-button is activated by pressing it through the prosthetic glove. A time lag prevents unintentional activation. In order to activate the chosen function keep the button pressed down for approx. 1 second. (Fig. 8).

#### 6.2 Size 7¼, 7¾, 8¼

The ON/OFF switch integrated in the black housing can be activated by pressing it through the prosthetic glove.

- |                      |              |          |
|----------------------|--------------|----------|
| Back of Hand region: | Function ON  | (Fig. 8) |
| Thumb region:        | Function OFF | (Fig. 9) |

## 7 Control modes

### Notice!

The photos mentioned are valid for all hand sizes.

For optimal patient fitting, the prosthetist may choose from two different control modes by inserting or removing the 13E185 Function Plug:

### 7.1 Function Plug inserted (as delivered)

The Function is identical to that of the current DMC control (Fig. 1, Fig. 2). Grip speed and grip force are exactly determined by the strength of the muscle signal. If the muscle signal strength varies, grip speed and grip force proportionally adapt to the new muscle signal.

### 7.2 Function Plug removed (Fig. 3, Fig. 4)

Grip speed and grip force are exactly determined by the strength of the muscle signal. If the muscle signal strength varies, grip speed and grip force proportionally adapt to the new muscle signal but after gripping once with maximum force the EMG signal required to open the Hand will be slightly higher. This reduces the risk of the Hand opening due to unintentional muscle contraction. Gripping security for situations such as when eating will be increased.

To change the control mode, remove the lock ring of inner hand. Pull the inner hand up enough to remove the Function Plug (Fig. 5, Fig. 6). Pull the inner hand back in place and re-install the lock ring. To reset the system so that it will recognise the new control mode, remove the battery, then re-set it.

## 8 Adjusting the electrodes

To gain the best possible function from the DMC plus Control the electrodes must be optimally adjusted. This is easy with the 757M11=X-CHANGE.

### Notice!

The dial on the electronic control is to be used exclusively for service work. Only certified Ottobock Myo-Service Facilities may change the position of the dial.

Procedure:

For detailed information please refer to the 647H265 Instructions for Use of the 757M11=X-CHANGE. In the MyoSoft/ PAULA program, select DMC plus Control. Each electrode should be adjusted so that the patient can maintain the corresponding muscle signal at a point barely above the value HIGH for about 2 seconds (Fig. 7)

### Attention!

Allow the patient to rest during the setting up of the electrodes to avoid fatigue, and thus incorrect electrode settings (i.e. excess sensitivity).

Please take care to place the electrodes so that the entire contact surface lies directly on intact skin. If strong interference by electrical appliances still occurs, please check the positioning of the electrode and – if necessary – reposition them. If the interference cannot be eliminated, please contact an Ottobock Myo-Service Facility.

## 9 Battery management

Only use fully charged batteries to operate the System Electric Hand with DMC plus Control. An intelligent “battery manager” alerts the patient to low voltage conditions: The System Electric

Hand DMC plus will become slower and the gripping force will decrease. At the same time, the battery is protected against harmful total discharge.

The System Electric Hand with DMC plus Control can be operated with the 757B20/757B21 EnergyPack Li-Ion Battery, the 757B35=\* MyoEnergy Integral or the 757B15 X-ChangePack. For more detailed information regarding the batteries, please refer to the individual Battery Manual.

 **Attention!**

As soon as the DMC plus Control recognises a fully charged 757B20/757B21 EnergyPack or 757B35=\* MyoEnergy Integral, it automatically switches to Li-Ion battery technology. Once this has happened, only this technology should be used to operate the system. If a 757B15 X-ChangePack is used, it is not possible to use the battery's full capacity. A reset to X-ChangePack technology can only be done with the Ottobock MyoSelect 757T13. For operational security and reliability, **only** 757B20, 757B21, 757B35=\* or 757B15 Batteries are to be used.

The default setting at delivery is optimised for NiMH Batteries

## 10 Cleaning and care

Clean the product with a damp, soft cloth and mild soap (e.g. Ottobock 453H10=1 Derma Clean) when needed. Make sure the product is not damaged and that no liquid can penetrate into system components.

Then, dry the prosthesis component with a soft cloth.

## 11 Technical data

Static Current (all sizes)	mA	1
Operating Temperature (all sizes)	°C	0 – 70
Opening Width Size 7	mm	79
Opening Width Size 7¼, 7¾, 8¼	mm	100
Proportional Speed (all sizes)	mm/s	15 – 130
Proportional Grip Force (all sizes)	N	0 – 90
Power Source		
EnergyPack	757B20,	7.2 V
EnergyPack	757B21,	7.2 V
MyoEnergy Integral	757B35=* ,	7.4 V
X-ChangePack	757B15,	6 V
Service life of the terminal device	5 years	
Battery service life	2 years	

### Ambient conditions

Storage (with and without packaging)	+5 °C/+41 °F to +40 °C/+104 °F Max. 85%relative humidity, non-condensing
Transport (with and without packaging)	-20 °C/-4 °F to +60 °C/+140 °F Max. 90% relative humidity, non-condensing
Operation	-5 °C/+23 °F to +45 °C/+113 °F Max. 95% relative humidity, non-condensing
Charging the battery	+5 °C/+41 °F to +40 °C/+104 °F Max. 85%relative humidity, non-condensing

## 12 Symbols used

**MD** Medical device

## 13 Liability

Otto Bock Healthcare Products GmbH, hereafter referred to as manufacturer, assumes liability only if the user complies with the processing, operating and maintenance instructions as well as the service intervals. The manufacturer explicitly states that this device may only be used in combination with components that were authorised by the manufacturer (see instructions for use and catalogues). The manufacturer assumes no liability for damage caused by component combinations (use of third-party products) and applications which it did not authorise.

The device may only be opened and repaired by authorised Ottobock technicians.

## 14 Trademarks

All denotations within this accompanying document are subject to the provisions of the respective-applicable trademark laws and the rights of the respective owners, with no restrictions.

All brands, trade names or company names may be registered trademarks and are subject to the rights of the respective owners.

Should trademarks in this accompanying document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

## 15 CE conformity

Otto Bock Healthcare Products GmbH hereby declares that the product is in compliance with applicable European requirements for medical devices. The product meets the requirements of the RoHS Directive 2011/65/EU on the restriction of the use of certain hazardous substances in electrical and electronic devices. The full text of the regulations and requirements is available at the following Internet address: <http://www.ottobock.com/conformity>

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