
ottobock.

Agilium Freestep 3.0


DE	Kurzinformation für Anwender	4
EN	Brief information for users.....	6
FR	Brèves informations destinées à l'utilisateur	8
IT	Brevi informazioni per l'utilizzatore	10
ES	Información breve para usuarios.....	12
PT	Nota informativa para o utilizador.....	15
NL	Beknopte handleiding voor gebruikers	17
SV	Snabbguide för användare.....	19
DA	Kort information til brugere	21
NO	Kort informasjon til bruker	23
FI	Pikatieote käyttäjille	25
PL	Skrócona informacja dla użytkowników	27
HU	Rövid tájékoztató a felhasználók számára	29
CS	Stručné informace pro uživatele	31
SL	Povzetek informacij za uporabnika	33
SK	Stručná informácia pre používateľa	35
HR	Kratke informacije za korisnika	37
BG	Кратка информация за потребителя.....	39
TR	Kullanıcı için kısa bilgi	41
EL	Σύντομη ενημέρωση για χρήστες.....	43
RU	Краткая памятка для пользователей	45
JA	使用者のための情報	48
ZH	致用户的简要信息	50


1

2

3

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Dear User,

EN

You have obtained the Agilium Freestep 3.0 lower leg orthosis, a medical device that is designed to treat, for example, unicompartmental osteoarthritis of the knee joint. The special feature of the Agilium Freestep 3.0 is that it works indirectly. Rather than being applied to the knee joint itself, it is worn in the shoe and on the lower leg. This means the mobility of the knee joint is not restricted by the orthosis, which is a crucial advantage.

What is osteoarthritis of the knee? 1

Osteoarthritis is the most common of all joint diseases. Approximately 60 per cent of all osteoarthritis patients are affected by osteoarthritis of the knee (also known as gonarthrosis). Osteoarthritis can be caused by excessive strain, joint misalignment (also known as knock-knees or bow legs) or injury to one or more joints. Mechanical stress – caused by inappropriate or excessive strain – wears down the cartilage in the joint, resulting in progressive damage to the joint.

Treatment for osteoarthritis of the knee begins with the foot 2

The Agilium Freestep 3.0 alters the load line of the knee, thus providing targeted relief for the damaged cartilage in the knee joint. The load line describes the point where the body weight primarily impacts the cartilage. Relief is provided via a rigid bridge over the ankle joint and the gentle lateral pressure the orthosis applies to the lower leg while walking or standing. Pain is effectively alleviated.

When a patient is diagnosed with medial (inner side) osteoarthritis of the knee, for example, the medial contact surfaces (inside of the knee) are especially sensitive to pain (algnesia). Relief therefore has to be provided by transferring the load to the central region of the knee while walking or standing.

We have compiled a few questions on the proper use and handling of the orthosis that are answered below. Your O&P professional will be happy to answer questions about your specific situation and the fitting of your orthosis.

How is the Agilium Freestep 3.0 adapted?

An expert O&P professional adapts the orthosis. They select the appropriate model and the correct size, adapting it once to your lower leg and your footwear. This adaptation is permanent and can only be readjusted by your O&P professional.

How do I put on the Agilium Freestep 3.0? 3

First take the insole (if removable or present) out of your shoe. Now position the Agilium Freestep 3.0 in the shoe. Then put the insole back into the shoe, on top of the orthosis sole plate. Now you can slide your foot into the shoe. Your O&P professional has adjusted the orthosis so the orthosis joint is at ankle height after putting on the shoe. Bring the upper section of the orthosis up along your lower leg and pass the calf shell around your calf. Make sure that the orthosis fits securely and comfortably. The strap should not be too tight and must not cause constriction. Fasten the quick-release closure on the strap. It should engage with an audible click.

What shoes should be worn with the Agilium Freestep 3.0?

The Agilium Freestep 3.0 always has to be worn in a shoe. Wearing closed loafers or low trainers is recommended to ensure the best possible functioning of the Agilium Freestep 3.0. The Agilium Freestep 3.0 cannot be worn in combination with sandals, boots, ballet flats or stilettos.

Are you also looking for a solution that works with any footwear, or for walking barefoot? Your O&P professional will be happy to advise you regarding suitable products that match your activities and indications.

How should I care for and clean the Agilium Freestep 3.0?

The orthosis should be wiped clean using water and/or a mild cleaning agent if it gets dirty. You can remove the strap and the pads that are attached to the joint and in the calf shell with hook-and-loop closures. These components can be hand washed at max. 30° Celsius with a standard mild detergent. Please do not use fabric softener. Allow the strap and pads to air dry.

What spare parts are available?

The strap, pads and the required screws can be reordered. Just speak to your O&P professional.

Can I use the Agilium Freestep 3.0 during sports?

Yes, sports activities are possible and in fact recommended.

Notice: Our experts have compiled special exercises for osteoarthritis of the knee for you on YouTube under "Agilium Move Program".

The Agilium Freestep can be used for sports such as walking, jogging or playing tennis without any problems. You should, however, avoid extreme strain. Sports that require special shoes, such as alpine skiing, are not possible with the Agilium Freestep 3.0. We recommend using our "Agilium Select" orthosis finder at www.ottobock.com to find a suitable alternative orthosis for your sports activities.

When should I contact my specialist orthopaedic dealer?

Check the orthosis regularly for any damage or wear and tear. Also have your orthosis inspected promptly if it feels different or gets uncomfortable to wear. Just speak to your O&P professional.

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