



8E38=7/8E39=7/8E41=7

DE	Gebrauchsanweisung	4
EN	Instructions for Use	10
FR	Instructions d'utilisation	15
IT	Istruzioni per l'uso	22
ES	Instrucciones de uso	28
PT	Manual de utilização	34
NL	Gebruiksaanwijzing	40
SV	Bruksanvisning	46
NO	Bruksanvisning	52
TR	Kullanma talimatı	57
EL	Οδηγίες χρήσης	63
JA	取扱説明書	69
ZH	使用说明书	74

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Date of the last update: 2021-04-19

- Please read this document carefully before using the product and observe the safety notices.
- Instruct the user in the safe use of the product.
- Please contact the manufacturer if you have questions about the product or in case of problems.
- Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- Please keep this document for your records.

1 Application

English

The 8E38=7, 8E39=7, 8E41=7 System Electric Hands are to be used **solely** for the exo-prosthetic fitting of upper limb prostheses.

The product is designed **exclusively** for use on one patient. Use of the product by other persons is not approved by the manufacturer.

1.1 Qualification

Patients may be fitted with the product only by qualified personnel authorised by Ottobock after completing the corresponding training.

2 Safety information

Failure to follow the safety instructions given below can lead to damage to or malfunction of the product. Follow the safety instructions and the precautions specified in this accompanying document.

- Do not use silicone spray when putting on the prosthetic glove; silicone spray may compromise the secure fit of the glove. 633S2 Ottobock Procomfort Gel is recommended as a donning aid.
- In order to reduce the malfunctioning risk of the System Electric Hand and the connected components due to strong electromagnetic radiation (e.g. when near high-tension power lines, transmitters, transformers or security systems for goods in department stores), set the electrode gains as low as possible.
- Before interrupting or establishing electrical connections (such as when removing the hand from the prosthesis), the system must be disconnected from the power supply. For this purpose, either remove the battery from the battery receptacle or turn the prosthesis off by pressing the button in the charging receptacle.
- If you do not achieve the expected results with the adjustments or with the chosen program, please contact an Ottobock Myo-Service Facility.
- Instruct the patient in the appropriate handling of the System Electric Hand as described in Chapter 3 "Patient Information".

3 Patient information

- Do not let foreign particles or liquids enter the Hand. The Hand should not be subjected to intense smoke, dust, mechanical vibrations, shocks or high temperatures.
- Avoid staying near high-tension power lines, transmitters, transformers or other sources of strong electromagnetic radiation (such as security systems for goods in department stores), as this can lead to malfunction of the Electric Hand.
- The System Electric Hand as well as damaged components of the Hand may be opened or repaired only by the certified Ottobock Myo-Service.

- The System Electric Hand was developed for every-day activities. It should not be used for exceptional activities such as extreme sports (free-climbing, paragliding, etc.). Careful handling of the prosthesis and its components increases not only the life-cycle of the hand but also ensures the personal safety of the amputee! If the System Electric Hand was exposed to extreme loads (e.g. during a fall,...) then have your Orthopaedic Technician check the hand for damages. If necessary the hand should be sent to an Ottobock Myo-Service facility.
- **If the hand is not in use, it should be left in the open position to protect the mechanism and prosthetic glove.**
- Please provide your patient with the 646D44 Patient Information.
- If the distance to HF communication devices (e.g. mobile phones, Bluetooth devices, WiFi devices) is too small, interference with internal data communications can result in unexpected prosthesis behaviour. Therefore, keeping the following minimum distances to these HF communication devices is recommended:
 - Mobile phone GSM 850/GSM 900: 0.99 m
 - Mobile phone GSM 1800/GSM 1900/UMTS: 0.7 m
 - DECT cordless phones incl. base station: 0.35 m
 - WiFi (routers, access points,...): 0.22 m
 - Bluetooth devices (third-party products not approved by Ottobock): 0.22 m
- When operating the product in the immediate vicinity of active, implantable systems (e.g. pacemakers, defibrillators, etc.), ensure that the minimum distances stipulated by the manufacturer of the implant are observed. Interference with active, implantable systems may occur due to the electromagnetic radiation emitted by the product. Furthermore, make sure to observe any operating conditions and safety instructions stipulated by the manufacturer of the implant.

3.1 Driving Vehicles

An upper limb amputee's ability to drive a vehicle is determined on a case-by-case basis. Factors include the type of fitting (amputation level, unilateral or bilateral, residual limb conditions, design of the prosthesis) and the amputee's abilities.

All persons are required to observe their country's national and state driving laws when operating vehicles. For insurance purposes, drivers should have their driving ability examined and confirmed by an authorized test center.

For maximum safety and convenience, Ottobock recommends that the vehicle should be adapted in a specialised garage (such as by installing a steering fork). It is imperative that the driver is able to operate the vehicle without any risk with the System Electric Hand turned off. Driving with the System Electric Hand switched ON may endanger other traffic participants should it malfunction.

Caution! If the electric hand is equipped with a quick-disconnect wrist, the gripping must be positioned prior to use in such a way that a slight turn cannot disconnect the hand from the prosthesis!

3.2 Maintenance

Performing regular maintenance (service inspections) every 24 months is recommended to prevent injuries and maintain the quality of the product. The grace period is no more than one month before or three months after maintenance is due. In general, all products are subject to compliance with the maintenance intervals during the warranty period. This is the only way to maintain full warranty cover. Additional services such as repairs may be provided in the course of maintenance. These additional services may be provided free of charge or can be billable according to an advance cost estimate, depending on the extent and validity of the warranty. The following components must always be sent in for maintenance and repairs: The product, battery charger and power

supply. The shipping container for the loaner unit you receive must be reused for sending back the components requiring inspection.

This product may not be disposed of with regular domestic waste. Disposal that is not in accordance with the regulations of your country may have a detrimental impact on health and the environment. Please observe the information provided by the responsible authorities in your country regarding return and collection processes.

4 Scope of delivery

1 pc. 8E**=7 System Electric Hand Digital Twin

1 pce. Operating Manual

5 Description and function

The 8E38=7, 8E39=7, 8E41=7 System Electric Hands are equipped with Digital Twin Control. The Digital Twin Control system provides both digital control mode (dual site) and a double channel control mode (single site) in the same terminal device. The desired control mode can be selected by inserting or removing the 13E185 Function Plug. The Digital Twin Control is located within a **black housing** with two **white dials**. The hand's coaxial bushing or hand cable is marked with a **white ring** for identification purposes.

The digital mode (dual site) is controlled by two electrodes and the double channel mode (single site) is controlled by only one electrode - either 13E200 or 13E202.

The System Electric Hand is delivered in digital mode (dual site) with the function plug inserted.

6 Switching the System Electric Hand ON and OFF

6.1 Size 7

The on / off button is activated with slight finger pressure on the prosthesis glove. A delay prevents accidental activation. In order to activate the desired function, the button must be pressed for approximately 1 second (Figure 11).

6.2 Size 7¼, 7¾, 8¼

The on / off button integrated into the black housing can be activated by pushing on the prosthesis glove.

Back of Hand region: Function ON (Fig. 11)

Thumb region: Funcion OFF (Fig. 12)

7 Control modes

Notice: The photos mentioned are valid for all hand sizes.

For optimal patient fitting, the prosthesis may choose from two different control modes by inserting or removing the Function Plug

7.1 7.1 Function Plug inserted (as delivered)

Digital Mode - Dual site (Fig.1, Fig. 2).

After the corresponding electrode signal has exceeded a certain threshold value, the hand opens and closes at a constant speed.

7.2 Function Plug removed

Double channel mode - single site (Fig. 3, Fig. 4).

The hand is opened and closed by one electrode only.

A slow and gentle muscle signal will close the hand.

A quick and strong muscle signal will open the hand.

8 Adjusting the electrode(s)

To gain the best possible function from the Digital Twin Control the electrodes must be optimally adjusted. This is easy with the 757M11 MyoBoy.

Notice: Dial 1 (Fig. 5, Fig 6, red arrow) is to be used exclusively for service work. Only certified Ottobock Service facilities may change the position of the dial.

8.1 Procedure with digital mode (with Function Plug)

For detailed information please refer to the 647G265 MyoBoy Instructions for Use.

In the MyoSoft/PAULA program, select „Digital Twin with Function Plug“. Adjust each electrode so that the muscle signal can be maintained at a point barely above the value ON (Fig. 7) for approximately two seconds. Dial 2 on the electronic control has no function in the digital mode.

8.2 Procedure with double channel mode - single site (without Function Plug)

For detailed information please refer to the 647G265 MyoBoy Instructions for Use.

Before changing the digital mode (dual site) to the double channel mode (single site), remove the inner hand lock ring. Then remove the inner hand and function plug (Fig. 9, Fig. 10). To reset the system so that it will recognize the new control mode, remove the battery, then re-insert it.

In the MyoSoft/PAULA program, select „Digital Twin without Function Plug“. Adjust the slider on the scale to 3.5. Set Dial 2 on the electronic control also to 3.5 (Fig. 5, Fig. 6).

Adjust the electrode so that the muscle signal reaches a value just below the second threshold marked in MyoSoft/PAULA by the broken line: the hand closes (Fig. 8).

Overcome the threshold with a quick, strong muscle signal to open the hand.

If the muscle signal is too slow or weak, set the threshold in MyoSoft/PAULA to a lower value (slider). Set Dial 2 (Fig. 5, Fig. 6) on the electronic control to the same level.

Pull on the inner hand and seal it with 633F11 Silicone Grease. Remount the lock ring.

Attention! Allow the patient to rest during the adjustment of the electrodes to avoid fatigue, and thus incorrect electrode settings (i.e. excess sensitivity). The electrodes should be placed on intact skin only and with as much electrode-skin contact as possible. In the case of strong interferences from electric devices, the position of the electrodes should be verified and, if need be, repositioned. If the interference cannot be eliminated, please contact the Ottobock Myo-Service.

9 Battery management

Only use fully charged Ottobock batteries to operate the System Electric Hand with Digital Twin Control.

The System Electric Hand with Digital Twin Control can be operated with the 757B20 / 757B21 Li-Ion battery EnergyPack, 757B35=* MyoEnergy Integral or 757B15 NiMH battery. For more information regarding the batteries please consult the Operating Manual accompanying the batteries.

Attention! As soon as the Digital Twin control unit **recognises a fully charged** 757B20/757B21 EnergyPack or 757B35=* MyoEnergy Integral, it automatically switches to the new Li-Ion battery technology. Subsequently, only this technology should be used to operate the system. If an Ottobock 757B15 X-ChangePack is used, its capacity can no longer be fully utilised. However, the 757T13 Ottobock MyoSelect can be used to reset the system for operation with X-ChangePack. For operational security and reliability, only 757B20, 757B21, 757B35=* or 757B15 Ottobock batteries are to be used.

The default setting at delivery is optimized for NiMH Batteries.

10 Technical data

Static Current (all sizes)	mA	1
Operating Temperature (all sizes)	°C	0–70
Opening Width Size 7	mm	79
Opening Width Size 7¼, 7¾, 8¼	mm	100
Grip Speed	mm/s	110
Grip Force (all sizes)	N	0–90
Power Source	Ottobock EnergyPack	757B20 – 7,2 V
	Ottobock EnergyPack	757B21 – 7,2 V
	MyoEnergy Integral	757B35=* – 7,4 V
	Ottobock X-ChangePack	757B15 – 6 V
Service life of the terminal device	5 years	
Battery service life	2 years	

Ambient conditions

Storage (with and without packaging)	+5 °C/+41 °F to +40 °C/+104 °F Max. 85%relative humidity, non-condensing
Transport (with and without packaging)	-20 °C/-4 °F to +60 °C/+140 °F Max. 90% relative humidity, non-condensing
Operation	-5 °C/+23 °F to +45 °C/+113 °F Max. 95% relative humidity, non-condensing
Charging the battery	+5 °C/+41 °F to +40 °C/+104 °F Max. 85%relative humidity, non-condensing

11 Symbols used

MD Medical device

12 Liability

Otto Bock Healthcare Products GmbH, hereafter referred to as manufacturer, assumes liability only if the user complies with the processing, operating and maintenance instructions as well as the service intervals. The manufacturer explicitly states that this product may only be used in combination with components that were authorized by the manufacturer (see instructions for use and catalogs). The manufacturer assumes no liability for damage caused by component combinations (use of third-party products) and applications which it did not authorise.

The device may only be opened and repaired by authorised Ottobock technicians.

13 Trademarks

All denotations within this accompanying document are subject to the provisions of the respective-applicable trademark laws and the rights of the respective owners, with no restrictions.

All brands, trade names or company names may be registered trademarks and are subject to the rights of the respective owners.

Should trademarks in this accompanying document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

14 CE conformity

Otto Bock Healthcare Products GmbH hereby declares that the product is in compliance with applicable European requirements for medical devices. The product meets the requirements of the RoHS Directive 2011/65/EU on the restriction of the use of certain hazardous substances in electrical and electronic devices. The full text of the regulations and requirements is available at the following Internet address: <http://www.ottobock.com/conformity>.

15 Appendices

15.1 Operating states/error signals

Battery charge level

An intelligent „battery manager“ alerts the patient to low voltage conditions: The System Electric Hand will become slower and the gripping force will decrease. At the same time, the battery is protected against harmful total discharge.

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