






4R120, 4R121=30, 4R121=34

DE	Gebrauchsanweisung	4
EN	Instructions for use	12
FR	Instructions d'utilisation	20
IT	Istruzioni per l'uso	28
ES	Instrucciones de uso	36
PT	Manual de utilização	45
NL	Gebruiksaanwijzing	53
SV	Bruksanvisning	62
DA	Brugsanvisning	70
NO	Bruksanvisning	78
PL	Instrukcja użytkowania	85
HU	Használati utasítás	94
HR	Upute za uporabu	102
EL	Οδηγίες χρήσης	110
RU	Руководство по применению	120
JA	取扱説明書	129
ZH	使用说明书	136

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11 Angewandte Symbole

Date of the last update: 2021-03-23

- Please read this document carefully before using the product and observe the safety notices.
- Instruct the user in the safe use of the product.
- Please contact the manufacturer if you have questions about the product or in case of problems. Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- Please keep this document for your records..

Please read the following instructions carefully to achieve the best possible fitting of the DeltaTwist. By understanding the function and fitting possibilities of the DeltaTwist you can better serve your patient.

1 Application

The Ottobock DeltaTwist 4R120, 4R121=30 and 4R121=34 with Rotational Function are to be used **solely** for the exoprosthetic fitting of the lower limbs. The Ottobock DeltaTwist can be used for transfemoral, transtibial, as well as distal and proximal applications. The DeltaTwist is not intended for high performance or athletic use. The DeltaTwist should not be used in bathing prostheses.

1.1 The DeltaTwist is available in the following variations:



- 4R120 with Ø 30 mm tube clamp
 4R121=30 with Ø 30 mm tube
 4R121=34 with Ø 34 mm tube

1.2 Field of application according to MOBIS

Field of application according to the **Ottobock Mobility System MOBIS**:

Recommendation for **mobility degrees 2 to 4** (walker confined to a limited radius outdoors, non-restricted outdoor walker and non-restricted outdoor walker with especially rigorous demands).

Approved for a patient with a weight of up to **100 kg** (4R120, 4R121=30) / **125 kg** (4R121=34).

1.3 Contraindications

All conditions which contradict or go beyond the specifications listed in the section on "Safety" and "Indications for use".

The Ottobock Mobis is not intended to be used as a guide to obtaining reimbursement for prosthetic components in the USA.

2 Safety information

- The DeltaTwist should be examined visually for any damage or corrosion during routine check-ups

by qualified personnel.

- Avoid contact with the skin during application of the 633F30 lubricant.
- Do not use the DeltaTwist in a bathing prosthesis.
- Any changes in technical parts are forbidden without the express written consent of the manufacturer.
- The DeltaTwist may only be repaired by an authorized Ottobock technician.
- Use the „Patient Information“ section 3 of this manual to instruct the patient in both the appropriate use and application of the DeltaTwist.

3 Patient information

- Handle the DeltaTwist carefully and examine regularly for any signs of damage.
- Prevent any unnecessary strain on, or impact to, the housing.
- Avoid letting the temperature of the DeltaTwist rise above 70 °C.
- Do not allow debris or liquids to get into the DeltaTwist. The DeltaTwist should not be subjected to intense smoke, dust, mechanical vibrations, shocks or high temperatures.
- The DeltaTwist is not intended for use in high performance or endurance sports.
- If you notice any changes in function, please contact your orthopedic technician. Patient Information

4 Contents of delivery

Every DeltaTwist model is delivered as a set (Fig. 1 shows the 4R120 as an example). The set includes:

Pos. 1	1 pc. DeltaTwist equipped with:		
	1 pc. Elastomer Rod durometer	<i>MEDIUM</i>	YELLOW
	2 pc. Elastomer Plates durometer	<i>MEDIUM</i>	YELLOW
Pos. 2	1 pc. Elastomer Rod Durometer	<i>SOFT</i>	RED
Pos. 3	1 pc. Elastomer Rod Durometer	<i>HARD</i>	GREEN
Pos. 4	2 pc. Elastomer Plates Durometer	<i>SOFT</i>	RED
Pos. 5	2 pc. Elastomer Plates Durometer	<i>HARD</i>	GREEN
Pos. 6	1 pc. Hook Wrench		
Pos. 7	1 pc. Tube Lubricant		
	1 pc. Manual		

4.1 Replacement parts and accessories

DeltaTwist components can be ordered individually from Ottobock. The following component amounts refer to a single package.

Pc.	Description		Color	Article #
1	Elastomer Rod Durometer	<i>SOFT</i>	red	709H6=1
1	Elastomer Rod Durometer	<i>MEDIUM</i>	yellow	709H6=2
1	Elastomer Rod Durometer	<i>HARD</i>	green	709H6=3
1	Elastomer Plate Durometer	<i>SOFT</i>	red	709H5=1
1	Elastomer Plate Durometer	<i>MEDIUM</i>	yellow	709H5=2
1	Elastomer Plate Durometer	<i>HARD</i>	green	709H5=3
1	Rotational Locking Plate	<i>EXTRA HARD</i>	black	709H5=4
1	Hook Wrench			709H4
1	Tube Lubricant			633F30
1	Connection Adapter, rotatable, with adjustment screw			4R72=D
1	Connection Adapter, rotatable, with adjustment screw			4R72=D-62

1	Connection Adapter, rotatable, with adjustment pyramid	4R84=D
1	Socket Wrench Size 11, ½",	709H7
1	Extension for Socket Wrench ½"	709H8
1	T-Handle for Extension ½"	709H9

5 Description

The Ottobock DeltaTwist combines the function of a shock absorber and the function of a torsional adapter in one unit. It allows for a dampening motion of up to 8 mm up and down as well as about +/- 20 degrees inner and outer rotation. Modular Elastomer components (plates and rods of varying durometers and a continuously adjustable tension) allow you to individualize the spring and dampening characteristics.

6 Function

6.1 Vertical compression

A short-term load occurs during initial contact while walking, which leads to force and moment impact in the extremities. The involved leg takes on body weight shortly after initial contact and the related increase in force may also feel like an impact to the patient. The DeltaTwist vertical compression characteristic dampens the dispersion of shock waves during initial contact while also reducing the rate of force during loading-response. In this way the locomotor system can better react to the loading-response: an improvement for both the comfort and loading-situation for the prosthesis wearer. The Shock Absorber resistance, with its exchangeable Elastomer bars, is continuously adjustable.

6.2 Torsion

During the stance phase of the gait cycle an internal rotation of the foot (relative to the hip) occurs. In order to have a natural, healthy gait it is important to obtain internal rotation for prosthetic limbs, too.

The high torsional stiffness of exoprosthetic systems and prevention of any rotation of the joints through the socket can increase the load on the residual limb. Furthermore, the patient might perform disadvantageous compensatory movements.

The rotational function of the shock-absorber makes it possible to compensate for the torsional stiffness of an exoprosthetic system by allowing for individual settings for both inner and outer rotation. A softer setting for the internal rotation insures the desired mobility while a harder setting for the external rotation offers stability and an increased sense of security for the patient. The resistance offered by the Elastomer plates gives the patient an increased sense of stability and more security in the starting position of the prosthetic foot.

7 Installing the DeltaTwist

7.1 General Information

To decrease inertial force and inertial moments, we recommend installing the DeltaTwist in the prosthesis as proximally as possible. The DeltaTwist can be installed both above and below the knee. To ensure the function and safety of the DeltaTwist connection it is necessary to observe the recommended torque values.

Proper use of the tube clamp and the socket adapter in connection with all Ottobock adapters fulfills the load capacity requirements of ISO 10328. If the DeltaTwist is used in combination with parts from other manufacturers', it is possible that the load capacity requirements of ISO 10328 cannot be fulfilled.

7.2 Using the socket adapter and adjustment screws

The contact surfaces of the adjustment pyramid and socket adapter must be clean and smooth. It is possible to use the four adjustment screws to make static corrections during the installation, trial fitting and even after the final fitting.

During disassembly or while changing the DeltaTwist the previously adjusted position can be maintained if the two adjustment screws that are set the deepest, and which lie next to each other, are removed. Adjustment screws in their final adjustment that stick out too far (and can possibly damage the foam cover) or which are too short (and might endanger stability) should be replaced with better fitting adjustment screws.

Adjustment screws are available in the following sizes.

Order number	Length
506G3=M8x10	10 mm
506G3=M8x12-V	12 mm
506G3=M8x14	14 mm
506G3=M8x16	16 mm

Adjustment screws should only be lightly tightened (max. 1 Nm) for the initial adjustment of the adjustment pyramid.

To tighten the adjustment screws, use a torque wrench such as the Ottobock 710 D4 and tighten to 15 Nm. Apply Loctite 636K13 to the adjustment screws when final adjustments have been made.

7.3 Model 4R120 with tube clamp

The DeltaTwist with Tube clamp 4R120 can be installed both proximally and distally. If the Tube clamp is fitted distally, make sure that the clamp screw lies at least 150 mm above the ankle axis. The Tube clamp is oriented posteriorly and tightened to 10 Nm.

The Tube clamp acts a receptacle for either a tube, a tube adapter, or a connection adapter, which can be inserted into the DeltaTwist up to the stop.

If a very low height is required we recommend using the 4R72=D or the 4R72=D-68 connection adapter.

7.4 Model 4R121=30 with tube

For transtibial fittings the DeltaTwist 4R121=30 Pyramid Receiver must be fitted distally. Please be sure that the tube section has the necessary length required by the respective connector, and also make sure that the tube is inserted deeply enough into the connection components.

7.5 Model 4R121=34 with tube

For transtibial fittings the DeltaTwist 4R121=34 Pyramid Receiver must be fitted distally. Please pay attention to the markings for minimal length to avoid damaging the DeltaTwist! Please also make sure that the tube is inserted deeply enough into the connection components.

7.6 Cutting the tube

Do not clamp tube in vise! The cut edge must be even and make a right angle with the length of the tube. To cut the DeltaTwist 4R121=30 and 4R121=34, please use a tube cutter (such as the Ottobock Tube Cutter 719R2) and apply at a right angle. The tube must be deburred inside and out (as with the Ottobock Pipe Reamer 718R1) and should fit evenly against the inside wall of the Tube clamp.

7.7 Maintenance instructions

Notice:

As a basic principle, all Ottobock modular adapters are subjected to tests involving three million load cycles. Depending on the amputee's activity this corresponds to a service life of three to five years.

We recommend carrying out regular safety checks once a year.

8 Exchanging the elastomer elements and individualizing the adjustment

8.1 General information

Never place the DeltaTwist in a vise, clamping chuck or similar tool to loosen the screws. Avoid placing any unusual strain on the DeltaTwist body. To loosen or tighten the Threaded Bushing with Inner Friction Bearing follow directions as shown analogously in Fig. 11. Notice that the Hook Wrench has a sharp tip at the end for removing the Elastomer Plates.

After disassembly the distal and proximal bearings as well as the Elastomer elements and the plastic rotational sections should be wiped clean of all lubricants. Apply liberal amounts of new lubricant before reassembly. Use the included 633F30 Lubricant for this purpose.

8.2 Customizing the vertical compression characteristics

The DeltaTwist comes with the MEDIUM Elastomer Rod pre-installed. The vertical compression can be individually adjusted for the patient by interchanging the Elastomer Rods.

8.2.1 Disassembly

Notice: To avoid damage to the DeltaTwist, tension on the DeltaTwist should be reduced to a minimum.

- 1) To reduce the tension to a minimum, turn the red tensioning screw (counter-clockwise) all the way to the stop (Fig. 2) with an Allen Wrench (Size 6).
- 2) Remove the green locking screw by unscrewing counter-clockwise with an Allen Wrench (Size 6) or a Socket Wrench (Size 11) (Fig. 3). For DeltaTwist models 4R121=30 and 4R121=34 use a Socket Wrench with Extension.
- 3) To remove the Elastomer Rod, shake it out of the DeltaTwist (Fig. 4).
- 4) Use the table below to select the durometer of the desired Elastomer Rod.:

Table 1: Selecting the Durometer

	SOFT (Red)	MEDIUM (Yellow)	HARD (Green)
< 75 kg	recommended	acceptable	not recommended
75 kg - 100 kg	not recommended	recommended	recommended
100 kg - 125 kg	not recommended	acceptable	recommended

Example 1: The patient weight is 72 kg. An Elastomer Rod with a SOFT durometer is recommended, although the harder option, MEDIUM is also acceptable. A HARD durometer Elastomer Rod will suppress the vertical compression too much.

Example 2: The patient weight is 83 kg. An Elastomer Rod with a MEDIUM or HARD Durometer is recommended. The SOFT variation is not recommended.

8.2.2 Assembly

- 1) Lubricate the selected Elastomer Rod with 633F30 Lubricant.
- 2) Place the Elastomer Rod in the DeltaTwist.
- 3) Turn the green locking screw into the DeltaTwist and tighten until it is hand tight (approx. 4 Nm).

8.2.3 Setting the tension

The tension setting prevents the compression of the DeltaTwist when standing on both legs, and also guarantees that the length of the unloaded prosthesis does not increase (i.e., in swing-phase). The patient can select more or less tension. Make sure that the installed Elastomer Rod has no play in it.

- 1) Use an Allen Wrench (Size 6) to unscrew the red tension screw (counter-clockwise) all the way to the stop. When you reach the stop you will find a noticeable resistance, which can be loaded with a maximum of 4 Nm (hand tight). This is the starting point for setting the tension.
- 2) Use Table 2 below to determine the number of turns needed to set the red tension screw.

Table 2: Selection of Tension (in Tension Screw Rotations)

	SOFT (Red)	MEDIUM (Yellow)	HARD (Green)
< 75 kg	max. 8	max. 5	—
75 kg - 100 kg	—	2 - 8	2 - 6
100 kg - 125 kg	—	8	6 - 8

- 3) Turn the red tension screw equal to the number of turns recommended in the table. The tension screw must not be tightened more than 8 turns. After 8 turns you will reach a stop, which will offer a noticeable resistance. This stop can be loaded with a maximum of 4 Nm (hand tight).

8.2.4 Setting the Amount of vertical compression


The pre-tension for the DeltaTwist should be set so that it compresses no more than 1 mm when standing on both legs, and no more than 4 mm while walking.

For ease of use, a black Indicator Ring has been included (Fig. 5a and 5b) that can stay on the DeltaTwist.

Ask the patient to stand equally on both legs (to check this, the Ottobock L.A.S.A.R. Posture can be employed). Press the Indicator Ring to the blue Wiper Ring (Fig. 5a). Finally, ask the patient to remove all weight from the DeltaTwist. The Indicator Ring should not rise completely above the edge of the Housing (Fig. 5b).

8.3 Adjusting the rotational function

The torsional stiffness of the DeltaTwist can be individually adapted to the patient by exchanging the Elastomer Plates. The DeltaTwist comes with MEDIUM (yellow) Elastomer Plates preinstalled.

The DeltaTwist rotational function can be set in a variety of ways (for example, the inner rotation can be set softer than the outer rotation if indicated).

Very stiff EXTRA-HARD durometer Elastomer Plates (Black), which suppress the rotational function, are available as an accessory (article number 709H5=4). The DeltaTwist must be disassembled in order to exchange the Elastomer Plates.

8.3.1 Disassembly

Notice: To avoid damage to the DeltaTwist, the Elastomer Rod must not be held under tension by the tensioning screw.

- 1) To set the tension to a minimum, use a Allen Head Wrench (Size 6) to turn the red tension screw counter-clockwise (Fig. 2).
- 2) Remove the Threaded Bushing with Inner Friction Bearing from the DeltaTwist housing with the hook wrench, as shown analogously in Fig. 11.
- 3) Turn the blue retaining screw counter-clockwise with a Allen Wrench (Size 10). Hold the DeltaTwist horizontally and make sure no inner parts fall out.
- 4) Turn the blue retaining screw until the inner parts can be taken out.
- 5) Lever the Elastomer Plate out with the tip located on the Combination Wrench (Fig. 6).
- 6) If the patient wants only the vertical compression or only one rotational direction, you can suppress the function by installing the EXTRA HARD Elastomer Plate (articlenumber 709H5=4). If the patient wants full rotational function, fit the Elastomer Plates for inner rotation and the Elastomer Plates for outer rotation based on the information in Table 3 below.

Notice: Remember that you can select a durometer for the inner rotation that is softer than the durometer for the outer rotation. (See section 6.2).

Table 3: Recommendations for Selecting Elastomer Plates

	Internal Rotation	External Rotation
< 75 kg	SOFT	MEDIUM/SOFT
75 kg - 100 kg	MEDIUM	HARD/MEDIUM
100 kg - 125 kg	HARD	X-HARD/HARD

Example: The patient weight is 83 kg. For the internal rotation a durometer of MEDIUM (yellow) can be employed. For the external rotation the HARD durometer can be employed, as well as the MEDIUM durometer.

8.3.2 Assembly

Positioning of Elastomer plates

- Fig. 7 (item 1) Elastomer plates, inner rotation, left leg
- Fig. 7 (item 2) Elastomer plates, outer rotation, left leg
- Fig. 7 (item 3) Elastomer plates, inner rotation, right leg
- Fig. 7 (item 4) Elastomer plates, outer rotation, right leg

- 1) Grease the selected Elastomer Plate.
- 2) Arch the Elastomer Plate slightly (Fig. 8) and place it between the rotational section (Fig. 7, Pos. 5) and the stop (Fig. 7, Pos. 6).
- 3) Press the Elastomer Plate against the inner surface of the housing until the Elastomer Plate lines up with the rotational section (Fig. 9).

NOTE: Insertion of the Rotational Locking Plate, durometer EXTRA-HARD, is made easier during assembly if the black Elastomer Insert is inserted last between the white rotational section and the black locking section. Disassembly is made easier if the black Elastomer Insert is removed first.

- 4) Press the inner part of the housing up to the stop
- 5) Set the DeltaTwist on a clean surface and hold the housing with one hand. While that hand presses the housing downwards so that the blue retaining screw can be reached, turn the blue retaining screw clockwise with the other hand (Fig. 10).
- 6) Tighten the blue retaining screw hand tight (approx. 4 Nm).

- 7) Insert the Threaded Bushing with Inner Friction Bearing and tighten it to between 15-20 Nm with the 709H4 Hook Wrench. To apply the correct torque use the Pyramid Receiver, as shown in Fig. 11, where a prosthetic foot is used as an example.
- 8) Reset the tension of the Elastomer Rod by turning the red tension screw (clockwise). Make sure not to loosen the blue locking screw while doing this.

9 Technical data

Article number	4R120	4R121=30	4R121=34
Proximal connection	Tube clamp Ø 30 mm	Tube Ø 30 mm	Tube Ø 34 mm
Distal connection	Pyramid receiver	Pyramid receiver	Pyramid receiver
Outer diameter	47 mm	47 mm	47 mm
Weight	~ 340 g	~ 530 g *	~ 585 g *
System height	117 mm	–	–
Min. system height	–	177 mm	218 mm
Max. system height	–	553 mm	578 mm
Material	Aluminum	Aluminum	Aluminum
Max. torsion inner/ outer	± 20°	± 20°	± 20°
Max. dampening	~ 8 mm	~ 8 mm	~ 8 mm
Mobility Grade	2 - 4		
Max. body weight	100 kg		125 kg

10 Legal information

10.1 Liability

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the Instructions for Use. The device must only be operated with tested modular components in accordance with the Ottobock Mobility System. The manufacturer is not liable for damage caused by component combinations that were not authorised by the manufacturer.

10.2 Trademarks

All denotations within this accompanying document are subject to the provisions of the respective applicable trademark laws and the rights of the respective owners, with no restrictions.

All brands, trade names or company names may be registered trademarks and are subject to the rights of the respective owners.

Should trademarks in this accompanying document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

10.3 CE conformity

Otto Bock Healthcare Products GmbH hereby declares that the product is in compliance with applicable European requirements for medical devices.

The full text of the regulations and requirements is available at the following Internet address: <http://www.ottobock.com/conformity>

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