1B1-2=* Meridium

DE	Gebrauchsanweisung (Fachpersonal)	3
EN	Instructions for use (qualified personnel)	43
FR	Instructions d'utilisation (Personnel spécialisé)	83
IT	Istruzioni per l'uso (Personale tecnico specializzato)	125
ES	Instrucciones de uso (Personal técnico especializado)	165
PT	Manual de utilização (Pessoal técnico)	205
NL	Gebruiksaanwijzing (Vakmensen)	245
SV	Bruksanvisning (Fackpersonal)	285
DA	Brugsanvisning (Faguddannet personale)	321
NO	Bruksanvisning (Fagpersonell)	359
FI	Käyttöohje (Ammattihenkilöstö)	397
PL	Instrukcja użytkowania (Personel fachowy)	435
HU	Használati utasítás (szakszemélyzet)	475
CS	Návod k použití (Odborný personál)	513
HR	Upute za uporabu (Stručno osoblje)	551
SK	Návod na používanie (Odborný personál)	589
TR	Kullanma talimatı (Uzman personel)	627
RU	Руководство по применению (Квалифицированный персонал)	665
JA	取扱説明書 (有資格担当者)	707
ZH	使用说明书 (专业人员)	743

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DE

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Table of contents

EN

1	Foreword	45
2	Product description	45
2.1	Design	45
2.2	Function	45
2.3	Combination possibilities	46
3	Intended use	46
3.1	Indications for use	46
3.2	Conditions of use	46
3.3	Indications	47
3.4	Contraindications	47
3.5	Qualification	47
4	Safety	47
4.1	Explanation of warning symbols	47
4.2	Structure of the safety instructions	47
4.3	General safety instructions	47
4.4	Information on the Power Supply/Battery Charging	49
4.5	Battery charger information	50
4.6	Information on Alignment/Adjustment	50
4.7	Information on Proximity to Certain Areas	51
4.8	Information on Use	52
4.9	Notes on the safety modes	53
4.10	Information on the use of a mobile device with the cockpit app	53
5	Scope of Delivery and Accessories	54
5.1	Scope of delivery	54
5.2	Accessories	54
6	Charging the battery	54
6.1	Connecting the power supply and battery charger	55
6.2	Charging the prosthesis battery	55
6.3	Display of the current charge level	55
7	Preparing the product for use	56
7.1	Applying/removing footshell	56
7.2	Alignment	57
7.2.1	Settings with the "M-Soft" adjustment software	57
7.2.1.1	Introduction	57
7.2.1.2	Data transfer between the product and the PC	57
7.2.1.3	Preparing the product to connect to the adjustment software	57
7.2.2	Bench alignment in alignment apparatus	57
7.2.3	Static alignment optimisation	58
7.2.4	Dynamic alignment optimisation	58
7.2.5	Mounting cover plate/connecting plate/cover cap	58
7.2.5.1	Mounting the cover cap	58
7.2.5.2	Mounting the cover plate	58
7.2.5.3	Mounting the Protector connecting plate	59
7.2.5.4	Mounting the connecting plate for the cosmetic foam cover	59
8	Cockpit app	59
8.1	Initial connection between cockpit app and prosthesis	59
8.1.1	Starting the cockpit app for the first time	60
8.2	Control elements for cockpit app	60
8.2.1	Cockpit app navigation menu	61
8.3	Managing components	61
8.3.1	Adding component	61
8.3.2	Deleting a component	62
8.3.3	Connecting component with multiple mobile devices	62

9	Use	62
9.1	Setting the heel height.....	62
9.1.1	Setting heel height using motion pattern.....	62
9.1.2	Setting the heel height with the cockpit app.....	63
9.2	Movement patterns in basic mode (mode 1).....	63
9.2.1	Standing.....	63
9.2.2	Walking.....	63
9.2.3	Sitting down/sitting.....	64
9.2.4	Standing up.....	64
9.2.5	Walking up stairs.....	64
9.2.6	Walking down stairs.....	64
9.2.6.1	Stair function.....	64
9.2.7	Walking up a ramp.....	65
9.2.8	Walking down a ramp.....	65
9.2.9	Walking backwards.....	65
9.2.10	Standing on inclined surface.....	65
9.2.11	Kneeling.....	66
9.2.12	Relief function.....	66
9.3	MyModes.....	66
9.3.1	Switching MyModes with the cockpit app.....	66
9.3.2	Switching MyModes using motion patterns.....	67
9.3.3	Switching on ankle lock.....	67
9.3.4	Switching from a MyMode back to basic mode.....	68
9.4	Changing prosthesis settings.....	68
9.4.1	Changing the prosthesis setting using the cockpit app.....	69
9.4.2	Overview of adjustment parameters in basic mode.....	69
9.4.3	Overview of adjustment parameters in MyModes.....	69
9.5	Switching off the product.....	70
9.6	Turning Bluetooth on the prosthesis on/off.....	70
9.6.1	Switching Bluetooth off/on using the cockpit app.....	70
9.7	Querying the prosthesis status.....	71
9.7.1	Query status through cockpit app.....	71
9.7.2	Status display in the cockpit app.....	71
10	Additional operating states (modes)	71
10.1	Empty battery mode.....	71
10.2	Mode for charging the prosthesis.....	71
10.3	Safety mode.....	71
10.4	Overheating mode.....	71
11	Cleaning	71
12	Maintenance	72
12.1	Identification of the product by the Service Center.....	72
13	Legal information	72
13.1	Liability.....	72
13.2	Trademarks.....	72
13.3	CE conformity.....	72
13.4	Local Legal Information.....	73
14	Technical data	74
15	Appendices	75
15.1	Symbols Used.....	75
15.2	Operating states/error signals.....	76
15.2.1	Signals for operating states.....	76
15.2.2	Warnings/error signals.....	77
15.2.3	Error messages while establishing a connection with the cockpit app.....	78
15.2.4	Status signals.....	79
15.3	Directives and manufacturer's declaration.....	79
15.3.1	Electromagnetic environment.....	79

1 Foreword

INFORMATION

Date of last update: 2021-05-25

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- ▶ Please keep this document for your records.

The product "1B1-2=* Meridium" is referred to as the product/component/prosthesis/prosthetic foot below.

These instructions for use provide you with important information on the use, adaptation and handling of the product.


Only put the product into use in accordance with the information contained in the accompanying documents supplied.

According to the manufacturer (Otto Bock Healthcare Products GmbH), the patient is the operator of the product according to the IEC 60601-1:2005/A1:2012 standard.

2 Product description

2.1 Design

The product consists of the following components:

1. Cover/connection plate
2. Cover cap with charging receptacle
3. Battery
4. Ankle spring
5. Heel bracket
6. Ankle axis
7. Heel spring
8. Hydraulic unit
9. Toe axis
10. Toe plate
11. Carbon fibre frame
12. Main electronics
13. Spherical cap with pyramid

2.2 Function

This product features microprocessor-controlled damping of plantar flexion (motion of the foot in the ankle joint in the direction of the sole of the foot) and dorsiflexion (motion of the foot in the ankle joint in the direction of the instep).

The microprocessor uses the measurements of an integrated sensor system as a basis to control a hydraulic unit that influences the damping behaviour of the product.

These sensor data are updated and evaluated 100 times per second. As a result, the behaviour of the product is adapted to the current motion situation (gait phase) dynamically and in real time.

Thanks to the microprocessor-controlled damping of the plantar flexion and dorsiflexion, the product can be individually adapted to the needs of the patient.

For this purpose, the product is configured with the "M-Soft" adjustment software.

The product features MyModes for special motion types (e.g. golf ...). These have default values configured using the adjustment software and can be accessed via the Cockpit app or using a special motion pattern (see page 66).

Additionally, provided it is configured using the adjustment software, a so-called lock mode (additional mode "Ankle lock") can also be selected, locking the ankle joint of the prosthetic foot in its current position.

In case of a product malfunction, safety mode makes restricted operation possible. Resistance parameters that are predefined by the product are configured for this purpose (see page 71).

Empty battery mode permits safe walking when the battery is drained. Resistance parameters that are predefined by the product are configured for this purpose (see page 71).

The microprocessor-controlled hydraulic unit offers the following advantages

- Approximation of the physiological gait pattern
- Stable standing on even and inclined surfaces
- Adaptation of product characteristics to various surfaces, inclines, gait situations, walking speeds and heel heights

2.3 Combination possibilities

This product can be combined with the following Ottobock components:

Knee joints

- Genium: 3B1, 3B1=ST, 3B1-2, 3B1-2=ST, 3B1-3, 3B1-3=ST
- Genium X3: 3B5-X3, 3B5-X3=ST, 3B5-2, 3B5-2=ST, 3B5-3, 3B5-3=ST
- 3C98-2, 3C88-2 C-Leg knee joint
- 3C98-3, 3C88-3 C-Leg

Cosmetic cover/protector

- 4X860=* C-Leg Protective Cover (w/o shield)
- 4X880=* Genium Protective Cover
- 3S26 cosmetic foam cover

Footshells

- Translucent: 2C7=[prosthesis side][foot size]/1
- Beige: 2C7=[prosthesis side][foot size]/4
- Brown: 2C7=[prosthesis side][foot size]/15

[prosthesis side]: L=left, R=right

[foot size]: 24, 25, 26, 27, 28, 29

Order example: Footshell for left prosthetic foot in foot size 25 and beige colour

Article number: 2C7=L25/4

3 Intended use

3.1 Indications for use

The product is to be used **solely** for lower limb exoprosthetic fittings.


3.2 Conditions of use

The product was developed for everyday use and must not be used for unusual activities. These unusual activities include, for example, sports with excessive impact load (tennis, basketball, running, ...) or extreme sports (free climbing, paragliding, etc.).

Permissible ambient conditions are described in the technical data (see page 74).

The product is intended **exclusively** for use on **one** patient. Use of the product by another person is not approved by the manufacturer.

Our components perform optimally when paired with appropriate components based upon weight and mobility grades identifiable by our MOBIS classification information and which have appropriate modular connectors.

The product is recommended for mobility grade 2 (restricted outdoor walker) and mobility grade 3 (unrestricted outdoor walker).

Foot size [cm]	24 to 25	26 to 29
----------------	----------	----------

Max. body weight [kg]	100	125
-----------------------	-----	-----

3.3 Indications

- For users with unilateral knee disarticulation and unilateral transfemoral amputation
- For users with unilateral or bilateral transtibial amputation
- The user must fulfil the physical and mental requirements for perceiving optical/acoustic signals and/or mechanical vibrations

3.4 Contraindications

- All conditions which contradict or go beyond the specifications listed in the section on "Safety" and "Indications for use".

3.5 Qualification

The product may be fitted only by qualified personnel authorised by Ottobock after completing the corresponding training.

4 Safety

4.1 Explanation of warning symbols

 WARNING	Warning regarding possible serious risks of accident or injury.
 CAUTION	Warning regarding possible risks of accident or injury.
 NOTICE	Warning regarding possible technical damage.

4.2 Structure of the safety instructions

 WARNING
The heading describes the source and/or the type of hazard
The introduction describes the consequences in case of failure to observe the safety instructions. Consequences are presented as follows if more than one consequence is possible:
> E.g.: Consequence 1 in the event of failure to observe the hazard
> E.g.: Consequence 2 in the event of failure to observe the hazard
▶ This symbol identifies activities/actions that must be observed/carried out in order to avert the hazard.

4.3 General safety instructions

 WARNING
Non-observance of safety notices
Personal injury/damage to the product due to using the product in certain situations.
▶ Observe the safety notices and the stated precautions in this accompanying document.

 WARNING
Use of damaged power supply unit, adapter plug or battery charger
Risk of electric shock due to contact with exposed, live components.
▶ Do not open the power supply unit, adapter plug or battery charger.
▶ Do not expose the power supply unit, adapter plug or battery charger to extreme loading conditions.
▶ Immediately replace damaged power supply units, adapter plugs or battery chargers.

 CAUTION
Failure to observe warnings/error signals
Falling due to unexpected behaviour of the prosthesis because of changed damping behaviour.
▶ Observe the warnings/error signals and corresponding change in damping settings.

⚠ CAUTION

Independent manipulation of the product and the components

Falling due to breakage of load-bearing components or malfunction of the product.

- ▶ Manipulations to the product other than the tasks described in these instructions for use are not permitted.
- ▶ The battery may only be handled by authorised, qualified Ottobock personnel (no replacement by the user).
- ▶ The product and any damaged components may only be opened and repaired by authorised, qualified Ottobock personnel.

⚠ CAUTION

Mechanical stress on the product

- > Falling due to unexpected product behaviour as the result of a malfunction.
- > Falling due to breakage of load-bearing components.
- > Skin irritation due to defects on the hydraulic unit with leakage of liquid.
- ▶ Do not subject the product to mechanical vibrations or impacts.
- ▶ Check the product for visible damage before each use.

⚠ CAUTION

Use of the product when battery charge level is too low

Falling due to unexpected behaviour of the prosthesis because of changed damping behaviour.

- ▶ Check the current charge level before use and charge the prosthesis if required.
- ▶ Note that the operating time of the product may be reduced at low ambient temperatures or due to ageing of the battery.

⚠ CAUTION

Penetration of liquids into the product

Falling due to unexpected product behaviour as the result of a malfunction.

- ▶ When the footshell is intact, the product is protected against splashed water from all directions. However, it is not protected against submersion, jets of water and steam.
- ▶ If water has penetrated the product, remove the footshell and allow both components to dry. The product must be inspected by an authorised Ottobock Service Center.
- ▶ If salt water has penetrated the prosthesis, the footshell must be removed immediately. The product must be inspected by an authorised Ottobock Service Center.
- ▶ Do not use the product for bathing prostheses.

⚠ CAUTION

Overloading due to unusual activities

- > Falling due to unexpected product behaviour as the result of malfunction.
- > Falling due to breakage of load-bearing components.
- > Skin irritation due to defects on the hydraulic unit with leakage of liquid.
- ▶ The product was developed for everyday use and must not be used for unusual activities. These unusual activities include, for example, sports with excessive impact load (tennis, basketball, running, ...) or extreme sports (free climbing, paragliding, etc.).
- ▶ Careful handling of the product and its components not only increases their service life but, above all, ensures the patient's personal safety!
- ▶ If the product and its components have been subjected to extreme loads (e.g. due to a fall, etc.), then the product must be inspected for damage immediately. If necessary, forward the product to an authorised Ottobock Service Centre.

⚠ CAUTION

Mechanical stress during transport

- > Falling due to unexpected product behaviour as a result of a malfunction.
- > Falling due to breakage of load-bearing components.
- > Skin irritation due to defects on the hydraulic unit with leakage of liquid.
- ▶ Only use the transport packaging for transportation.

⚠ CAUTION**Signs of wear and tear on the product components**

Falling due to damage or malfunction of the product.

- ▶ Regular service inspections (maintenance) are mandatory in the interest of patient safety and in order to maintain operating reliability and protect the warranty.

⚠ CAUTION**Use of unapproved accessories**

> Falling due to product malfunction as a result of reduced interference resistance.

> Interference of other electronic devices due to increased emissions.

- ▶ Use the product only in combination with the accessories, signal converters and cables listed in the sections "Scope of delivery" (see page 54) and "Accessories" (see page 54).

NOTICE**Improper product care**

Damage to the product due to the use of incorrect cleaning agents.

- ▶ Only clean the product with a damp cloth and mild soap (e.g. 453H10=1-N Ottobock DermaClean).

NOTICE**Mechanical damage to the product**

Change in or loss of functionality due to damage.

- ▶ Use caution when working with the product.
- ▶ If the product is damaged, check it for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, do not continue using the product (see "Signs of changes in or loss of functionality during use" in this section).
- ▶ Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

INFORMATION

When exoprosthetic components are used, the use of control functions which are performed hydraulically or movement of the component in the footshell can cause movement noise. This kind of noise is normal and unavoidable. It generally does not cause any problems. If movement noise increases noticeably during the lifecycle of the component, the component should be inspected by an authorised Ottobock Service Centre immediately.

Signs of changes in or loss of functionality during use

Decreased forefoot resistance or changes in roll-over behaviour are noticeable indications of loss of functionality.

4.4 Information on the Power Supply/Battery Charging**⚠ CAUTION****Charging the prosthesis without taking it off**

> Falling due to walking and getting caught on the connected battery charger.

> Falling due to unexpected behaviour of the prosthesis because of changed damping behaviour.

- ▶ Instruct the patient that the prosthesis must be taken off before it is charged.

⚠ CAUTION**Charging the product with damaged power supply unit/charger/charger cable**

Falling due to unexpected behaviour of the product caused by insufficient charging.

- ▶ Check the power supply unit, charger and charger cable for damage before use.
- ▶ Replace any damaged power supply unit, charger or charger cable.

NOTICE

Use of incorrect power supply unit/battery charger

Damage to product due to incorrect voltage, current or polarity.

- ▶ Use only power supply units/battery chargers approved for this product by Ottobock (see instructions for use and catalogues).

4.5 Battery charger information

NOTICE

Penetration of dirt and humidity into the product

Lack of proper charging functionality due to malfunction.

- ▶ Ensure that neither solid particles nor liquids can penetrate into the product.

NOTICE

Mechanical stress on the power supply/battery charger

Lack of proper charging functionality due to malfunction.

- ▶ Do not subject the power supply/battery charger to mechanical vibrations or impacts.
- ▶ Check the power supply/battery charger for visible damage before each use.

NOTICE

Operating the power supply unit/charger outside of the permissible temperature range

Lack of proper charging functionality due to malfunction.

- ▶ Only use the power supply unit/charger for charging within the allowable temperature range. The section "Technical data" contains information on the allowable temperature range (see page 74).

NOTICE

Independent changes or modifications carried out to the battery charger

Lack of proper charging functionality due to malfunction.

- ▶ Have any changes or modifications carried out only by Ottobock authorised, qualified personnel.

4.6 Information on Alignment/Adjustment

⚠ CAUTION

Use of unsuitable prosthesis components

Falling due to unexpected behaviour of the product or breakage of load-bearing components.

- ▶ Use the product only in combination with components listed in the section "Combination possibilities" (see page 46).

⚠ CAUTION

Use of unapproved accessories

- > Falling due to product malfunction as a result of reduced resistance to interference.
- > Interference of other electronic devices due to increased emissions.

- ▶ Use the product only in combination with the accessories, signal converters and cables listed in the sections "Combination possibilities" (see page 46), "Scope of delivery" (see page 54) and "Accessories" (see page 54).

⚠ CAUTION

Incorrect alignment or assembly

Falling due to damage to the prosthesis components.

- ▶ Observe the alignment and assembly instructions.

⚠ CAUTION**Disconnecting/establishing the connection during the adjustment process with the adjustment software**

Falling due to unexpected behaviour of the product.

- ▶ When wearing the product, the patient must not remain unattended during the configuration process while connected to the adjustment software.
- ▶ Observe the maximum range of the Bluetooth connection.
- ▶ During the data transfer (PC to product), the patient must stand or sit without moving.
- ▶ If the connection fails while making adjustments, the prosthetist must immediately warn and secure the patient.
- ▶ The connection to the product must always be disconnected after adjustments have been completed.

⚠ CAUTION**Operator errors during the adjustment process with the adjustment software**

Falling due to unexpected behaviour of the product.

- ▶ Participation in an Ottobock product training course is mandatory prior to using the product. Additional product training courses may become necessary to qualify for software updates.
- ▶ Correctly entering the patient data (e.g. body weight) is an important criterion for the quality of the fitting.
- ▶ Always enter the weight in the specified unit.
- ▶ Ensure that the required dimensions are selected and entered correctly.
- ▶ If the patient uses walking aids (e.g. crutches or walking canes) during the adjustment process, readjustment is required as soon as the patient stops using these walking aids.

4.7 Information on Proximity to Certain Areas**⚠ CAUTION****Insufficient distance to HF communication devices (e.g. mobile phones, Bluetooth devices, WiFi devices)**

Falling due to unexpected behaviour of the product caused by interference with internal data communication.

- ▶ Therefore, keeping a minimum distance of 30 cm to HF communication devices is recommended.

⚠ CAUTION**Operating the product in very close proximity to other electronic devices**

Falling due to unexpected behaviour of the product caused by interference with internal data communication.

- ▶ Do not operate the product in the immediate vicinity of other electronic devices.
- ▶ Do not stack the product with other electronic devices during operation.
- ▶ If simultaneous operation cannot be avoided, monitor the product and verify proper use in the existing setup.

⚠ CAUTION**Proximity to sources of strong magnetic or electrical interference (e.g. theft prevention systems, metal detectors)**

Falling due to unexpected behaviour of the product caused by interference with internal data communication.

- ▶ Ensure that the patient is not in the vicinity of sources of strong magnetic and electrical interference during trial fitting (such as theft prevention systems, metal detectors...).
If this cannot be avoided, ensure at least that the patient has a safeguard when walking or standing (e.g. a handrail or the support of another person).
- ▶ In general, monitor the product for unexpected changes in the damping behaviour when electronic or magnetic devices are in the immediate vicinity.

⚠ CAUTION

Entering a room or area with strong magnetic fields (e.g. magnetic resonance tomographs, MRT (MRI) equipment...)

- > Falling due to unexpected restriction of the product's range of motion caused by metallic objects adhering to the magnetised components.
- > Irreparable damage to the product due to the effect of strong magnetic fields.
- ▶ Make sure that the patient takes off the product before entering the room or area and stores the product outside this room or area.
- ▶ Damage to the product caused by exposure to strong magnetic fields cannot be repaired.

⚠ CAUTION

Remaining in areas outside the allowable temperature range

Falling due to malfunction or the breakage of load-bearing product components.

- ▶ Ensure that the patient is not in areas outside the permissible temperature range (see page 74) during trial fitting.

4.8 Information on Use

⚠ CAUTION

Walking up stairs

Falling due to foot placed incorrectly on stair.

- ▶ Ensure that the patient uses the handrail when walking up stairs and places most of the sole of the foot on the stair surface. The toe plate may fold away if only the front part of the foot is placed on the edge of the stair.
- ▶ Inform the patient that particular caution is required when carrying children while walking up stairs.

⚠ CAUTION

Walking down stairs

Falling due to foot being placed incorrectly on stair as a result of changed damping behaviour.

- ▶ Ensure that the patient uses the handrail when walking down stairs and places most of the sole of the foot on the stair surface. Rolling the foot over the edge of the stair is not necessary.
- ▶ Inform the patient that particular caution is required when carrying children while walking down stairs.

⚠ CAUTION

High temperature of the hydraulic unit due to uninterrupted, increased activity (e.g. extended walking downhill)

Falling due to unexpected behaviour of the product caused by use in overheating mode.

- ▶ Be sure to pay attention when pulsating vibration signals start. They indicate the risk of overheating.
- ▶ As soon as these pulsating vibration signals begin, you must reduce your level of activity so the hydraulic unit can cool down.
- ▶ Please note that the range of motion in the ankle joint is reduced as the temperature increases, and this may ultimately result in the ankle joint being locked completely. Special caution is therefore required, especially when walking down stairs.
- ▶ Full activity may be resumed after the pulsating vibration signals stop.

⚠ CAUTION

Improper mode switching

Falling due to unexpected behaviour of the product because of changed damping behaviour.

- ▶ Ensure that the patient stands securely during all switching processes.
- ▶ Inform the patient that the changed damping characteristics have to be verified after switching and feedback from the acoustic signal emitter must be observed.
- ▶ Switching back to basic mode is mandatory once the activities in MyMode have been completed.
- ▶ If required, take the weight off the product and correct the switching.

⚠ CAUTION**Using prosthetic foot without footshell**

Falling due to slipping when walking on smooth floor surfaces (tiles).

- ▶ Do not use the prosthetic foot without the accompanying footshell.

⚠ CAUTION**Using prosthetic foot with a damaged foot shell**

- > Falling due to unexpected product behaviour as the result of malfunction.
- > Falling due to breakage of load-bearing components.

- ▶ Do not use the prosthetic foot when the foot shell is damaged. Replace damaged foot shells immediately, before using the product again.

4.9 Notes on the safety modes**⚠ CAUTION****Using the product in safety mode**

Falling due to unexpected product behaviour because of changed damping behaviour.

- ▶ The warnings/error signals (see page 77) have to be observed.

⚠ CAUTION**Safety mode cannot be activated due to malfunction caused by water penetration or mechanical damage**

Falling due to unexpected product behaviour because of changed damping behaviour.

- ▶ Using the product when it is defective is prohibited.
- ▶ The product must be inspected by an authorised Ottobock Service Centre.

⚠ CAUTION**Safety mode cannot be deactivated**

Falling due to unexpected product behaviour because of changed damping behaviour.

- ▶ If safety mode cannot be deactivated by recharging the battery, a permanent error has occurred.
- ▶ Using the product when it is defective is prohibited.
- ▶ The product must be inspected by an authorised Ottobock Service Centre.

⚠ CAUTION**Safety signal occurs (ongoing vibration)**

Falling due to unexpected product behaviour because of changed damping behaviour.

- ▶ The warnings/error signals (see page 77) have to be observed.
- ▶ After the safety signal has been emitted, further use of the product is prohibited.
- ▶ The product must be inspected by an authorised Ottobock Service Centre.

4.10 Information on the use of a mobile device with the cockpit app**⚠ CAUTION****Improper use of the mobile device**

Falling due to altered damping behaviour as a result of unexpected switching into a MyMode.

- ▶ Use the instructions for use (user) to instruct the patient on the proper handling of the mobile device with the Cockpit app.

⚠ CAUTION**Independently applied changes or modifications made to the mobile device**

Falling due to altered damping behaviour as a result of unexpected switching to a MyMode.

- ▶ Do not make any independent changes to the hardware of the mobile device on which the app is installed.
- ▶ Do not make any independent changes to the software/firmware of the mobile device that are not included in the update function of the software/firmware.

⚠ CAUTION

Improper mode switching with the mobile device

Falling due to unexpected product behaviour because of changed damping behaviour.

- ▶ Ensure that the patient stands securely during all switching processes.
- ▶ Inform the patient that the changed damping characteristics have to be verified after switching, and feedback from the acoustic signal emitter and the mobile device display must be observed.
- ▶ Switching back to basic mode is mandatory once the activities in the MyMode have been completed.

NOTICE

Failure to observe the system requirements for the installation of the Cockpit app

Mobile device malfunction.

- ▶ The Cockpit App should only be installed on mobile devices and versions which comply with the specifications in the respective online stores (e.g. Apple App Store, Google Play Store, ...)

INFORMATION

The illustrations in these instructions for use are only examples and may deviate from the respective mobile device being used and the version.

5 Scope of Delivery and Accessories

5.1 Scope of delivery

- 1 pc. 1B1-2 Meridium
- 1 pc. 757L16-4 power supply
- 1 pc. 4E50-2 C-Leg battery charger
- 1 pc. 646C107 Bluetooth PIN card
- 1 pc. 647F542 prosthesis passport
- 1 pc. Instructions for use (qualified personnel)
- 1 pc. Instructions for use (user)
- 1 pc. 2C7 footshell inc. instructions for use
- 1 pc. 4G872=* set of cover caps
- 1 pc. 2C101 footshell replacement tool

Cockpit app for download from the website: <https://www.ottobock.com/cockpitapp>

- "Cockpit 4X441-IOS=V*" app for iOS
- "Cockpit 4X441-ANDR=V*" app for Android

5.2 Accessories

The following components are not included in the scope of delivery and may be ordered separately:

- 60X5 BionicLink PC
- 1 pc. 757P48 Y adapter cable
This is for charging more than one product at the same time (e.g. 1B1-2; 1B1; 3B1/3B1=ST; 3B1-2/3B1-2=ST; 3B5-X3/3B5-X3=ST; 3C98-2/3C88-2; 3C98-3/3C88-3; 3C96-1/3C86-1) with the 757L16-4 power supply.
- 704G30 clamping tool


6 Charging the battery

The following points must be observed when charging the battery:

- Use the 757L16-4 power supply and 4E50-2 battery charger to charge the battery.
- The capacity of a fully charged battery is sufficient for one full day.
- We recommend charging the product once a day when used by the patient on a daily basis.
- For the maximum operating time with one battery charge, disconnecting the battery charger from the product only immediately before using the product is recommended.

- The battery should be charged until the yellow LED on the battery charger turns off prior to initial use, and for at least 4 hours. This calibrates the charge level indicator via the Cockpit app and by turning over the prosthesis. If the battery charger is disconnected from the prosthesis too soon, the charge level indicator via the Cockpit app and by turning over the prosthesis may not correspond to the actual charge level.
- The ankle joint of the prosthetic foot is locked during the charging process.
- The battery may discharge while the product is not being used.

6.1 Connecting the power supply and battery charger

- 1) Slide the country-specific plug adapter onto the power supply until it locks into place (see fig. 1).
- 2) Connect the round, **four-pin** plug of the charging cable to the **OUT** receptacle on the battery charger so that the plug locks into place (see fig. 2).
INFORMATION: Ensure correct polarity (guide lug). Do not use force when connecting the cable plug to the battery charger.
- 3) Connect the round, **three-pin** plug of the power supply to the **12 V** receptacle on the battery charger so that the plug locks into place (see fig. 2).
INFORMATION: Ensure correct polarity (guide lug). Do not use force when connecting the cable plug to the battery charger.
- 4) Plug the power supply unit into the wall socket.
→ The green LED on the back of the power supply and the green LED on the battery charger light up (see fig. 3).
→ If the green LED on the power supply and the green LED on the battery charger do not light up, there is an error (see page 77).

6.2 Charging the prosthesis battery

- 1) Open the charging receptacle cover.
- 2) Connect the charging plug to the charging receptacle of the product.
INFORMATION: Take note of the plug direction!
→ The correct connection of the battery charger to the product is indicated by feedback (see page 79).
- 3) The charging process starts.
→ Once the product battery is fully charged, the yellow LED on the battery charger turns off.
- 4) Disconnect the product after the charging process is complete.
→ This is followed by an electronics self-test which is confirmed by feedback (see page 79).
- 5) Close the charging receptacle cover.

6.3 Display of the current charge level

INFORMATION

The charge level cannot be displayed during the charging process.

- 1) Turn the prosthesis by 180° (the sole of the foot must face up).
- 2) Hold the prosthesis still and wait for beep signals.

Prosthetic foot with knee joint:

The beep signal for the knee joint sounds after approx. 2 seconds.

The beep signal for the prosthetic foot sounds after approx. 4 seconds.

Prosthetic foot without knee joint:

The beep signal for the prosthetic foot sounds after approx. 2 seconds.

Beep signal	Vibration signal	Battery charge level
5x short	–	More than 80%
4x short	–	66% to 80%
3x short	–	51% to 65%
2x short	–	36% to 50%
1x short	3x long	20% to 35%
1x short	5x long	Less than 20%

INFORMATION
 If the **Volume** parameter is set to '0' in the Cockpit app, there are no beep signals (see page 68).

Display of the current charge level using the cockpit app:

Once the Cockpit app has been started, the current charge level is displayed in the bottom line of the screen:

1. 38% – Charge level of battery for currently connected component

7 Preparing the product for use

7.1 Applying/removing footshell

For applying/removing the footshell, please consult the instructions for use included with the footshell.

CAUTION
Incorrect assembly/disassembly and improper use of footshell

- > Falling due to unexpected product behaviour as the result of malfunction.
- > Falling due to breakage of load-bearing components.
- ▶ Use the prosthetic foot only in combination with the footshell.
- ▶ Remove the footshell from the prosthetic foot only when necessary.
- ▶ Only use the 2C101 changing device and the 704G30 clamping device for mounting/removing.
- ▶ Replace worn footshells. Do not use the prosthetic foot when the footshell is damaged.

INFORMATION
 Before mounting the component on the tube adapter, make a note of the serial number of the component. It is found next to the pyramid on the spherical cap (see page 45).
 The serial number is required to connect to the adjustment software and the Cockpit app and must be entered in the prosthesis passport.

7.2 Alignment

7.2.1 Settings with the "M-Soft" adjustment software

7.2.1.1 Introduction

The "M-Soft" adjustment software makes it possible to optimise the product settings for a patient. The adjustment software provides step-by-step guidance through the adjustment process. After the settings are configured, the data for them can be saved and printed for documentation. These data can be retrieved if required and imported into the product.

Please consult the integrated online help in the adjustment software for further information.

Updating the M-Soft adjustment software

- 1) Click "**Help > About**" in the menu bar of the Data Station when you are connected to the Internet.
 - The window opens with the versions of the previously installed programs and the manufacturer's address.
- 2) Click the "**Check for updates**" button in this window.
 - A search for updates of previously installed software products and components is performed via the Internet.
- 3) If updates are available, click "**Download**" in the column on the right in order to download and save the update.
- 4) Extract the "ZIP file" and execute it.

INFORMATION

Cybersecurity

- ▶ Keep your operating system up to date and always install any available security updates.
- ▶ Protect your computer from unauthorised access (e.g., by using virus scans, password protection etc.).
- ▶ Do not use unsecured networks.
- ▶ Please contact the manufacturer if you suspect cybersecurity problems.

7.2.1.2 Data transfer between the product and the PC

Product settings using the adjustment software can only be made via Bluetooth data transfer. For this purpose, a Bluetooth wireless connection must be established between the product and the PC using the "60X5 BionicLink PC" Bluetooth adapter. The installation and use of the "60X5 BionicLink PC" adapter are described in the instructions for use supplied with the adapter.

7.2.1.3 Preparing the product to connect to the adjustment software

If the product does not emit any signals when querying the charge level (Display of battery charge level without additional devices), the battery is drained or the product is switched off.

Switching on the product

- 1) Connect the power supply and battery charger to the wall socket.
 - 2) Connect the battery charger to the product.
 - 3) Wait for feedback signals.
 - 4) Disconnect the battery charger from the product.
- After feedback signals are emitted (self-test), the product is switched on.

Switching on Bluetooth

Upon delivery, the Bluetooth function of the prosthesis is switched on.

The Bluetooth function can be switched off using the Cockpit app or the adjustment software. When the Bluetooth function is switched off, it is only turned on temporarily for 2 minutes after connecting/disconnecting the battery charger and is then turned off again automatically. When a connection with the PC is active (the  icon is lit up), the Bluetooth function is not switched off automatically.

7.2.2 Bench alignment in alignment apparatus

With correct bench alignment, e.g. in the PROS.A. assembly alignment apparatus (743A200), the benefits of the product are realised to best advantage. If the L.A.S.A.R. assembly alignment apparatus (743L200) is available, it can be used as well.

Alignment can also be carried out using LaserLine/plumb line.

The following points must be observed during alignment:

- The alignment recommendations provided must be observed in order for the prosthetic foot to function correctly.

- Static alignment in the **alignment apparatus** must always be carried out **without shoes**, otherwise a correct adjustment will not be possible.
- A marking is located on the exterior distal side of the footshell. The marking serves as orientation for the alignment point on the foot.
- Before beginning the alignment, the prosthetic foot must be put into alignment mode using the adjustment software ("Alignment" tab, "Alignment Recommendation" tab). Only in alignment mode will the prosthetic foot be locked in neutral position and allow static alignment to be carried out correctly.

Foot size (cm)	Mid-foot anterior to the alignment reference line	Heel height
24 – 29	30 mm	0 mm

Bench alignment for transtibial (TT) modular lower limb prostheses

The adjustment software must be used to provide information for alignment of TT modular lower limb prostheses (also see TT modular lower limb prostheses: 646F336).

Bench alignment for transfemoral (TF) modular lower limb prostheses

The alignment recommendations for the appropriate Ottobock knee joint must be observed for TF modular lower limb prostheses (also see TF modular prostheses: 646F219).

7.2.3 Static alignment optimisation

- Use the L.A.S.A.R. Posture for static alignment optimisation of the prosthesis, if available.
- Use the alignment recommendations of the manufacturer (TF modular lower limb prostheses: **646F219**, TT modular lower limb prostheses: **646F336**).

7.2.4 Dynamic alignment optimisation

Make sure the heel height is correctly set before beginning the dynamic alignment optimisation!

- Dynamic alignment optimisation is not necessary if the static situation is satisfactory. The adjustment software is used to make adjustments to ensure correct heel contact, smooth rollover and optimum weight transfer on the contralateral side.
- Adjust the prosthesis in the frontal plane (ML) by using the medial-lateral screws to change its angle or shift it (see 646F336) so as to minimise any axial run-out in the knee.

7.2.5 Mounting cover plate/connecting plate/cover cap

Depending on the cosmetic that is being used (cosmetic foam cover, Protector), the suitable element must be selected from the set of cover caps and mounted on the product.

7.2.5.1 Mounting the cover cap

The cover cap forms the top of the footshell.


- ▶ Mount the cover cap according to the instructions for use of the footshell.


7.2.5.2 Mounting the cover plate

If the cover plate is being used, there should be no cosmetic (Protector, cosmetic foam cover) on the product.

- ▶ Fit the cover plate on the cover cap which has already been mounted.

7.2.5.3 Mounting the Protector connecting plate


The Protective Cover connecting plate acts as the connection between the prosthetic foot and the Protective Cover.

- 1) Apply adhesive to the cover cap in 4 places (adhesive areas about 15 mm long).
- 2) Fit the Protective Cover connecting plate to the cover cap which has already been mounted.
- 3) Fasten the clips of the foot cuff to the connecting plate of the Protective Cover.

INFORMATION: Please note the 647G1113/647G942 instructions for use.

7.2.5.4 Mounting the connecting plate for the cosmetic foam cover

The connecting plate for the cosmetic foam cover acts as the connection between the prosthetic foot and the cosmetic foam cover.

- 1) Fit the connecting plate for the cosmetic foam cover to the cover cap which has already been mounted.
- 2) Trim the cosmetic foam cover to fit the connecting plate for the cosmetic foam cover.
- 3) Apply a layer of 636N9 contact adhesive to the contact surfaces of the cosmetic foam cover.
- 4) Place the cosmetic foam cover on the connecting plate for the cosmetic foam cover.

8 Cockpit app

The cockpit app enables switching from basic mode into the pre-configured MyModes. In addition, information about the product (step counter, battery charge level, etc.) can be called up.

The behaviour of the product can be changed to a certain extent on a day-to-day basis using the app (e.g. while becoming accustomed to the product). The adjustment software can be used to trace the change at the patient's next appointment.

Information on the Cockpit app

- The Cockpit app can be downloaded free of charge from the respective online store. For more information, please visit the following website: <https://www.ottobock.com/cockpitapp>. To download the Cockpit app, the QR code on the supplied Bluetooth PIN card can also be read with the mobile device (requirement: QR code reader and camera).
- The language of the Cockpit app user interface can only be changed using the adjustment software.
- Depending on the version of the Cockpit app being used, the language of the user interface in the Cockpit app corresponds to the language of the mobile device on which the Cockpit app is being used.
- The serial number of the component to be connected has to be registered with Ottobock the first time it is connected. If the registration is not accepted, use of the Cockpit app for this component will be limited.
- Bluetooth on the prosthesis must be turned on in order to use the Cockpit app.
If Bluetooth is switched off, it can be turned on by turning the prosthesis upside-down (sole of the foot must point up) or by connecting/disconnecting the battery charger. Bluetooth is then turned on for approx. 2 minutes. During this time, the app must be started and used to establish a connection. If required, Bluetooth on the prosthesis can be switched on permanently afterwards (see page 70).
- Keep the mobile app up to date at all times.
- Please contact the manufacturer if you suspect cybersecurity problems.

8.1 Initial connection between cockpit app and prosthesis

The following points need to be observed before establishing the connection:

- Bluetooth of the component must be switched on (see page 70).

- Bluetooth on the mobile device must be switched on.
- The mobile device must not be in "flight mode" (offline mode), otherwise all wireless connections are turned off.
- **The mobile device must be connected to the Internet.**
- The serial number and Bluetooth PIN of the component being connected must be known. They are found on the enclosed Bluetooth PIN card. The serial number begins with the letters "SN".


INFORMATION
 If the Bluetooth PIN card with the Bluetooth PIN and serial number of the component is lost, the Bluetooth PIN can be determined using the adjustment software.


8.1.1 Starting the cockpit app for the first time

- 1) Tap the symbol of the Cockpit app ().
 → The end user license agreement (EULA) is displayed.
 - 2) Accept the end user license agreement (EULA) by tapping the **Accept** button. If the end user license agreement (EULA) is not accepted, the Cockpit app cannot be used.
 → The welcome screen appears.
 - 3) Hold the prosthesis with the sole of the foot facing up, or connect and then disconnect the battery charger, in order to activate recognition (visibility) of the Bluetooth connection for 2 minutes.
 - 4) Tap the **Add component** button.
 → The Connection Wizard opens and guides you through the process of establishing a connection.
 - 5) Follow the subsequent instructions on the screen.
 - 6) After the Bluetooth PIN is entered, a connection to the component is established.
 → While the connection is being established, 3 beep signals sound and the  symbol appears.
 The  symbol is displayed when the connection has been established.
- Once the connection has been established, the data are read from the component. This process may take up to a minute.
 Then the main menu appears with the name of the connected component.


INFORMATION
 After the initial connection to the component has been established successfully, the app will connect automatically each time it is started. No further steps are required.

8.2 Control elements for cockpit app

1.  Access the navigation menu (see page 61)
2. Product
 The component name can only be changed with the adjustment software.
3. If connections to more than one component have been saved, you can switch between the saved components by tapping the **change** item (see page 61).
4. MyModes configured with the adjustment software.
 Switch the mode by tapping the corresponding symbol and confirming by tapping "OK".
5. Currently selected mode
6. Charge level of the component.
 Component battery fully charged
 Component battery empty
 Component battery charging
 The current charge level is also displayed in %.
7. Display and name of the currently selected mode (e.g. **1. Basic Mode**)
8.  Connection to component has been established
 Connection to component has been interrupted. The app is attempting to re-establish the connection automatically.
 No existing connection to the component.

8.2.1 Cockpit app navigation menu

Tap the ☰ symbol in the menus to display the navigation menu. Additional settings for the connected component can be configured in this menu.

Product

Name of the connected component

MyModes

Return to the main menu to switch MyModes

Heel Height

Setting the heel height (see page 63)

Functions

Access additional functions of the component (e.g. turn off Bluetooth) (see page 70)

Settings

Change settings of the selected mode (see page 68)

Status

Query status of the connected component (see page 71)

Manage components

Add or delete components (see page 61)

Imprint/Info

Display information/legal notices for the Cockpit app

8.3 Managing components

Connections with up to four different components can be stored in the app. However, a component can only be connected to one mobile device at a time.

INFORMATION

Before establishing the connection, observe the points in the section "Initial connection between Cockpit app and component" (see page 59).

8.3.1 Adding component

- 1) Tap the ☰ symbol in the main menu.
 - The navigation menu opens.
- 2) In the navigation menu, tap the "**Manage components**" item.
- 3) Hold the prosthesis with the sole of the foot facing up, or connect and then disconnect the battery charger, in order to activate recognition (visibility) of the Bluetooth connection for 2 minutes.
- 4) Tap the "+" button.
 - The Connection Wizard opens and guides you through the process of establishing a connection.
- 5) Follow the subsequent instructions on the screen.
- 6) After the Bluetooth PIN is entered, a connection to the component is established.
 - While the connection is being established, three beep signals sound and the (📶) symbol appears. The (📶) symbol is displayed when the connection has been established.
 - Once the connection has been established, the data are read from the component. This process may take up to a minute. The main menu will then appear with the name of the connected component.

INFORMATION

If establishing a connection to a component is not possible, perform the following steps:

- ▶ Delete the component from the Cockpit app if applicable (see the section "Deleting a component")
- ▶ Add the component again in the Cockpit app (see the section "Adding a component")

INFORMATION

After activating the "visibility" of the component (holding the component with the sole of the foot facing up, or connecting and then disconnecting the battery charger), the component can be recognised by another device (e.g. smartphone) within 2 minutes. If registration or establishing the connection takes too long, the process of establishing a connection is cancelled. In this case, hold the component with the sole of the foot facing up again, or connect and then disconnect the battery charger.


8.3.2 Deleting a component

- 1) Tap the ☰ symbol in the main menu.
→ The navigation menu opens.
- 2) In the navigation menu, tap the entry "**Manage components**".
- 3) Tap the "**Edit**" button.
- 4) Tap the 🗑 symbol under the component you want to delete.
→ The component is deleted.

8.3.3 Connecting component with multiple mobile devices

The connection for a component can be stored on more than one mobile device. However, only one mobile device can be connected to the component at one time.

If there is an existing connection between the component and a different mobile device, the following information appears while the connection is being established with the current mobile device:

- ▶ Tap the "**OK**" button.
- The connection to the last connected mobile device is broken off and established with the current mobile device.

9 Use**9.1 Setting the heel height**


The heel height of the prosthesis must be set on a level surface. A sloped ground surface distorts the heel height measurement and leads to improper control of damping behaviour.

A heel height that is too high may cause control of the prosthetic foot to function incorrectly because of too little motion in the ankle joint. This occurs especially in case of small feet, heels that are moved forward, when walking down stairs and ramps, and standing on surfaces that are sloping downwards. The maximum heel height in the section "Technical data" must therefore be taken into account (see page 74).

9.1.1 Setting heel height using motion pattern

- 1) Put on shoes with the new heel height.
 - 2) Stretch out foot with the prosthetic foot sideways.
 - 3) Swing the foot to the side three times.
→ A beep signal will sound to confirm that the motion pattern has been recognised.
 - 4) Put feet at the same level and make sure the heel and toes are resting on the ground surface.
 - 5) Distribute weight on feet evenly.
→ A confirmation signal will sound to indicate that the new heel height was successfully stored.
- INFORMATION: If there is no feedback (e.g. beep signal), this indicates that the new heel height could not be stored. Repeat the heel height measurement.**


9.1.2 Setting the heel height with the cockpit app

- 1) Once the component is connected and in the desired mode, tap the ☰ symbol in the main menu.
→ The navigation menu opens.
- 2) Tap the "**Heel Height**" menu option.
- 3) Follow the on-screen instructions.
- 4) Tap the "**Set the heel height**" option.
- 5) Follow the subsequent instructions on the screen.

9.2 Movement patterns in basic mode (mode 1)

9.2.1 Standing


The intuitive stance recognises any situation where the prosthesis is kept still while standing. The standing position is stabilised through high dorsiflexion damping with vertical lower leg. Plantar flexion is only dampened to a minor extent to allow the lower leg to be retracted to the neutral position for correction of the standing position.

This function is deactivated during rollover to the front or by lifting the prosthesis off the floor.

Coming to a stop with the prosthesis side after walking can lead to the knee joint sinking in as a result of the ankle position during the rollover. To restore a stable standing position, place the leg back under the body and stretch the leg or put weight on the heel.

The relief function can be used when standing (see page 66).

9.2.2 Walking

Initial attempts at walking with the prosthetic foot always require instruction by trained, qualified personnel.

When walking, the damping of the dorsiflexion and plantar flexion is adapted to the current walking phase, making a physiological gait possible.


Plantar flexion damping is increased during heel strike in order to provide support for stance phase flexion in the knee.

By increasing dorsiflexion damping in the stance phase, the lower leg is led into an extended rollover angle. The rollover behaviour automatically adjusts itself to the walking speed.

Plantar flexion damping is increased during the transition to the swing phase in order to prevent the toes from dropping and to maintain ground clearance.

During the swing phase, plantar flexion damping is adjusted to the current position of the lower leg at all times. This ensures a comfortable tread, with suitable heel leverage for the respective step length.

9.2.3 Sitting down/sitting


Sitting down

- 1) Place both feet side by side at the same level.
- 2) While sitting down, distribute weight evenly between both legs and use armrests, if available.
- 3) Move the buttocks in the direction of the backrest and lean the upper body forward.

Sitting

The relief function can be used when sitting; the tip of the foot lowers to reach a more natural foot position. (see page 66).

9.2.4 Standing up


- 1) Place the feet at the same level. Make sure the foot is positioned perpendicularly under the knee or shifted further forward, and that weight is distributed evenly over the feet.

INFORMATION: The ankle joint may lock if the prosthetic foot is set further back than at a perpendicular angle under the knee.

- 2) Lean the upper body forward.
- 3) Put the hands on armrests, if available.
- 4) Stand up with support from the hands, while keeping weight evenly distributed on the feet.


9.2.5 Walking up stairs

The position is stabilised through high dorsiflexion damping with vertical lower leg. Walking up stairs step-over-step is possible, depending upon the type of fitting.

Always hold on to the handrail with one hand when walking up steps.

9.2.6 Walking down stairs

This function must be practised and executed consciously. Only by properly stepping down with the sole can the system switch correctly and permit controlled rollover. The motion must be carried out in a continuous pattern in order to allow the motion to proceed in a fluid manner.

A stair function can be activated by using the adjustment software. See the following section for further information regarding the stair function.

- 1) Hold on to the handrail with one hand.
- 2) Position the leg with the prosthetic foot on the step so that as much of the sole of the foot as possible is on the step.


INFORMATION: Rolling the foot over the edge of the stair is not necessary.

- 3) Place the contralateral side onto the next step.
While doing so, check whether the knee joint and prosthetic foot will allow this motion.
- 4) Place the leg with the prosthetic foot on one step beyond the next step.
- 5) At the end of the stairs, take a bigger step when transitioning to level ground, so that the prosthetic foot switches over correctly from walking down the stairs to normal walking phase.

9.2.6.1 Stair function

The stair function increases the rollover angle when walking down stairs. This function should be switched on when walking down stairs step-over-step. If walking down stairs step-over-step is not required, this function can be switched off. Further information on switching on/off see page 69.


9.2.7 Walking up a ramp

With the first step, the foot is already adjusted to the inclination of the ramp, making a rollover possible if the step is on the heel or middle of the foot. In order to do so, the lower leg should be positioned almost perpendicular to the surface of the ramp and the full area of the foot should be placed on the surface.

When stepping on the forefoot with the lower leg at a steep angle (e.g. on very steep ramps), the foot secures dorsiflexion, which enables the body to be lifted in a stable manner.

9.2.8 Walking down a ramp


With the first step, the foot is already adjusted to the inclination of the ramp, making it possible for the heel to strike with extended plantar flexion so the complete area of the foot is on the surface during the rollover.

After the foot is put down on the ramp, you should not use the knee to work against it, but should allow flexion in the knee joint for the heel strike (yielding). Only in this way will the foot be able to recognise the motion as walking and allow an extended rollover. This allows for a guided lowering of the body's centre of gravity.


Plantar flexion is limited for walking with a prosthetic knee joint (higher amputation level than transtibial amputation) in order to support flexion in the knee joint during heel strike (yielding).

9.2.9 Walking backwards

The foot enables plantar flexion from the stance phase when walking backwards. When the toe is subsequently put down, the ankle joint drops to the neutral position in the dorsiflexion direction.

9.2.10 Standing on inclined surface

Standing on an inclined ground surface is no different to standing on a level surface. The foot stabilises in dorsiflexion when the lower leg is perpendicular. Put weight on the heel in order to lower the forefoot (e.g. when standing in a downward direction).

Perform one of the following motions in order to continue to walk on a downward-sloping ground surface from the standing position:

- Start the first step with the prosthesis side.
- Deliberately produce a rollover motion with the prosthesis side.

The prosthetic foot then yields in dorsiflexion, in order to allow the body's centre of gravity to drop before the heel of the other leg strikes.

The relief function can be used when standing on an inclined surface (see page 66).


Wearing shoes with heels limits the range of incline, so it may not be possible to achieve a perpendicular lower leg.

9.2.11 Kneeling

If the leg with the joint is tilted backwards, plantar flexion damping is reduced, making it possible for the foot to bend so the lower leg can lie flatter to the ground surface.

9.2.12 Relief function





When weight is evenly placed on the heel without movement for more than 2 seconds, the tip of the foot lowers to reach a more natural foot position.

Possible applications include sitting with the heel in front of the knee axis, standing while leaning against something, and standing on a downward sloping surface.

9.3 MyModes

In addition to basic mode (mode 1), MyModes can be activated and configured with the adjustment software. They can be called up by the patient using the Cockpit app or movement patterns. Switching by using movement patterns has to be activated in the adjustment software.

These modes are intended for specific motion patterns or postures (e.g. golfing, ...). Default settings for these movement patterns and postures can be accessed and individually adapted using the adjustment software.

Settings can also be adjusted by the patient using the Cockpit app (see page 69).

9.3.1 Switching MyModes with the cockpit app

INFORMATION


Bluetooth on the prosthesis must be turned on in order to use the Cockpit app.

If Bluetooth is switched off, it can be turned on by turning the prosthesis upside-down (function only available in basic mode) or by connecting/disconnecting the battery charger. Bluetooth is then turned on for approx. 2 minutes. During this time, the app must be started and used to establish a connection. If required, Bluetooth on the prosthesis can be switched on permanently afterwards (see page 70).

INFORMATION

If the **Volume** parameter is set to '0' in the Cockpit app, there are no beep signals (see page 68).

Once a connection to a prosthesis has been established, the cockpit app can be used to switch between the MyModes.

- 1) Tap the symbol of the MyMode (1) you want in the main menu of the app.
→ A security question for changing the MyMode appears.
- 2) If you want to change the mode, tap the "OK" button.
→ A beep signal sounds to confirm the switch.
- 3) After switching, a symbol (2) is displayed to identify the active mode.
→ The current mode is also indicated by the name on the lower edge of the screen (3).

9.3.2 Switching MyModes using motion patterns

Information on switching

- Switching and the number of motion patterns must be activated in the adjustment software.
- Before further activities, always check whether the selected mode corresponds to the required motion type.
- If the **Volume** parameter is set to "0" in the Cockpit app, there are no beep signals (see page 68).

Switching process

- 1) Hold the prosthetic leg under the body.
- 2) Using the heel of the prosthetic foot, knock on a solid obstacle behind you (e.g. a wall) according to the number of the configured MyMode (MyMode 1 = three times, MyMode 2 = four times, MyMode 3 = five times). You can also knock against the tip of the shoe on the contralateral leg.
→ A beep and vibration signal will occur to confirm that the movement pattern has been recognised.
INFORMATION: If this beep and vibration signal is not emitted, the knocking was not recognised.
- 3) Lean the prosthetic foot slightly to the back and place weight on the forefoot.
INFORMATION: If the prosthetic foot is strongly dorsiflexed, weight can be placed on the heel.
→ A confirmation signal will sound to indicate that the foot successfully switched to the corresponding mode (two times = MyMode 1, three times = MyMode 2, four times = MyMode 3).
INFORMATION: If this confirmation signal does not occur, weight was not placed on the prosthetic foot correctly or not for enough time. Repeat the process for correct switching.
- 4) Remove the load from the prosthetic leg.
→ The mode has been changed.

9.3.3 Switching on ankle lock

Information on switching

- The ankle lock must be activated as MyMode "Ankle lock" and the number of motion patterns used to switch it on must be activated in the adjustment software.
- Before further activities, always check whether the selected mode corresponds to the required motion type.
- If the **Volume** parameter is set to "0" in the Cockpit app, there are no beep signals (see page 68).

Switching process

- 1) Hold the prosthetic leg under the body.
- 2) Using the heel of the prosthetic foot, knock on a solid obstacle behind you (e.g. a wall) according to the number of the configured MyMode (MyMode 1 = 3 times, MyMode 2 = 4 times, MyMode 3 = 5 times). It is also possible to knock against the tip of the shoe of the contralateral leg.
→ A beep and vibration signal will occur to confirm that the movement pattern has been recognised.
- 3) Lean the prosthetic foot slightly to the back and place weight on the forefoot.
INFORMATION: If the prosthetic foot is strongly dorsiflexed, weight can be placed on the heel.

→ A confirmation signal will sound to indicate that the prosthesis has successfully switched to the corresponding mode (2 times = MyMode 1, 3 times = MyMode 2, 4 times = MyMode 3).

INFORMATION: If this confirmation signal does not occur, weight was not placed on the prosthetic foot correctly or not for enough time. Repeat the process to correctly switch to the required mode.

- 4) Remove the load from the prosthetic leg.
 - The mode has been changed.
- 5) Within 2 seconds, lower the prosthetic leg and assume the required position of the ankle angle.
 - At the end of this time, a signal will occur to indicate the locking of the ankle joint.

9.3.4 Switching from a MyMode back to basic mode

Information on switching

- Regardless of the configuration of additional MyModes in the adjustment software, it is always possible to switch back to basic mode (mode 1) with a motion pattern.
- It is always possible to switch back to basic mode (mode 1) by connecting/disconnecting the battery charger.
- Before further activities, always check whether the selected mode corresponds to the required motion type.
- If the **Volume** parameter is set to "0" in the Cockpit app, there are no beep signals (see page 68).

Switching process

- 1) Hold the prosthetic leg under the body.
- 2) Using the heel of the prosthetic foot, knock on a solid obstacle behind you at least 3 times but no more than 5 times.
 - A beep and vibration signal will occur to confirm that the motion pattern has been recognised.
- 3) Lean the prosthetic foot slightly to the back and place weight on the forefoot.

INFORMATION: If the prosthetic foot is strongly dorsiflexed, weight can be placed on the heel.

 - A confirmation signal will sound to indicate that the prosthesis has successfully switched over to basic mode.
 - INFORMATION: If this confirmation signal does not occur, weight was not placed on the prosthetic foot correctly or not for enough time. Repeat the process to correctly switch to the required mode.**
- 4) Remove the load from the prosthetic leg.
 - The mode has been changed.
 - Before further activities, always check whether the selected mode corresponds to the required motion type.

9.4 Changing prosthesis settings

Once an active connection to a component has been established, the settings **of the respective active mode** can be changed using the Cockpit app.


INFORMATION

Bluetooth on the prosthesis must be switched on to change the prosthesis settings. If Bluetooth is switched off, it can be turned on by turning the prosthesis upside-down or by connecting/disconnecting the battery charger. Bluetooth is then turned on for approx. 2 minutes. The connection must be established during this period.

Information for changing the prosthesis settings

- Before changing settings, always check the main menu of the Cockpit app to make sure the correct component has been selected. Otherwise parameters could be changed for the wrong component.
- It is not possible to change prosthesis settings nor to switch to a different mode while the prosthesis battery is being charged. Only the status of the prosthesis can be called up. Instead of the  symbol, the  symbol appears in the bottom row of the screen in the cockpit app.
- The O&P professional's setting is in the middle of the scale. After making adjustments, this setting can be restored by tapping the "**Standard**" button in the Cockpit app.
- Prosthesis settings should be optimised using the adjustment software. The Cockpit app is not intended for use by the O&P professional to set up the prosthesis. The patient can use the app to change the behaviour of the prosthesis to a certain extent during everyday use (e.g. while becoming accustomed to the prosthesis). The O&P professional can use the adjustment software to track these changes at the patient's next appointment.
- If the settings of a MyMode are to be modified, one must first switch to this MyMode.

9.4.1 Changing the prosthesis setting using the cockpit app

- 1) Once the component is connected and in the desired mode, tap the ☰ symbol in the main menu.
→ The navigation menu opens.
- 2) Tap the "**Settings**" menu option.
→ A list appears with the parameters for the currently selected mode.
- 3) Change the setting of the desired parameter by tapping the "<", ">" symbols.

INFORMATION: The O&P professional's setting is marked and, after the setting has been changed, can be restored by tapping the "Standard" button.

9.4.2 Overview of adjustment parameters in basic mode

The parameters in basic mode describe the dynamic behaviour of the prosthesis in a normal gait cycle. These parameters act as basic settings for automatically adjusting the damping behaviour to the current motion situation (e.g. ramps, slow walking speed, etc.).

The following parameters can be modified:

Parameter	Adjustment software range	Setting range, app	Meaning
Pitch	1000 Hz–4000 Hz	1000 Hz–4000 Hz	Pitch (frequency) of beep signal for confirmation tones
Volume	0–4	0–4	Volume of beep signal for confirmation tones (e.g. when checking the charge level, switching MyModes). The "0" setting deactivates the audible feedback signals. However, warning signals are still generated when errors occur (see page 77).
Heel Resist.	10–60	± 20	Plantar flexion damping. How fast the forefoot is lowered when weight is placed on the heel.
Rollover Resist.	110–170	± 10	This parameter defines how easy the rollover is.
Stair Function	ON–OFF	ON–OFF	Switching on this function increases the rollover angle when walking down stairs. In order to do so, this function needs to be enabled in the adjustment software.

9.4.3 Overview of adjustment parameters in MyModes

The parameters in the MyModes describe the static behaviour of the prosthesis for a specific motion pattern such as golf. Damping behaviour is not automatically controlled and adjusted in MyModes.

The following parameters can be modified in MyModes:

Parameter	Adjustment software range	Setting range, app	Meaning
Heel Resist.	0–195	± 20	Plantar flexion damping. How fast the forefoot is lowered when weight is placed on the heel.
Rollover Resist.	0–195	± 10	Dorsiflexion damping. How easily the value of the "Stop angle" parameter can be reached, in other words how strong the resistance is to reaching the value of the "Stop angle" parameter.
Stop angle	-200–200	± 10 displayed in 0.1°	Ankle angle from which motion in the rollover direction (in direction of dorsiflexion) is locked.

9.5 Switching off the product

⚠ CAUTION

Using the product while switched off

Falling due to unexpected behaviour of the product because of changed damping behaviour.

- ▶ Before using the product, switch it on by connecting the power supply and battery charger.

If the prosthesis is not being worn, it switches to energy saving mode after some time (15 minutes). All sensors are turned off. Energy saving mode is deactivated again as soon as the prosthesis is moved.

In certain cases, e.g. for storage or transportation, the prosthesis can be purposely switched off. It can only be switched on by connecting to a live outlet, a power supply and a battery charger.

Switching off

- ▶ Connect/disconnect the battery charger to/from the product 3x. A wait time less than 3 seconds is required before disconnecting the battery charger again.
- After disconnecting 3 times, a descending sequence of 5 beeps is emitted and the product is then switched off.

Switching on

- 1) Connect the power supply with battery charger to the socket.
- 2) Connect the battery charger to the product.
 - The correct connection of the battery charger to the product is indicated by feedback (see page 76 and see page 79).

9.6 Turning Bluetooth on the prosthesis on/off

INFORMATION

Bluetooth on the prosthesis must be turned on in order to use the Cockpit app.

If Bluetooth is switched off, it can be turned on by turning the prosthesis upside-down (function only available in basic mode) or by connecting/disconnecting the battery charger. Bluetooth is then turned on for approx. 2 minutes. During this time, the app must be started and used to establish a connection. If required, Bluetooth on the prosthesis can be switched on permanently afterwards (see page 70).

9.6.1 Switching Bluetooth off/on using the cockpit app

Switching off Bluetooth

- 1) When the component is connected, tap the ☰ symbol in the main menu.
 - The navigation menu opens.
- 2) In the navigation menu, tap the entry "**Functions**".
- 3) Tap the entry "**Deactivate Bluetooth**".
- 4) Follow the on-screen instructions.

Switching on Bluetooth

- 1) Rotate the component or connect/disconnect the battery charger.

- Bluetooth is switched on for approx. 2 minutes. During this time, the app must be started in order to establish a connection to the component.
- 2) Follow the on-screen instructions.
- If Bluetooth is switched on, the  symbol appears on the screen.

9.7 Querying the prosthesis status

9.7.1 Query status through cockpit app

- 1) When the component is connected, tap the  symbol in the main menu.
- 2) In the navigation menu, tap the entry "**Status**".

9.7.2 Status display in the cockpit app

Menu option	Description	Possible actions
Trip: 1747	Daily step counter	Reset the counter by tapping the " Reset " button.
Step: 1747	Total step counter	For informational purposes only
Batt.: 68	Current prosthesis charge level, as a percentage	For informational purposes only

10 Additional operating states (modes)

10.1 Empty battery mode

Beeps and vibration signals are emitted if the available battery charge level drops to 0% (see page 77). During this time, damping settings are set to their safety mode values. The prosthesis is then switched off. You can switch back to basic mode (mode 1) from empty battery mode by charging the product.

10.2 Mode for charging the prosthesis

The ankle joint of the prosthetic foot is locked during the charging process.

10.3 Safety mode

The product automatically switches to safety mode if a critical error occurs (e.g. failure of a sensor signal) or if the battery is drained. Safety mode remains in effect until the error has been rectified.

Default damping values are activated in safety mode. This makes limited walking possible for the user even though the product is not active.

The switch to safety mode is indicated by beeps and vibration signals immediately prior to switching (see page 77).

Safety mode can be disabled by connecting then disconnecting the battery charger. If the product switches into safety mode again, this means a permanent error exists. The product must be inspected by an authorised Ottobock Service Centre.

10.4 Overheating mode

To prevent overheating of the hydraulic unit due to uninterrupted, increased activity (e.g. extended walking downhill), the range of motion in the ankle joint is restricted as the temperature rises. Depending on the temperature, this restriction may also result in the ankle joint being locked completely. When the hydraulic unit cools down, the product switches back to the settings that existed before the overheating mode.

Overheating mode is indicated by a brief vibration every 5 seconds.

11 Cleaning

- 1) Switch off the product before cleaning.
- 2) Clean the product with a damp cloth and mild soap when needed.
Make sure that no liquids get into the product and product components.
- 3) Dry the product with a lint-free cloth and allow it to air dry fully.

12 Maintenance

INFORMATION

The foot shell for the prosthetic foot has an expected working life of approx. one year when professionally mounted and properly used. Damaged foot shells must be replaced immediately, before the prosthetic foot is used again.

Regular maintenance (service inspections) at 24-month intervals is mandatory in the interest of patient safety and in order to maintain operating reliability and protect the warranty, maintain basic safety and the essential performance characteristics, and ensure safety in regards to EMC.

When maintenance is due, this is indicated by feedback after disconnecting the battery charger (see the section "Operating states/error signals", see page 76). The manufacturer grants a grace period of no more than one month before, or two months after, the due date.

Additional services such as repairs may be provided in the course of maintenance. These additional services may be provided free of charge or can be billable according to an advance cost estimate, depending on the extent and validity of the warranty.

The following components must always be sent in for maintenance and repairs:

The product, battery charger and power supply. The shipping container for the loaner unit you receive must be reused for sending back the components requiring inspection.

12.1 Identification of the product by the Service Center

The product may have been identified by an authorised Ottobock Service Center:

Factory setting

The patient-specific product settings have been reset to the state at delivery (factory setting).

User setting

The settings already configured using the adjustment software were not changed.

CAUTION

Use of the prosthesis with incorrect setting data

Falling due to unexpected prosthesis behaviour caused by triggering the swing phase at the wrong time.

- ▶ The prosthesis settings (parameters) have to be checked using the corresponding adjustment software and changed as needed.

13 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

13.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

13.2 Trademarks

All product names mentioned in this document are subject without restriction to the respective applicable trademark laws and are the property of the respective owners.

All brands, trade names or company names may be registered trademarks and are the property of the respective owners.

Should trademarks used in this document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

13.3 CE conformity

Otto Bock Healthcare Products GmbH hereby declares that the product is in compliance with European requirements for medical devices.

This product meets the requirements of the 2014/53/EU directive.

The full text of the regulations and requirements is available at the following Internet address:
<http://www.ottobock.com/conformity>

13.4 Local Legal Information

Legal information that applies **exclusively** to specific countries is written in the official language of the respective country of use in this chapter.

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:

- 1) This device may not cause harmful interference, and
- 2) This device must accept any interference received, including interference that may cause undesired operation.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/ TV technician for help.

Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

Caution: Exposure to Radio Frequency Radiation.

This device must not be co-located or operating in conjunction with any other antenna or transmitter.

Responsible party:

Otto Bock Health Care, LP
 3820 West Great Lakes Drive
 Salt Lake City, Utah 84120-7205 USA
 Phone + 1-801-956-2400
 Fax + 1-801-956-2401

This device complies with RSS 210 of Industry Canada.

Operation is subject to the following two conditions:

- (1) this device may not cause interference, and
- (2) this device must accept any interference, including interference that may cause undesired operation of this device.

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Caution: Exposure to Radio Frequency Radiation.

The installer of this radio equipment must ensure that the antenna is located or pointed such that it does not emit RF field in excess of Health Canada limits for the general population; consult Safety Code 6, obtainable from Health Canada's website

<http://www.hc-sc.gc.ca/rpb>.

Responsible party:

Otto Bock Healthcare Canada Ltd.
 5470 Harvester Road
 L7L 5N5 Burlington, Ontario
 Canada

Caution: Federal law (USA) restricts this device to sale by or on the order of a practitioner licensed by law of the State in which he/she practices to use or order the use of the device.

14 Technical data

Environmental conditions						
Storage and transport in original packaging (≤3 months)	-20 °C/-4 °F to +40 °C/+104 °F					
Storage and transport without packaging (<48 hours)	-25 °C/-13 °F to +70 °C/+122 °F Max. 93% relative humidity, non-condensing					
Long-term storage (>3 months)	-20°C/-4°F to +20°C/+68°F Max. 93% relative humidity, non-condensing					
Operation	-10°C/+14°F to +40°C/+104°F Max. 93% relative humidity, non-condensing					
Charging the battery	+10 °C/+50 °F to +45 °C/+113 °F					
Product						
Reference number	1B1-2					
Maximum adjustable heel height	50 mm/2 inches					
Dorsiflexion with heel height of 1 cm/0.39 inches	14.5°					
Plantar flexion with heel height of 1 cm/0.39 inches	22°					
Mobility grade (MOBIS)	2–3					
Footshell colours	Translucent, beige, brown					
Max. system height with heel height of 2 cm/0.79 inches	18.5 cm/7.28 inches					
Protection rating	IP54					
Water resistance	Weatherproof but not corrosion-resistant Not designed for underwater use or submersion					
Range of Bluetooth connection	Max. 10 m					
Information on the product's ruleset and firmware version	Accessible via the Cockpit app navigation menu and the menu item "Imprint/Info"					
Expected lifetime given compliance with the prescribed maintenance intervals	6 years					
Test procedure (foot sizes 24 and 25)	ISO 22675-P5-100 kg/2 million load cycles					
Test procedure (foot sizes 26 through 29)	ISO 22675-P6-125 kg/2 million load cycles					
Foot size [cm]	24	25	26	27	28	29
Max. body weight	100 kg/220 lbs	125 kg/275 lbs	125 kg/275 lbs			
Max. weight incl. footshell	approx. 1275 g/45 oz	approx. 1485 g/52 oz	approx. 1555 g/55 oz			
Data transfer						
Wireless technology	Bluetooth Smart Ready					
Range	approx. 10 m / 32.8 ft					
Frequency range	2402 MHz to 2480 MHz					
Modulation	GFSK, π/4 DQPSK, 8DPSK					
Data rate (over the air)	2178 kbps (asymmetrical)					
Maximum output power (EIRP):	+8.5 dBm					
Prosthesis battery						
Battery type	Li-Ion					
Charging cycles (charging and discharging cycles) after which at least 80% of the original battery capacity remains available	500					

Prosthesis battery	
Charging time until battery is fully charged	8 hours
Behaviour of prosthetic foot while being charged	The ankle joint of the prosthetic foot is locked
Operating time of prosthesis with fully charged battery	1 day with average use
Power supply unit	
Reference number	757L16-4
Type	FW8001M/12
Storage and transport in original packaging	-40 °C/-40 °F to +70 °C/+158 °F 10% to 95% relative humidity, non-condensing
Storage and transport without packaging	-40 °C/-40 °F to +70 °C/+158 °F 10% to 95% relative humidity, non-condensing
Operation	0 °C/+32 °F to +50 °C/+122 °F Max. 95% relative humidity Air pressure: 70–106 kPa (up to 3,000 m without pressure equalisation)
Input voltage	100 V~ to 240 V~
Mains frequency	50 Hz to 60 Hz
Output voltage	12 V ==
Battery charger	
Reference number	4E50-2
Storage and transport in original packaging	-25 °C/-13 °F to +70 °C/+158 °F
Storage and transport without packaging	-25 °C/-13 °F to +70 °C/+158 °F Max. 93% relative humidity, non-condensing
Operation	0 °C/+32 °F to +40 °C/+104 °F Max. 93% relative humidity, non-condensing
Input voltage	12 V ==
Lifetime	8 years
Cockpit app	
Reference number	4X441-IOS=*/4X441-Andr=V* Cockpit
Supported operating system	See the information in the respective online store (e.g. Apple App Store, Google Play Store, etc.) regarding compatibility with mobile devices and versions.
Website for download	https://www.ottobock.com/cockpitapp

15 Appendices

15.1 Symbols Used

Manufacturer

Compliance with the requirements according to "FCC Part 15" (USA)

Compliance with the requirements under the "Radiocommunications Act" (AUS)

Non-ionising radiation

IP54

Protected against dust, protected against splashing water

In some jurisdictions it is not permissible to dispose of these products with unsorted household waste. Disposal that is not in accordance with the regulations of your country may have a detrimental impact on health and the environment. Please observe the instructions of your national authority pertaining to return and collection.

DUAL

The product's Bluetooth wireless module can establish a connection to mobile devices with the following operating systems: iOS (iPhone, iPad, iPod...) and Android

Declaration of conformity according to the applicable European directives

Serial number (YYYY WW NNN)
 YYYY – year of manufacture
 WW – week of manufacture
 NNN – sequential number

Lot number (PPPP YYYY WW)
 PPPP – plant
 YYYY – year of manufacture
 WW – week of manufacture

Article number

Medical device

Please note the instructions for use

Check the product settings using the corresponding adjustment software of the Ottobock Data Station.

15.2 Operating states/error signals

The prosthesis indicates operating states and error messages through beeps and vibration signals.

15.2.1 Signals for operating states

Battery charger connected/disconnected

Beep signal	Vibration signal	Event
1x short		Battery charger is connected or battery charger already disconnected prior to start of charging mode
	3x short	Charging mode started (3 sec. after connection of battery charger)
1x short	1x before beep signal	Battery charger disconnected after start of charging mode

Mode switching

INFORMATION

If the **Volume** parameter is set to '0' in the Cockpit app, there are no beep signals (see page 68).

Beep signal	Vibration signal	Additional action performed	Event
1 x short	1 x short	Mode switching using the Cockpit app	Mode switching is performed using the Cockpit app
1 x short	1 x short	Tapping with heel to switch modes or	Motion pattern recognised

Beep signal	Vibration signal	Additional action performed	Event
		swing sideways three times to set heel height	
1 x short	1 x short	Weight placed on prosthetic leg and leg kept still for 1 second to switch mode or feet placed on same level and weight distributed evenly on both feet to set heel height	Switching to basic mode (mode 1) carried out.
2 x short	2 x short	Weight placed on prosthetic leg and leg kept still for 1 second	Switching to MyMode 1 (mode 2) carried out.
3 x short	3 x short	Weight placed on prosthetic leg and leg kept still for 1 second	Switching to MyMode 2 (mode 3) carried out.

15.2.2 Warnings/error signals

Error during use

Beep signal	Vibration signal	Event	Required action
–	1 x long at interval of approx. 5 seconds	Overheated hydraulic unit	Reduce activity.
–	3 x long	Battery charge level under 25%	Charge battery soon.
–	5 x long	Battery charge level under 15%	Charge battery immediately; the product will be switched off after the next warning signal.
10 x short	10 x long	Charge level 0% After the beep and vibration signals, the product switches to empty battery mode and then switches off.	Charge the battery.
30 x long	1 x long, 1 x short repeated every 3 seconds	Severe error/indication that safety mode is activated e.g. one or more sensors are not operational.	Walking possible with restrictions. Please note the possible change in flexion/extension resistance. Attempt to reset this error by connecting/disconnecting the battery charger. The battery charger must be connected for at least 5 seconds before it is disconnected. If the error persists, use of the product is prohibited. The product must be inspected by an authorised Ottobock Service Centre.
–	Continuous	Total failure Electronic control no longer possible. Safety mode active or undetermined valve state. Undetermined product behaviour.	Attempt to reset this error by connecting/disconnecting the battery charger. If the error persists, use of the product is prohibited. The product must be inspected by an authorised Ottobock Service Centre.

Error while charging the product

LED on power supply	LED on battery charger	Error	Resolution
○	🔌 ○ ○ ①	Country-specific plug adapter not fully engaged on power supply	Check whether the country-specific plug adapter is fully engaged on the power supply.
		Non-functional socket	Check socket with another electric device.
		Defective power supply	The battery charger and power supply must be inspected by an authorised Ottobock Service Centre.
●	🔌 ○ ○ ①	No connection between battery charger and power supply	Check whether the charging cable plug is fully engaged on the battery charger.
		Defective battery charger	The battery charger and power supply must be inspected by an authorised Ottobock Service Centre.
●	🔌 ○ ● ①	Battery is fully charged (or connection with product is interrupted).	Take note of the confirmation signal for differentiation. When the battery charger is connected or disconnected, a self-test is conducted and confirmed by a beep and vibration signal. The battery is fully charged if this signal is heard. If no signal is emitted, the connection to the product is interrupted.
			If the connection to the product is interrupted, an authorised Ottobock Service Centre must inspect the product, battery charger and power supply.
Beep signal	Error	Resolution	
4 x short at intervals of approx. 20 sec. (continuously)	Charging the battery outside the allowable temperature range	Check whether the specified ambient conditions for charging the battery are met (see page 74).	

15.2.3 Error messages while establishing a connection with the cockpit app

Error message	Cause	Correction
Component was connected to another device. Establish connection?	The component was connected to another mobile device	To disconnect the original connection, tap the "OK" button. If the original connection is not to be disconnected, tap the "Cancel" button.
Mode change failed	An attempt was made to switch to a different MyMode while the component was in motion (e.g. while walking)	For safety reasons, switching MyModes is only permitted when components are at rest, e.g. while standing or sitting.
	A current connection to the prosthesis was interrupted	Check the following points: <ul style="list-style-type: none"> • Distance from the prosthesis to the mobile device • Charge level of the prosthesis • Bluetooth on the prosthesis switched on? (see page 70) • Hold the component with the sole of the foot facing up to make the component "visible" for 2 minutes. • Prosthesis switched on? (see page 70) • If multiple prostheses were stored, was the correct prosthesis selected?

15.2.4 Status signals

Battery charger is connected

LED on power supply	LED on battery charger	Event
●	🔌 ○ ● ①	Power supply and battery charger operational

Battery charger disconnected

Beep signal	Vibration signal	Result
1 x short	1 x short	Self-test completed successfully. Product is operational.
3 x short	–	Maintenance note Conduct the self-test again by connecting/disconnecting the battery charger. If the beep signal sounds again, product maintenance should be carried out by an authorised Ottobock Service Center. The product can be used without restrictions. However, vibration signals may not be generated.
–	–	Conduct the self-test again by connecting/disconnecting the battery charger. If no beep/vibration signal is emitted after connecting/disconnecting the battery charger again, the product must be inspected by an authorised Ottobock Service Center.

Battery charge level

Battery charger	
🔌 ● ● ①	Battery is being charged, battery charge level is less than 50%
🔌 ☀️ ● ①	Battery is being charged, battery charge level is over 50%
🔌 ○ ● ①	Battery is fully charged (or connection with product is interrupted). Take note of the confirmation signal for differentiation. When the battery charger is connected or disconnected, a self-test is conducted and confirmed by a beep and vibration signal. The battery is fully charged if this signal is heard. If no signal is emitted, the connection to the product is interrupted.

15.3 Directives and manufacturer's declaration

15.3.1 Electromagnetic environment

This product is designed for operation in the following electromagnetic environments:

- Operation in a professional healthcare facility (e.g. hospital, etc.)
- Operation in areas of home healthcare (e.g. use at home, use outdoors)

Observe the safety notices in the section "Information on proximity to certain areas" (see page 51).

Electromagnetic emissions

Interference measurements	Compliance	Electromagnetic environment directive
HF emissions according to CISPR 11	Group 1/class B	The product uses HF energy exclusively for its internal functioning. Its HF emissions are therefore very low, and interference with neighbouring electronic devices is unlikely.
Harmonics according to IEC 61000-3-2	Not applicable – power below 75 W	–
Voltage fluctuations/flicker according to IEC 61000-3-3	Product meets the requirements of the standard.	–

Electromagnetic interference immunity

Phenomenon	EMC basic standard or test procedure	Interference immunity test level
Electrostatic discharge	IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air,
High-frequency electro-magnetic fields	IEC 61000-4-3	10 V/m 80 MHz to 2.7 GHz 80% AM at 1 kHz
Magnetic fields with rated power frequencies	IEC 61000-4-8	30 A/m 50 Hz or 60 Hz
Electrical fast transients/bursts	IEC 61000-4-4	± 2 kV 100 kHz repetition rate
Surges Line against line	IEC 61000-4-5	± 0.5 kV, ± 1 kV
Conducted interference induced by high-frequency fields	IEC 61000-4-6	3 V 0.15 MHz to 80 MHz 6 V in ISM and amateur frequency bands between 0.15 MHz and 80 MHz 80% AM at 1 kHz
Voltage drops	IEC 61000-4-11	0% U _T ; 1/2 period At 0, 45, 90, 135, 180, 225, 270 and 315 degrees
		0% U _T ; 1 period and 70% U _T ; 25/30 periods Single phase: at 0 degrees
Voltage interruptions	IEC 61000-4-11	0% U _T ; 250/300 periods

Interference resistance against wireless communication devices

Test frequency [MHz]	Frequency band [MHz]	Radio service	Modulation	Maximum power [W]	Distance [m]	Interference immunity test level [V/m]
385	380 to 390	TETRA 400	Pulse modulation 18 Hz	1.8	0.3	27
450	430 to 470	GMRS 460, FRS 460	FM ± 5 kHz deviation 1 kHz sine	1.8	0.3	28
710	704 to 787	LTE band 13, 17	Pulse modulation 217 Hz	0.2	0.3	9
745						
780						
810	800 to 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, GSM 800/900, LTE band 5	Pulse modulation 18 Hz	2	0.3	28
870						
930						
1,720	1,700 to 1,990	GSM 1800; CDMA 1900;		2	0.3	28
1,845						
1,970	1,700 to 1,990	GSM 1900; DECT; LTE band 1, 3, 4, 25;	Pulse modulation 217 Hz	2	0.3	28

Test frequency [MHz]	Frequency band [MHz]	Radio service	Modulation	Maximum power [W]	Distance [m]	Interference immunity test level [V/m]
		UMTS				
2,450	2,400 to 2,570	Bluetooth WLAN 802.11-b/g/n, RFID 2450 LTE band 7	Pulse modulation 217 Hz	2	0.3	28
5,240	5,100 to 5,800	WLAN 802.11-a/n	Pulse modulation 217 Hz	0.2	0.3	9
5,500						
5,785						

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