

## 17AD1\*, 17AF2\*



<b>DE</b>	Gebrauchsanweisung .....	5
<b>EN</b>	Instructions for use .....	12
<b>FR</b>	Instructions d'utilisation .....	19
<b>IT</b>	Istruzioni per l'uso .....	26
<b>ES</b>	Instrucciones de uso .....	34
<b>PT</b>	Manual de utilização .....	41
<b>NL</b>	Gebruiksaanwijzing .....	48
<b>SV</b>	Bruksanvisning .....	56
<b>DA</b>	Brugsanvisning .....	62
<b>NO</b>	Bruksanvisning .....	69
<b>PL</b>	Instrukcja użytkowania .....	75
<b>HU</b>	Használati utasítás .....	83
<b>CS</b>	Návod k použití .....	90
<b>TR</b>	Kullanma talimatı .....	96
<b>EL</b>	Οδηγίες χρήσης .....	103
<b>RU</b>	Руководство по применению .....	111


1

2

## INFORMATION

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### INFORMATION

Date of last update: 2021-05-07

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

These instructions for use provide you with important information on the processing of the 17AD1\* and 17AF2\* orthotic ankle joints.

## 2 Product description

### 2.1 Available sizes

The orthotic ankle joints are available in two versions and are delivered in pairs.

Size	Weight (per pair)	Max. body weight up to	Lower leg length (MTP – floor)	Overall length upper section – foot stirrup	Material
17AD1=120	160 g	100 kg	55 cm	120 mm	Stainless steel
17AD1=93	130 g	35 kg	35 cm	93 mm	

Size	Weight (per pair)	Max. body weight up to	Lower leg length (MTP – floor)	Overall length upper section – foot stirrup	Material
17AF2=77	60 g	100 kg	55 cm	77.5 mm	Stainless steel
17AF2=63	30 g	35 kg	35 cm	63.5 mm	

## 2.2 Components

### 17AD1\*

Scope of delivery (see Fig. 1)			
Item	Designation	Spare part	Article number
1	Oval head screw		501S89=M5x6 501S89=M5x8 (for heavier plastic)
2	Set screw		506G21=M6x6
3	Compression spring	X	513D83= 1.1x3.7
4	Ball		509Y1=5.0
5	Shoulder screw for parallel retainer		501S98=M4x30
6	Parallel retainer		29PK5
7	Joint screw	X	501S32=M4X7.5X10
not illustrated	Bearing nut	X	17Y93=M4X6X6.1
not illustrated	Red washer	X	507U5=9.8X7.4X1.5
8	Stop pin		506A8=2x18 506A8=3x18 <b>506A8=5x18</b> 506A8=5x20

### Service set

Article number	Consisting of
17AS=19	Bearing nut Compression spring Joint screw Red washer

## Accessories

Article number	Designation	Size
743Y57	Alignment insert	5

### 17AF2\*

Scope of delivery (see Fig. 2)		
Item	Component	Article number
1	Oval head screw	501S89=M5x6 501S89=M5x8 (for heavier plastic)
2	Shoulder screw for parallel retainer	501S98=M4x30
3	Parallel retainer	29PK5

### 2.3 Design

The orthotic ankle joints are made of stainless steel.

To limit movement, e.g. for dorsiflexion, the 17AD1\* can be equipped with a stop pin or compression spring.

The 17AF2\* is a freely moving orthotic joint with no movement limitation.

## 3 Intended use

### 3.1 Indications for use

The orthotic ankle joints are to be used **exclusively** in pairs in thermoplastic ankle-foot orthoses (AFOs) for treatment of the lower limbs.

### 3.2 Indications

- Partial or total paralysis of the leg muscles.
  - Orthopaedic diseases of the lower limbs
- Indications must be determined by the physician.

### 3.3 Lifetime

The product is designed for a lifetime of **3 years** when used as intended and assembled professionally.

### 3.4 Qualification

Patients may be fitted with the product only by trained qualified personnel. The qualified personnel must be familiar with the handling of the various techniques, materials, machines and tools.

## 4 Safety

### 4.1 Explanation of warning symbols

 **CAUTION** Warning regarding possible risks of accident or injury.

 **NOTICE** Warning regarding possible technical damage.

### 4.2 General safety instructions

 **CAUTION!**

#### Risk of injury and risk of product damage

- ▶ Use the product as intended and only for the defined area of application.
- ▶ If the product has been exposed to extreme strain (e.g. due to falling), take appropriate measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).
- ▶ Use caution when working with the product.
- ▶ Check the product it for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, do not continue using the product and have it checked by authorised, qualified personnel.
- ▶ Use the product on one patient only.
- ▶ Observe the maintenance recommendations.

**NOTICE!**

#### Risk of product damage and limited functionality

- ▶ Avoid storage in condensing ambient humidity.
- ▶ Avoid contact with abrasive substances (e.g. sand, dust).
- ▶ Do not expose the product to temperatures below **-10 °C** or above **+60 °C** (e.g. sauna, excessive sunlight, drying on a radiator).

## 5 Preparing the product for use

 **CAUTION**

#### **Incorrect alignment, assembly or adjustment**

Risk of injury due to change in or loss of functionality

- ▶ Assembly, adjustment and maintenance operations may only be completed by qualified personnel.
- ▶ Observe the alignment, assembly and adjustment instructions.

## NOTICE

### Thermal overloading of the product

Damage due to improper thermal treatment

- ▶ Do not carry out any heat treatment at temperatures above **300 °C (570 °F)**.
- ▶ Prior to thermal treatment, remove all temperature-critical components (such as plastic parts).

### Vacuum forming the orthosis

- > **Required materials:** 743A7 pivot point adjustment aid, 623T3 Perlon stockinette, Loctite 241, plaster negative, metal rod, parallel retainer, Plasticine or 636K8 Plastaband, thermoplastic material
  - > **Prerequisite:** A patient-specific plaster model has been fabricated.
- 1) Mark the compromise axis of rotation of the anatomical ankle joint on the plaster negative. Use the 743A7 pivot point adjustment aid for this purpose (see fig. 3).
  - 2) Open the plaster negative at the marked compromise pivot points (see fig. 4).
  - 3) Insert the parallel retainer. Ensure that the parallel retainer does not bend (see fig. 5). **Optional:** Use the 743Y57 alignment insert, size 5.
  - 4) Insert a metal rod into the plaster negative. Ensure that the metal rod and the parallel retainer do not touch (see fig. 6).
  - 5) Fabricate the plaster positive.
  - 6) Remove the plaster negative.
  - 7) Shorten the parallel retainer. The parallel retainer has to be flush with the contour of the plaster positive (see fig. 7).
  - 8) Mark the alignment reference lines of the orthotic joints on the plaster positive. Use the 743A7 pivot point adjustment aid for this purpose (see fig. 8).
  - 9) Adjust the orthotic joints and the foot stirrups. Use bending irons (see fig. 9, see fig. 10).
  - 10) Attach the orthotic joints to the parallel retainer with the shoulder screws.
  - 11) Attach the orthotic joints to the plaster positive with Plasticine. The Plasticine has to be flush with the plaster positive (see fig. 11).
  - 12) Pull two layers of Perlon stockinette over the plaster positive.
  - 13) Vacuum form the heated thermoplastic material over the plaster positive (see fig. 12).
  - 14) Mark the cut edges of the orthosis (see fig. 13).

- 15) Cut the orthosis free and deburr the edges. Ensure that the cut edge at the foot stirrup runs directly below the ankle joint. This results in the highest possible supporting surface of the foot stirrup.
- 16) Transfer the bore holes of the bar and the foot stirrup to the orthosis shell (see fig. 14).
- 17) Screw the orthotic joints to the orthosis shell.
- 18) Secure screw connections with Loctite 241 (see fig. 15).

### 17AD1\*

The 17AD1\* orthotic ankle joints are adapted to the needs of the patient through the individual combination of the ball, stop pin and compression spring with the set screws.

The following combinations are possible:

Combination		Intended use
Installation position Anterior	Installation position Posterior	
Stop pin	Stop pin	e.g. ICP or spina bifida
Compression spring	Stop pin	Support for dorsal extension, e.g. in case of lower leg muscle weakness, to prevent hyper-extension of the knee joint or for energy return during toe-off
Compression spring	Compression spring	Dorsal and plantar support, e.g. in case of lower leg muscle weakness
Stop pin	Compression spring	Dorsiflexion assist with adjustable dorsal stop, e.g. in case of peroneal paralysis with or without own muscular knee joint protection

- ▶ Assemble the threaded holes with ball, stop pin and/or compression spring in the required combination and close the hole with the set screw.
- ▶ For final assembly, secure the set screws with Loctite 241.

## 6 Cleaning

Promptly clean the product after contact with water containing salt, chlorine or soap, or if it gets dirty.

- 1) Rinse the product with clear fresh water.
- 2) Dry the product with a cloth or allow it to air dry. Avoid exposure to direct heat (e. g. from an oven or radiator).

## 7 Maintenance

### INFORMATION

**The product may be exposed to increased loads by the patient.**

- ▶ Shorten the maintenance intervals according to the expected loads.

The manufacturer requires an inspection of the product every **6 months** to verify functionality and check for wear at a minimum.

Spare parts are listed in the section "Components/design".

#### 17AD1\*

A service kit can be ordered for product maintenance.

#### Replacing the compression spring

- ▶ Replace the compression spring with a new spring when there are signs of wear and tear, but no later than after **12 months**.
- ▶ Lubricate the orthotic ankle joints, using only 633F7 special lubricant.

Before handing the orthosis over to the patient, secure all screw connections with Loctite 241.

#### 17AF2\*

- ▶ The orthotic ankle joints are maintenance-free.
- If you have any doubts, please contact the manufacturer.

## 8 Disposal

Dispose of the product in accordance with national regulations.

## 9 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

### 9.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

### 9.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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