21Y14




21Y15

21Y14, 21Y15

DE	Gebrauchsanweisung	5
EN	Instructions for use	11
FR	Instructions d'utilisation	18
IT	Istruzioni per l'uso	25
ES	Instrucciones de uso	31
PT	Manual de utilização	38
NL	Gebruiksaanwijzing	45
SV	Bruksanvisning	51
DA	Brugsanvisning	58
NO	Bruksanvisning	64
PL	Instrukcja użytkowania	70
HU	Használati utasítás	77
CS	Návod k použití	83
EL	Οδηγίες χρήσης	90
RU	Руководство по применению	97
JA	取扱説明書	104
ZH	使用说明书	110

INFORMATION

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INFORMATION

Last update: 2021-06-10

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- ▶ Please keep this document for your records.

Explanation of symbols

 Warnings regarding possible risks of accident or injury.

 Warnings regarding possible technical damage.

 Additional information on the fitting / use.

1 Description

1.1 Indications for use

The Ottobock PushValves und MagValves are intended **exclusively** for the prosthetic fitting of lower limb amputees.

1.2 Field of application

1.2.1 PushValve

The PushValve is suitable for amputees with the following requirements:

- Easy handling (threadless).
- Extended lever arms facilitating valve opening in case of restricted finger functionality.
- The valve opening is large enough to allow easy removal of donning assists.
- Easier position control of the residual limb.
- Audible feedback when the valve has been positioned correctly.

1.2.2 MagValve

The MagValve is suitable for amputees with the following requirements:

- Easy handling (threadless), especially for patients with restricted finger functionality.
- Smaller in size and lower structural height than the PushValve.
- Easier position control of the residual limb.

2 Safety instructions

Please forward the following safety instructions to your patients:

 **CAUTION**

Improper combination of prosthetic components

Risk of injury due to breakage or deformation of the product

- ▶ Only combine the product with prosthetic components that are approved for that purpose.
- ▶ Based on the instructions for use of the prosthetic components, verify that they may be combined with each other.

CAUTION

Reuse on other persons and improper cleaning

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ The product may be used by one person only.
- ▶ Clean the product regularly.

NOTICE

Use under unallowable environmental conditions

Damage to product due to unallowable environmental conditions

- ▶ Do not expose the product to unallowable environmental conditions.
- ▶ If the product has been exposed to unallowable environmental conditions, check it for damage.
- ▶ If damage is apparent or in case of doubt, do not continue using the product.
- ▶ Take suitable measures if required (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop, etc.).

NOTICE

Use of powder (e.g. baby powder, talcum)

Loss of function of prosthetic components due to clogging with particles or deprivation of lubricant

- ▶ Keep the product away from powders.

NOTICE

Incorrect alignment, assembly or adjustment

Product damage due to incorrectly installed or adjusted prosthetic components

- ▶ Assembly, adjustment and maintenance operations may only be completed by qualified personnel.
- ▶ Observe the alignment, assembly and adjustment instructions.

3 Function

3.1 PushValve

The PushValve is a threadless one-way valve. It is opened and closed by pressing together two levers. It is used to evacuate air from the prosthetic socket.

3.2 MagValve

The MagValve is a threadless one-way valve with low structural height. It is closed with magnetic force. It is used to evacuate air from the prosthetic socket.

4 Components

21Y14 PushValve (Fig. 1)	21Y15 MagValve (Fig. 2)
(1) Lower screw ring	(1) Lower screw ring
(2) Upper screw ring	(2) Upper screw ring
(3) Valve upper part	(3) Valve upper part
(4) Air release button	(4) Air release button
(5) Vacuum forming dummy	(5) Vacuum forming dummy
(6) Lamination dummy	(6) Lamination dummy
(7) 501S86=M5x12 Countersunk Head Screw	(7) 501T19=M5x10 Countersunk Head Screw
(8) 503S3 Screw	(8) 503S3 Screw

21Y14 PushValve (Fig. 1)	21Y15 MagValve (Fig. 2)
(9) Pin wrench	(9) Pin wrench

5 Handling

The dummy provided is used during vacuum forming or lamination and allows for precise fitting into the socket made of different materials. Please always make sure that the full dummy surface rests on the plaster positive.

5.1 Setup

NOTICE

Risk of damage to the valve. The valve is not intended for the fixation of inner socket and frame socket. It has no holding function. Only attach the valve to the test socket / to the flexible inner socket.

5.1.1 Installation preparation for Sockets Made of ThermoLyn®

The enclosed vacuum forming dummy (5) provides for that the lower screw ring is flush on the inner side of the socket.

- Soak 99B25 Nylon Stockinette with 85F1 Plaster Parting Agent and pull it over the plaster positive.
- Secure the vacuum forming dummy to the plaster positive with the enclosed 503S3 Screw (Fig. 3).

5.1.2 Completion of a Test Socket

- After vacuum forming of the test socket of ThermoLyn®, grind away the material covering the vacuum forming dummy until it is flush (Fig. 4).
- Remove the 503S3 Screw and vacuum forming dummy.
- Disassemble the valve and apply 617H46 Bonding Agent to the contact surface of the lower screw ring if need be (Fig. 5). Position the lower screw ring in the socket (Fig. 6). Tighten the upper screw ring. The pin wrench can be used for tightening the lower and upper screw rings (Fig. 7).

5.1.3 Completion of a Flexible Inner Socket with Carbon Frame

The enclosed lamination dummy (6) provides for the space required for the upper screw ring in the laminated socket.

- After vacuum forming of the definitive inner socket of ThermoLyn®, grind away the material covering the vacuum forming dummy until it is flush (Fig.4).
- Remove the 503S3 Screw.
- Put the lamination dummy onto the vacuum forming dummy and secure with the enclosed countersunk head screw (Fig. 8). The 5Z1 Lamination Kit with the aid of the Technical Information 2.3.5 or the 5Z4 Lamination Kit for CAT-CAM can be used for reinforcement of the frame.
- After laminating, grind away the material covering the lamination dummy until it is flush and remove the countersunk head screw.
- Remove the carbon frame with inner socket from the plaster positive.
- Disassemble the valve and apply 617H46 Bonding Agent to the contact surface of the lower screw ring if need be (Fig. 5). Position the lower screw ring in the socket (Fig. 6). Tighten the upper screw ring. The pin wrench can be used for tightening the lower and upper screw rings (Fig. 7).
- Position the flexible inner socket in the carbon frame. Make sure that the valve upper part can be easily installed into and removed from the valve base part. Remove more excessive material if necessary.

5.1.4 Installation Preparation for Laminated Sockets

The enclosed vacuum forming dummy provides for that the lower screw ring is flush with the inner side of the socket.

- Pull 99B25 Nylon Stockinette over the plaster positive.
- Secure the vacuum forming dummy to the plaster positive with the enclosed 503S3 Screw (Fig.3).
- Pull 99B81 PVA Bag over the plaster positive and reinforce the laminated socket adequately.

5.1.5 Completion of a Laminated Socket without Flexible Inner Socket

- After laminating the lamination resin socket, grind away the material covering the vacuum forming dummy until it is flush.
- Remove the 503S3 Screw and vacuum forming dummy.
- Disassemble the valve and apply 617H46 Bonding Agent to the contact surface of the lower screw ring if need be (Fig. 5). Position the lower screw ring in the socket (Fig. 6). Tighten the upper screw ring. The pin wrench can be used for tightening the lower and upper screw rings (Fig. 7).
- The lamination dummy is not needed here.

5.2 Instructions for use

5.2.1 PushValve

Place the PushValve upper part into the valve's base part by simply pressing. An audible feedback is emitted when the valve upper part has been pressed in properly. To remove the valve upper part, press the air release button, lift both wings with the thumb and index finger and pull the valve upper part out of the base part.

5.2.2 MagValve

Place the MagValve upper part into the valve's base part by simply pressing. An audible feedback is emitted when the contact between the magnet and metal base has been established. To remove the valve upper part, press the air release button and pull the valve upper part out of the base part.

5.3 Cleaning

The valve should be cleaned regularly with clear, warm water.

6 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

6.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

6.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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