






50K14 Genu Arexa PCL

DE Gebrauchsanweisung	4
EN Instructions for use	8
FR Instructions d'utilisation	13
IT Istruzioni per l'uso	18
ES Instrucciones de uso	23
PT Manual de utilização	28
NL Gebruiksaanwijzing	33
SV Bruksanvisning	38
DA Brugsanvisning	43
FI Käyttöohje	47
PL Instrukcja użytkowania	52
SK Návod na používanie	57
RU Руководство по применению	61
ZH 使用说明书	67

7	Seite/Side	Größe/Size	Umfang/Circumference			
			A (cm)	A (inch)	B (cm)	B (inch)
	L/R	S	32-37	12.6-14.6	39-46	15.4-18.1
	L/R	M	37-41	14.6-16.1	46-53	18.1-20.9
	L/R	L	41-46	16.1-18.1	53-61	20.9-24.0
	L/R	XL	46-50	18.1-19.7	61-69	24.0-27.2
	L/R	XXL	50-56	19.7-22.1	69-79	27.2-31.1

Material	AL, SS, klettbarer Velours aus PA/velastic velour of PA, PU-Schaum/PU-foam, EVA-Schaum/EVA-foam, PA, POM, TPU, PP
-----------------	---

INFORMATION

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INFORMATION

Date of last update: 2021-08-21

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.

- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

These instructions for use provide important information on the fitting and application of the 50K14 Genu Arexa PCL knee joint orthosis.

2 Product description

Scope of delivery (see fig. 1)			
Item	Component	Item	Component
1	Strap guide loop	6	Pad retainer
2	Pad button	7	Extension stop (angular hand grip)
3	Tibia pad	8	Flexion stop (round hand grip)
4	Clip rivet	-	Closures 1 through 6
5	Condyle pad		

Scope of delivery (see fig. 2)			
Item	Component	Item	Component
1	Calf pad	5	Screw connection
2	Dial	6	Fastening strap
3	Scale	7	Quick-release closure
4	Closure holder		

3 Intended use

3.1 Indications for use

The orthosis is intended **exclusively** for orthotic fittings of the lower limbs and **exclusively** for contact with intact skin.

The orthosis must be used in accordance with the indications.

3.2 Indications

- Rupture of the posterior cruciate ligament (acute treatment, conservative therapy, postoperative treatment)

Indications must be determined by the physician.

3.3 Contraindications

3.3.1 Absolute Contraindications

Not known.

3.3.2 Relative Contraindications

- Sensory disorders of the lower limb (e.g. polyneuropathy)
- Circulatory disorders
- Skin changes in the treatment area (e.g. skin diseases, skin injuries, scars)
- Inflammation, reddening and hyperthermia in the treated body region
- Circulatory disorders, lymphatic and venous insufficiencies, severe varicose veins in the lower limbs

3.4 Modes of action

The orthosis relieves and stabilises the knee joint.

3.5 Lifetime

The product is designed for a maximum lifetime of **2 years**.

3.6 Qualification

Patients may be fitted with the product only by trained qualified personnel. The qualified personnel must be familiar with the handling of the various techniques, materials, machines and tools.

4 Safety

4.1 Explanation of warning symbols

CAUTION Warning regarding possible risks of accident or injury.

NOTICE Warning regarding possible technical damage.

4.2 General safety instructions

CAUTION

Use of the product on damaged skin

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ Only use the product on healthy skin.

CAUTION

Reuse on other persons and improper cleaning

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ The product may be used by one person only.
- ▶ Clean the product regularly.

CAUTION

Improper moulding or application

Injuries or damage to the orthosis due to overloading of the material and improper fit of the orthosis due to breakage of load-bearing components

- ▶ Do not make any improper changes to the orthosis.
- ▶ Always apply the orthosis according to the information in the instructions.

CAUTION

Contact with heat, embers or fire

Risk of injury (such as burns) and risk of product damage.

- ▶ Keep the product away from open flames, embers and other sources of heat.

NOTICE

Contact with oils, salves, lotions or other products that contain oils or acids

Insufficient stabilisation due to loss of material functionality

- ▶ Do not expose the product to oils, salves, lotions or other products that contain oils or acids.

NOTICE

Incorrect setting

Damage to the orthosis due to overloading of the material and improper fit of the orthosis due to breakage of load-bearing components

- ▶ The orthosis may only be adjusted by qualified personnel.
- ▶ Do not make any improper changes to the settings.

NOTICE

Use of a worn or damaged product

Limited effectiveness

- ▶ Before each use, check the product for functional reliability and for possible wear or damage.
- ▶ Do not continue using a product that is no longer functional, or that is worn or damaged.

5 Handling

INFORMATION

- ▶ The daily duration of use and period of application are generally determined by the physician.
- ▶ The initial fitting and application of the product must be carried out by qualified personnel.
- ▶ Instruct the patient in the handling and care of the product.
- ▶ Instruct the patient to see a physician immediately if any exceptional changes are noted (e.g. worsening of the complaint).

5.1 Size Selection

- 1) Measure the circumference of the knee joint (circumference A).
- 2) Measure the circumference of the thigh (circumference B) 15 cm above the centre of the patella.
- 3) Determine the size of the orthosis (see page 3, Fig. 7).

5.2 Adaptation

Adjusting the range of motion

INFORMATION

The joint stops (medial and lateral) must have the same degree settings.

Available extension stops (square grip):	0°, 10°, 20°, 30°, 45°
Available flexion stops (round grip):	0°, 10°, 20°, 30°, 45°, 60°, 75°, 90°

- 1) Remove the pad retainer (and the spacer rings if applicable) (see fig. 4).
- 2) Remove the pinch protection stop (see fig. 1, item 8).
- 3) Insert the desired extension and flexion stops into the joint as specified by the doctor (see fig. 5).
→ The stops engage with slight resistance.
- 4) **Optional:** Snap the required number of spacer rings onto the joint (see fig. 6).
- 5) Attach the pad retainer.

Adapting the orthosis

- > **Prerequisite:** The patient is sitting on the front edge of a chair and bends the knee at **70°**. All straps are opened.
- 1) Set the maximum width of the top lower leg closure on dial 5 (see fig. 2, item 2).
 - 2) Place the orthosis on the leg.

- 3) Align the centre of the orthotic joints according to the compromise pivot point according to Nietert (see fig. 3).

Adjusting the strap length

- > The orthosis has been taken off.

 - 1) Remove the pad.
 - 2) Adjust the strap length on the hook-and-loop strap.
 - 3) Attach the pad.

Adjusting the knee width

The knee width of the orthosis can be adjusted with the included spacer rings and the condyle pads.

- 1) Remove the condyle pads and the pad retainers (see fig. 4).
- 2) Attach the required number of spacer rings to the orthotic joint (see fig. 6).
- 3) Attach the pad retainers.
- 4) Attach the condyle pads to the pad retainers.

Adjusting the top lower leg closure

- 1) Fasten closure 1 of the top lower leg closure (see fig. 2, item 7).

INFORMATION: The position of the calf pad is specified by the mounting hole in the orthosis frame.

- 2) **Optional:** Move the top lower leg closure down by one hole.
- 3) Pre-tension the calf pad using the dial.

NOTICE! To avoid a risk of breakage, do not over-tighten the dial.

- 4) Check the pre-tension with the patient standing up and correct as needed.
- 5) Mark the setting for the patient with a pen on the scale of the fastening strap.

Repositioning the top lower leg closure

- 1) Remove the screw connections of the closure holder and the closure on both sides.
- 2) Reposition the closure holder and then insert the screws and tighten them.

5.3 Application

⚠ CAUTION

Incorrect or excessively tight application

Risk of local pressure and constriction of blood vessels and nerves due to improper or excessively tight application

- ▶ Ensure that the product is applied properly and fits correctly.

> **Prerequisite:** All straps have been opened, and the orthosis has been applied.

- 1) Fasten top lower leg closure 1 (see fig. 3).
- 2) Fasten closures 2, 3 and 4 one after the other.
- 3) Turn dial 5 of the calf pad until the pad fits firmly and comfortably.

INFORMATION: Note the position of the fastening strap that was determined by the qualified personnel during the initial fitting and marked on the scale.

- 4) Fasten closure 6.
- 5) Tighten the closures so the orthosis fits firmly and comfortably.
- 6) After **15 min** of active movement, retighten the closures in the specified sequence.

5.4 Removal

- 1) Loosen the pre-tension of the calf pad with dial 5.
- 2) Remove quick-release closure 1 from the pad button.
- 3) Open closures 2, 3 and 4.
- 4) Remove the orthosis from the leg.

6 Cleaning

NOTICE

Use of improper cleaning agents

Damage to the product due to use of improper cleaning agents

- ▶ Only clean the product with the approved cleaning agents.

Clean the orthosis regularly:

- 1) Remove the straps and pads from the orthosis.
- 2) Fasten all hook-and-loop closures.
- 3) Hand-wash the straps and pads in warm water at **30 °C (86 °F)** using neutral detergent. Rinse thoroughly.
- 4) Allow to air dry. Do not expose to direct heat sources (e.g. sunlight, stove or radiator).
- 5) If the orthosis comes into contact with salt water or dirt, rinse with clear water and allow to air dry.

7 Disposal

Dispose of the product in accordance with national regulations.

8 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

8.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

8.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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