1K10, 1K30, 1S30

DE Gebrauchsanweisung	3	CS Návod k použití.....	44
EN Instructions for use	6	RO Instrucțiuni de utilizare.....	47
FR Instructions d'utilisation.....	9	HR Upute za uporabu	50
IT Istruzioni per l'uso	12	SL Navodila za uporabo	53
ES Instrucciones de uso	15	SK Návod na používanie	56
PT Manual de utilização.....	18	BG Инструкция за употреба	59
NL Gebruiksaanwijzing.....	22	TR Kullanma talimatı	62
SV Bruksanvisning	25	EL Οδηγίες χρήσης	65
DA Brugsanvisning	28	RU Руководство по применению	68
NO Bruksanvisning	31	JA 取扱説明書.....	72
FI Käyttöohje.....	34	ZH 使用说明书.....	75
PL Instrukcja użytkowania.....	37	KO 사용 설명서.....	77
HU Használati utasítás	40		

1

1K10

1K30

1S30

2

1

2

1 Produktbeschreibung

Deutsch

INFORMATION

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INFORMATION

Date of last update: 2021-09-14

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

1.1 Construction and Function

The 1K10, 1K30 and 1S30 prosthetic feet are suitable for use in both modular and exoskeletal prostheses. The functional properties are achieved through the combination of a contoured core and functional foam.

1.2 Combination possibilities

This prosthetic component is compatible with Ottobock's system of modular connectors for children. Functionality with components of other manufacturers that have compatible modular connectors has not been tested.

2 Intended use

2.1 Indications for use

The product is intended exclusively for lower limb exoprosthetic fittings.

2.2 Area of application

- The maximum approved body weight is specified in the technical data (see page 8).

2.3 Environmental conditions

Storage and transport

Temperature range -20 °C to +60 °C (-4 °F to +140 °F), relative humidity 20 % to 90 %, no mechanical vibrations or impacts

Allowable environmental conditions

Temperature range: -10 °C to +45 °C (14 °F to 113 °F)

Moisture: relative humidity: 20% to 90%, non-condensing

Unacceptable environmental conditions

Chemicals/liquids: fresh water, salt water, perspiration, urine, acids, soapsuds, chlorine water

Solids: dust, sand, highly hygroscopic particles (e. g. talcum)

Other environmental conditions

UV light Not resistant

Hydrolysis Not resistant

2.4 Lifetime

This product has been load-tested by the manufacturer. The maximum lifetime is 2 years.

3 Safety

3.1 Explanation of warning symbols

 **CAUTION** Warning regarding possible risks of accident or injury.

NOTICE

Warning regarding possible technical damage.

3.2 General safety instructions

CAUTION!

Risk of injury and risk of product damage

- ▶ Comply with the product's field of application and do not expose it to excessive strain (see page 6).
- ▶ The product may only be used by one person; reuse on another person is prohibited.
- ▶ Note the combination possibilities/comboination exclusions in the instructions for use of the products.

NOTICE!

Risk of product damage and limited functionality

- ▶ To prevent mechanical damage, use caution when working with the product.
- ▶ If you suspect the product is damaged, check it for proper function and readiness for use.
- ▶ Do not use the product if its functionality is restricted. Take suitable measures (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop).
- ▶ Do not expose the product to prohibited environmental conditions.
- ▶ Check the product for damage if it has been exposed to prohibited environmental conditions.
- ▶ Do not use the product if it is damaged or in a questionable condition. Take suitable measures (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop).

Signs of changes in or loss of functionality during use

Decreased forefoot resistance or changes in roll-over behaviour are noticeable indications of loss of functionality.

4 Scope of delivery

Quantity	Designation
1	Instructions for use
1	Prosthetic foot

Spare parts/accessories (not included in the scope of delivery)

Fig.	Item	Designation	Reference number
For 1K10 prosthetic feet size 14 – 17, 1K30 size 14 – 17, 1S30 size 12 – 13			
2	1	Foot adapter with screw connection	2R40=2
For 1K10 prosthetic feet size 18 – 21, 1K30 size 18 – 21			
2	2	Foot adapter with screw connection	2R40=1

5 Preparing the product for use

CAUTION!

Incorrect alignment or assembly

- Risk of injury due to damaged prosthetic components
- ▶ Observe the alignment and assembly instructions.

5.1 Mounting the adapter

- ▶ Install the foot adapter on the prosthetic foot according to its instructions for use.

5.2 Alignment

NOTICE

Grinding the prosthetic foot

Premature wear resulting from damage to the prosthetic foot

- ▶ Do not grind the prosthetic foot.

5.2.1 Bench Alignment

Foot size	Mid-foot anterior to the alignment reference line	Heel height
12 – 16 cm	15 mm	see Technical Data
17 – 21 cm	20 mm	

5.2.2 Static Alignment

- Ottobock recommends checking the alignment of the prosthesis using the L.A.S.A.R. Posture and adapting it as needed.
- If necessary, the alignment recommendations (TF modular leg prostheses: **646F219***, TT modular leg prostheses: **646F336***) may be requested from Ottobock.

5.2.3 Dynamic Trial Fitting

- Adapt the alignment of the prosthesis in the frontal plane and the sagittal plane (e.g. by making angle or slide adjustments) to ensure an optimum gait pattern.
- **TT fittings:** Make sure that physiological knee movement in the sagittal and frontal plane is achieved when the leg begins to bear weight after the heel strike. Avoid medial movement of the knee joint. If the knee joint moves in the medial direction in the first half of the stance phase, move the prosthetic foot in the medial direction. If the medial movement occurs in the second half of the stance phase, reduce the exterior rotation of the prosthetic foot.

6 Cleaning and Care

- 1) Clean the product with a damp, soft cloth.
- 2) Dry the product with a soft cloth.
- 3) Allow to air dry in order to remove residual moisture.

7 Maintenance

- ▶ A visual inspection and functional test of the prosthetic components should be performed after the first 30 days of use.
- ▶ Inspect the entire prosthesis for wear during normal consultations.
- ▶ Have the prosthetic components inspected every 3 months.

8 Disposal

In some jurisdictions it is not permissible to dispose of the product with unsorted household waste. Improper disposal can be harmful to health and the environment. Observe the information provided by the responsible authorities in your country regarding return, collection and disposal procedures.

9 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

9.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

9.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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