




Harmony P3 4R147

DE Gebrauchsanweisung	5
EN Instructions for use	13
FR Instructions d'utilisation	21
IT Istruzioni per l'uso	30
ES Instrucciones de uso	39
PT Manual de utilização	48
NL Gebruiksaanwijzing	56
SV Bruksanvisning	65
DA Brugsanvisning	73
NO Bruksanvisning	81
PL Instrukcja użytkowania	89
CS Návod k použití	97
TR Kullanma talimatı	105
EL Οδηγίες χρήσης	114
RU Руководство по применению	123
JA 取扱説明書	132

INFORMATION

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INFORMATION

Date of last update: 2021-08-23

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

1.1 Construction and Function

The 4R147 Harmony P3 mechanical vacuum pump generates a vacuum in the prosthetic socket that improves the connection with the residual limb. The vacuum pump is usually installed directly below the prosthetic socket.

While walking, the vacuum pump is alternately compressed and released again. This generates the vacuum and simultaneously dampens vertical impact loads. The design also permits slight torsion of the prosthetic socket relative to the foot.

The Harmony system consists of a vacuum pump, a total surface weight-bearing socket, a polyurethane liner and a sealing sleeve.

1.2 Combination possibilities

This prosthetic component is compatible with Ottobock's system of modular connectors. Functionality with components of other manufacturers that have compatible modular connectors has not been tested.

2 Intended use

2.1 Indications for use

The product is intended exclusively for lower limb exoprosthetic fittings.

2.2 Area of application

Our components perform optimally when paired with appropriate components based upon weight and mobility grades identifiable by our MOBIS classification information and which have appropriate modular connectors.

The product is recommended for mobility grade 2 (restricted outdoor walker), mobility grade 3 (unrestricted outdoor walker) and mobility grade 4 (unrestricted outdoor walker with particularly high demands).

- The maximum approved body weight is specified in the technical data (see page 21).

2.3 Environmental conditions

Storage and transport
Temperature range -20 °C to +60 °C (-4 °F to +140 °F), relative humidity 20 % to 90 %, no mechanical vibrations or impacts
Allowable environmental conditions
Temperature range: -10 °C to +45 °C (14 °F to 113 °F)
Moisture: relative humidity: 20% to 90%, non-condensing
Chemicals/liquids: fresh water as dripping water, occasional contact with salty air (e.g. near the ocean)
Solids: dust

Prohibited environmental conditions
Chemicals/moisture: salt water, perspiration, urine, acids, soapsuds, chlorine water
Solids: dust in high concentrations (e.g. construction site), sand, highly hygroscopic particles (e.g. talcum)

2.4 Lifetime

This product was tested by the manufacturer with 3 million load cycles. Depending on the user's activity level, this corresponds to a maximum lifetime of 5 years.

2.5 Qualification

- **Fabricating the total surface weight-bearing socket:** Fabricating a total surface weight-bearing socket is required for the use of the product.

2.6 Absolute Contraindications

- Effectively non-weight-bearing residual limb end
- Neuromas in the area of the residual limb

- User limitations which prevent the safe use of the product

3 Safety

3.1 Explanation of warning symbols

 CAUTION	Warning regarding possible risks of accident or injury.
 NOTICE	Warning regarding possible technical damage.

3.2 General safety instructions

Risk of injury and risk of product damage

- ▶ Comply with the product's field of application and do not expose it to excessive strain (see page 14).
- ▶ Note the combination possibilities/combination exclusions in the instructions for use of the products.
- ▶ Do not expose the product to prohibited environmental conditions.
- ▶ Check the product for damage if it has been exposed to prohibited environmental conditions.
- ▶ Do not use the product if it is damaged or in a questionable condition. Take suitable measures (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop).
- ▶ To avoid the risk of injury and product damage, do not use the product beyond the tested lifetime.
- ▶ To avoid the risk of injury and product damage, only use the product for a single patient.
- ▶ To prevent mechanical damage, use caution when working with the product.
- ▶ If you suspect the product is damaged, check it for proper function and readiness for use.
- ▶ Do not use the product if its functionality is restricted. Take suitable measures (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop).
- ▶ Only use suitable liners (without a textile coating or with a partial textile coating) for prostheses with Harmony vacuum pumps. This prevents skin irritation due to vacuum on the skin.

Signs of changes in or loss of functionality during use

Among other factors, changes in functionality can be indicated by reduced prosthetic suspension.

4 Scope of delivery

Quantity	Designation	Reference number
1	Instructions for use	–
1	Harmony P3 with functional ring	–
2	Socket connector	2R134
1	Sound absorber	4Y360
1	Flange kit	4Y383

Spare parts	
Designation	Reference number
Functional ring	4X147=*
Service set for Harmony P3	4X148
Socket connector	2R117=0
Filter unit	4X632

5 Preparing the product for use

CAUTION

Incorrect alignment or assembly

Risk of injury due to damaged prosthetic components

- ▶ Observe the alignment and assembly instructions.

CAUTION

Improper assembly of the screw connections

Risk of injury due to breakage or loosening of the screw connections

- ▶ Clean the threads before every installation.
- ▶ Apply the specified torque values.
- ▶ Follow the instructions regarding the length of the screws and about how to secure the screws.

5.1 Pre-compressing the functional ring

In order to achieve its full flexibility, the functional ring has to be compressed before the first use.

> **Required tools:** 4X247 Pretensioning tool

- 1) Using an Allen key, unscrew the mounting screw in the distal end of the base body (see fig. 2).
- 2) Pull the adapter with pyramid up and out of the base body (see fig. 3).
- 3) Compress the functional ring for **15 minutes** using the pretensioning tool (see fig. 4, see fig. 5).
- 4) Disassemble the pretensioning tool.

5.2 Installing the socket connector

The vacuum pump is connected to the prosthetic socket by a socket connector. The hose can either be considerably shortened or wrapped around the adapter to secure it closely on the prosthesis.

A hole in the prosthetic socket is required for assembly of the socket connector. It is either drilled or an existing hole is used (from an installed check valve, for example). The socket connector is inserted through the hole and glued to the prosthetic socket.

> **Required materials:** Sandpaper, adhesive (Ottobock recommends Fabtech +PLUSeries Composite 1 Minute Adhesive or Ottobock SuperGlue), 627B40 polyethylene adhesive tape, rounded knife

- 1) **If not already present:** Drill a **Ø 6 mm** hole in the distal, posterior area of the prosthetic socket (see fig. 6).
- 2) Deburr the hole.
- 3) Roughen the adhesive surface on the prosthetic socket with sandpaper (see fig. 7).
- 4) Clean the roughened surface and the flange of the socket connector with a degreasing cleaner.
- 5) Apply the adhesive to the flange of the socket connector (see fig. 8). In doing so, make sure that no adhesive gets onto the end of the flange that tapers to a point.
- 6) Insert the tapered end of the flange into the bore hole, aligning the hose in the direction of the vacuum pump.
- 7) Wrap the prosthetic socket with polyethylene adhesive tape at the level of the socket connector to generate contact pressure (see fig. 9).
- 8) Allow the adhesive to cure and then remove the polyethylene adhesive tape.
- 9) Using a suitable knife, cut off the protruding end of the flange inside the prosthetic socket.

5.3 Installing the vacuum pump

The vacuum pump has a pyramid receiver on the proximal end and a receiver for a 34 mm tube adapter on the distal end. The receiver for the tube adapter has an insertion depth of **42 mm**.

Connecting the pyramid adapter and pyramid receiver

The pyramid adapter is fixed with the set screws of the pyramid receiver.

> **Required materials:** torque wrench (e.g. 710D20), 636K13 Loctite 241

1) **Fitting:**

Screw in the set screws.

Use the torque wrench to tighten the set screws (**10 Nm**).

2) **Definitive assembly:**

Secure the set screws with Loctite.

Screw in the set screws.

Pre-tighten the set screws with the torque wrench (**10 Nm**) and then tighten them (**15 Nm**).

3) Replace any set screws that are protruding or are recessed too much with appropriate set screws (see selection chart).

Selection table for set screws	
Reference number	Length (mm)
506G3=M8X12-V	12
506G3=M8X14	14
506G3=M8X16	16

Connecting the tube adapter

> **Required materials:** torque wrench (e.g. 710D20), degreasing cleaner (e.g. acetone), 636K13 Loctite 241

1) Loosen the tube clamp cap screw by **two turns**.

2) **For final assembly:** Degrease the insertion range of the receiver and tube adapter.

3) Insert the tube adapter into the receiver up to the stop.

4) **For final assembly:** Apply Loctite to the tube clamp cap screw.

5) Tighten the cap screw (tightening torque: **10 Nm**).

Connecting the socket connector

1) Shorten the hose of the socket connector or wrap it around the prosthesis so it does not stick out.

2) Connect the hose to the intake valve of the vacuum pump.

5.4 Installing the outlet hose

The outlet hose is attached to the outlet valve of the vacuum pump. Air and fluids (such as perspiration) are drained through the outlet hose. Fluids must not be drained onto metal components of the prosthesis. The longer the outlet hose, the quieter the vacuum pump.

► Connect the outlet hose to the outlet valve of the vacuum pump.

5.5 Checking the stiffness of the functional ring

The stiffness of the functional ring influences the level of comfort while walking. With the Harmony vacuum pump, the functional ring also determines the efficiency of vacuum generation. The stiffness is specified on the functional ring (the higher the number, the greater the stiffness).

Functional ring while walking	Cause	Correction/solution
Soft damping for the patient. The functional ring is not compressed to the stop.	The functional ring has the appropriate stiffness	–
Patient hits the stop during damping. The functional ring is fully compressed.	Functional ring too soft	Install a harder functional ring (see page 20)
No damping for the patient. The functional ring is not compressed.	Functional ring too hard	Install a softer functional ring (see page 20)

Selection of the functional ring	
Body weight	Stiffness
40 to 47	0
48 to 55	1
56 to 65	2
66 to 75	3
76 to 87	4
88 to 100	5
101 to 112	6
113 to 125	7

5.6 Checking vacuum generation

Vacuum generation is tested during the dynamic fitting.

> **Recommended materials:** 755Z37 pressure gauge (with T-piece)

- 1) Remove the hose from the intake valve of the vacuum pump.
- 2) Connect the hose to the T-piece of the pressure gauge.
- 3) Connect the pressure gauge with the T-piece to the intake valve on the vacuum pump.
- 4) Attach the pressure gauge to the prosthesis so that the patient is able to move freely.
- 5) Have the patient walk **50** steps to activate the vacuum pump.
- 6) Check whether the vacuum is between **340 hPa** and **680 hPa** and is maintained for **20 seconds**.

If a sufficient vacuum is built up and the vacuum does not decrease within **20 seconds**, the vacuum system is functioning correctly.

If a sufficient vacuum is built up but the vacuum decreases within **20 seconds**, determine the cause (see page 19).

- 7) Remove the pressure gauge after completing the test and reinstall the hose on the intake valve of the vacuum pump.

6 Use

CAUTION

Use of powder (e.g. baby powder or talcum)

Skin irritation on the residual limb and loss of function of prosthetic components due to clogging with particles or deprivation of lubricant

- Keep the product away from powders.

6.1 Donning the prosthesis

- 1) Roll the liner over the residual limb. Avoid wrinkles, air pockets and soft tissue displacement as you do so.
- 2) **For liners without a textile coating:** Put on a suitable stocking (e.g. 451F21).
- 3) Slide into the prosthetic socket and unroll the sealing sleeve over the prosthetic socket up to the thigh.

6.2 Doffing the prosthesis

- 1) Unroll the sealing sleeve from the thigh onto the prosthetic socket.
- 2) Doff the residual limb and liner out of the prosthetic socket.

7 Cleaning

7.1 Cleaning the prosthetic socket

- Clean the product with a damp, soft cloth and then let it dry.

7.2 Cleaning the liner

- ▶ Clean the liner as described in the liner instructions for use.

7.3 Cleaning the vacuum pump

> **Required materials:** 30 ml to 60 ml distilled water, vessel

- 1) Fill distilled water into the prosthetic socket.

Alternatively, a residual limb sock can be soaked with distilled water to flush the vacuum pump while walking.

- 2) Position the vessel under the outlet of the vacuum pump.
- 3) Alternately load and unload the prosthesis to pump all the water through the vacuum system.
- 4) Clean and dry the interior surface of the prosthetic socket.

8 Maintenance

⚠ CAUTION

Failure to follow the maintenance instructions


Risk of injuries due to changes in or loss of functionality and damage to the product

- ▶ Observe the following maintenance instructions.

- ▶ A visual inspection and functional test of the prosthetic components should be performed after the first 30 days of use.
- ▶ Inspect the entire prosthesis for wear during normal consultations.
- ▶ Conduct annual safety inspections.
- ▶ Check the condition of the functional ring as part of each inspection. If the user frequently compresses the functional ring to the stop while walking, or if the functional ring exhibits visible damage, replace the functional ring (see page 20).

8.1 Checking for leaks

This section gives a methodical troubleshooting procedure for when the vacuum system leaks.

To check for leaks, the prosthesis is divided into areas.

- **Area 1:** Prosthetic socket and sealing sleeve
- **Area 2:** Socket connector and connection to the vacuum pump
- **Area 3:** Vacuum pump

Checking area 1

First check the donning process and the sealing sleeve. No pressure gauge is needed. The vacuum suction the sealing sleeve against the socket brim. The socket brim is clearly visible under the sealing sleeve during the stance phase. If no vacuum is generated or the vacuum decreases, this can be identified on the sealing sleeve.

- Ask the user to put on the prosthesis. Check the process to exclude errors during donning.
- Get the user to walk in order to build up a vacuum. Get the user to stand and check whether the vacuum is maintained.
- **If the vacuum is not maintained:** Replace the sealing sleeve with a new one.

- Get the user to walk in order to build up a vacuum. Get the user to stand and check whether the vacuum is maintained.
- **If the vacuum is still not maintained:** The fault is in area 2 or 3.

Checking area 3

- Connect the pressure gauge directly to the vacuum pump.
- Activate the vacuum pump two to four times (vacuum of **340 hPa to 680 hPa**).
- **If the vacuum is maintained:** The vacuum pump is OK. The fault is in area 2.
- **If the vacuum is not maintained:** Replace the valves of the vacuum pump. Then check whether this measure was successful.
- **If this measure was not successful:** Replace the functional ring. Then check whether this measure was successful.
- **If this measure was not successful:** Send the vacuum pump to the manufacturer for inspection.

Checking area 2

- Connect the pressure gauge with a T-piece between the socket connector and the vacuum pump.
- Seal the opening to the socket connector on the inside of the prosthetic socket with air-tight adhesive tape (e.g. PVC tape).
- Activate the vacuum pump until a vacuum of **340 hPa to 680 hPa** is reached.
- **If the vacuum is not maintained:** The socket connector is defective. Replace the socket connector. Then check whether this measure was successful.
- **If the vacuum is maintained:** The glued joint of the socket connector is likely defective.
- Seal the opening to the socket connector on the inside of the prosthetic socket with a small piece of paper and air-tight adhesive tape (e.g. PVC tape).
- Activate the vacuum pump until a vacuum of **340 hPa to 680 hPa** is reached.
- **If the vacuum is not maintained:** Glue the socket connector again. Then check whether this measure was successful.

8.2 Replacing the valves

- 1) Pull the old valves off the functional ring, noting the flow direction (arrows on the valves).
- 2) Clean the functional ring with compressed air.
- 3) Plug the new valves into the openings of the functional ring in the correct flow direction.

8.3 Replacing the functional ring

To replace the functional ring, the adapter with pyramid is removed from the base body. The O-rings on the adapter and the mounting screw are also replaced together with the functional ring.

- > **Required materials:** 4X247 pretensioning tool, torque wrench (e.g. 710D20), degreasing cleaner (e.g. 634A58 isopropyl alcohol), Loctite 243 or Loctite 246
- 1) Remove the vacuum pump from the prosthesis.
 - 2) Insert an Allen key (**6 mm**) into the hole in the pyramid and unscrew the mounting screw in the distal end of the base body (see fig. 2).
 - 3) Pull the adapter with pyramid up and out of the base body (see fig. 3).
 - 4) Remove the functional ring.
 - 5) **If the O-rings are worn:** Remove the O-rings from the adapter (see fig. 10).
 - 6) **If the O-ring is worn:** Remove the O-ring from the mounting screw (see fig. 10).
 - 7) Clean the mounting screw and adapter including thread with degreasing cleaner (see fig. 10).
 - 8) Slide new O-rings onto the adapter and apply silicone grease. Make sure that no silicone grease gets into the threads (see fig. 10).
 - 9) Slide a new O-ring onto the mounting screw (see fig. 10).
 - 10) Apply Loctite to the thread of the mounting screw.
 - 11) Slide the new functional ring onto the adapter. The structure of the functional ring permits installation only with the correct alignment.

- 12) Insert the adapter into the base body and attach it with the mounting screw (tightening torque: **15 Nm**). To do so, insert the torque wrench through the hole in the pyramid.
- 13) Pre-compress the new functional ring (see page 16).
- 14) Reinstall the vacuum pump in the prosthesis.

9 Disposal

In some jurisdictions it is not permissible to dispose of the product with unsorted household waste. Improper disposal can be harmful to health and the environment. Observe the information provided by the responsible authorities in your country regarding return, collection and disposal procedures.

10 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

10.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

10.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

10.3 Warranty

The manufacturer warrants this device from the date of purchase. The warranty covers defects that can be proven to be a direct result of flaws in the material, production or construction and that are reported to the manufacturer within the warranty period.

Further information on the warranty terms and conditions can be obtained from the competent manufacturer distribution company.

11 Technical data

Reference number	4R147
Weight [g]	400
System height [mm]	175
Build height [mm]	117
Body weight [kg]	40 to 125
Mobility grade	2 to 4

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