

## 3R80, 3R80=ST

<b>DE</b> Gebrauchsanweisung .....	18
<b>EN</b> Instructions for use .....	28
<b>FR</b> Instructions d'utilisation .....	38
<b>IT</b> Istruzioni per l'uso .....	48
<b>ES</b> Instrucciones de uso .....	59
<b>PT</b> Manual de utilização .....	69
<b>NL</b> Gebruiksaanwijzing .....	80
<b>SV</b> Bruksanvisning .....	90
<b>PL</b> Instrukcja użytkowania .....	100
<b>HU</b> Használati utasítás .....	110
<b>CS</b> Návod k použití .....	120
<b>RO</b> Instrucțiuni de utilizare .....	129
<b>HR</b> Upute za uporabu .....	140
<b>TR</b> Kullanma talimatı .....	150
<b>RU</b> Руководство по применению .....	159
<b>JA</b> 取扱説明書 .....	171
<b>ZH</b> 使用说明书 .....	180


1 **1.1 Kenndaten (Produkt, Dokument)**  
**1.1 Identification data (Product, Document)**


2 **3.3.1 Empfohlene Kombinationen | Unzulässige Kombinationen**  
**3.3.1 Recommended combinations | Unallowable combinations**


**3** **i** **5** Lieferumfang  
**5** Scope of delivery

**4** **i** **6** Technische Daten  
**6** Technical data


5 **i** 6 Technische Daten  
6 Technical data

235 mm

235 mm


150°


Ø 34 mm

①

163 mm

179 mm

②


28 mm


44 mm

③


135 mm




6 **i** 7.2 Grundaufbau der Prothese  
**i** 7.2 Bench alignment of the prosthesis








7 **i** 7.2 Grundaufbau der Prothese  
**i** 7.2 Bench alignment of the prosthesis

8 7.2 Grundaufbau der Prothese  
 1 7.2 Bench alignment of the prosthesis

9 7.2 Grundaufbau der Prothese  
 1 7.2 Bench alignment of the prosthesis

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### INFORMATION

Date of last update: 2021-09-09

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.


This document and the product are intended exclusively for the fabrication of a prosthesis by qualified personnel with technical knowledge of lower limb prosthetics.

### 1.1 Characteristics (product, document)

The following characteristics are important for identification: (1) – page: 2)


- ① 3R80, 3R80=ST (product reference number)
- ② YYYYYWNNN (product serial number)  
YYYY (year of manufacture)  
WW (calendar week)  
NNN (number)
- ③ 647G403=all\_INT-VV-YYMM (standard line of the document)  
647G403=all\_INT (document reference number)  
VV (version number)  
YYMM (publication date) – YY (year); MM (month)

### 1.2 Explanation of warning symbols

 **CAUTION** Warning regarding possible risks of accident or injury.

 **NOTICE** Warning regarding possible technical damage.

### 1.3 Meanings of pictograms in the illustration

	Numbering for the illustrations		Reference to the number of the corresponding section
	Numbering for a defined sequence		Numbering for the parts of an illustration
	Right		Wrong
	Observe the safety notice in the section		Risk of falling
	Replace		

## 2 Product description

The product (3R80\*) has the following key features:

- Monocentric prosthetic knee joint with rotation hydraulics
- Product components for stance phase stability:
  - Adjustable stance phase flexion resistance (hydraulic damping)
  - Adjustable trigger threshold for switching from swing phase extension resistance to stance phase flexion resistance
  - Lock (activated and deactivated by the patient)
- Product components for swing phase control:
  - Adjustable swing phase flexion resistance (hydraulic damping)
  - Adjustable swing phase extension resistance (hydraulic damping)
  - Hydraulic unit with extension assist function (the extension assist function can be enhanced by installing an additional spring)

## 3 Intended use

### 3.1 Indications for use

The product is intended exclusively for lower limb exoprosthetic fittings.

### 3.2 Area of application

 <b>CAUTION</b>
<b>Excessive strain on the product</b>
Fall due to breakage of load-bearing components
▶ Only use the product according to its allowable field of application.

Allowable field of application (3R80, 3R80=ST)	
Recommended mobility grade: <b>3 + 4</b>	Everyday prosthesis
Allowable body weight: <b>≤ 150 kg</b>	

### 3.3 Combination possibilities

 <b>CAUTION</b>
<b>Improper combination of prosthetic components</b>
Injuries, malfunctions or product damage due to unallowable combination of prosthesis components
▶ Based on the instructions for use of all prosthetic components used, verify that they may be combined with each other and are approved for the patient's area of application.

<b>INFORMATION</b>
In a prosthesis, all prosthetic components have to meet the patient's requirements regarding the amputation level, body weight, activity level, environmental conditions and field of application.

This prosthetic component is compatible with Ottobock's system of modular connectors. Functionality with components of other manufacturers that have compatible modular connectors has not been tested.

Our components perform optimally when paired with appropriate components based upon weight and mobility grades identifiable by our MOBIS classification information and which have appropriate modular connectors.

### 3.3.1 Recommended combinations | Unallowable combinations

- ▶ Note the illustration regarding the recommended combinations and impermissible combinations. (2) – page: 2)

### 3.4 Environmental conditions

#### **⚠ CAUTION**

##### **Use under unallowable environmental conditions**

Fall due to damaged product

- ▶ Do not expose the product to unallowable environmental conditions (see the table "Unallowable environmental conditions" in this section).
- ▶ If the product was exposed to unallowable environmental conditions, take suitable steps (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop etc.).

#### **Storage and transport**

Temperature range -20 °C to +60 °C (-4 °F to +140 °F), relative humidity 20 % to 90 %, no mechanical vibrations or impacts

#### **Allowable environmental conditions**

**Temperature range:** -10 °C to +45 °C (14 °F to 113 °F)

**Chemicals/liquids:** fresh water, soapsuds, chlorine water

**Moisture:** Submersion: max. 1 h in 2 m depth, relative humidity: no restrictions

**Solids:** Dust, occasional contact with sand

**Clean the product after contact with humidity/chemicals/solids, in order to avoid increased wear and damage** (see page 36).

#### **Prohibited environmental conditions**

**Solids:** highly hygroscopic particles (e.g. talcum), dust in high concentrations (e.g. construction site), intensive contact with sand

**Chemicals/liquids:** salt water, perspiration, urine, acids, continuous use in liquid media

### 3.5 Reuse and lifetime

#### **⚠ CAUTION**

##### **Reuse on another patient**

Fall due to loss of functionality as well as damage to the product

- ▶ Only use the product for a single patient.

#### **⚠ CAUTION**

##### **Exceeding the lifetime**

Fall due to change in or loss of functionality and damage to the product

- ▶ Ensure that the maximum lifetime defined in this section is not exceeded.

This prosthetic component has been load tested by the manufacturer according to ISO 10328. The maximum lifetime is 5 years.

## 4 General safety instructions

### CAUTION

#### **Reaching into the area of the joint mechanism**

Pinching of limbs (e.g. fingers) and the skin due to uncontrolled joint movement

- ▶ Do not reach into the joint mechanism during use.
- ▶ Pay special attention during assembly and adjustment tasks.

### CAUTION

#### **Mechanical damage to the product**


Risk of injury due to change in or loss of functionality

- ▶ Use caution when working with the product.
- ▶ If the product is damaged, check it for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, do not continue using the product (see "Signs of changes in or loss of functionality during use" in this section).
- ▶ Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

#### **Signs of changes in or loss of functionality during use**

Changes in functionality can manifest themselves, for example, through a changed gait (swing phase, stance phase), incomplete extension, stiffness and the development of noise.

## 5 Scope of delivery

The scope of delivery is shown in the following figures:  – page: 3

Only product components with reference numbers in the illustration can be reordered separately.

In the illustration, a reference number shown in bold stands for a package in which the product components shown in the frame can be reordered.

## 6 Technical data

The technical data are shown in the following figures:  – page: 3,  – page: 4

①: system height; ②: proximal system height; ③: distal system height

## 7 Preparing the product for use

### 7.1 Information on fabrication of a prosthesis

#### CAUTION

#### **Incorrect alignment, assembly or adjustment**

Injury due to incorrectly installed or adjusted as well as damaged prosthetic components

- ▶ Observe the alignment, assembly and adjustment instructions.

#### CAUTION

#### **Initial use of the prosthesis by the patient**

Fall due to lack of patient experience, incorrect alignment or incorrect adjustment of the prosthesis

- ▶ For the safety of the patient, use a suitable device (e.g. parallel bars, handrail, wheeled walker) during initial standing and walking.

## 7.2 Bench alignment

- ▶ **CAUTION!** To enable the patient to stand safely, align the prosthesis using the alignment images as well as the instructions for use of all prosthetic components being used. (6 – page: 5, 7 – page: 5)

### Use 2Z11=KIT

**INFORMATION:** The connection area of the prosthetic knee joint can be protected against scratches with the protective film of the 2Z11=KIT during alignment in the workshop and testing in the trial fitting area.

- ▶ Use the protective film as illustrated in the accompanying document for the 2Z11=KIT.
- ▶ Remove the protective film before the patient leaves the fitting area.

### Adapter inserts

- ▶ Mount the adapter inserts shown in the illustration when using an alignment apparatus. (8 – page: 6)
  - ①: The prosthetic knee joint reference numbers shown in square brackets are marked on the adapter inserts. The adapter insert reference numbers and mm dimensions in the round brackets are not found on the adapter inserts.

## Installing the tube adapter

### ⚠ CAUTION

#### Improper assembly of the screw connections

Risk of injury due to breakage or loosening of the screw connections

- ▶ Clean the threads before every installation.
- ▶ Apply the specified torque values.
- ▶ Follow the instructions for thread lock.

### ⚠ CAUTION

#### Incorrect processing of tube

Fall due to damage to the tube

- ▶ Do not clamp the tube into a vice.
- ▶ To shorten the tube, use only a tube cutter or a cutting device.

### ⚠ CAUTION

#### Incorrect mounting of the tube

Risk of injury due to breakage of load-bearing components

- ▶ Clean the contact surfaces of the tube and the tube adapter using a degreasing cleaning agent.
- ▶ Position the tube in the tube adapter so the minimum insertion depth is met.

- ▶ **CAUTION!** Position the tube adapter only in the permitted insertion range shown in the illustration. (9 – page: 6)
  - ①: Permitted insertion range: 0 mm
  - ②: Impermissible insertion range: > 0 mm

## 7.3 Optimising the static alignment

**INFORMATION:** Additional technical knowledge is provided by the instructions for use of the measuring apparatus being used, the TF alignment poster and the Ottobock seminars.

The 743L500 3D L.A.S.A.R. Posture measuring apparatus is needed to optimise the static alignment.

The following optimisation methods are available:

- ①: Optimisation with 3D mode activated
- ②: Optimisation with 3D mode deactivated

The optimisation method ② can be used as an optional second step. With the 743L100 L.A.S.A.R. Posture, optimisation is possible only based on the values of the optimisation method ②.

**10 (Page: 7): ①: Optimisation with 3D mode activated**

- ▶ Optimise the static alignment so the patient stands in a relaxed position, the values for the reference points shown in the illustration are met and the following points are observed:
  - ①: 3D mode is activated (3D symbol – colour: light grey).
    - Sequence for optimisation in reference to the planes:
      - ①③: Sagittal plane – ①②: Frontal plane
  - ②: The prosthetic knee joint is in extension. The lock, if present, is activated.
  - ③: The reference line lies on the sagittal reference point of the prosthetic knee joint (pivot point).
  - ④: The reference points are in the centre (50 : 50: pyramids of the prosthetic knee joint and the prosthetic foot, shoes).
  - ⑤: The extension of the reference line points to the reference point (anterior superior iliac spine).
  - ⑥: The reference line lies on the frontal reference point of the prosthetic knee joint (pyramid).

**11 (Page: 7): ②: Optimisation with 3D mode deactivated**

- ▶ Optimise the static alignment so the values for the reference points shown in the illustration are met and the following points are observed:
  - ①: 3D mode is deactivated (3D symbol – colour: dark grey).
    - Sequence for optimisation in reference to the planes:
      - ②③: Sagittal plane – ②②: Frontal plane
  - ②: The prosthetic knee joint is in extension. The lock, if present, is activated.
  - ③: The reference line lies on the sagittal reference point of the prosthetic knee joint (pivot point).
  - ④: The reference points are in the centre (50 : 50: pyramids of the prosthetic knee joint and the prosthetic foot, shoes).
  - ⑤: The extension of the reference line points to the reference point (anterior superior iliac spine).
  - ⑥: The reference line lies on the frontal reference point of the prosthetic knee joint (pyramid).
  - ⑦: The extension of the reference line points to the reference point (pyramid).

## 7.4 Optimising during dynamic fitting

### 7.4.1 Notices regarding dynamic fitting

**⚠ CAUTION**

#### Adjusting the Settings

Fall due to incorrect or unfamiliar settings

- ▶ Only adapt the settings to the patient gradually.
- ▶ Explain the effects of the adjustments on the use of the prosthesis to the patient.

### 7.4.2 Overview of adjustment possibilities

#### Stance phase flexion resistance (12 – page: 8 – ④)

①: Adjustment range; ②: Decrease stance phase flexion resistance; ③: Increase stance phase flexion resistance

#### Switching threshold (13 – page: 8 – ④)

①: Adjustment range; ②: Decrease switching threshold (lower weight required to initiate the stance phase); ③: Increase switching threshold (higher weight required to initiate the stance phase)

**Swing phase flexion resistance** (14 – page: 9 – 4)

①: Adjustment range; ②: Decrease swing phase flexion resistance; ③: Increase swing phase flexion resistance

**Swing phase extension resistance** (15 – page: 9 – 4)

①: Adjustment range; ②: Decrease swing phase extension resistance; ③: Increase swing phase extension resistance

**Lock** (16 – page: 10 – 3)

**CAUTION**

**Using the lock**

Falling due to non-activation or unintentional deactivation of the lock

- ▶ Each time after activating and deactivating the lock, carefully check the proper functioning of the prosthesis.
  - ▶ Lock the prosthetic knee joint before using the prosthesis in wet areas.
  - ▶ Adjust the walking speed to the ambient conditions.
  - ▶ Make sure the lock is not deactivated or activated through unintentional operation.
  - ▶ On wet surfaces, only use a prosthetic foot with a slip-resistant sole.
- ▶ Activate the lock as shown in the illustration for use in wet areas (e.g. showering and swimming) and relaxed standing. Deactivate the lock for walking. Put weight on the prosthetic knee joint and push against the extension stop during activation and deactivation.
- ①: Activate lock
  - ②: Deactivate lock

**7.4.3 Checking settings and restoring factory settings**

▶ Prior to optimisation during the dynamic fitting, check the settings with the patient and restore the factory settings.

17 – Page: 10 – ① Stance phase flexion resistance; ② Switching threshold to initiate the stance phase; ③ Bumper without play, pre-stressed

18 – Page: 11 – ① Swing phase flexion resistance; ② Swing phase extension resistance

19 – Page: 11 – ① Lock (deactivated)

Function	Means of adjustment	Factory setting	Meaning
Stance phase flexion resistance	Adjustment ring "S"	Stop "+"	Maximum resistance
Switching threshold to initiate the stance phase	Adjustment ring "I"	Bumper without play, pre-stressed	Low trigger threshold (Stance phase flexion resistance is activated with low load)
Swing phase flexion resistance	Adjustment valve "F"	Stop "-"	Minimum resistance
Swing phase extension resistance	Adjustment valve "E"	Stop "-"	Minimum resistance

**7.4.4 Switching between the swing phase and stance phase**

20 – Page: 12 – ①: Stance phase flexion resistance activated at heel strike; ②: Stance phase flexion resistance deactivated when heel is relieved

- ① The stance phase flexion resistance is activated at heel strike, at the beginning of the stance phase.
- ② The stance phase flexion resistance is deactivated when the heel is relieved during the stance phase.

## 7.4.5 Exercises and adjustments

### Recommended sequence of the exercises and settings (21 – page: 12)

- ▶ **CAUTION!** If the settings cannot be fully adjusted to the comfort needs of the patient, the settings must be established primarily based on safety aspects. The extension position has to be reached with every step and at any walking speed.
- ▶ Adapt the settings of the prosthetic knee joint to the patient by means of fine-tuning and exercises.
- ▶ Check the adjustment of the prosthesis settings during normal consultations and the annual safety inspection.  
Advise the patient to check the prosthesis if changes in function occur.
- ▶ Observe the recommended sequence of the exercises and settings.

### Sitting down (22 – page: 13)

- ▶ Practise sitting down with the patient.

### Walking (23 – page: 13)

- ▶ Adjust all settings in small increments.
- ▶ Check the effect on the gait pattern after each change.
- ▶ Adjust the settings first at normal walking speed, then with short, quick steps and finally with long, fast steps.
- ▶ **CAUTION!** Allow the patient to walk only with safeguards in place, as the prosthetic knee joint may switch to swing phase flexion resistance too late with the factory settings.  
If stance phase resistance cannot be activated at the beginning of the stance phase, decrease the switching threshold in small increments until it can be activated even with different walking speeds and stride lengths.
- ▶ Adjust the swing phase flexion resistance so the lower leg of the prosthesis does not swing through too far in the dorsal direction but reaches full extension in time for the next heel strike.
- ▶ Adjust the swing phase extension resistance so the prosthetic knee joint does not swing too hard against the extension stop but reaches full extension in time for the next heel strike.
- ▶ If necessary, install the additional supplied spring according to figure 24 on page: 14 to enhance the extension assist function so the prosthetic knee joint reaches full extension at the next heel strike.

### Walking down ramps (25 – page: 14)

- ▶ Adjust the stance phase flexion resistance so it provides the patient with adequate safety yet does not generate excessive resistance.

### Walking down stairs (26 – page: 15)

- ▶ Adjust the stance phase flexion resistance so it provides the patient with adequate safety yet does not generate excessive resistance.

### Walking (repeat for readjustment) (27 – page: 15)

- ▶ Check the swing phase flexion resistance again and, if necessary, adjust it in small increments so the lower leg of the prosthesis does not swing through too far in the dorsal direction and reaches full extension in time for the next heel strike.

### Using the lock (28 – page: 16)

- ▶ Practise using the lock.

### Cycling (29 – page: 16)

- ▶ If the patient wants to participate in cycling or similar sports, perform corresponding exercises (e.g. on an exercise machine) with them. The lock must be deactivated for these sports.

## 7.5 Attaching the cosmetic cover

The illustration shows the recommended foam cover. (30 – page: 17)

- ▶ Note the following points for attaching the cosmetic cover:
  - Check the settings from the dynamic fitting.

### **⚠ CAUTION**

#### **Use of highly hygroscopic particles (grease-absorbent substances such as talcum)**

Risk of injury, damage to the product due to lack of lubrication

- ▶ Do not allow the product to come into contact with highly hygroscopic particles.

- ▶ To reduce friction and to eliminate noise, apply 519L5 Silicone Spray directly onto the contact surfaces of the cosmetic foam cover.

## **7.6 Finishing the prosthesis**

- ▶ **CAUTION!** [31] – page: 17 – To avoid product damage and the risk of falling, finish the prosthesis by replacing set screws that are too short or too long and by tightening all screw connections. In doing so, observe the instructions for use for all prosthetic components regarding assembly tightening torques and thread lock.

## **8 Use**

### **8.1 Information for use**

### **⚠ CAUTION**

#### **Overheating of the hydraulics due to excessive strain**

Burns, injuries caused by falling due to changes in functionality and damage to prosthetic components

- ▶ Do not touch overheated prosthetic components.
- ▶ Reduce all activities if there are changes in functionality, in order to allow the overheated prosthetic components to cool down.
- ▶ If overheating or changes in functionality occur, have the prosthesis checked by authorised qualified personnel.

### **NOTICE**

#### **Mechanical overload**

Impaired functionality due to mechanical damage

- ▶ Check the product for damage prior to each use.
- ▶ Do not use the product if its functionality has been impaired.
- ▶ Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

### **INFORMATION**

#### **Leak in the product's hydraulic system**

Environmental damage or skin irritation due to leaking hydraulic oil

- ▶ Before each use, check the product for leaking hydraulic oil.
- ▶ In case of leaks, do not continue using the product and have it repaired promptly.

## **8.2 Cleaning**

### **⚠ CAUTION**

#### **Use of unsuitable cleaning agents or disinfectants**

Impairment of functionality and damage due to incorrect cleaning agents or disinfectants

- ▶ Only clean the product with the approved cleaning agents.
- ▶ Only disinfect the product with the approved disinfectants.
- ▶ Observe the instructions for cleaning and care.

## Permissible cleaning agents and disinfectants

### Requirements:

Free of solvents, chlorine and phosphates

### INFORMATION

- ▶ Clean the product if it gets dirty.
- ▶ Observe the instructions for cleaning and drying in the section “Environmental conditions” (see page 30).
- ▶ Only use cleaning agents and disinfectants that do not affect the materials of the product. Test the chosen agent for material compatibility in an inconspicuous location.
- ▶ Do not spray the product directly with a pump spray or propellant spray cleaner and disinfectant.
- ▶ Only apply cleaners and disinfectants using a soft, lint-free cloth.
- ▶ Dry with a soft, lint-free cloth. Allow to air dry in order to remove residual moisture.
- ▶ Observe the cleaning instructions for all prosthetic components.

## 9 Maintenance

### CAUTION

#### Failure to follow the maintenance instructions

Risk of injuries due to changes in or loss of functionality and damage to the product

- ▶ Observe the following maintenance instructions.

- ▶ **NOTICE! Do not lubricate and grease the prosthetic joint.**
- ▶ **NOTICE! Repair work must be performed exclusively by manufacturer service.**
- ▶ Arrange regular maintenance intervals with the patient depending on the level of use.
- ▶ Following an individual period for the patient to get accustomed to the prosthesis, check the settings of the prosthetic joint and adapt them to the patient's requirements again as needed.
- ▶ The prosthetic components should be inspected after the first 30 days of use.
- ▶ Inspect the entire prosthesis for wear during normal consultations.
- ▶ Conduct annual safety inspections.
- ▶ As part of the safety inspections, inspect the prosthetic joint for wear and proper functionality. Special attention should be paid to movement resistance, bearings and abnormal noises. Full flexion and extension must be ensured at all times. Readjust as required.

## 10 Disposal

In some jurisdictions it is not permissible to dispose of the product with unsorted household waste. Improper disposal can be harmful to health and the environment. Observe the information provided by the responsible authorities in your country regarding return, collection and disposal procedures.

## 11 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

### 11.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

## 11.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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