





50K13 Genu Arexa



DE Gebrauchsanweisung	5
EN Instructions for use	10
FR Instructions d'utilisation	15
IT Istruzioni per l'uso	21
ES Instrucciones de uso	26
PT Manual de utilização	32
NL Gebruiksaanwijzing	37
SV Bruksanvisning	43
DA Brugsanvisning	48
FI Käyttöohje	53
PL Instrukcja użytkowania	59
HU Használati utasítás	64
CS Návod k použití	70
SK Návod na používanie	75
BG Инструкция за употреба	80
TR Kullanma talimatı	86
EL Οδηγίες χρήσης	91
RU Руководство по применению	97
JA 取扱説明書	103
ZH 使用说明书	108

	Seite/ Side	Größe/ Size	Umfang/Circumference			
	A (cm)	A (inch)	B (cm)	B (inch)		
	L/R	S	32-37	12,6-14,6	39-46	15,4-18,1
	L/R	M	37-41	14,6-16,1	46-53	18,1-20,9
	L/R	L	41-46	16,1-18,1	53-61	20,9-24,0
	L/R	XL	46-50	18,1-19,7	61-69	24,0-27,2
L/R	XXL	50-56	19,7-22,1	69-79	27,2-31,1	

Material	AL, SS, klettbarer Velours aus PA/velastic velour of PA, PU-Schaum/PU-foam, EVA-Schaum/EVA-foam, PA, POM, TPU, PP
-----------------	---

INFORMATION

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INFORMATION

Date of last update: 2021-09-22

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

These instructions for use provide important information on the fitting and application of the 50K13 Genu Arexa knee joint brace.

2 Product description

see fig. 1

Item	Component	Item	Component
1	Strap guide loop	6	Pad retainer
2	Pad button	7	Extension stop (square grip)
3	Clip rivet	8	Flexion stop (round grip)
4	Tibia pad	-	Straps 1 through 5
5	Condyle pad	-	Strap 6 (CI version, optional ACL/PCL version)

3 Intended use

3.1 Indications for use

The brace is intended **exclusively** for orthotic fittings of the lower limbs and **exclusively** for contact with intact skin.

The brace must be used in accordance with the indications.

3.2 Indications

- Conservative treatment of injuries (e.g. after distension, partial ruptures, ruptures) of the anterior and/or posterior cruciate ligament and/or the medial and/or lateral collateral ligament of the knee joint
- Functional follow-up treatment after ligament surgery of the anterior and/or posterior cruciate ligament and/or the medial and/or lateral collateral ligament of the knee joint
- Conservative therapy of chronic complex knee joint instability

The indication must be determined by the physician.

3.3 Contraindications

3.3.1 Absolute Contraindications

None known.

3.3.2 Relative Contraindications

The following indications require consultation with a physician: skin diseases/injuries; inflammation; prominent, swollen scars; reddening and hyperthermia of the treated limb; pronounced varicose veins, especially with impaired return flow; lymphatic flow disorders, including unclear soft tissue swelling distal to the body area where the device will be applied; sensory and circulatory disorders in the legs, e. g. associated with diabetic neuropathy.

3.4 Mechanism of Action

The brace relieves and stabilises the knee joint.

3.5 Lifetime

The product is designed for a maximum lifetime of **2 years**.

3.6 Qualification

Patients may be fitted with the product only by trained qualified personnel. The qualified personnel must be familiar with the handling of the various techniques, materials, machines and tools.

4 Safety

4.1 Explanation of warning symbols

Warning regarding possible risks of accident or injury.

Warning regarding possible technical damage.

4.2 General safety instructions

Use of the product on damaged skin

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ Only use the product on healthy skin.

Reuse on other persons and improper cleaning

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ The product may be used by one person only.
- ▶ Clean the product regularly.

⚠ CAUTION

Exceeding the lifetime and reuse

Risk of injury due to loss of functionality as well as damage to the product

- ▶ Ensure that the maximum lifetime is not exceeded.
- ▶ Only use the product for one person.

⚠ CAUTION

Contact with heat, embers or fire

Risk of injury (such as burns) and risk of product damage

- ▶ Keep the product away from open flames, embers and other sources of heat.

NOTICE

Contact with oils, salves, lotions or other products that contain oils or acids

Insufficient stabilization due to loss of material functionality

- ▶ Do not expose the product to oils, salves, lotions or other products that contain oils or acids.

NOTICE

Incorrect setting

Damage to the brace due to overloading of the material and improper fit of the brace due to breakage of load-bearing components

- ▶ The brace may only be adjusted by qualified personnel.
- ▶ Do not make any improper changes to the settings.

5 Use

INFORMATION

- ▶ The daily duration of use and period of application are generally determined by the physician.
- ▶ The initial fitting and application of the product must be carried out by qualified personnel according to the instructions of the treating physician.
- ▶ Consult a physician if any exceptional changes are noted (such as worsening of the complaint).

5.1 Selecting the Size

1. Measure the circumference of the knee joint (circumference A).
2. Measure the circumference of the thigh (circumference B) **15 cm** above the center of the patella.
3. Determine the size of the brace (see size chart).

5.2 Fitting and Application

⚠ CAUTION

Incorrect or excessively tight application

Risk of local pressure and constriction of blood vessels and nerves due to improper or excessively tight application

- ▶ Ensure that the product is applied properly and fits correctly.

NOTICE

Use of a worn or damaged product

Limited effectiveness

- ▶ Before each use, check the product for functional reliability and for possible wear or damage.
- ▶ Do not continue using a product that is no longer functional, or that is worn or damaged.

5.2.1 Adjusting the Range of Motion (Flexion/Extension)

Available extension stops (square grip):	0°, 10°, 20°, 30°, 45°
Available flexion stops (round grip):	0°, 10°, 20°, 30°, 45°, 60°, 75°, 90°

INFORMATION

The joint stops (medial and lateral) must have identical degree settings.

- 1) Remove the pad retainer (and the spacer rings if applicable) (see fig. 5).
- 2) Insert the desired extension and flexion stops into the joint (see fig. 6).
→ The stops will engage with slight resistance.
- 3) **Optional:** Snap the required number of spacer rings onto the joint (see fig. 7).
- 4) Attach the pad retainer.

5.2.2 Adjusting the strap height

The height of straps 1 and 5 can be adjusted for knee joint relief in case of post-operative therapy and/or to adapt the position of the straps to the anatomical shape of the leg.

- 1) Remove the clip rivet of the strap holder for straps 1 and 5 on both sides.
- 2) Reposition the strap holder.
- 3) Attach the strap holder to the brace frame with the supplied new clip rivets.

5.2.3 Applying the ACL Version

- > **Prerequisite:** The brace is in the delivery condition (ACL version) (see fig. 2).
 - > **Prerequisite:** The patient is sitting on the front edge of a chair and bends the knee at **45°**.
- 1) Position the brace with the middle of the joints approximately even with the top of the patella.
→ Although this puts the brace higher than its final position, it settles into the correct position on its own.
 - 2) Fasten strap 1 (see fig. 2).
 - 3) Fasten strap 2.
 - 4) Engage strap 3 on the pad button with the strap guide loop and tighten the strap.
→ The strap guide loop engages audibly.
 - 5) Tighten strap 4.
 - 6) Push the brace towards the thigh.
 - 7) Fasten strap 5.
 - 8) Tighten the brace straps so that the brace fits firmly but is comfortable to wear.
 - 9) After **15 min** of active movement (e.g. slight knee bends): if necessary, slightly retighten the straps in the specified sequence (see fig. 2, straps 1-5).

Optional: Use the included strap 6 if the knee joint requires more stabilization. Open the hook-and-loop closure of the strap holder for strap 6 on both sides and mount the strap on the front.

5.2.4 Conversion to and application of the PCL version

- > **Prerequisite:** The brace is in the delivery condition (ACL version) and has to be converted to the PCL version (see fig. 3).
- 1) Turn the strap holder for strap 2 from the rear to the front of the brace and align the strap (see fig. 3).

- 2) **Optional:** Use the included strap 6 if the knee joint requires more stabilization. In doing so, strap 2 remains on the front of the brace and strap 6 is attached with hook-and-loop to the strap holder on both sides on the back of the brace.
- 3) Put on the brace and fasten the numbered straps in sequence (see fig. 3, straps 1-5, **optional** strap 6).

5.2.5 Conversion to and application of the CI version

- > **Prerequisite:** The brace is in the delivery condition (ACL version) and has to be converted to the CI version (see fig. 4).
- 1) Attach the supplied strap 6 to the strap holder on both sides on the front of the brace with hook-and-loop (see fig. 4).
 - 2) Put on the brace and fasten the numbered straps 1-6 in sequence (see fig. 4).

5.3 Delivery

NOTICE

Incorrect setting

Damage to the brace due to overloading of the material and improper fit of the brace due to breakage of load-bearing components

- ▶ The brace may only be adjusted by qualified personnel.
 - ▶ Do not make any improper changes to the settings.
- ▶ Check that the brace fits the leg properly at the time of delivery.
→ The centre of the patella and the centre of the joints are approximately aligned.

6 Cleaning

NOTICE

Use of improper cleaning agents

Damage to the product due to use of improper cleaning agents

- ▶ Only clean the product with the approved cleaning agents.

INFORMATION

- ▶ Ensure proper alignment when inserting the tibia pad. The wide side of the pad must face up.
 - ▶ **Allowable cleaning agent:** pH neutral soap (e.g. 453H10 Derma Clean)
- 1) Remove the straps and pads from the brace.
 - 2) Fasten all hook-and-loop closures.
 - 3) Hand-wash the straps and pads in warm water at **30 °C (86 °F)** using neutral detergent. Rinse thoroughly.
 - 4) Allow to air dry. Do not expose to direct heat sources (e.g. sunlight, stove, or radiator).
 - 5) Attach the straps and pads to the brace.
- ▶ **If the brace comes into contact with salt water or dirt:** Rinse the brace with clear water and allow to air dry.

7 Information on re-use (applies for Germany only)

The product is approved for re-use **exclusively** in Germany. Products that are being re-used are subject to increased strain. Re-use is only permitted up to a maximum total period of use of **2 years**.

Reconditioning

Decontamination of the product (according to the guidelines of the Robert Koch Institute) is required prior to reuse.

All components subject to heavy wear and tear or that come into direct contact with the skin are not suitable for reuse. These components are contained in the 29K110-2 replacement set.

All reusable components have to be inspected for function and damage prior to reuse.

7.1 Replacing the Straps and Pads

Straps:

- 1) Release straps 3 and 4 from the frame with the quick-release closures and replace them.
- 2) Loosen straps 1, 2, 5 and 6 from the strap guide loops and replace them. Note section 5.2.2 for the adjustment of straps 1 and 5.

Strap and Frame Pads

The strap and frame pads are attached with hook-and-loop closures and are easy to replace.

8 Disposal

Dispose of the product in accordance with national regulations.

9 Legal Information

All legal conditions are subject to the respective national law of the country of use and may vary accordingly.

9.1 Local Legal Information

Legal information that applies **exclusively** to specific countries is written in the official language of the respective country of use after this chapter.

9.2 Liability

The manufacturer shall be liable in the event that the product is used in accordance with the descriptions and instructions in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorized modification of the product.

9.3 CE Conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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