## 50K15 Genu Sensa

<b>DE</b>	Gebrauchsanweisung .....	4
<b>EN</b>	Instructions for use .....	8
<b>FR</b>	Instructions d'utilisation .....	12
<b>IT</b>	Istruzioni per l'uso .....	16
<b>ES</b>	Instrucciones de uso .....	21
<b>PT</b>	Manual de utilização .....	25
<b>NL</b>	Gebruiksaanwijzing .....	29
<b>SV</b>	Bruksanvisning .....	33
<b>DA</b>	Brugsanvisning .....	37
<b>NO</b>	Bruksanvisning .....	41
<b>FI</b>	Käyttöohje .....	45
<b>PL</b>	Instrukcja użytkowania .....	50
<b>HU</b>	Használati utasítás .....	54
<b>CS</b>	Návod k použití .....	58
<b>HR</b>	Upute za uporabu .....	62
<b>SL</b>	Navodila za uporabo .....	66
<b>SK</b>	Návod na používanie .....	69
<b>BG</b>	Инструкция за употреба .....	73
<b>TR</b>	Kullanma talimatı .....	78
<b>EL</b>	Οδηγίες χρήσης .....	81
<b>RU</b>	Руководство по применению .....	86
<b>ZH</b>	使用说明书 .....	90

	Größe Size	Umfang Circumference					
		A		B		B V-MAX	
		cm	inch	cm	inch	cm	inch
	XXS	25-28	9.1-11.0	35-38	13.8-15.0	-	-
	XS	28-31	11.0-12.2	38-41	15.0-16.1	-	-
	S	31-34	12.2-13.4	41-44	16.1-17.3	-	-
	M	34-37	13.4-14.6	44-47	17.3-18.5	47-50	18.5-19.7
	L	37-40	14.6-15.8	47-50	18.5-19.7	50-53	19.7-20.9
	XL	40-43	15.8-16.9	50-53	19.7-20.9	53-56	20.9-22.0
	XXL	43-46	16.9-18.1	53-56	20.9-22.0	56-59	22.0-23.2

<b>Material</b>	Gestrick / Knit Fabric: PA, EL
	Frotteebezüge / Terry Cloth Covers: Baumwolle / Cotton, PES
	Spiralfedern / Spiral Springs: verzinkter Federstahl / zinc-coated Spring Steel, PE
	Pelotte / Pad: viskoelastisches Technogel® / viscoelastic Technogel®

## INFORMATION

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### INFORMATION

Date of last update: 2021-12-08

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

These instructions for use provide important information on the fitting and application of the 50K15 Genu Sensa knee support.

## 2 Intended use

### 2.1 Indications for use

The support is intended **exclusively** for fittings of the lower limbs and **exclusively** for contact with intact skin.

The support must be used according to the indication.

### 2.2 Indications

- Feeling of instability
- Irritation (postoperative, posttraumatic and chronic)
- Knee pain due to stress (tendomyopathy, insertion ligamentosis)
- Joint effusion
- Swelling

Indications must be determined by the physician.

## 2.3 Contraindications

### 2.3.1 Absolute Contraindications

Not known.

### 2.3.2 Relative Contraindications

The following indications require consultation with a physician: skin diseases/injuries, inflammation, prominent scars that are swollen, reddening and hyperthermia of the fitted limb; pronounced varicose veins, especially with return flow impairment, lymphatic flow disorders, including unclear soft tissue swelling distant to the body area to which the support will be applied; sensory and circulatory disorders in the area of the lower leg and the foot, e.g. in case of diabetic neuropathy.

## 2.4 Lifetime

The product is designed for a maximum lifetime of 6 months.

## 2.5 Effects

The form-knitted support promotes sensorimotor functions and thereby activates muscular stabilisation of the knee joint. The support with integrated shape-retaining pad offers relief for the kneecap.

## 3 Safety

### 3.1 Explanation of warning symbols

 **CAUTION** Warning regarding possible risks of accident or injury.

 **NOTICE** Warning regarding possible technical damage.

### 3.2 General safety instructions

 **CAUTION**

#### **Contact with heat, embers or fire**

Risk of injury (such as burns) and risk of product damage.

- ▶ Keep the product away from open flames, embers and other sources of heat.

 **CAUTION**

#### **Reuse on other persons and improper cleaning**

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ The product may be used by one person only.

- ▶ Clean the product regularly.

#### NOTICE

#### Contact with oils, salves, lotions or other products that contain oils or acids

Insufficient stabilisation due to loss of material functionality

- ▶ Do not expose the product to oils, salves, lotions or other products that contain oils or acids.

## 4 Handling

#### INFORMATION

- ▶ The daily duration of use and period of application are generally determined by the physician.
- ▶ The initial fitting and application of the product must be carried out by qualified personnel.
- ▶ Instruct the patient in the handling and care of the product.
- ▶ Instruct the patient to see a physician immediately if any exceptional changes are noted (e.g. worsening of the complaint).

### 4.1 Size Selection

- 1) Measure the circumference of the thigh 14 cm above the centre of the patella.
- 2) Measure the circumference of the lower leg 12 cm below the centre of the patella.
- 3) Determine the size of the support (see size chart).

### 4.2 Application

#### CAUTION

#### Incorrect or excessively tight application

Skin irritation due to overheating, local pressure points due to tight application to the body

- ▶ Ensure that the product is applied properly and fits correctly.
- ▶ When in doubt, do not continue to use the product if skin irritation occurs.

## NOTICE

### Use of a worn or damaged product

Limited effectiveness

- ▶ Before each use, check the product for functional reliability and for possible wear or damage.
- ▶ Do not continue using a product that is no longer functional, or that is worn or damaged.

- 1) Grasp the support with both hands on the side and on the upper end of the spiral springs.
- 2) Pull the support over the knee (see fig. 1).
- 3) Position the knee pad so that the kneecap is surrounded without pressure (see fig. 2).
- 4) Correct the position of the knee pad if required (see fig. 3).
- 5) When providing the support to the patient, the proper fit and positioning of the support must be verified.

### 4.3 Removal

- 1) Grasp the support with both hands on the side and on the lower end of the spiral springs.
- 2) Pull the support down over the foot.

### 4.4 Cleaning

## NOTICE

### Use of improper cleaning agents

Damage to the product due to use of improper cleaning agents

- ▶ Only clean the product with the approved cleaning agents.

Clean the support regularly:

- 1) **Recommendation:** Use a laundry bag or net.
- 2) Wash the support in warm water at 30 °C with standard mild detergent. Do not use fabric softener. Rinse thoroughly.
- 3) Allow to air dry. Do not expose to direct heat (e.g. sunshine, stove or radiator).

### 5 Disposal

Dispose of the product in accordance with national regulations.

## 6 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

### 6.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

### 6.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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