

## 50R236 Smartspine SI Support

<b>DE</b> Gebrauchsanweisung .....	4
<b>EN</b> Instructions for use .....	7
<b>FR</b> Instructions d'utilisation .....	10
<b>IT</b> Istruzioni per l'uso .....	14
<b>ES</b> Instrucciones de uso .....	17
<b>PT</b> Manual de utilização .....	20
<b>NL</b> Gebruiksaanwijzing .....	24
<b>SV</b> Bruksanvisning .....	27
<b>DA</b> Brugsanvisning .....	30
<b>NO</b> Bruksanvisning .....	33
<b>FI</b> Käyttöohje .....	36
<b>IS</b> Notkunarleiðbeiningar .....	40
<b>PL</b> Instrukcja użytkowania .....	43
<b>HU</b> Használati utasítás .....	46
<b>CS</b> Návod k použití .....	50
<b>RO</b> Instrucțiuni de utilizare .....	53
<b>HR</b> Upute za uporabu .....	56
<b>SK</b> Návod na používanie .....	59
<b>BG</b> Инструкция за употреба .....	62
<b>TR</b> Kullanma talimatı .....	66
<b>RU</b> Руководство по применению .....	69
<b>JA</b> 取扱説明書 .....	73
<b>ZH</b> 使用说明书 .....	76
<b>KO</b> 사용 설명서 .....	78

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### INFORMATION

Date of last update: 2021-12-08

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

These instructions for use provide important information on the fitting and application of the 50R236 Smartspine SI sacroiliac support.

## 2 Intended use

### 2.1 Indications for use

The support is intended **exclusively** for orthotic fittings of the pelvis and **exclusively** for contact with intact skin.

The support must be used according to the indication.

### 2.2 Indications

- Sciatica
- Irritation of the sacroiliac joint
- Instability of the sacroiliac joint
- Symphysiolysis

The indication must be determined by the physician.

## 2.3 Contraindications

### 2.3.1 Absolute Contraindications

None known.

### 2.3.2 Relative Contraindications

The following indications require consultation with a physician: skin diseases/injuries, inflammation, raised scars with swelling, reddening, and warming in the area where the support is fitted, tumors, lymphatic congestion – including unclear soft tissue swelling distal to the support, and dysesthesia in the pelvis.

## 2.4 Lifetime

The product is designed for a maximum lifetime of **1 year**.

## 2.5 Mechanism of Action

The support stabilizes the pelvic girdle by immobilizing the sacroiliac joints and symphysis, and can alleviate pain.

Pads that can be individually positioned promote blood flow and reduce cramping.

The applied compression can be individually controlled with the Mechanical Advantage Pulley System.

## 3 Safety

### 3.1 Explanation of warning symbols

**CAUTION** Warning regarding possible risks of accident or injury.

**NOTICE** Warning regarding possible technical damage.

### 3.2 General safety instructions

**CAUTION**

#### Contact with heat, embers or fire

Risk of injury (such as burns) and risk of product damage

- ▶ Keep the product away from open flames, embers and other sources of heat.

**CAUTION**

#### Reuse on other persons and improper cleaning

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ The product may be used by one person only.
- ▶ Clean the product regularly.

**NOTICE**

#### Improper use and changes

Change in or loss of functionality as well as damage to the product

- ▶ Use this product with care and only for its intended purpose.
- ▶ Do not make any improper changes to the product.

**NOTICE**

#### Contact with oils, salves, lotions or other products that contain oils or acids

Insufficient stabilization due to loss of material functionality

- ▶ Do not expose the product to oils, salves, lotions or other products that contain oils or acids.

## 4 Handling

### INFORMATION

- ▶ The daily duration of use and period of application are generally determined by the physician.
- ▶ The initial fitting and application of the product must be carried out by qualified personnel.
- ▶ Instruct the patient in the handling and care of the product.
- ▶ Instruct the patient to see a physician immediately if any exceptional changes are noted (e.g. worsening of the complaint).

### 4.1 Selecting the Size

- ▶ Select the support size based on the circumference of the pelvis (see size chart).

### 4.2 Fitting and Application

#### ⚠ CAUTION

##### **Incorrect or excessively tight application**

Skin irritation due to overheating, local pressure points due to tight application to the body

- ▶ Ensure that the product is applied properly and fits correctly.
- ▶ When in doubt, do not continue to use the product if skin irritation occurs.

#### ⚠ CAUTION

##### **Direct skin contact with the product**

Skin irritation due to chafing or sweat formation

- ▶ Do not wear this product directly on the skin.

### NOTICE

##### **Use of a worn or damaged product**

Limited effectiveness

- ▶ Before each use, check the product for functional reliability and for possible wear or damage.
  - ▶ Do not continue using a product that is no longer functional, or that is worn or damaged.
- > The patient is standing.
  - > The hook-and-loop closure of the abdominal strap is open.
  - > The pull cable is completely loosened.
- 1) **Optional – inserting the pad:** Use the hook-and-loop closure to attach the pad to the inside of the support (see fig. 1).
  - 2) Guide the support around the torso and apply equal force to both sides to pull it forward (see fig. 2).
  - 3) Attach the right side to the left side of the abdominal strap (see fig. 3). Ensure that the hook-and-loop sections are lying flat against each other and are not in contact with clothing or skin.
  - 4) Pull the end of the pull cable forward and attach it to the abdominal strap with the hook-and-loop fastener (see fig. 4).
  - 5) **Optional – shortening the pull cable:** Pull the winding card out of the pocket (see fig. 5), wrap the pull cable around the winding card (see fig. 6) and put the winding card back into the pocket.
  - 6) **Optional – changing the direction of the pull cable from left to right:** The pull cable can be changed from the left to the right side by rotating the brace by 180° when putting it on.

## 4.3 Cleaning

### NOTICE

#### Use of improper cleaning agents

Damage to the product due to use of improper cleaning agents

- ▶ Only clean the product with the approved cleaning agents.

Clean the support regularly:

- 1) Fasten all hook-and-loop closures.
- 2) Hand wash the support in warm water at **30 °C** (86 °F) with standard mild detergent.
- 3) Rinse the support well.
- 4) Allow to air dry. Do not expose to direct heat sources (e.g. sunlight, stove or radiator).

## 5 Disposal

Dispose of the product in accordance with national regulations.

## 6 Legal Information

All legal conditions are subject to the respective national law of the country of use and may vary accordingly.

### 6.1 Local Legal Information

Legal information that applies **exclusively** to specific countries is written in the official language of the respective country of use after this chapter.

### 6.2 Liability

The manufacturer shall be liable in the event that the product is used in accordance with the descriptions and instructions in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorized modification of the product.

### 6.3 CE Conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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