
50R231 Smartspine Extension Brace

DE Gebrauchsanweisung	5
EN Instructions for use	8
FR Instructions d'utilisation	12
IT Istruzioni per l'uso	16
ES Instrucciones de uso	19
PT Manual de utilização	23
NL Gebruiksaanwijzing	27
SV Bruksanvisning	31
DA Brugsanvisning	34
NO Bruksanvisning	38
FI Käyttöohje	41
IS Notkunarleiðbeiningar	45
PL Instrukcja użytkowania	48
HU Használati utasítás	52
CS Návod k použití	56
RO Instrucțiuni de utilizare	59
HR Upute za uporabu	63
SK Návod na používanie	67
BG Инструкция за употреба	70
TR Kullanma talimatı	74
RU Руководство по применению	78
JA 取扱説明書	82
ZH 使用说明书	85
KO 사용 설명서	88

	Größe / Size	Taillenumfang* / Waist circumference	
		cm	inch
	S	76,2–88,9	30–35
	M	88,9–101,6	35–40
	L	101,6–114,3	40–45
XL	114,3–127	45–50	

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INFORMATION

Date of last update: 2021-12-08

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

These instructions for use provide important information on the fitting and application of the 50R231 Smartspine Extension thoracolumbosacral brace.

2 Intended use

2.1 Indications for use

The brace is intended **exclusively** for orthotic fittings of the back and **exclusively** for contact with intact skin.

The brace must be used in accordance with the indications.

2.2 Indications

- Thoracic spine pain (mechanical)
- Thoracic kyphosis (hunchback) (without fused vertebrae, mechanical)
- Compression fractures of the thoracic spine (with no neurological deficits)
- Degenerative changes of the thoracic spine
- Facet syndrome of the thoracic spine
- Osteochondrosis of the thoracic spine
- Osteoporosis
- Osteoporotic fractures of the thoracic and/or lumbar spine

The indication must be determined by the physician.

2.3 Contraindications

2.3.1 Absolute Contraindications

None known.

2.3.2 Relative Contraindications

The following conditions require consultation with a physician: skin diseases/injuries, inflammation, raised scars with swelling, reddening, and warming in the area where the device is fitted, tumours, lymphatic congestion – including unclear soft tissue swelling distal to the device, and dysaesthesia in the trunk.

2.4 Lifetime

The product is designed for a maximum lifetime of **1 year**.

2.5 Mechanism of Action

The brace supports muscular stabilization of the thoracic spine and straightens it actively. It promotes sensorimotor function, helps alleviate pain and relieves the thoracic spine through active straightening and muscular stabilization and the lumbar spine by increasing intra-abdominal pressure and providing external stabilization.

The applied compression can be individually controlled with the Mechanical Advantage Pulley System.

3 Safety

3.1 Explanation of warning symbols

 CAUTION	Warning regarding possible risks of accident or injury.
---	---

 NOTICE	Warning regarding possible technical damage.
--	--

3.2 General safety instructions

 CAUTION

Contact with heat, embers or fire

Risk of injury (such as burns) and risk of product damage

- ▶ Keep the product away from open flames, embers and other sources of heat.

⚠ CAUTION

Reuse on other persons and improper cleaning

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ The product may be used by one person only.
- ▶ Clean the product regularly.

NOTICE

Improper use and changes

Change in or loss of functionality as well as damage to the product

- ▶ Use this product with care and only for its intended purpose.
- ▶ Do not make any improper changes to the product.

NOTICE

Contact with oils, salves, lotions or other products that contain oils or acids

Insufficient stabilization due to loss of material functionality

- ▶ Do not expose the product to oils, salves, lotions or other products that contain oils or acids.

4 Handling

INFORMATION

- ▶ The daily duration of use and period of application are generally determined by the physician.
- ▶ The initial fitting and application of the product must be carried out by qualified personnel.
- ▶ Instruct the patient in the handling and care of the product.
- ▶ Instruct the patient to see a physician immediately if any exceptional changes are noted (e.g. worsening of the complaint).

4.1 Selecting the Size

- 1) Measure the waist circumference approx. **2 cm** above the iliac crest.
- 2) Determine the brace size (see size chart).

4.2 Adaptation

- > The patient is standing.
 - > The hook-and-loop closures on the abdominal strap are open.
 - > The pull cable is fully loosened.
- 1) Measure the upper part of the brace over the shoulders like a backpack.
 - 2) Center the back plate of the brace (see fig. 1).
 - 3) **Optional:** Adjust the height position of the abdominal strap. In order to do so, detach the Mechanical Advantage Pulley System from the hook-and-loop fastener on the back plate and reposition it (see fig. 2).

INFORMATION: Ensure that the position of the hook-and-loop surface is the same on both sides.

- 4) Guide the brace around the torso and apply equal force to both sides to pull it forward.
- 5) Attach the right side to the left side of the abdominal strap (see fig. 3). Ensure that the hook-and-loop sections are lying flat against each other and are not in contact with clothing or skin.
- 6) Pull the end of the pull cable forward and attach it to the abdominal strap with the hook-and-loop fastener.
- 7) Pull the ends of the straps forward (see fig. 4).
- 8) **Optional:** Shorten the straps to the desired length (see fig. 5) and attach the Y-hook-and-loop to the end of the straps (see fig. 6).

- 9) Attach the ends of the straps to the abdominal strap with the hook-and-loop fasteners (see fig. 7).

4.3 Application

CAUTION

Direct skin contact with the product

Skin irritation due to chafing or sweat formation

- ▶ Do not wear this product directly on the skin.

CAUTION

Incorrect or excessively tight application

Risk of local pressure and constriction of blood vessels and nerves due to improper or excessively tight application

- ▶ Ensure that the product is applied properly and fits correctly.

NOTICE

Use of a worn or damaged product

Limited effectiveness

- ▶ Before each use, check the product for functional reliability and for possible wear or damage.
- ▶ Do not continue using a product that is no longer functional, or that is worn or damaged.

- > The patient is standing.
 - > The hook-and-loop closures on the abdominal strap are open.
 - > The pull cable is fully loosened.
- 1) Place the upper part of the brace over the shoulders like a backpack (see fig. 8).
 - 2) Attach the right side to the left side of the abdominal strap (see fig. 9). Ensure that the hook-and-loop sections are lying flat against each other and are not in contact with clothing or skin.
 - 3) Pull the ends of the straps forward and attach them to the abdominal strap with the hook-and-loop fasteners (see fig. 10).
 - 4) Pull the end of the pull cable forward and attach it to the abdominal strap with the hook-and-loop fastener (see fig. 11).

4.4 Cleaning

NOTICE

Use of improper cleaning agents

Damage to the product due to use of improper cleaning agents

- ▶ Only clean the product with the approved cleaning agents.

Clean the brace regularly:

- 1) Remove all boards/plates.
- 2) Fasten all hook-and-loop closures.
- 3) Hand wash the brace in warm water at **30 °C** with standard mild detergent.
- 4) Rinse well.
- 5) Allow to air dry. Do not expose to direct heat sources (e.g. sunlight, stove or radiator).

5 Disposal

Dispose of the product in accordance with national regulations.

6 Legal Information

All legal conditions are subject to the respective national law of the country of use and may vary accordingly.

6.1 Local Legal Information

Legal information that applies **exclusively** to specific countries is written in the official language of the respective country of use after this chapter.

6.2 Liability

The manufacturer shall be liable in the event that the product is used in accordance with the descriptions and instructions in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorized modification of the product.

6.3 CE Conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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