## 1C20 ProSymes

<b>DE</b> Gebrauchsanweisung .....	5
<b>EN</b> Instructions for use .....	12
<b>FR</b> Instructions d'utilisation .....	19
<b>IT</b> Istruzioni per l'uso .....	26
<b>ES</b> Instrucciones de uso .....	33
<b>PT</b> Manual de utilização .....	40
<b>NL</b> Gebruiksaanwijzing .....	47
<b>SV</b> Bruksanvisning .....	54
<b>DA</b> Brugsanvisning .....	61
<b>NO</b> Bruksanvisning .....	67
<b>FI</b> Käyttöohje .....	74
<b>PL</b> Instrukcja użytkowania .....	81
<b>HU</b> Használati utasítás .....	88
<b>CS</b> Návod k použití .....	95
<b>RO</b> Instrucțiuni de utilizare .....	101
<b>HR</b> Upute za uporabu .....	108
<b>SL</b> Navodila za uporabo .....	115
<b>SK</b> Návod na používanie .....	122
<b>BG</b> Инструкция за употреба .....	128
<b>TR</b> Kullanma talimatı .....	135
<b>EL</b> Οδηγίες χρήσης .....	142
<b>RU</b> Руководство по применению .....	150
<b>JA</b> 取扱説明書 .....	157
<b>ZH</b> 使用说明书 .....	164
<b>KO</b> 사용 설명서 .....	169

1


2


3

4


11

12

13

## INFORMATION

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### INFORMATION

Date of last update: 2022-01-24

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

### 1.1 Construction and Function

The 1C20 ProSymes prosthetic foot is used for the fitting of Syme ankle amputations, under certain conditions for Pirogoff amputations, and for modified Chopart residual limbs. The prosthetic foot features a lamination anchor and is connected directly to the prosthetic socket. Thanks to its low structural height, it is particularly suited to long residual limbs with low ground clearance. The functional properties of the prosthetic foot are achieved by a design that features double carbon spring elements with a PU elastomer layer in between.

Characteristics of the prosthetic foot	
Structural height (with lamination anchor and footshell)	43 mm or 52 mm (with spacer plate)
Plantar flexion and dorsiflexion	±6 °
Abduction and adduction	±6 °
Rotation	±30 °

### 1.2 Combination possibilities

This prosthetic component is compatible with Ottobock's system of modular connectors. Functionality with components of other manufacturers that have compatible modular connectors has not been tested.

## 2 Intended use

### 2.1 Area of application

Our components perform optimally when paired with appropriate components based upon weight and mobility grades identifiable by our MOBIS classification information and which have appropriate modular connectors.

### Size 25 stiffness 1, size 26 to 28 stiffness 2

The product is recommended for mobility grade 2 (restricted outdoor walker) and mobility grade 3 (unrestricted outdoor walker). Approved for a body weight of up to **100 kg (220 lbs)**.

### Size 26 to 28 stiffness 3

The product is recommended for mobility grade 2 (restricted outdoor walker) and mobility grade 3 (unrestricted outdoor walker). Approved for a body weight of **up to 125 kg (275 lbs)**.

## 2.2 Environmental conditions

### Allowable environmental conditions

Temperature range for use: -10 °C to +60°C

Allowable relative humidity 0 % to 90 %, non-condensing

### Unallowable environmental conditions

Mechanical vibrations or impacts

Perspiration, urine, fresh water, salt water, acids

Dust, sand, highly hygroscopic particles (e.g. talcum)

## 2.3 Lifetime

The product was tested by the manufacturer with 2 million load cycles. Depending on the user's activity level, this corresponds to a maximum lifetime of 3 years.

## 3 Safety

### 3.1 Explanation of warning symbols

**CAUTION** Warning regarding possible risks of accident or injury.

**NOTICE** Warning regarding possible technical damage.

### 3.2 General safety instructions

#### Risk of injury and risk of product damage

- ▶ Comply with the product's field of application and do not expose it to excessive strain (see page 12).
- ▶ Note the combination possibilities/combination exclusions in the instructions for use of the products.
- ▶ Do not expose the product to prohibited environmental conditions.
- ▶ Check the product for damage if it has been exposed to prohibited environmental conditions.
- ▶ Do not use the product if it is damaged or in a questionable condition. Take suitable measures (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop).
- ▶ To avoid the risk of injury and product damage, do not use the product beyond the tested lifetime.

- ▶ To avoid the risk of injury and product damage, only use the product for a single patient.
- ▶ To prevent mechanical damage, use caution when working with the product.
- ▶ If you suspect the product is damaged, check it for proper function and readiness for use.
- ▶ Do not use the product if its functionality is restricted. Take suitable measures (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop).

### Signs of changes in or loss of functionality during use

Decreased forefoot resistance or changes in roll-over behaviour are noticeable indications of loss of functionality.

## 4 Scope of delivery

1C20 ProSymes			
Fig.	Quantity	Designation	Reference number
–	1	Instructions for use	–
1	1	Prosthetic foot	–
1	1	Lamination anchor with lamination cover	2G120
–	1	Protective sock	SL=Spectra-Sock

Spare parts/accessories (not included in the scope of delivery)	
Designation	Reference number
Footshell (including connection cap)	2C2
Setting aid with screw	2Z328
Screw set	2Z120

## 5 Preparing the product for use

### **⚠ CAUTION**

#### **Incorrect alignment or assembly**

Risk of injury due to damaged prosthetic components

- ▶ Observe the alignment and assembly instructions.

### **NOTICE**

#### **Grinding the prosthetic foot**

Premature wear resulting from damage to the prosthetic foot

- ▶ Do not grind the prosthetic foot.

The three-dimensional alignment of the prosthetic socket affects the function of the prosthesis. The patient can only benefit from the special characteristics of the prosthetic foot when the product is correctly aligned.

The residual limb position, e.g. abduction, adduction and flexion, must be considered when positioning the lamination anchor. This guarantees the maximum freedom of play for adjustments.

The alignment of the prosthesis is performed in the following steps:

- Preparing the prosthetic socket
- Bench alignment
- Dynamic trial fitting
- Finishing the prosthesis

## 5.1 Preparing the prosthetic socket

### 5.1.1 Fabricating the prosthetic socket

#### > Required materials:

- 99B81 PVA Bag, 623T9 Nyloglas Stockinette, 616G12 Carbon Fibre Cloth, 616G13 Fibreglass Woven Hose, 617H55 C-Orthocryl
- 1) Pull on a soaked PVA bag that is heat welded at the distal end over the model.
  - 2) Pull over two layers of Nyloglas stockinette.
  - 3) Place a strip of carbon fibre cloth on the medial and lateral sides from the distal end of the residual limb up as far as the condyles.
  - 4) Place a strip of carbon fibre cloth (width: approx. **8 cm**) circularly around the MPT (mid-patella tendon) point.
  - 5) Cover the surfaces free of carbon around the attachment of the lamination anchor using carbon fibre cloth.
  - 6) Cut a piece of fibreglass woven hose (**1.5 times** the length of the plaster positive).
  - 7) Pull the fibreglass woven hose over the model to the edge.
  - 8) Tie off the fibreglass woven hose at the distal end and pull it over the model.
  - 9) Perform the casting process with C-Orthocryl.

### 5.1.2 Trial fitting of the prosthetic socket

#### > Required tools and materials:

- 743L30 LaserLine or a plumb line
- > The patient should be standing at the height of the heel without shoes on. He/she should be standing comfortably with the knees in a normal position.
  - > Compensate for the distance from the socket to the floor.
- 1) Check the upright pelvis position and make corrections if necessary.
  - 2) Measure and note down the distance from the end of the socket to the floor.
  - 3) Mark up the position of the residual limb on the socket in the frontal and sagittal planes using the LaserLine or a plumb line (see fig. 3, see fig. 4). When doing this, draw on frontally from the middle of the patella down to the end of the residual limb.

## 5.2 Bench Alignment

### 5.2.1 Optional: Inserting the spacer plate

Inserting the spacer plate increases the structural height of the prosthetic foot from 43 mm to 52 mm.

#### > Required tools and materials:

- 710D4 Torque Wrench
- 1) **If necessary:** Remove the lamination anchor from the concave plate (see fig. 5).
  - 2) Remove the 2 cap screws from the sides of the concave plate (see fig. 6).
  - 3) Loosen the countersunk head screw on the surface of the concave plate (see fig. 7).
  - 4) Lift off the concave plate.
  - 5) Remove the short fine-threaded screw from the counter plate and replace it with the long fine-threaded screw (see fig. 8).
  - 6) Remove the washer from the short fine-threaded screw and place it on the long fine-threaded screw with the concave side facing upwards.
  - 7) Put on the concave plate (see fig. 9).
  - 8) **INFORMATION: Screw on the screws without tightening.**  
Screw in the two cap screws and the countersunk head screw fully.
  - 9) First tighten the countersunk head screw (installation torque: **46 Nm**)
  - 10) Tighten the two cap screws (installation torque: **6 Nm**)
  - 11) Place the spacer plate on the concave plate so that the dowel pin is positioned in the head of the countersunk head screw.

### 5.2.2 Installing the lamination anchor

> **Required tools and materials:**

710D4 Torque Wrench

- 1) Position the lamination anchor with the two anchor arms close to each other towards the posterior and carefully place it on the fine-threaded screw.
- 2) Allow the setting aid for the lamination anchor to click in to the guiding groove on the concave plate.
- 3) Tighten the fine-threaded screw (installation torque: **46 Nm**)

### 5.2.3 Installing the prosthetic foot

> **Required tools and materials:**

743L200 or 743L300 L.A.S.A.R. Assembly, 711S4\* Bending Iron, degreasing cleaner (e.g. 634A3 Acetone), 627B2=19 Monofilament Adhesive Tape (alternatively: 699G30=\* Cellacast)

- 1) Set the prosthetic foot in the alignment apparatus to the heel height.
- 2) Bring the prosthetic foot into outer rotation and clamp it in the alignment apparatus.
- 3) Place the socket in the lamination anchor and align it.  
Position the drawn on line in the medial third of the prosthetic foot from the frontal perspective.  
Align the posterior end of the socket in a natural line with the heel from the sagittal perspective.
- 4) Adapt the arms of the lamination anchor to the shape of the residual limb using the bending iron (see fig. 10).
- 5) Roughen the distal socket end.
- 6) Clean the lamination anchor using a degreasing cleaner.
- 7) Protect the prosthetic foot from dirt by covering it.
- 8) Create a thick mixture of sealing resin and talcum.
- 9) Apply the sealing resin and talcum mixture to the lamination anchor and the socket, evening out any small differences in height.
- 10) Position the socket on the lamination anchor, align it and leave the sealing resin and talcum mixture to harden.
- 11) Take the prosthesis out of the alignment apparatus.
- 12) Create a thick mixture of sealing resin and talcum.
- 13) Apply the resin-talcum mixture in the region of the anchor arms for support and smooth it (see fig. 11).
- 14) Leave the resin-talcum mixture to harden.
- 15) Wrap glass-fibre reinforced adhesive tape or Cellacast circularly around the area of the anchor arms and secure it.

### 5.2.4 Installing the footshell

> **Required tools and materials:**

Footshell, foam for inserting into cavities in the forefoot, sanding tool

- 1) Sand down the foam for fitting in cavities in the forefoot region.
- 2) **Optional:** Cut the footshell in the proximal area and sand it to fit the shape of the socket.
- 3) Put the footshell on the prosthetic foot until it clicks onto the heel spring and slide the foam into the forefoot area.

### 5.3 Dynamic Trial Fitting

The alignment of the prosthetic foot is checked and corrected during the dynamic trial fitting. If additional lateralisation of the prosthetic foot is necessary, the lamination anchor must be repositioned on the prosthetic socket.

> Inform the patient that the foot is only provisionally secured to the prosthetic socket.

- 1) Check the length of the prosthesis and the alignment of the prosthetic foot.
- 2) **If necessary:** Adjust the alignment of the prosthetic foot (see page 17).

### 5.3.1 Aligning the prosthetic foot

The positioning of the foot relative to the prosthetic socket can be adjusted by moving the lamination anchor. The following adjustments can be made:

- Plantar flexion/dorsiflexion:  $\pm 6^\circ$
- Rotation:  $\pm 30^\circ$
- Adduction/abduction:  $\pm 6^\circ$

> **Required tools and materials:**

710D4 Torque Wrench

- 1) Loosen the fine-threaded screw until the lamination anchor can be moved.
- 2) Adjust the setting aid towards the anterior for maximum plantar flexion or towards the posterior for maximum dorsiflexion (see fig. 2).
- 3) Click the setting aid in to the guiding groove.
- 4) Adjust for the desired adduction or abduction.
- 5) Adjust for the desired rotation.
- 6) Tighten the fine-threaded screw (installation torque: **46 Nm**)

## 5.4 Finishing the Prosthesis

### 5.4.1 Finishing the prosthetic socket

After the trial fitting of the prosthesis and the alignment of the prosthetic foot, the prosthetic socket is finished. Before finishing, the position of the prosthetic foot is marked and the foot is removed.

#### Marking the Foot Position

> **Required tools and materials:**

Pen, scribe, Allen key (2.5 mm and 3 mm), sanding tool, 638K8 Plastaband

- 1) Mark the position of the lamination anchor using a scribe next to the slots in the concave plate.
- 2) Mark the position of the setting aid with a pen.
- 3) Mark the sliding block position on the lamination anchor with a pen.
- 4) Loosen the fine-threaded screw and remove the foot from the socket.
- 5) Trace the markings clearly on the surfaces of the lamination anchor.
- 6) Undo the countersunk head screw on the setting aid using an Allen key and remove the setting aid.
- 7) Sand down the excess resin-talcum mixture on the socket as far as the anchor arms.
- 8) Seal the openings on the sliding block flush with the surface of the lamination anchor using Plastaband.
- 9) Place the lamination dummy on the lamination anchor and screw it on with the screws provided.
- 10) Seal the screw holes on the lamination anchor with Plastaband.

#### Finishing the Prosthetic Socket

> **Required tools and materials:**

Soaked 99B81 PVA Bag, 623T9 Nylglas Stockinette, 616G12 Carbon Fibre Cloth, 99B25 Tube Sock, 616G15 Carbon Fibre Woven Hose, 616G13 Fibreglass Woven Hose, 617H55 C-Orthocryl

- 1) Place two layers of carbon fibre cloth over the anchor arms (see fig. 12).
- 2) Pull over a tube sock.
- 3) Pull a carbon fibre woven hose over the distal **third** and tie it off in the groove of the lamination anchor (see fig. 13).
- 4) Pull the carbon fibre woven hose over the distal **quarter**.
- 5) Cut a Nylglas stockinette to twice the length of the socket.
- 6) Pull the Nylglas stockinette once over the socket and tie it off in the groove of the lamination anchor.
- 7) Place a layer of carbon fibre cloth over the anchor arms.

- 8) If a patient with a mobility grade of 3 weighs over **100 kg (220 lbs)**, place additional carbon fibre cloth over the anchor arms.
- 9) Fold over the Nylglas stockinette.
- 10) Perform the lamination process with C-Orthocryl.
- 11) **Optional:** Leave the casting resin to harden and perform a cast with socket décor fabric.
- 12) **NOTICE! Do not sand the lamination anchor sliding surface or groove.**  
After the lamination resin has cured and cooled, remove the lamination dummy.

#### 5.4.2 Installing the prosthetic foot

##### > Required materials and tools:

Allen key (2.5 mm), 710D4 Torque Wrench

- 1) Install the setting aid on the lamination anchor.
- 2) Place the prosthetic foot on the lamination anchor and bring it to the marked position.
- 3) Screw on and tighten the fine-threaded screw (installation torque: **46 Nm**).
- 4) **Optional:** Cut the footshell in the proximal area and sand it to fit the shape of the socket.
- 5) Put the footshell on the prosthetic foot until it clicks onto the heel spring and slide the foam into the forefoot area.

## 6 Maintenance

- ▶ A visual inspection and functional test of the prosthetic components should be performed after the first 30 days of use.
- ▶ Inspect the entire prosthesis for wear during normal consultations.
- ▶ Conduct annual safety inspections.

## 7 Disposal

In some jurisdictions it is not permissible to dispose of the product with unsorted household waste. Improper disposal can be harmful to health and the environment. Observe the information provided by the responsible authorities in your country regarding return, collection and disposal procedures.

## 8 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

### 8.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

### 8.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

## 9 Technical data

Size [cm]	25	26	27	28
Heel height [mm]	10 ±5			
Structural height without spacer plate [mm]	43			
Structural height with spacer plate [mm]	52			
Product weight with lamination anchor, without foot-shell [g]	465	475	490	495

Size [cm]	25	26	27	28
Maximum body weight [kg]	100	100 (stiffness 2) 125 (stiffness 3)		
Mobility Grade	2 + 3			

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