
50R232 Smartspine LSO high


50R233 Smartspine LSO low

50R235 Smartspine LSO standard

DE Gebrauchsanweisung	4
EN Instructions for use	8
FR Instructions d'utilisation	12
IT Istruzioni per l'uso	16
ES Instrucciones de uso	20
PT Manual de utilização	24
NL Gebruiksaanwijzing	28
SV Bruksanvisning	33
DA Brugsanvisning	37
NO Bruksanvisning	40
FI Käyttöohje	44
IS Notkunarleiðbeiningar	48
PL Instrukcja użytkowania	52
HU Használati utasítás	56
CS Návod k použití	61
RO Instrucțiuni de utilizare	65
HR Upute za uporabu	69
SK Návod na používanie	73
BG Инструкция за употреба	77
TR Kullanma talimatı	81
RU Руководство по применению	85
JA 取扱説明書	90
ZH 使用说明书	94
KO 사용 설명서	97

	Größe / Size	Tailenumfang* / Waist circumference	
		cm	inch
	S	63,5–76,2	25–30
	M	76,2–88,9	30–35
	L	88,9–101,6	35–40
	XL	101,6–114,3	40–45
	XXL	114,3–127	45–50
	XXXL	127–139,7	50–55

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INFORMATION

Date of last update: 2021-12-08

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

These instructions for use provide important information on the fitting and application of the 50R232 Smartspine LSO high, 50R233 Smartspine LSO low and 50R235 Smartspine LSO standard lumbosacral braces.

2 Intended use

2.1 Indications for use

The brace is intended **exclusively** for orthotic fittings of the back and **exclusively** for contact with intact skin.

The brace must be used in accordance with the indications.

2.2 Indications

50R232 Smartspine LSO high

- Stenosis of the lumbar spine
- Degeneration of the intervertebral discs in the lumbar spine
- Spondylolisthesis of the lumbar spine
- Facet syndrome of the lumbar spine
- Spondylolysis of the lumbar spine
- Stable vertebral body fractures of the lumbar spine (without neurological deficits)
- Postoperative stabilization after spinal fusion
- Low back/thoracic spine pain (mechanical)
- Compression fractures of the thoracic and lumbar spine (with no neurological deficits)
- Intervertebral disc prolapse in the thoracic and lumbar spine

50R233 Smartspine LSO low

- Stenosis of the lumbar spine
- Degeneration of the intervertebral discs in the lumbar spine
- Spondylolisthesis of the lumbar spine
- Facet syndrome of the lumbar spine

50R235 Smartspine LSO standard

- Stenosis of the lumbar spine
- Degeneration of the intervertebral discs in the lumbar spine

- Spondylolisthesis of the lumbar spine
- Facet syndrome of the lumbar spine
- Spondylolysis of the lumbar spine
- Stable vertebral body fractures of the lumbar spine (without neurological deficits)

The indication must be determined by the physician.

2.3 Contraindications

2.3.1 Absolute Contraindications

None known.

2.3.2 Relative Contraindications

The following conditions require consultation with a physician: skin diseases/injuries, inflammation, raised scars with swelling, reddening, and warming in the area where the device is fitted, tumours, lymphatic congestion – including unclear soft tissue swelling distal to the device, and dysaesthesia in the trunk.

2.4 Lifetime

The product is designed for a maximum lifetime of **1 year**.

2.5 Mechanism of Action

The brace promotes sensorimotor function and activates muscular stabilization of the lumbar spine. It supports pain relief and relieves the lumbar spine by increasing intra-abdominal pressure and by providing external stabilization.

The applied compression can be individually controlled with the Mechanical Advantage Pulley System.

3 Safety

3.1 Explanation of warning symbols

Warning regarding possible risks of accident or injury.

Warning regarding possible technical damage.

3.2 General safety instructions

Contact with heat, embers or fire

Risk of injury (such as burns) and risk of product damage

- ▶ Keep the product away from open flames, embers and other sources of heat.

Reuse on other persons and improper cleaning

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ The product may be used by one person only.
- ▶ Clean the product regularly.

Improper use and changes

Change in or loss of functionality as well as damage to the product

- ▶ Use this product with care and only for its intended purpose.
- ▶ Do not make any improper changes to the product.

NOTICE

Contact with oils, salves, lotions or other products that contain oils or acids

Insufficient stabilization due to loss of material functionality

- ▶ Do not expose the product to oils, salves, lotions or other products that contain oils or acids.

4 Handling

INFORMATION

- ▶ The daily duration of use and period of application are generally determined by the physician.
- ▶ The initial fitting and application of the product must be carried out by qualified personnel.
- ▶ Instruct the patient in the handling and care of the product.
- ▶ Instruct the patient to see a physician immediately if any exceptional changes are noted (e.g. worsening of the complaint).

4.1 Selecting the Size

- 1) Measure the waist circumference approx. **2 cm** above the iliac crest.
- 2) Determine the brace size (see size chart).

4.2 Adapt to patient

The height of the back plate can be customized to the needs of the patient.

- 1) Loosen the hook-and-loop closure of the back pad.
- 2) Mount the elastic bands of the abdominal strap at the desired position (see fig. 1).
- 3) Close the hook-and-loop closure of the back pad.

4.3 Application

⚠ CAUTION

Direct skin contact with the product

Skin irritation due to chafing or sweat formation

- ▶ Do not wear this product directly on the skin.

⚠ CAUTION

Incorrect or excessively tight application

Risk of local pressure and constriction of blood vessels and nerves due to improper or excessively tight application

- ▶ Ensure that the product is applied properly and fits correctly.

NOTICE

Use of a worn or damaged product

Limited effectiveness

- ▶ Before each use, check the product for functional reliability and for possible wear or damage.
 - ▶ Do not continue using a product that is no longer functional, or that is worn or damaged.
- > The patient is standing.
 - > The hook-and-loop closure of the abdominal strap is open.
 - > The pull cable is completely loosened.
- 1) Center the back plate of the brace (see fig. 2).
 - 2) Guide the brace around the torso and apply equal force to both sides to pull it forward.

- 3) Attach the right side to the left side of the abdominal strap (see fig. 3). Ensure that the hook-and-loop sections are lying flat against each other and are not in contact with clothing or skin.
- 4) Pull the end of the pull cable forward and attach it to the abdominal strap with the hook-and-loop fastener (see fig. 4).
- 5) **Optional – shortening the pull cable:** Pull the winding card out of the pocket (see fig. 5), wrap the pull cable around the winding card (see fig. 6) and put the winding card back into the pocket.
- 6) **Optional – changing the direction of the pull cable from left to right:** The pull cable can be changed from the left to the right side by rotating the brace by 180°. The back pad has to be turned around in the process.

4.4 Optional – only 50R232 Smartspine LSO high: Remove components

The side plates, abdominal plate and back plate can be removed. The removal of the components is directed by the physician.

Side plates

- 1) Open the hook-and-loop closures between the side plates and back plate.
- 2) Remove the side plates. To do so, pull the side plates over the abdominal plate (see fig. 7).

Abdominal plate

- 1) Open the hook-and-loop closures on the abdominal plate (see fig. 8).
- 2) Remove the abdominal plate.

Back plate

- 1) Open the hook-and-loop closure on the back plate.
- 2) Remove the elastic bands from the back plate (see fig. 1).
- 3) Close the hook-and-loop closures on the back plate.
- 4) Remove the back plate. To do so, pull the abdominal strap out of the back plate.

4.5 Cleaning

NOTICE

Use of improper cleaning agents

Damage to the product due to use of improper cleaning agents

- ▶ Only clean the product with the approved cleaning agents.

Clean the brace regularly:

- 1) Remove all boards/plates.
- 2) Fasten all hook-and-loop closures.
- 3) Hand wash the brace in warm water at **30 °C** with standard mild detergent.
- 4) Rinse well.
- 5) Allow to air dry. Do not expose to direct heat sources (e.g. sunlight, stove or radiator).

5 Disposal

Dispose of the product in accordance with national regulations.

6 Legal Information

All legal conditions are subject to the respective national law of the country of use and may vary accordingly.

6.1 Local Legal Information

Legal information that applies **exclusively** to specific countries is written in the official language of the respective country of use after this chapter.

6.2 Liability

The manufacturer shall be liable in the event that the product is used in accordance with the descriptions and instructions in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorized modification of the product.

6.3 CE Conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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