



## 50C91 Smartspine Universal Collar

<b>DE</b>	Gebrauchsanweisung .....	4
<b>EN</b>	Instructions for use .....	9
<b>FR</b>	Instructions d'utilisation .....	13
<b>IT</b>	Istruzioni per l'uso .....	18
<b>ES</b>	Instrucciones de uso .....	23
<b>PT</b>	Manual de utilização .....	28
<b>NL</b>	Gebruiksaanwijzing .....	33
<b>SV</b>	Bruksanvisning .....	38
<b>DA</b>	Brugsanvisning .....	42
<b>NO</b>	Bruksanvisning .....	47
<b>FI</b>	Käyttöohje .....	51
<b>IS</b>	Notkunarleiðbeiningar .....	56
<b>PL</b>	Instrukcja użytkowania .....	60
<b>HU</b>	Használati utasítás .....	65
<b>CS</b>	Návod k použití .....	70
<b>RO</b>	Instrucțiuni de utilizare .....	75
<b>HR</b>	Upute za uporabu .....	80
<b>SK</b>	Návod na používanie .....	84
<b>BG</b>	Инструкция за употреба .....	89
<b>TR</b>	Kullanma talimatı .....	94
<b>JA</b>	取扱説明書 .....	98
<b>ZH</b>	使用说明书 .....	103
<b>KO</b>	사용 설명서 .....	107

		Größe / Size	Abstand* / Distance	
			cm	inch
	50C91	universal	6,4–12,7	2,5–5

\* von Kinnspeitze zu Oberkante Sternum / from the tip of the chin to the top edge of the sternum

<b>Material</b>	50C91: HDPE, PU Schaum/Pu foam, PA, POM
-----------------	---

## INFORMATION

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## INFORMATION

Date of last update: 2021-12-08

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

These instructions for use provide important information on the fitting and application of the 50C91 Smartspine Universal Collar cervical brace.

## 2 Intended use

### 2.1 Indications for use

The brace is intended **exclusively** for orthotic fittings of the cervical spine and **exclusively** for contact with intact skin.

The brace must be used in accordance with the indications.

### 2.2 Indications

#### 50C91 Smartspine Universal Collar

- Acute pain syndromes of the cervical spine
- Stenosis of the cervical spine
- Whiplash injuries of the cervical spine
- Fractures in the cervical spine (simple and stable)
- Dislocation & dislocation fractures of the cervical spine (postoperative)
- Traumatic injury to the cervical spine
- Radiculopathies of the cervical spine

The indication must be determined by the physician.

### 2.3 Contraindications

#### 2.3.1 Absolute Contraindications

None known.

#### 2.3.2 Relative Contraindications

Consult the physician if any of the following symptoms occur: skin diseases/injuries, inflammation, raised scars with swelling, reddening, and

warming in the area where the brace is fitted, tumors, lymphatic congestion and dysaesthesia in the cervical spine.

## 2.4 Lifetime

The product is designed for a maximum lifetime of **1 year**.

## 2.5 Mechanism of Action

The brace immobilizes the cervical spine.

# 3 Safety

## 3.1 Explanation of warning symbols

**CAUTION** Warning regarding possible risks of accident or injury.

**NOTICE** Warning regarding possible technical damage.

## 3.2 General safety instructions

### **CAUTION**

#### **Contact with heat, embers or fire**

Risk of injury (such as burns) and risk of product damage

- ▶ Keep the product away from open flames, embers and other sources of heat.

### **CAUTION**

#### **Reuse on other persons and improper cleaning**

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ The product may be used by one person only.
- ▶ Clean the product regularly.

### **NOTICE**

#### **Improper use and changes**

Change in or loss of functionality as well as damage to the product

- ▶ Use this product with care and only for its intended purpose.
- ▶ Do not make any improper changes to the product.

#### NOTICE

#### Contact with oils, salves, lotions or other products that contain oils or acids

Insufficient stabilization due to loss of material functionality

- ▶ Do not expose the product to oils, salves, lotions or other products that contain oils or acids.

## 4 Handling

#### INFORMATION

- ▶ The daily duration of use and period of application are generally determined by the physician.
- ▶ The initial fitting and application of the product must be carried out by qualified personnel.
- ▶ Instruct the patient in the handling and care of the product.
- ▶ Instruct the patient to see a physician immediately if any exceptional changes are noted (e.g. worsening of the complaint).

### 4.1 Selecting the Size

The brace is available in a universal size. Distances between the tip of the chin and the top of the sternum in the range of 6.4 cm to 12.7 cm can be fitted.

### 4.2 Adaptation

The brace is delivered ready for trial fitting. Adjustment is not required.

- 1) Turn the brace lock on both sides counter-clockwise to the "open" position (arrowheads in vertical position) (see fig. 1).
- 2) Engage the brace on both sides in the corresponding position (1-9) (see fig. 2).

**CAUTION! Ensure that the brace setting is the same on both sides.**

- 3) Turn the brace lock clockwise to the "closed" position (arrowheads in horizontal position) (see fig. 3).
- 4) **Adjusting the inner part of the brace for small neck diameters:** Hold the neckpiece of the inner part of the brace with both hands and press it together. This will allow you to bend the brace round.
- 5) **Optional:** Shorten the hook-and-loop closure to the desired length (see fig. 6).

## 4.3 Application

### CAUTION

#### **Use with heat-sensitive skin**

Skin irritation due to overheating

- ▶ Do not use the product in case of a known heat allergy.
- ▶ Do not continue to use the product if skin irritation occurs.

### CAUTION

#### **Incorrect or excessively tight application**

Risk of local pressure and constriction of blood vessels and nerves due to improper or excessively tight application

- ▶ Ensure that the product is applied properly and fits correctly.

### NOTICE

#### **Use of a worn or damaged product**

Limited effectiveness

- ▶ Before each use, check the product for functional reliability and for possible wear or damage.
- ▶ Do not continue using a product that is no longer functional, or that is worn or damaged.

- > The patient should be standing or sitting down.
- > Open the hook-and-loop closure on one side of the brace.
  - 1) Place the inside of the brace on the neck from the front (see fig. 4).
  - 2) Place the outside of the brace on the neck from behind.
    - The chin rests on the recess.
  - 3) Fasten the hook-and-loop closure again (see fig. 5).

## 5 Cleaning

### NOTICE

#### **Use of improper cleaning agents**

Damage to the product due to use of improper cleaning agents

- ▶ Only clean the product with the approved cleaning agents.

Clean the brace regularly:

- 1) Fasten all hook-and-loop closures.
- 2) Hand wash the brace in warm water at **30 °C** (86 °F) with standard mild detergent.

- 3) Rinse well.
- 4) Allow to air dry. Do not expose to direct heat sources (e.g. sunlight, stove or radiator).

## 6 Disposal

Dispose of the product in accordance with national regulations.

## 7 Legal Information

All legal conditions are subject to the respective national law of the country of use and may vary accordingly.

### 7.1 Local Legal Information

Legal information that applies **exclusively** to specific countries is written in the official language of the respective country of use after this chapter.

### 7.2 Liability

The manufacturer shall be liable in the event that the product is used in accordance with the descriptions and instructions in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorized modification of the product.

### 7.3 CE Conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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