## bebalance+ 1.4 560X12=V1.4

<b>DE</b> Gebrauchsanweisung [kurz] (Fachpersonal) .....	3
<b>EN</b> Instructions for use [short] (qualified personnel) .....	5

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**INFORMATION**

Date of last update: 2022-02-09

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

## Note the instructions for use

This document contains the sections "Intended use" and "Safety", and does not replace the instructions for use. The instructions for use are available in digital form (PDF format) (see page 2 of this document). To view PDF files, you need an Internet browser or other program capable of opening this file format.

You may request a printed copy of the instructions for use under the link on page 2. This will be delivered to the specified address within 7 days after the order is received.

The "560X12=V.14 bebalance+ 1.4" adjustment software is referred to simply as bebalance+, software or adjustment software below.

Basic UDI-DI: 4064411000560X12=V1.4TP

## 2 Intended use

### 2.1 Indications for use

The software allows for patient-specific adaptation of the bebionic hand control device.

### 2.2 Conditions of use

The software is intended exclusively for adjusting the bebionic hand in the course of unilateral or bilateral fittings of the upper limbs.

The B33061 Bluetooth Long Range Dongle is the only permissible method of data transfer between the system components. It can be purchased from Ottobock.

An existing RF dongle can be used only for older bebionic hands of the M and L models.

## 3 Safety

### 3.1 Explanation of warning symbols

**⚠ CAUTION** Warning regarding possible risks of accident or injury.

### 3.2 Structure of the safety instructions

**⚠ CAUTION**

#### **The heading describes the source and/or the type of hazard**

The introduction describes the consequences in case of failure to observe the safety instructions. Consequences are presented as follows if more than one consequence is possible:

- > E.g.: Consequence 1 in case of failure to observe the hazard
- > E.g.: Consequence 2 in case of failure to observe the hazard
- ▶ This symbol identifies activities/actions that must be observed/carried out in order to avert the hazard.

### 3.3 General safety instructions

#### CAUTION

##### **Non-observance of safety instructions**

Personal injury/damage to the product due to using the product in certain situations.

- ▶ Observe the safety instructions and the stated precautions in this accompanying document.

#### CAUTION

##### **Operator errors during the adjustment process**

Injury due to unexpected prosthesis behaviour.

- ▶ Do not charge the prosthesis battery during the adjustment process since the prosthesis is not functional while the battery is being charged.
- ▶ The prosthesis must not remain unattended during the adjustment process while connected to the adjustment software and when being worn by the patient.
- ▶ Observe the maximum range of the Bluetooth connection and note that obstacles may limit this range.
- ▶ The Bluetooth dongle must not be unplugged from the computer during data transfer.
- ▶ Inform the patient immediately if the data connection is accidentally interrupted during the adjustment process.
- ▶ The connection to the prosthesis must always be disconnected after adjustments have been completed.
- ▶ Participation in an Ottobock product training course is recommended.
- ▶ Use the online help which is integrated into the software.

#### CAUTION

##### **Errors while establishing a connection with the prosthesis**

Injury resulting from connection with more than one or the wrong prosthesis.

- ▶ When establishing a connection using the RF dongle, make sure only one RF device is activated in the immediate vicinity.
-

Otto Bock Healthcare Products GmbH  
Brehmstraße 16 · 1110 Wien · Austria  
T +43-1 523 37 86 · F +43-1 523 22 64  
info.austria@ottobock.com · www.ottobock.com