
DMC plus 8E38=6*, 8E39=6*, 8E44=6*

Digital Twin 8E38=7*, 8E39=7*

DE Gebrauchsanweisung (Benutzer)	3
EN Instructions for use (user)	14
FR Instructions d'utilisation (Utilisateur)	26
IT Istruzioni per l'uso (Utilizzatore)	38
ES Instrucciones de uso (Usuario)	50
PT Manual de utilização (Usuário)	62
NL Gebruiksaanwijzing (Gebruiker)	74
SV Bruksanvisning (Användare)	86
DA Brugsanvisning (Bruger)	97
PL Instrukcja użytkownika (Użytkownik)	108
HU Használati utasítás (felhasználó)	120
CS Návod k použití (Uživatel)	132
RO Instrucțiuni de utilizare (Utilizator)	143
TR Kullanma talimatı (Kullanıcı)	155
EL Οδηγίες χρήσης (Χρήστης)	166

INFORMATION

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INFORMATION

Date of last update: 2022-01-19

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Obtain instruction from the qualified personnel in the safe use of the product.
- ▶ Please contact the qualified personnel if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

The 8E38=6*, 8E38=7*, 8E39=6*, 8E39=7* and 8E44=6* System Electric Hands are referred to hereafter as the product/terminal device.

These instructions for use provide you with important information on the use, adaptation and handling of the product.

Only put the product into use in accordance with the information contained in the accompanying documents supplied.

2 Product description

2.1 System Electric Hand DMC plus

The DMC plus control device encompasses the conventional DMC control device and a DMC control device with a higher activation threshold (DMC plus). The desired control variant is selected with the 13E185 function plug. The DMC plus control device is installed in a black housing with a black adjusting ring. The coaxial bushing and hand cable are marked with a brown ring.

Two different operating modes can be selected with the function plug of the product with the DMC plus control device.

2.1.1 Operating mode DMC control device: Function plug connected

The gripping speed and gripping force are determined by the strength of the muscle signal. When the strength of the muscle signal changes, the gripping speed and gripping force immediately adjust proportionally to the changed muscle signal.

The gripping process is significantly more physiological and fragile objects can be gently grasped with ease.

2.1.2 Operating mode DMC plus control device: Function plug not connected

The gripping speed and gripping force are determined by the strength of the muscle signal. When the strength of the muscle signal changes, the gripping speed and gripping force immediately adjust proportionally to the changed muscle signal.

After gripping with maximum force, the activation threshold in the OPEN direction is raised to a higher value. This increase reduces the risk of opening the hand with unintended muscle signals. The reliability of the grip is improved, for example, for eating.

2.2 System Electric Hand Digital Twin

The Digital Twin control device features a digital and a dual-channel control mode. The desired control variant is selected with the 13E185 function plug. The Digital Twin control device is installed in a black housing with two white adjusters. The coaxial bushing and hand cable are marked with a white ring.

Digital mode is controlled with two electrodes and dual-channel mode with just one electrode.

Two different operating modes can be selected with the function plug of the product with the Digital Twin control device.

2.2.1 Operating mode of digital control device: Function plug connected

The hand opens and closes at a constant speed after the respective electrode signal exceeds a threshold.

2.2.2 Operating mode of dual-channel control device: Function plug not connected

The hand is opened or closed with one electrode.

Slow, weak muscle signal: The hand closes.

Fast, strong muscle signal: The hand opens.

3 Intended use

3.1 Indications for use

The product is to be used exclusively for exoprosthetic fittings of the upper limbs.

3.2 Conditions of use

The product is intended **exclusively** for use on **one** user. Use of the product by another person is not approved by the manufacturer.

The product was developed for everyday use and must not be used for unusual activities. These unusual activities include, for example, sports which involve excessive strain and/or shocks to the wrist joint (pushups, downhill racing, mountain biking, ...) or extreme sports (free climbing, paragliding, etc.). Furthermore, the product should not be used to operate motor vehicles, heavy equipment (e.g. construction machines), industrial machines or motor-driven equipment.

Permissible ambient conditions are described in the technical data (see page 24).

3.3 Indications

- Amputation level transcarpal, transradial, transhumeral and shoulder disarticulation
- For unilateral or bilateral amputation
- Dismelia of the forearm or upper arm
- The user must be able to understand usage instructions and safety notices and put them into practice.
- The user must fulfil the physical and mental requirements for perceiving optical/acoustic signals and/or mechanical vibrations

3.4 Contraindications

- All conditions which contradict or go beyond the specifications listed in the section on "Safety" and "Indications for use".

3.5 Qualification

The fitting of a patient with the product may only be carried out by O&P professionals who have been authorised with the corresponding Ottobock training.

4 Safety

4.1 Explanation of warning symbols

 WARNING	Warning regarding possible serious risks of accident or injury.
 CAUTION	Warning regarding possible risks of accident or injury.
 NOTICE	Warning regarding possible technical damage.

4.2 Structure of the safety instructions

 WARNING
The heading describes the source and/or the type of hazard
The introduction describes the consequences in case of failure to observe the safety instructions. Consequences are presented as follows if more than one consequence is possible:
> E.g.: Consequence 1 in the event of failure to observe the hazard
> E.g.: Consequence 2 in the event of failure to observe the hazard
▶ This symbol identifies activities/actions that must be observed/carried out in order to avert the hazard.

4.3 General safety instructions

 WARNING
Non-observance of safety notices
Personal injury/damage to the product due to using the product in certain situations.
▶ Observe the safety notices and the stated precautions in this accompanying document.

⚠ WARNING

Operating a vehicle with the prosthesis

Accidents due to unexpected actions of the prosthesis.

- ▶ The prosthesis should not be used for the operation of motor vehicles and heavy equipment (e.g. construction machines).

⚠ WARNING

Using the prosthesis while operating machines

Injury due to unexpected prosthesis behaviour.

- ▶ The prosthesis should not be used to operate industrial machines or motor-driven equipment.

⚠ WARNING

Operating the prosthesis near active implanted systems

Interference with active implantable systems (e.g., pacemaker, defibrillator, etc.) due to electromagnetic emissions of the prosthesis.

- ▶ When operating the prosthesis in the immediate vicinity of active implantable systems, ensure that the minimum distances stipulated by the manufacturer of the implant are observed.
- ▶ Make sure to observe any operating conditions and safety instructions stipulated by the manufacturer of the implant.

⚠ WARNING

Skin contact with leaking lubricants due to mechanical defects

Injury due to skin irritation.

- ▶ Avoid contact between leaking lubricants and the mouth, nose and eyes.
- ▶ The product must be inspected by an authorised Ottobock Service Center.

⚠ CAUTION

Penetration of product with dirt and humidity

Injury due to unexpected product behaviour or malfunction.

- ▶ Ensure that neither solid particles nor liquids can penetrate into the product.

⚠ CAUTION

Independent manipulation of the product

Injury due to malfunction and resulting unexpected prosthesis actions.

- ▶ Manipulations to the product other than the tasks described in these instructions for use are not permitted.
- ▶ The product and any damaged components may only be opened and repaired by authorised, qualified Ottobock personnel.

4.4 Information on the Power Supply/Battery Charging

⚠ WARNING

Use of damaged power supply unit, adapter plug or battery charger

Risk of electric shock due to contact with exposed, live components.

- ▶ Do not open the power supply unit, adapter plug or battery charger.

- ▶ Do not expose the power supply unit, adapter plug or battery charger to extreme loading conditions.
- ▶ Immediately replace damaged power supply units, adapter plugs or battery chargers.

⚠ WARNING

Charging the prosthesis without taking it off

Risk of electric shock due to defects in the power supply unit or in the battery charger.

- ▶ For safety reasons, remove the prosthesis prior to charging it.

⚠ CAUTION

Prosthesis use when battery charge level is too low

Injury due to unexpected prosthesis behaviour

- ▶ Check the current charge level before use and charge the prosthesis if required.
- ▶ Note that the operating time of the prosthesis may be reduced at low ambient temperatures or due to ageing of the battery.
- ▶ Note that the actions/reactions of the terminal device become slower when the battery voltage gets very low.
- ▶ Note that only a few grips or actions will still be possible with the terminal device when the battery voltage gets very low.
- ▶ A small opening width may indicate a low battery voltage.

4.5 Information on Proximity to Certain Areas

⚠ CAUTION

Insufficient distance to HF communication devices (e.g. mobile phones, Bluetooth devices, WiFi devices)

Injury due to unexpected behaviour of the product caused by interference with internal data communication.

- ▶ Therefore, keeping a minimum distance of 30 cm from HF communication devices is recommended.

⚠ CAUTION

Operating the product in very close proximity to other electronic devices

Injury due to unexpected behaviour of the product caused by interference with internal data communication.

- ▶ Do not operate the product in the immediate vicinity of other electronic devices.
- ▶ Do not stack the product with other electronic devices during operation.
- ▶ If simultaneous operation cannot be avoided, monitor the product and verify proper use in the existing setup.

⚠ CAUTION

Proximity to sources of strong magnetic or electrical interference (e.g. theft prevention systems, metal detectors)

Injury due to unexpected behaviour of the product caused by interference with internal data communication.

- ▶ Avoid remaining in the vicinity of visible or concealed theft prevention systems at the entrance/exit of stores, metal detectors/body scanners for persons (e.g. in airports) or other sources of strong magnetic and electrical interference (e.g. high-voltage lines, transmitters, transformer stations, computer tomographs, magnetic resonance tomographs, etc.).

- ▶ When walking through theft prevention systems, body scanners or metal detectors, watch for unexpected behaviour of the product.

⚠ CAUTION

Wearing in extreme ambient temperatures

Injury due to faulty control or malfunction of the Axon-Bus prosthetic system.

- ▶ Avoid wearing the prosthesis in areas with temperatures outside of the permissible range (see page 24).

4.6 Information on Use

⚠ CAUTION

Mechanical stress on the product

Injury due to faulty control or malfunction of the product.

- ▶ Do not subject the product to mechanical vibrations or impacts.
- ▶ Check the product for visible damage before each use.

⚠ CAUTION

Improper handling

Injury due to faulty operation or malfunction of the product.

- ▶ Seek expert guidance in the proper use of the product.

⚠ CAUTION

Improper product care

- > Injuries due to faulty control/malfunction of the product or damage to the mechanical components
- > Damage or breakage due to brittleness of plastics caused by the use of acetone, petrol or similar solvents.
- ▶ Clean the product only as described in the section "Cleaning and care" (see page 24).
- ▶ Do not clean the product under running water.
- ▶ When using a prosthetic glove, also note the instructions for use of the glove.

⚠ CAUTION

Grasping objects with incorrect gripping forces

Injury due to unexpected product behaviour.

- ▶ Note that the gripping has to be controlled manually depending on the consistency (soft/hard) of the object being grasped.

⚠ CAUTION

Overloading due to unusual activities

Injury due to unexpected product behaviour as a result of a malfunction.

- ▶ The product was developed for everyday use and must not be used for unusual activities. These unusual activities include, for example, sports with excessive strain and/or shocks to the wrist joint (pushups, downhill mountain biking, ...) or extreme sports (free climbing, paragliding, etc.).
- ▶ Careful handling of the product and its components not only increases their service life but, above all, ensures the patient's personal safety!

- ▶ If the product and its components have been subjected to extreme loads (e.g. due to a fall, etc.), then the product must be inspected for damage immediately. If necessary, forward the product to an authorised Ottobock Service Centre.

⚠ CAUTION

Risk of pinching between the fingertips

Injury due to pinching parts of the body.

- ▶ Ensure that no body parts are between the fingertips when using the product.
- ▶ Ensure that no body parts are between the fingertips when closing the hand.
- ▶ Make sure the product is switched off for cleaning.

⚠ CAUTION

Insufficient distance from sources of high heat

Inflammation of the product.

- ▶ Do not expose the product to sources of high heat (fire, stovetop, space heater, radiator etc.).
- ▶ Do not grasp or hold any red-hot objects with the product.

⚠ CAUTION

Unintentional unlocking of the terminal device

Injury due to releasing the terminal device from the forearm (e.g. while carrying objects).

- ▶ When connecting the hand to the prosthetic socket or components, ensure that the connection is carried out correctly.

⚠ CAUTION

Insufficient skin contact of the electrodes

Injury due to unexpected product behaviour as a result of uncontrolled activation of the prosthesis component.

- ▶ The electrodes are to be placed on intact skin only and with as much electrode-skin contact as possible.
- ▶ Make sure that skin contact of the electrodes continues to be given when carrying heavy loads.
- ▶ If the product cannot be controlled correctly using muscle signals, turn off the entire prosthesis and contact your O&P professional.

⚠ CAUTION

Using the prosthesis with pointed or sharp objects (e.g. knives in the kitchen)

Injury due to unintentional movements.

- ▶ Use extra caution when using the prosthesis for handling pointed or sharp objects.

INFORMATION

Avoid exposing the product without a prosthetic glove to direct sunlight or UV light (solarium) for extended periods.

5 Scope of delivery

- 1 pc. 8E38 System Electric Hand
- or
- 1 pc. 8E39 System Electric Hand

or

- 1 pc. 8E44 Transcarpal Hand
- 1 pc. Instructions for use (user)

6 Use

6.1 Terminal device on/off

The ON/OFF button is located directly under the palm of the hand.

System Electric Hand size 7¼, 7¾, 8¼

The ON/OFF button can be operated by pressing on the prosthetic glove.

Area	Function
Back of the hand	ON (see fig. 1)
Thumb	OFF (see fig. 2)

System Electric Hand size 7, Transcarpal Hand

To avoid unintentional activation, the ON/OFF button has to be pressed for about 1 second to activate and deactivate the prosthesis. (see fig. 3)

Press the ON/OFF button only when you want to maintain a grip for a longer period of time, while holding cutlery, writing or during longer passive pauses, for example, while travelling by rail or air. This prevents unintentional opening of the hand caused by inadvertent tensing of the muscles or extreme electrical interference. The operating time of the prosthesis battery is also extended.

The hand switch is easy to operate with your other hand or by pressing against the leg or the edge of a table.

The ON/OFF button can be operated by pressing on the prosthetic glove.

Area	Function
Back of the hand	ON/OFF (see fig. 3)

6.2 Quick-disconnect wrist

The quick-disconnect wrist makes switching between the product and the System Electric Greifer quick and easy.

If the product is equipped with a quick-disconnect wrist, the following options are available:

- Bring the terminal device into the most favourable gripping position by rotating it (passive pronation and supination).
- Easily disconnect the terminal device from the prosthetic socket to replace it with a System Electric Greifer.

To do so, turn the terminal device once around its own axis (360°), clockwise or anticlockwise, until a slight resistance is felt. The terminal device can be removed after overcoming this resistance. (see fig. 4)

Engaging the quick-disconnect wrist

- 1) Insert the quick-disconnect wrist of the prosthesis into the lamination ring and press it in firmly.
- 2) Turn the prosthesis slightly to the left or right.
- 3) Pressure triggers the locking mechanism so the hand is firmly anchored in the lamination ring.

INFORMATION: After it engages, check that the prosthesis or System Electric Greifer is firmly attached.

6.3 Charging the battery

The following points must be observed when charging the battery:

- The capacity of a fully charged battery is sufficient for one full day.
- We recommend charging the product every day when used on a daily basis.
- The battery should be charged for at least 3 hours prior to initial use.

6.3.1 757L20 battery charger and 757B2* EnergyPack

NOTICE

Deep discharge of the battery

Irreparable damage to the product and resulting loss of functionality in case of irregular and incomplete charging of the battery.

- ▶ In case of long-term storage, fully charge the battery every 4 to 6 months.

- 1) Plug the charging station with power supply unit into the wall socket.
→ The green LED on the power supply unit and the green LED in the middle of the LED row light up.
- 2) Insert one or two rechargeable batteries into the charging station.
- 3) The green LED for the respective charging socket lights up ● and the charging process starts.
- 4) When the rechargeable battery is fully charged, the green LED for the respective charging socket flashes :◉:
- 5) Remove the rechargeable battery after the charging process is complete.

6.3.1.1 Display of the current charge level

When inserting/engaging the rechargeable battery in the prosthesis, a charge level indicator is activated on the battery for a few seconds.

LED indicator	Event
●	Charge level above 50% (illuminated in green)
● and ●	Charge level below 50% (flashing alternating green and orange)
●	Charge level below 5% (illuminated in orange)

6.3.1.2 Troubleshooting

Charging socket	LED for the charging socket	Event
Empty	● or ● or ●	The battery charger is defective The battery charger and power supply must be inspected by an authorised Ottobock Service Centre.
Empty or rechargeable battery inserted	☀	
Rechargeable battery inserted	●	The rechargeable battery is charging
	☀	The rechargeable battery is fully charged
	●	The temperature of the rechargeable battery is too high Remove the rechargeable battery and let it cool down
	●	The inserted rechargeable battery is defective The rechargeable battery must be inspected by an authorised Ottobock Service Centre.

6.3.2 757L35 battery charger and 757B35=* rechargeable battery

6.3.2.1 Connect battery charger to the product

- 1) Connect the charging plug to the charging receptacle on the prosthetic socket.
→ A correct connection between the battery charger and the prosthesis is indicated by feedback (Status signals).
- 2) The charging process starts.
→ The prosthesis is switched off automatically.
- 3) Disconnect the prosthesis when the charging process is complete.

6.3.2.2 Safety Shutoff

The purpose of the prosthesis safety shutoff is to protect the battery; it is triggered in case of:

- Excessively high or low temperature
- Overvoltage or undervoltage
- Short circuit

After a short circuit, the charging plug must be plugged into the charging receptacle and then unplugged in order for the electronics to be activated.

6.3.2.3 Display of the current charge level

The battery charge level can be queried at any time.

- 1) With the prosthesis switched on, press the charging receptacle button and hold for less than one second.
- 2) The LED display on the charging receptacle provides information on the current battery charge level.

Charging receptacle	Event
●	Battery fully charged (illuminated in green)
●	Battery 50% charged (illuminated in yellow)
●	Battery drained (illuminated in orange)

6.3.2.4 Beep signals

Beep signal	Additional display	Event
1 x long	-	<ul style="list-style-type: none"> • Turn off prosthesis on charging receptacle • Charging begins (charging plug connected to charging receptacle) • Charging ends (charging plug disconnected from charging receptacle)
2 x short	LED on charging receptacle lights up briefly	Switch on prosthesis
3 x short	-	Battery voltage too low, prosthesis shuts off automatically

7 Storage

If the System Electric Hand will not be used, store the System Electric Hand in the open position to protect the sensors and mechanism.

8 Cleaning and Care

1) Clean the product with a damp cloth and mild soap (e.g. Ottobock 453H10=1-N Derma Clean) when needed.

Ensure that no liquid penetrates into the system component(s).

2) Dry the product with a lint-free cloth and allow it to air dry fully.

Please see the information included with the glove for prosthetic glove care instructions. Further information is available from your O&P professional.

9 Maintenance

Performing regular maintenance (service inspections) every 24 months is recommended to prevent injuries and maintain the quality of the product.

In general, all products are subject to compliance with the maintenance intervals during the warranty period. This is the only way to maintain full warranty cover.

Additional services such as repairs may be provided in the course of maintenance. These additional services may be provided free of charge or can be billable according to an advance cost estimate, depending on the extent and validity of the warranty.

10 Legal information

10.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

10.2 Trademarks

All product names mentioned in this document are subject without restriction to the respective applicable trademark laws and are the property of the respective owners.

All brands, trade names or company names may be registered trademarks and are the property of the respective owners.

Should trademarks used in this document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

10.3 CE conformity

Otto Bock Healthcare Products GmbH hereby declares that the product is in compliance with applicable European requirements for medical devices.

The product meets the requirements of the RoHS Directive 2011/65/EU on the restriction of the use of certain hazardous substances in electrical and electronic devices.

The full text of the regulations and requirements is available at the following Internet address: <http://www.ottobock.com/conformity>

11 Technical data

Ambient conditions	
Storage (with and without packaging)	+5 °C/+41 °F to +40 °C/+104 °F Max. 85% relative humidity, non-condensing
Transport (with and without packaging)	-20 °C/-4 °F to +60 °C/+140 °F Max. 90% relative humidity, non-condensing
Operation	-5 °C/+23 °F to +45 °C/+113 °F Max. 95% relative humidity, non-condensing

Reference number	8E38=6*, 8E39=6*,	8E44=6*	8E38=7*, 8E39=7*
Opening width	79–100 mm	79–100 mm	100 mm
Speed	15–130 mm/s	15–160 mm/s	110 mm/s (constant)
Proportional gripping force	0–90 N		
Lifetime	5 years		

Rechargeable battery reference number	757B20	757B21	757B35 =1	757B35 =3	757B35 =5
Temperature range for charging [°C]	+5 to +40				
Capacity [mAh]	900 ¹ / 950 ¹	680 ¹ / 800 ¹	600	1150	3450
Rechargeable battery lifetime [years]	2				
Product behaviour during the charging process	The product is non-functional				
Operating time of the product with fully charged battery [grip cycles]	Approx. 2500– 3000	Approx. 2000– 2500 (depend- ing on the bat- tery capacity)	Approx. 1000– 2000	Approx. 4000	Approx. 10000
Charging times (rechargeable battery fully dis- charged) [hours]	Approx. 3.5	Approx. 3	Approx. 2.5		Approx. 3
Nominal voltage [V]	Approx. 7.4				
Battery technology	Li-Ion		Li-Po	Li-Ion	
Approved battery chargers	757L20		757L35		

¹ See labeling on battery

12 Appendix

12.1 Symbols Used

Manufacturer

In some jurisdictions it is not permissible to dispose of these products with unsorted household waste. Disposal that is not in accordance with the regulations of your country may have a detrimental impact on health and the environment. Please observe the instructions of your national authority pertaining to return and collection.

Declaration of conformity according to the applicable European directives

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