











50K4=3 Agilium Freestep 3.0




DE	Gebrauchsanweisung	7
EN	Instructions for use	14
FR	Instructions d'utilisation	21
IT	Istruzioni per l'uso	28
ES	Instrucciones de uso	35
PT	Manual de utilização	43
NL	Gebruiksaanwijzing	50
SV	Bruksanvisning	57
DA	Brugsanvisning	63
NO	Bruksanvisning	70
FI	Käyttöohje	76
PL	Instrukcja użytkowania	83
HU	Használati utasítás	90
CS	Návod k použití	97
HR	Upute za uporabu	104
SL	Navodila za uporabo	110
SK	Návod na používanie	117
BG	Инструкция за употреба	123
TR	Kullanma talimatı	131
EL	Οδηγίες χρήσης	137
RU	Руководство по применению	145
JA	取扱説明書	153
ZH	使用说明书	159

1

Größe / Size	Außengelenkhöhe / Lateral Joint Height	
XS	53 – 62 mm	2,1 – 2,4 inch
S	63 – 72 mm	2,5 – 2,8 inch
M	73 – 82 mm	2,9 – 3,2 inch
L	83 – 92 mm	3,3 – 3,6 inch
XL	93 – 103 mm	3,7 – 4,1 inch

Material	Stahl/Steel, Kunststoff/Plastic, TPE, Textilkaschierter Schaum/Fabric-lined foam
-----------------	---

INFORMATION

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INFORMATION

Date of last update: 2022-03-09

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

These instructions for use provide important information on the fitting and application of the 50K4=3 Agilium Freestep 3.0 brace for osteoarthritis of the knee.

2 Product description

2.1 Components (see Fig. 1)

Item:	Designation	Item:	Designation
1	Bars	2	Foot stirrup
3	Sole	4	Brace ankle joint
5	Allen head screw	6	Calf shell with strap

Components not illustrated	
Height adjuster	Allen wrench

2.2 Design

The brace consists of a foot component and a lower leg component joined by a brace ankle joint. The foot component consists of a sole that can be cut to size and a foot stirrup. The lower leg component consists of the brace ankle joint, the bars and the calf shell with the strap. Adjustments to the brace are made by angling the bars.

3 Intended use

3.1 Indications for use

The brace is intended **exclusively** for orthotic fittings of the lower limbs and **exclusively** for contact with intact skin.

The brace must be used in accordance with the indications.

3.2 Indications

- Unicompartmental knee osteoarthritis
- Conditions/injuries requiring unicompartmental relief (e.g. post-operative treatment following meniscus reconstruction or ligament injuries which require unilateral relief)

The indication must be determined by the physician.

3.3 Contraindications

3.3.1 Absolute Contraindications

None known.

3.3.2 Relative Contraindications

The following indications require consultation with a physician: skin diseases/injuries; inflammation; prominent, swollen scars; reddening and hyperthermia of the treated limb; pronounced varicose veins, especially with impaired return flow; lymphatic flow disorders, including unclear soft tissue swelling distal to the body area where the device will be applied; sensory and circulatory disorders in the legs, e. g. associated with diabetic neuropathy.

- Bicompartamental osteoarthritis of the knee

3.4 Lifetime

The product is designed for a maximum lifetime of **2 years**.

3.5 Mechanism of Action

The brace reduces the varus or valgus moment on the knee joint depending on the selected brace.

Two mechanisms influence the extent of the varus moment on the knee joint. The rigid bridging of the ankle joint in the frontal plane initially reduces ever-

sion in the lower ankle joint. Adjustable forces are simultaneously transferred to the lower leg region at the lateral contact point of the brace.

In the treatment of medial osteoarthritis of the knee (brace bar lateral), this effect shifts the vector of the ground reaction force laterally. This immediately results in the reduction of the varus moment in the knee joint and relief for the medial compartment. In the treatment of lateral osteoarthritis of the knee (brace bar medial), the force transmission point is shifted in the medial direction. This increases the varus moment and reduces the valgus moment in the knee joint, relieving the lateral compartment. Strain can be relieved on either the lateral or medial compartment, depending on the selected brace.

4 Safety

4.1 Explanation of warning symbols

CAUTION Warning regarding possible risks of accident or injury.

NOTICE Warning regarding possible technical damage.

4.2 General safety instructions

CAUTION

Reuse on other persons and improper cleaning

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ The product may be used by one person only.
- ▶ Clean the product regularly.

CAUTION

Contact with heat, embers or fire

Risk of injury (such as burns) and risk of product damage

- ▶ Keep the product away from open flames, embers and other sources of heat.

NOTICE

Contact with oils, salves, lotions or other products that contain oils or acids

Insufficient stabilization due to loss of material functionality

- ▶ Do not expose the product to oils, salves, lotions or other products that contain oils or acids.

5 Handling

INFORMATION

- ▶ The daily duration of use and period of application are generally determined by the physician.
- ▶ The initial fitting and application of the product must be carried out by qualified personnel.
- ▶ Instruct the patient in the handling and care of the product.
- ▶ Instruct the patient to see a physician immediately if any exceptional changes are noted (e.g. worsening of the complaint).

5.1 Selecting the Size

- 1) Determine the height of the lateral malleolus, including existing orthopedic insoles. Select the appropriate brace size (see sizing table, see page 5-6).
- 2) Verify that the calf shell is at least **3 cm** below the fibular head. If the calf shell extends closer to the fibular head, select the next smaller brace size.

5.2 Adaptation

⚠ CAUTION

Incorrect alignment, assembly or adjustment

Injury and damage to product components

- ▶ Observe the alignment, assembly and adjustment instructions.
- ▶ Work on the product may only be carried out by qualified personnel.

NOTICE

Improper molding or application

Damage to the product due to overloading of the material and improper fit of the product due to breakage of load-bearing components

- ▶ The product may only be molded by qualified personnel.
- ▶ Do not make any improper changes to the product.
- ▶ Always apply the product according to the information in the instructions.

Aligning and adapting the brace

The brace joint and ankle joint are aligned with each other and the brace sole is cut to size in the course of adapting the brace. The brace sole can be shortened if required. If the sole is shortened, the brace may also be worn in **approved safety shoes**.

- 1) **INFORMATION: If orthopedic insoles are used, position them on the brace sole.** Place the foot on the brace sole or on the sole with the orthopedic insole.
- 2) Shift the foot in the anterior-posterior direction so the axis of the brace joint lines up with the axis of the ankle joint (see fig. 2).
Adjust the medio-lateral distance so the brace joint does not touch the ankle (see fig. 3).
- 3) **Optional: Bend the foot stirrup with bending irons** (see fig. 4).
- 4) Mark the foot contour on the brace sole (see fig. 5).
- 5) Cut the brace sole along the contour (see fig. 6).
- 6) **Optional: Shorten the brace sole on one of the thinly tapered cut lines using scissors** (see fig. 6).
- 7) **Optional: Round the cut edges using a suitable tool.**

Optional: Gluing on the insole

- 1) Roughen the top of the brace sole.
- 2) Apply adhesive to the roughened sole and glue on the insole.

Adjusting the correction strength

- 1) Insert the brace into the shoe (see fig. 7).
- 2) Ottobock recommends using a height adjuster on the contralateral side to compensate for the sole thickness of up to **6 mm**.
- 3) Put the brace on together with the shoe.
- 4) Loosen the Allen head screw (see fig. 8).
- 5) Position the calf shell in front of the lower leg (see fig. 9).
- 6) The lower leg and foot are in the neutral-zero position in the frontal plane. Adjust the calf shell so that it is at approx. 1/3 of the width of the lower leg in the frontal plane (see fig. 9).
- 7) Tighten the Allen head screw (tightening torque **5 Nm**) (see fig. 10).

Optional: Shortening the strap

- 1) Pass the strap around the calf and fasten it on the closure in the strap holder (see fig. 11).
- 2) Determine the required strap length.
- 3) To make adjustments, open the strap at the strap holder and take off the brace.
- 4) Open the strap on the side opposite the closure (see fig. 12).

- 5) Pull off the hook-and-loop patch and place it in the desired position.
- 6) Trim the excess strap (see fig. 13).

Dynamic fitting

- ▶ Conduct the dynamic fitting. Note the position of the ankle joint in doing so. Reduce the correction force if the swing phase or heel strike is not physiological.
- **Lateral foot stirrup:** Excessive correction leads to pronounced supination
- **Medial foot stirrup :** Excessive correction leads to pronounced pronation

5.3 Application

CAUTION

Using the brace without padding

Risk of local pressure points, superficial injuries or skin irritation due to skin contact with rigid or sharp brace components

- ▶ Only use the brace with undamaged padding.

CAUTION

Brace is put on too tightly

Risk of local pressure points and restriction of adjacent blood vessels or nerves

- ▶ Check that the brace is correctly positioned and fitted.

NOTICE

Use of a worn or damaged product

Limited effectiveness

- ▶ Before each use, check the product for functional reliability and for possible wear or damage.
- ▶ Do not continue using a product that is no longer functional, or that is worn or damaged.

- 1) Open the strap closure.
- 2) Insert the brace into the shoe (see fig. 7).
- 3) Put the brace on together with the shoe.
- 4) Tighten the strap and fasten it on the closure on the strap holder (see fig. 14).
- 5) Listen for the "click" when the strap engages in the strap holder (see fig. 15).

5.4 Removal

- 1) Open the strap closure.
- 2) Remove the shoe together with the brace.

6 Cleaning

NOTICE

Use of improper cleaning agents

Damage to the product due to use of improper cleaning agents

- Only clean the product with the approved cleaning agents.

Clean the brace regularly:

Strap and pads

- 1) Remove the strap and pads.
- 2) Hand wash the strap and pads in **30 °C** warm water using a standard mild detergent. Do not use fabric softener. Rinse thoroughly.
- 3) Allow to air dry. Do not expose to direct heat (e.g. from the sun, an oven or radiator).

Brace frame

- 1) Remove soiling on the brace frame with a damp cloth.
- 2) Allow to air dry. Do not expose to direct heat (e.g. from the sun, an oven or radiator).

7 Disposal

Dispose of the product in accordance with national regulations.

8 Legal Information

All legal conditions are subject to the respective national law of the country of use and may vary accordingly.

8.1 Liability

The manufacturer shall be liable in the event that the product is used in accordance with the descriptions and instructions in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorized modification of the product.

8.2 CE Conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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