
Michelangelo Hand 8E500

DE Gebrauchsanweisung [kurz] (Fachpersonal)	3
EN Instructions for use [short] (qualified personnel)	7

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INFORMATION

Date of last update: 2022-03-02

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

Note the instructions for use

This document contains the sections "Intended use" and "Safety", and does not replace the instructions for use. The instructions for use are available in digital form (PDF format) (see page 2 of this document). To view PDF files, you need an Internet browser or other program capable of opening this file format.

You may request a printed copy of the instructions for use under the link on page 2. This will be delivered to the specified address within 7 days after the order is received.

The 8E500=R/L Michelangelo Hand is referred to as the product / Axon-Bus terminal device / Michelangelo Hand below.

Basic UDI-DI: 4064411000000008E500KX

2 Intended use

2.1 Indications for Use

The Michelangelo Hand 8E500 is to be used **exclusively** for exoprosthetic fittings of the upper limbs.

2.2 Conditions of use

The product was developed for everyday use and must not be used for unusual activities. These unusual activities include, for example, sports which involve excessive strain and/or shocks to the wrist joint (pushups, downhill racing, mountain biking, ...) or extreme sports (free climbing, paragliding, etc.). Furthermore, the product should not be used to operate motor vehicles, heavy equipment (e.g. construction machines), industrial machines or motor-driven equipment.

The product is intended **exclusively** for use on **one** patient. Use of the product by another person is not approved by the manufacturer.

Permissible ambient conditions are described in the technical data (Technical data).

2.3 Indications

- Amputation level transradial and transhumeral
- For unilateral or bilateral amputation
- Dysmelia of the forearm or upper arm
- The user must be able to understand usage instructions and safety notices and put them into practice.
- The patient must fulfil the physical and mental requirements for perceiving visual/acoustic signals and/or mechanical vibrations.

2.4 Contraindications

- All conditions which contradict or go beyond the specifications listed in the section on "Safety" and "Indications for use".

2.5 Qualification

The fitting of a patient with the product may only be carried out by O&P professionals who have been authorised with the corresponding Ottobock training.

3 Safety

3.1 Explanation of warning symbols

 WARNING	Warning regarding possible serious risks of accident or injury.
--	---

 CAUTION	Warning regarding possible risks of accident or injury.
--	---

3.2 Structure of the safety instructions

WARNING

The heading describes the source and/or the type of hazard

The introduction describes the consequences in case of failure to observe the safety instructions. Consequences are presented as follows if more than one consequence is possible:

- > E.g.: Consequence 1 in the event of failure to observe the hazard
- > E.g.: Consequence 2 in the event of failure to observe the hazard
- ▶ This symbol identifies activities/actions that must be observed/carried out in order to avert the hazard.

3.3 General safety instructions

WARNING

Non-observance of safety notices

Personal injury/damage to the product due to using the product in certain situations.

- ▶ Observe the safety notices and the stated precautions in this accompanying document.

WARNING

Use of damaged power supply unit, adapter plug or battery charger

Risk of electric shock due to contact with exposed, live components.

- ▶ Do not open the power supply unit, adapter plug or battery charger.
- ▶ Do not expose the power supply unit, adapter plug or battery charger to extreme loading conditions.
- ▶ Immediately replace damaged power supply units, adapter plugs or battery chargers.

WARNING

Operating the prosthetic system near active implanted systems

Interference with active implantable systems (e.g., pacemaker, defibrillator, etc.) due to electromagnetic interference of the prosthetic system.

- ▶ When operating the prosthetic system in the immediate vicinity of active implantable systems, ensure that the minimum distances stipulated by the manufacturer of the implant are observed.
- ▶ Make sure to observe any operating conditions and safety instructions stipulated by the manufacturer of the implant.

CAUTION

Independent manipulation of the product

Injury due to malfunction and resulting unexpected prosthesis actions.

- ▶ Manipulations to the product other than the tasks described in these instructions for use are not permitted.
- ▶ The product and any damaged components may only be opened and repaired by authorised, qualified Ottobock personnel.

CAUTION

Signs of wear on the product

Injury due to faulty control or malfunction of the product.

- ▶ In the interest of the patient's safety and in order to maintain operating reliability, the product should be serviced at regular intervals.

3.4 Information on Proximity to Certain Areas

CAUTION

Proximity to sources of strong magnetic or electrical interference (e.g. theft prevention systems, metal detectors)

Injury due to unexpected behaviour of the prosthesis system caused by interference with internal data communication.

- ▶ Avoid remaining in the vicinity of visible or concealed theft prevention systems at the entrance/exit of stores, metal detectors/body scanners for persons (e.g. in airports) or other sources of strong magnetic and electrical interference (e.g. high-voltage lines, transmitters, transformer stations, computer tomographs, magnetic resonance tomographs, etc.).
- ▶ When walking through theft prevention systems, body scanners or metal detectors, watch for unexpected behaviour of the prosthetic system.

CAUTION

Insufficient distance to HF communication devices (e.g. mobile phones, Bluetooth devices, WiFi devices)

Injury due to unexpected behaviour of the product caused by interference with internal data communication.

- ▶ Therefore, keeping a minimum distance of 30 cm from HF communication devices is recommended.

CAUTION

Storage of the product in the closed state

Injury due to faulty control or malfunction of the product as a result of damage to the sensors and mechanism.

- ▶ Store the product only in the neutral position or in the open state.

3.5 Information on Alignment/Adjustment

CAUTION

Operator errors during the adjustment process with the adjustment software

Injury due to unexpected product behaviour.

- ▶ Participation in an Ottobock product training course is mandatory prior to using the product. During this product training course you will receive a password giving you access to the adjustment software. Additional product training courses may become necessary to qualify for software updates.
- ▶ Do not share your unlock PIN.
- ▶ Use the online help which is integrated into the software.

CAUTION

Incorrect electrode settings

Injury due to unexpected product behaviour.

- ▶ The electrodes are to be placed on intact skin only and with as much electrode-skin contact as possible. In the case of strong interference from electronic devices, the position of the electrodes should be checked and changed if necessary. If the interference cannot be eliminated or if you do not achieve the expected results by adjustment or selection of the appropriate control programme, please contact the Ottobock branch responsible for your country.

- ▶ Set the electrode gain as low as possible in order to reduce interference from powerful electromagnetic radiation (e.g. visible or concealed theft prevention systems at the entrance/exit of stores), metal detectors/body scanners for persons (e.g. in airports) or other sources of strong electromagnetic interference (e.g. high-voltage lines, transmitters, transformer stations, computer tomographs, magnetic resonance tomographs, etc.).

3.6 Information on Use

CAUTION

Improper use

Injury due to malfunction and resulting unexpected prosthesis actions.

- ▶ Instruct the patient in the proper use of the product.

CAUTION

Mechanical stress on the product

Injury due to faulty control or malfunction of the product.

- ▶ Do not subject the product to mechanical vibrations or impacts.
- ▶ Check the product for visible damage before each use.

CAUTION

Operation outside of the permissible temperature range

Injury due to faulty control or malfunction of the product.

- ▶ Avoid operation in areas outside the permissible temperature range (Technical data).

CAUTION

Penetration of dirt and humidity into the prosthetic components

Injury due to faulty control or malfunction of the prosthetic components.

- ▶ Ensure that neither solid particles nor liquids can penetrate into the product or prosthetic components (e.g. the terminal device).

CAUTION

Changing Axon-Bus gripping components when switched on

Injury due to faulty control or malfunction of the Axon-Bus prosthetic system.

- ▶ Power down the Axon-Bus prosthetic system by pressing the button in the charging receptacle before changing any Axon-Bus components (e.g. Axon-Bus gripping component).

CAUTION

Unintentional unlocking of the Axon-Bus gripping component

Injury due to releasing the Axon-Bus gripping component from the forearm (e.g. while carrying objects).

- ▶ The two release buttons should only be activated in order to change the Axon-Bus gripping component, deliberately and with due consideration of the respective situation.

3.7 Information on the Power Supply/Battery Charging

CAUTION

Charging the product with soiled or damaged contacts

Injury due to unexpected behaviour of the product caused by insufficient charging.

- ▶ Make sure that the contacts are always clean and free of grease.

- ▶ Clean the electrical contacts of the charging plug and charging receptacle regularly using cotton swabs and a mild soap solution.
 - ▶ Take care to avoid damaging the contacts with pointed or sharp objects.
-

The 8E500 Michelangelo is covered by the following patents:

Australia:	AU 2006 332 253; AU 2006 332 292; AU 2006 332 315; AU 2006 332 317; AU 2006 332 318; AU 2006 332 316;
Canada:	CA 2 631 970; CA 2 631 982; CA 2 632 241; CA 2 632 551; CA 2 632 240; CA 2 676 196; CA 2 678 987; CA 2 631 966
China:	CN 101 340 864; CN 101 340 865; CN 101 340 866; CN 101 340 867; CN 101 346 107; CN 101 346 109; CN 101 681 704
Japan:	JP 4 843 055; JP 5 242 409; JP 5 123 206; JP 5 074 414; JP 5 155 183; JP 5 389 667; JP 5 074 415
Mexico	MX 285 670; MX 286 932; MX 291 871; MX 291 872; MX 292 784
Russia:	RU 2 387 412; RU 2 414 871; RU 2 416 379; RU 2 423 952; RU 2 427 348; RU 2 429 803; RU 2 469 429
South Korea	KR 101 131 692; KR 101 169 834; KR 101 178 679; KR 101 178 646; KR 101 509 265; KR 101 509 264; KR 101 353 867; KR 101 265 934
Taiwan:	R.O.C. Invention Patent No. I421884
USA:	US 7 867 287; US 8 016 969; US 8 257 446; US 8 188 835; US 8 579 991; US 8 663 339; US 8 690 963
European Patent	EP 1962732 in AT, DE, FR, GB, IT, NL, SE, TR; EP 1962734 in AT, DE, FR, GB, IT, NL, SE, TR; EP 1962737 in AT, DE, FR, GB, IT, NL, SE, TR; EP 1962738 in DE, FR, GB, TR; EP 1971297 in AT, DE, FR, GB, IT, NL, SE, TR; EP 2129340 in DE, FR, GB, IS, IT, NL, SE, TR; EP 2115752 in AT, CH, DE, ES, FR, GB, IT, NL, SE, TR; EP 2528549 in DE, FR, GB, IS, PL, TR

Patents pending in Brazil, EPA, Germany, India, Mexico and USA.

The 8E500 Michelangelo or parts thereof are covered by the following registered designs and design patents:

Australia:	317212; 317213; 317214; 317789; 317790; 317791; 317792; 317866; 317867; 317868; 317869; 317870; 317871
Canada:	©Ottobock No. 122162; 122163
China:	ZL 200730154423.X; ZL 200730154429.7; ZL 201130050582.1; ZL 201130050654.2; ZL 200730154424.4; ZL 200730154425.9
European Design:	No.000786421; No.000786694; No.001824004
Germany:	40701345.8; 40701357.1
India:	212295; 212296; 212297; 212298; 212299; 212300; 212301; 212302; 212303
Japan:	Registered Design No. 1 365 277; 1 365 278
Russia:	69 071; 70 542
Taiwan:	R.O.C. Design Patent D 128 170; D 128 171
USA:	Patent US D 595,854; D 597,672; D 694,189; D697 030

Michelangelo is a tradename of Otto Bock Healthcare GmbH.

Michelangelo is a registered trademark in many countries of the world, beside others registered at the US Patent and Trade-mark Office, Reg.–No. 4008171.

Otto Bock Healthcare Products GmbH
Brehmstraße 16 · 1110 Wien · Austria
T +43-1 523 37 86 · F +43-1 523 22 64
info.austria@ottobock.com · www.ottobock.com