3S80, 3S80=1 Sport





DE Gebrauchsanweisung	7
EN Instructions for use	20
FR Instructions d'utilisation	34
IT Istruzioni per l'uso	49
ES Instrucciones de uso	63
PT Manual de utilização	78
NL Gebruiksaanwijzing	92
SV Bruksanvisning	106
DA Brugsanvisning	120
NO Bruksanvisning	133
FI Käyttöohje	146
HU Használati utasítás	159
HR Upute za uporabu	173
JA 取扱説明書	186

Grundaufbau Bench Alignment 1E90 Sprinter

**Anfänger
Beginners**

**Fortgeschrittene
Advanced users**



5


Grundaufbau Bench Alignment

1E91 Runner

1E93 Runner junior

743A200
743A220

743Y580=

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INFORMATION

Date of last update: 2022-03-29

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.

- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

INFORMATION

These instructions for use are only valid starting with the following serial number: 201429001

1.1 Explanation of warning symbols

- | | |
|---|---|
|  WARNING | Warning regarding possible serious risks of accident or injury. |
|  CAUTION | Warning regarding possible risks of accident or injury. |
|  NOTICE | Warning regarding possible technical damage. |

1.2 Meanings of pictograms in the illustration

	Numbering for the illustrations		Reference to the number of the corresponding section
	Numbering for a defined sequence		Numbering for the parts of an illustration

2 Product description

2.1 Construction and Function

The 3S80* Sport is a sport knee joint designed for sports with no double support phase while running. These include in particular jogging, sprinting and related running disciplines in athletics that encompass corresponding movement sequences. For this purpose, the prosthetic knee joint has rotation hydraulics to control the swing phase.

The prosthetic knee joint is not suitable for an everyday prosthesis since it has no stance phase control and is designed for the considerably greater flexion angles that occur during sports.

It has adjustable extension and flexion damping for adaptation to the chosen sport. The prosthetic knee joint has an integrated lock.

The hydraulic oil used in the 3S80=1 has a considerably lower viscosity than the oil used in the 3S80. The less viscous oil facilitates smooth-running

movements, what can be advantageous e. g. in sport prostheses for children and especially in sprint prostheses.

3 Intended use

3.1 Indications for use

The product is intended exclusively for lower limb exoprosthetic fittings.

3.2 Area of application

CAUTION

Excessive strain on the product

Fall due to breakage of load-bearing components

- ▶ Only use the product according to its allowable field of application.

Approved for a body weight up to **max. 100 kg**.

Permitted sports

For use in sport prostheses (recreational sports, competitive sports) for jogging, sprinting and related running disciplines in athletics

Support for other sports

- ▶ Contact Ottobock if a prosthesis is needed for sports other than the permitted sports. Ottobock is available for enquiries regarding application in other sports.

3.3 Combination possibilities

CAUTION

Improper combination of prosthetic components

Injuries, malfunctions or product damage due to unallowable combination of prosthesis components

- ▶ Based on the instructions for use of all prosthetic components used, verify that they may be combined with each other and are approved for the patient's area of application.

INFORMATION

In a prosthesis, all prosthetic components have to meet the patient's requirements regarding the amputation level, body weight, activity level, environmental conditions and field of application.

INFORMATION

In the sport prosthesis, only use prosthetic components approved for the respective sport or a body weight of 150 kg.

This prosthetic component is compatible with Ottobock's system of modular connectors. Functionality with components of other manufacturers that have compatible modular connectors has not been tested.

Prosthetic components that are especially well suited for combination with the product are listed below.

Designation	Reference number
Lamination anchor	4R111, 4R116
Socket adapter	4R54, 4R77, 4R55, 4R51
Double adapter	4R72*
Foot adapter	1E90 Sprinter: 4R204, 4R206 1E91 Runner: 4R218, 4R216 1E93 Runner junior: 4R224
Prosthetic feet	1E90 Sprinter, 1E91 Runner, 1E93 Runner junior
Foot soles	1E90 Sprinter: 2Z500, 2Z501 1E91 Runner: 2Z540, 2Z541 1E93 Runner junior: 2Z543

3.4 Environmental conditions

CAUTION

Use under unallowable environmental conditions

Fall due to damaged product

- ▶ Do not expose the product to unallowable environmental conditions (see the table "Unallowable environmental conditions" in this section).
- ▶ If the product was exposed to unallowable environmental conditions, take suitable steps (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop etc.).

Storage and transport

Temperature range -20 °C to +60 °C (-4 °F to +140 °F), relative humidity 20 % to 90 %, no mechanical vibrations or impacts

Allowable environmental conditions

Temperature range: -10 °C to +45 °C (14 °F to 113 °F)

Chemicals/liquids: fresh water, soapsuds, chlorine water

Allowable environmental conditions

Moisture: Submersion: max. 1 h in 2 m depth, relative humidity: no restrictions

Solids: Dust, occasional contact with sand

Clean the product after contact with humidity/chemicals/solids, in order to avoid increased wear and damage (see page 31).

Prohibited environmental conditions

Solids: highly hygroscopic particles (e.g. talcum), dust in high concentrations (e.g. construction site), intensive contact with sand

Chemicals/liquids: salt water, perspiration, urine, acids, continuous use in liquid media

3.5 Reuse and lifetime

CAUTION

Reuse on other patients and improper maintenance

Fall due to loss of functionality or damage to the product

- ▶ Only use the product for a single patient.
- ▶ Maintain the product regularly to ensure a long service life (see page 32).

CAUTION

Exceeding the lifetime

Fall due to change in or loss of functionality and damage to the product

- ▶ Ensure that the maximum lifetime defined in this section is not exceeded.

This product was tested by the manufacturer with 3 million load cycles.

Recreational sports: The maximum lifetime is 2 years.

Intensive use and competitive sports: The lifetime is 1 year.

4 General safety instructions

WARNING

Using the product

Severe falls due to improper use of the product

- ▶ Only use the product after receiving full instruction by qualified personnel.
- ▶ Do not use the product as an everyday prosthesis.

⚠ CAUTION**Reaching into the area of the joint mechanism**

Pinching of limbs (e.g. fingers) and the skin due to uncontrolled joint movement

- ▶ Do not reach into the joint mechanism during daily use.
- ▶ Close attention is required during assembly and adjustment tasks.

⚠ CAUTION**Mechanical damage to the product**

Risk of injury due to change in or loss of functionality

- ▶ Use caution when working with the product.
- ▶ If the product is damaged, check it for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, do not continue using the product (see "Signs of changes in or loss of functionality during use" in this section).
- ▶ Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

Signs of changes in or loss of functionality during use

Changes in functionality can manifest themselves, for example, through a changed gait (swing phase, stance phase), incomplete extension, stiffness and the development of noise.

5 Scope of delivery

The scope of delivery for the product is illustrated on page 1 (title image).

The following single components and accessories are included in the scope of delivery in the specified quantities and are available for reordering as single components (■), single components with minimum order quantities (▲) or single-component packs (●):

Fig.		Item no.	Quantity	Designation	Reference number
–	■	–	1	Instructions for use	647G813=all_I-NT
–	■	–	1	Quick reference guide	646D1563=all_INT
–	–	–	1	Sports knee joint	–
–	■	–	1	Adjustment wrench	710H10=2x3

Fig.		Item no.	Quantity	Designation	Reference number
-	■	-	2	Pyramid Dome Protector Kit	2Z11=KIT

6 Preparing the product for use

⚠ CAUTION

Incorrect alignment, assembly or adjustment

Injury due to incorrectly installed or adjusted as well as damaged prosthetic components

- ▶ Observe the alignment, assembly and adjustment instructions.

⚠ CAUTION

Use of talcum

Risk of injury, damage to the product due to lack of lubrication

- ▶ Do not use talcum on the product or other prosthesis components.

The three-dimensional alignment of the prosthetic socket and the modular components affects the function of the prosthesis. The patient will only be able to benefit from the special characteristics of the prosthetic knee joint if the alignment is correct.

The position of the residual limb, e.g. abduction, adduction and flexion, must be taken into account for positioning the socket connector. Proceeding according to the manufacturer's alignment guidelines is recommended.

The alignment is performed in two steps:

- Bench alignment
- Static alignment

The individual adjustment of the prosthesis and the prosthetic knee joint is carried out during the dynamic fitting.

Use 2Z11=KIT

INFORMATION: The connection area of the prosthetic knee joint can be protected against scratches with the protective film of the 2Z11=KIT during alignment in the workshop and testing in the trial fitting area.

- ▶ Use the protective film as illustrated in the accompanying document for the 2Z11=KIT.
- ▶ Remove the protective film before the patient leaves the fitting area.

6.1 Bench Alignment

Bench alignment of the prosthesis: 1E90 Sprinter

Optimising the prosthetic foot height

In the initial bench alignment, the sports prosthesis is fabricated as a trial prosthesis with the test foot adapter. This is used to check whether the correct stiffness of the prosthetic foot was chosen. Then the height of the 1E90 prosthetic foot is optimally adapted to the patient's requirements. The movement of the hips and, in particular, the shoulders provides important information on the optimal prosthesis length. When the prosthesis is the optimal length, dropping the shoulder on the prosthesis side while running is prevented. See the 647H537 quick reference guide for further information and exercises.

The bench alignment is performed with the installed 2Z500 or 2Z501 outer sole (spikes not installed!).

Bench alignment process	
+ = anterior location / - = posterior location (in relation to the alignment reference line)	
Item	Neue Callout-Grafik
	Required tools and materials: 743A80 50:50 gauge, alignment apparatus (e.g. 743L200 L.A.S.A.R. Assembly or 743A200 PROS.A. Assembly) Prosthetic knee joint, locked (note the section "Use" – see page 30)
①	a-p position of the prosthetic foot's alignment reference point to the alignment reference line: +40 mm (beginner) to 0 mm (advanced)
②	External foot rotation: Approx. 5°
③	Build height of the prosthetic knee joint at the alignment reference point (rotation axis of the prosthetic knee joint): MTP to floor measurement + 20 mm
	a-p position of the alignment reference point to the alignment reference line: 0 mm
	External rotation of the prosthetic knee joint: Approx. 5°
④	Join the prosthetic foot with the prosthetic knee joint using the foot adapter and the socket adapter.

Bench alignment process	
+ = anterior location / - = posterior location (in relation to the alignment reference line)	
Item	Neue Callout-Grafik
	Be sure to follow the adapter instructions for use during adaptation and assembly.
5	Mark two centre points laterally on the prosthetic socket: 1. At the level of the greater trochanter 2. In the distal region
6	Join the two points with a line.
7	Join the prosthetic socket and the prosthetic knee joint using the lamination anchor. Be sure to follow the lamination anchor instructions for use during adjustment and assembly. Positioning of the prosthetic socket to the alignment reference line: <ul style="list-style-type: none"> • a-p position: 0 mm (proximal point of the prosthetic socket) • Socket flexion: Hip flexion contracture (determined with Thomas test) + at least 4° Take note of the patient's individual situation.

Bench alignment of the prosthesis: 1E91 Runner, 1E93 Runner junior

- see fig. 5

► Please read this document carefully: 647G1145

6.2 Static Alignment

Static alignment process	
+ = anterior location / - = posterior location (to the load line)	
Item	1E90 Sprinter: see fig. 6; 1E91 Runner/1E93 Runner junior: see fig. 7
	Required materials and tools: 743L100 L.A.S.A.R. Posture Prosthetic knee joint, locked (note the section "Use" – Use)
1	To determine the load line, position the patient on the L.A.S.A.R. Posture as follows: <ul style="list-style-type: none"> • Prosthetic foot on force plate (sufficient load: > 35% body weight) • Other foot (with shoe) on height compensation plate • The tip of the shoe and prosthetic foot are aligned with each other

Static alignment process	
+ = anterior location / - = posterior location (to the load line)	
Item	1E90 Sprinter: see fig. 6; 1E91 Runner/1E93 Runner junior: see fig. 7
2	<p>Optimise the static alignment solely by changing the plantar flexion on the socket adapter of the foot adapter.</p> <p>a-p position of the alignment reference point (rotation axis of the prosthetic knee joint) to the load line:</p> <ul style="list-style-type: none"> • 1E90 Sprinter: Beginner: > -80 mm Advanced: -80 mm • 1E91 Runner/1E93 Runner junior: ~ 40 mm

6.3 Dynamic Trial Fitting

During the dynamic fitting, the alignment of the prosthesis is checked, plantar flexion is adapted and flexion damping and extension damping are adjusted. The correct use of the sports prosthesis has to be practised to ensure safe running and walking. Basic information regarding the exercises is found in this section.

Walking exercises

Use parallel bars for safety while trying out the prosthesis. After each change to the prosthesis, begin the exercises between parallel bars. Always start with slow, short steps.

Perform the walking exercises with the prosthetic knee joint locked until the patient gets a feel for the prosthetic foot.

Then unlock the prosthetic joint to develop a feel for the functionality of the prosthetic joint.

Increase the stride length and walking speed as the patient gains experience.

Perform the walking exercises until the patient feels safe while walking.

Running exercises

Use suitable protective equipment (e.g. for inline skating) during the initial running exercises. Choose a surface that is not too hard (e.g. all-weather track, indoors or short, firm lawn), because a hard surface (e.g. asphalt or concrete) puts more strain on the patient and can cause pain.

Adapt the prosthesis to the sport and the requirements of the patient by making corresponding adjustments and changing settings.

Completing the prosthesis

After successful optimisation, the trial prosthesis is converted to the definitive prosthesis by replacing the test foot adapter with the sports foot adapter. Further information can be found in the instructions for use of the prosthetic foot and the foot adapters.

6.3.1 Adjusting Flexion and Extension Damping

Valve	Symbol	Factory setting	
Flexion valve	(F)	Stop +	High
Extension valve	(E)	Stop –	Low

> Required tools and materials:

710H10=2x3 Adjustment Wrench

- 1) Use the adjustment wrench to adapt the flexion resistance to the needs of the patient.
- 2) Use the adjustment wrench to adapt the extension resistance to the needs of the patient.
- 3) Check the adjustments in dynamic trial fitting and correct as required according to the table below.

Situation	Measure	Adjustment
Lower leg swings through too far	Increase damping	Turn flexion valve to the right (+)
Lower leg does not swing through sufficiently	Reduce damping	Turn flexion valve to the left (-)
Lower leg impact too hard on the extension stop	Increase damping	Turn extension valve to the right (+)
Lower leg does not reach full extension before heel strike	Reduce damping	Turn extension valve to the left (-)

7 Use

7.1 Information for use

CAUTION

Overheating of the hydraulics due to excessive strain

Burns, injuries caused by falling due to changes in functionality and damage to prosthetic components

- ▶ Do not touch overheated prosthetic components.
- ▶ Reduce all activities if there are changes in functionality, in order to allow the overheated prosthetic components to cool down.

- ▶ If overheating or changes in functionality occur, have the prosthesis checked by authorised qualified personnel.

NOTICE

Mechanical overload

Impaired functionality due to mechanical damage

- ▶ Check the product for damage prior to each use.
- ▶ Do not use the product if its functionality has been impaired.
- ▶ Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

INFORMATION

Leak in the product's hydraulic system

Environmental damage or skin irritation due to leaking hydraulic oil

- ▶ Before each use, check the product for leaking hydraulic oil.
- ▶ In case of leaks, do not continue using the product and have it repaired promptly.

INFORMATION

Safe standing

The prosthetic joint has a lever to lock it for standing (see fig. 3).

- ▶ While standing under load (full extension), push the lever down to activate the standing lock.
 - ▶ While standing under load (full extension), pull the lever up to release the standing lock.
- ▶ Inspect the prosthetic knee joint for wear and tear and to verify proper functionality after each use.

7.2 Cleaning

⚠ CAUTION

Use of unsuitable cleaning agents or disinfectants

Impairment of functionality and damage due to incorrect cleaning agents or disinfectants

- ▶ Only clean the product with the approved cleaning agents.
- ▶ Only disinfect the product with the approved disinfectants.
- ▶ Observe the instructions for cleaning and care.

Permissible cleaning agents and disinfectants

Requirements:

Free of solvents, chlorine and phosphates

INFORMATION

- ▶ Clean the product if it gets dirty.
- ▶ Observe the instructions for cleaning and drying in the section "Environmental conditions" (see page 23).
- ▶ Only use cleaning agents and disinfectants that do not affect the materials of the product. Test the chosen agent for material compatibility in an inconspicuous location.
- ▶ Do not spray the product directly with a pump spray or propellant spray cleaner and disinfectant.
- ▶ Only apply cleaners and disinfectants using a soft, lint-free cloth.
- ▶ Dry with a soft, lint-free cloth. Allow to air dry in order to remove residual moisture.
- ▶ Observe the cleaning instructions for all prosthetic components.

8 Maintenance

CAUTION

Failure to follow the maintenance instructions

Risk of injuries due to changes in or loss of functionality and damage to the product

- ▶ Observe the following maintenance instructions.
- ▶ **NOTICE! Repair work must be performed exclusively by manufacturer service.**
- ▶ The prosthetic components should be inspected after the first 30 days of use.
- ▶ Following an individual period for the patient to get accustomed to the prosthesis, check the settings of the prosthetic joint and adapt them to the patient's requirements again as needed.
- ▶ Arrange regular maintenance intervals with the patient depending on the level of use.
- ▶ Conduct annual safety inspections.
- ▶ Inspect the entire prosthesis for wear during normal consultations.

- ▶ As part of the safety inspections, inspect the prosthetic joint for wear and proper functionality. Special attention should be paid to movement resistance, bearings and abnormal noises. Full flexion and extension must be ensured at all times. Readjust as required.
- ▶ Direct the patient to observe the following maintenance instructions:
Check the entire prosthesis for damage after each use.
In case of changes in or loss of functionality and damage, do not continue using the prosthesis and have it checked by qualified personnel.
Comply with the maintenance intervals.

9 Disposal

In some jurisdictions it is not permissible to dispose of the product with unsorted household waste. Improper disposal can be harmful to health and the environment. Observe the information provided by the responsible authorities in your country regarding return, collection and disposal procedures.

10 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

10.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

10.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

11 Technical data

Reference number	3S80, 3S80=1 Sport
Weight (g)	682
System height (mm)	48
Proximal system height up to the alignment reference point (mm)	28

Reference number	3S80, 3S80=1 Sport
Distal system height up to the alignment reference point (mm)	20
Proximal connection	Pyramid
Distal connection	Pyramid
Max. flexion angle (°)	135
Max. body weight (kg)	100

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