



## 17AF10=\*

<b>DE</b> Gebrauchsanweisung .....	3	<b>CS</b> Návod k použití.....	28
<b>EN</b> Instructions for use .....	4	<b>RO</b> Instrucțiuni de utilizare.....	30
<b>FR</b> Instructions d'utilisation.....	6	<b>HR</b> Upute za uporabu .....	32
<b>IT</b> Istruzioni per l'uso .....	8	<b>SL</b> Navodila za uporabo .....	34
<b>ES</b> Instrucciones de uso .....	10	<b>SK</b> Návod na používanie .....	36
<b>PT</b> Manual de utilização.....	12	<b>BG</b> Инструкция за употреба .....	38
<b>NL</b> Gebruiksaanwijzing.....	14	<b>TR</b> Kullanma talimatı .....	40
<b>SV</b> Bruksanvisning .....	16	<b>EL</b> Οδηγίες χρήσης .....	41
<b>DA</b> Brugsanvisning .....	18	<b>RU</b> Руководство по применению .....	44
<b>NO</b> Bruksanvisning .....	20	<b>JA</b> 取扱説明書.....	46
<b>FI</b> Käyttöohje.....	22	<b>ZH</b> 使用说明书.....	48
<b>PL</b> Instrukcja użytkowania.....	24	<b>KO</b> 사용 설명서.....	49
<b>HU</b> Használati utasítás .....	26		

# 1 Vorwort

Deutsch

## INFORMATION

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### INFORMATION

Date of last update: 2022-03-28

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- ▶ Please keep this document for your records.

These instructions for use provide you with important information on the processing of the 17AF10=\* ankle joints.

## 2 Intended use

### 2.1 Indications for use

The product is intended **exclusively** for use in pairs, in orthotic fittings of the lower limbs and **exclusively** for contact with intact skin.

The orthosis must be used in accordance with the indications.

### 2.2 Indications

- Dorsiflexor weakness

Indications must be determined by the physician.

### 2.3 Lifetime

The product is designed for the following lifetime when used as intended and assembled properly:

Joint size S: **1 year**

Joint size L, M: **2 years**

### 2.4 Qualification

Patients may be fitted with the product only by trained qualified personnel. The qualified personnel must be familiar with the handling of the various techniques, materials, machines and tools.

## 3 Safety

### 3.1 Explanation of warning symbols

Warning regarding possible risks of accident or injury.

Warning regarding possible technical damage.

### 3.2 General safety instructions

#### Incorrect alignment, assembly or adjustment

Injury and damage to product components

- ▶ Observe the alignment, assembly and adjustment instructions.
- ▶ Work on the product may only be carried out by qualified personnel.

#### Use of the product in case of known allergies to chrome or nickel

Allergic reaction

- ▶ Avoid skin contact with the product in case of known allergies to chrome or nickel.

#### Mechanical damage to the product

Injuries due to changes in or loss of functionality

- ▶ Use caution when working with the product.
- ▶ Check the product for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, discontinue use of the product and have it checked by authorised, qualified personnel.

## 4 Preparing the product for use

### 4.1 Environmental conditions

**Permitted:** Temperature range for use  $-10^{\circ}\text{C}$  to  $+60^{\circ}\text{C}$ ; relative humidity 0 % to 90 %, non-condensing.

**Not permitted:** Mechanical vibrations, impacts, perspiration, urine, fresh water, salt water, acids, dust, sand, highly hygroscopic particles (such as talcum).

### 4.2 Processing

For further processing instructions, see the 646T523=\* technical information.

- 1) Place the dummy on the joint pivot point of the plaster model and secure it with the nails (see fig. 1).
- 2) Vacuum-form or laminate the orthosis.
- 3) Remove the orthosis and dummies from the plaster model after it has hardened/cooled.
- 4) Separate the foot and calf component with a saw and deburr the edges (see fig. 2).
- 5) For pre-flexed joints make a U-shaped cutout, for straight joints make a V-shaped cutout at the front sides of the joint adapters (see fig. 3).
- 6) Drill the mounting holes for the joints (diameter 4.5 mm).
- 7) Insert the joints and attach them to the orthosis with the screws (see fig. 4) (torque value: 0.5 Nm).
- 8) Cover both joints with the padding.

## 5 Cleaning

Promptly clean the product after contact with water containing salt, chlorine or soap, or if it gets dirty.

- 1) Rinse the product with clear fresh water.
- 2) Dry the product with a cloth or allow it to air dry. Avoid exposure to direct heat (e. g. from an oven or radiator).

## 6 Maintenance

### INFORMATION

**The product may be exposed to increased loads by the patient.**

- ▶ Shorten the maintenance intervals according to the expected loads.

## 7 Disposal

Dispose of the product in accordance with national regulations.

## 8 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

### 8.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

### 8.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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