

50K324 Agilium Reactive

DE	Gebrauchsanweisung	5
EN	Instructions for use	11
FR	Instructions d'utilisation	17
IT	Istruzioni per l'uso	24
ES	Instrucciones de uso	30
PT	Manual de utilização	36
NL	Gebruiksaanwijzing	43
SV	Bruksanvisning	49
DA	Brugsanvisning	55
NO	Bruksanvisning	61
FI	Käyttöohje	67
PL	Instrukcja użytkowania	73
HU	Használati utasítás	79
CS	Návod k použití	85
SK	Návod na používanie	91
TR	Kullanma talimatı	97
RU	Руководство по применению	103
JA	取扱説明書	110
ZH	使用说明书	116

Material

Al, PP, PA, PU, PES, EL, PE Schaum/PE foam, Stahl/steel

INFORMATION

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INFORMATION

Date of last update: 2022-03-09

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

These instructions for use provide important information on the fitting and application of the 50K324 Agilium Reactive knee joint brace.

2 Product description

Fig. 1/ Item:	Designation	Fig. 1/ Item:	Designation
1	Brace	7	Joint stop for flexion limitation
2	Strap pads	8	Tibia pad
3	Condyle pad	9	Strap guide loop
4	Allen wrench	10	Top lower leg strap closure clip
5	Screwdriver	11	BOA dial

Fig. 1/ Item:	Designation	Fig. 1/ Item:	Designation
6	Joint stop for extension limitation		

3 Intended use

3.1 Indications for use

The brace is intended **exclusively** for orthotic fittings of the lower limbs and **exclusively** for contact with intact skin.

The brace must be used in accordance with the indications.

3.2 Indications

- Unicompartmental knee osteoarthritis
- Post-operative fitting after meniscus reconstruction
- Ligament injuries requiring unilateral relief

The indication must be determined by the physician.

3.3 Contraindications

3.3.1 Absolute Contraindications

None known.

3.3.2 Relative Contraindications

The following indications require consultation with a physician: skin diseases/injuries; inflammation; prominent, swollen scars; reddening and hyperthermia of the treated limb; pronounced varicose veins, especially with impaired return flow; lymphatic flow disorders, including unclear soft tissue swelling distal to the body area where the device will be applied; sensory and circulatory disorders in the legs, e. g. associated with diabetic neuropathy.

- Bicompartamental osteoarthritis of the knee

3.4 Lifetime

The product is designed for a maximum lifetime of **2 years**.

3.5 Mechanism of Action

The brace relieves the medial or lateral compartment of the knee joint using the three-point principle. In the process, corrective force is applied to the contralateral condyle.

4 Safety

4.1 Explanation of warning symbols

CAUTION Warning regarding possible risks of accident or injury.

NOTICE Warning regarding possible technical damage.

4.2 General safety instructions

CAUTION

Reuse on other persons and improper cleaning

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ The product may be used by one person only.
- ▶ Clean the product regularly.

CAUTION

Breakage of the brace frame or joint due to improper shaping

Injuries caused by sharp edges, loss of function

- ▶ Avoid adapting the brace repeatedly.
- ▶ Do not make any improper changes to the brace.
- ▶ **Inform your patient.**

CAUTION

Incorrect or excessively tight application

Risk of local pressure and constriction of blood vessels and nerves due to improper or excessively tight application

- ▶ Ensure that the product is applied properly and fits correctly.

CAUTION

Contact with heat, embers or fire

Risk of injury (such as burns) and risk of product damage

- ▶ Keep the product away from open flames, embers and other sources of heat.

NOTICE

Contact with oils, salves, lotions or other products that contain oils or acids

Insufficient stabilization due to loss of material functionality

- ▶ Do not expose the product to oils, salves, lotions or other products that contain oils or acids.

NOTICE

Improper use and changes

Change in or loss of functionality as well as damage to the product

- ▶ Use this product with care and only for its intended purpose.
- ▶ Do not make any improper changes to the product.

NOTICE

Use of a worn or damaged product

Limited effectiveness

- ▶ Before each use, check the product for functional reliability and for possible wear or damage.
- ▶ Do not continue using a product that is no longer functional, or that is worn or damaged.

5 Handling

INFORMATION

- ▶ The daily duration of use and period of application are generally determined by the physician.
- ▶ The initial fitting and application of the product must be carried out by qualified personnel according to the instructions of the treating physician.
- ▶ Consult a physician if any exceptional changes are noted (such as worsening of the complaint).

5.1 Application

- > The patient should be seated.
 - > The knee joint should be slightly bent and relieved of pressure.
- 1) Apply the brace to the leg from the side. In doing so, the brace joint is on the lateral side of the affected compartment.
 - 2) Position the brace joint at the height of the middle of the patella (see fig. 2).
 - 3) Position the tibia pad on the lower leg and fasten the top lower leg strap. To do so, apply the snap-fit of the strap at an angle and push against the leg until the closure engages. In doing so, ensure that the tibia pad is under the closure (see fig. 3).

- 4) Push the dial down against the leg (see fig. 4).
- 5) Turn the dial clockwise until the lower leg strap is sitting firmly and comfortably (see fig. 5).
- 6) If the top lower leg strap is too tight, turn the dial anticlockwise (see fig. 5). This loosens the tieback in the lower leg strap. If the tieback is too loose, turn the dial clockwise to tighten the strap again.
- 7) Fasten the remaining straps. Do this by pulling the straps through the strap guide loops and fastening them. Do this in the following order: the bottom lower leg strap, the bottom thigh strap, the top thigh strap (see fig. 6, see fig. 7, see fig. 8, see fig. 9).
- 8) Ensure that the brace is in a firm and comfortable position.

5.2 Adaptation

Adjusting the corrective force

- > **Prerequisite:** The patient should be seated.
 - > **Prerequisite:** The knee joint should be relieved of pressure and slightly bent.
- 1) Adjust the corrective force of the brace using the Allen head screw above the brace joint (see fig. 10). This changes the varus/valgus angle of the brace.
 - 2) Check the corrective force while standing and walking.
 - 3) Adapt the corrective force until the desired adjustment has been achieved.

Optional: shorten the bottom lower leg and upper leg straps

- 1) Remove the Y-hook-and-loop.
- 2) Shorten the straps to the required length.
- 3) Stick the Y-hook-and-loop onto the strap.

Optional: shorten the top lower leg strap

- 1) Loosen the hook-and-loop closure on the strap guide loop.
- 2) Remove the hook-and-loop circle.
- 3) Shorten the strap to the desired length.
- 4) Fasten the hook-and-loop circle to the top lower leg strap.

Optional: limiting the knee flexion and extension

The extension limitation at delivery is 0°. The flexion is not limited.

- 1) Remove the condyle pad (see fig. 11).
- 2) To limit the extension, remove the anterior screw (see fig. 12). To limit the flexion, remove the posterior screw.
- 3) Choose the joint stops appropriate for the desired flexion and extension limitation (item 6, item 7, see fig. 1).

- 4) Insert the joint stops into the brace knee joint (see fig. 13). Insert the joint stop for extension limitation at the front and the joint stop for flexion limitation at the back. Ensure that the heads of the joint stops are pointing outwards. Secure the joint stop with the screw.
- 5) Attach the condyle pad.

5.3 Removal

- > The patient should be seated.
 - > The knee joint should be slightly bent and relieved of pressure.
- 1) Pull the dial out and turn it anticlockwise. This loosens the tieback in the top lower leg strap.
 - 2) Using your hand, partially pull out the top lower leg strap (see fig. 14).
 - 3) Pull the dial out of the closure clip at an angle (see fig. 15).
 - 4) Loosen all straps and take the brace off the leg.

6 Cleaning

NOTICE

Use of improper cleaning agents

Damage to the product due to use of improper cleaning agents

- Only clean the product with the approved cleaning agents.

- 1) Remove the straps and pads from the brace.
- 2) Fasten all hook-and-loop closures.
- 3) Hand-wash the straps and pads in warm water at **30 °C (86 °F)** using neutral detergent. Rinse thoroughly.
- 4) Allow to air dry. Do not expose to direct heat sources (e.g. sunlight, stove or radiator).
- 5) If the brace comes into contact with salt water or dirt, rinse with clear water and allow to air dry.

7 Disposal

Dispose of the product in accordance with national regulations.

8 Legal Information

All legal conditions are subject to the respective national law of the country of use and may vary accordingly.

8.1 Local Legal Information

Legal information that applies **exclusively** to specific countries is written in the official language of the respective country of use after this chapter.

8.2 Liability

The manufacturer shall be liable in the event that the product is used in accordance with the descriptions and instructions in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorized modification of the product.

8.3 CE Conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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