

50R59 Dorso Direxa Posture

| | | |
|-----------|----------------------------------|----|
| DE | Gebrauchsanweisung | 4 |
| EN | Instructions for use | 9 |
| FR | Instructions d'utilisation | 13 |
| IT | Istruzioni per l'uso | 18 |
| ES | Instrucciones de uso | 23 |
| NL | Gebruiksaanwijzing | 28 |
| SV | Bruksanvisning | 33 |
| DA | Brugsanvisning | 38 |
| NO | Bruksanvisning | 43 |
| FI | Käyttöohje | 47 |
| IS | Notkunarleiðbeiningar | 52 |
| PL | Instrukcja użytkowania | 56 |
| HU | Használati utasítás | 61 |
| CS | Návod k použití | 66 |
| SK | Návod na používanie | 71 |
| TR | Kullanma talimatı | 75 |
| EL | Οδηγίες χρήσης | 80 |
| RU | Руководство по применению | 85 |
| JA | 取扱説明書 | 91 |
| ZH | 使用说明书 | 95 |

9

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INFORMATION

Date of last update: 2022-03-31

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

These instructions for use provide important information on the fitting and application of the 50R59 Dorso Direxa Posture thoracolumbar brace.

2 Intended use

2.1 Indications for use

The brace is intended **exclusively** for orthotic fittings of the back and **exclusively** for contact with intact skin.

The brace must be used in accordance with the indications.

2.2 Indications

- Poor body posture (hunchback without fused vertebrae, Scheuermann's disease)
- Osteoporosis of the thoracic and lumbar spine
- Aftercare of osteoporotic fractures of the thoracic and/or lumbar spine
- Back pain (lumbar and lower thoracic spine), sciatica, muscular imbalance of the thoracic and lumbar spine, thoracic and lumbar spine ligament sprain
- Slight degenerative changes of the thoracic and lumbar spine (e.g. osteochondrosis, facet joint syndrome)

The indication must be determined by the physician.

2.3 Contraindications

2.3.1 Absolute Contraindications

None known.

2.3.2 Relative Contraindications

The following conditions require consultation with a physician: skin diseases/injuries, inflammation, raised scars with swelling, reddening, and warming in the area where the device is fitted, tumours, lymphatic congestion – including unclear soft tissue swelling distal to the device, and dysaesthesia in the trunk.

2.4 Lifetime

The product is designed for a maximum lifetime of **2 years**.

2.5 Mechanism of Action

The brace supports muscular stabilization of the thoracic spine and straightens it actively. It promotes sensorimotor function, helps alleviate pain and relieves the thoracic spine through active straightening and muscular stabilization and the lumbar spine by increasing intra-abdominal pressure and providing external stabilization.

3 Safety

3.1 Explanation of warning symbols

 **CAUTION** Warning regarding possible risks of accident or injury.

 **NOTICE** Warning regarding possible technical damage.

3.2 General safety instructions

 **CAUTION**

Contact with heat, embers or fire

Risk of injury (such as burns) and risk of product damage

- ▶ Keep the product away from open flames, embers and other sources of heat.

 **CAUTION**

Reuse on other persons and improper cleaning

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ The product may be used by one person only.
- ▶ Clean the product regularly.

NOTICE**Contact with oils, salves, lotions or other products that contain oils or acids**

Insufficient stabilization due to loss of material functionality

- ▶ Do not expose the product to oils, salves, lotions or other products that contain oils or acids.

4 Handling

INFORMATION

- ▶ The daily duration of use and period of application are generally determined by the physician.
- ▶ The initial fitting and application of the product must be carried out by qualified personnel.
- ▶ Instruct the patient in the handling and care of the product.
- ▶ Instruct the patient to see a physician immediately if any exceptional changes are noted (e.g. worsening of the complaint).

4.1 Selecting the Size

- 1) Measure the waist circumference approx. **2 cm** above the iliac crest.
- 2) Determine the brace size (see size chart).

4.2 Fitting and Application

⚠ CAUTION**Direct skin contact with the product**

Skin irritation due to chafing or sweat formation

- ▶ Do not wear this product directly on the skin.

⚠ CAUTION**Incorrect or excessively tight application**

Risk of local pressure and constriction of blood vessels and nerves due to improper or excessively tight application

- ▶ Ensure that the product is applied properly and fits correctly.

NOTICE**Use of a worn or damaged product**

Limited effectiveness

- ▶ Before each use, check the product for functional reliability and for possible wear or damage.
- ▶ Do not continue using a product that is no longer functional, or that is worn or damaged.

- 1) **Fitting the metal rods** (see fig. 1): Remove the 2 metal rods from the pockets and fit them to the spine individually. Put the metal rods back in the pockets.
- 2) **Optional – install the pad** (see fig. 2): Use the hook-and-loop closure to attach the pad to the inside of the brace so that it is at least 1 cm above the seat surface when seated.
- 3) **Application:** Place the upper part of the brace over the shoulders like a backpack (see fig. 3). Place the lower part of the brace around the trunk. Place the fingers in the hand straps of the abdominal closure. Apply equal force to both sides to pull the brace forward (see fig. 4).
- 4) **Hook-and-loop material:** Attach the right part to the top of the left part of the abdominal closure (see fig. 5). Remove the fingers from the left and then from the right hand strap. Ensure that the hook-and-loop sections lie flat against each other and are not in contact with clothing or skin. Route the ends of the shoulder strap to the front around the chest (see fig. 6) and secure the hook-and-loop closure below the chest (see fig. 7).
- 5) **Increasing compression:** Pull the ends of the elastic straps forward (see fig. 8) and attach them to the bottom closure (see fig. 9).
- 6) **Final check:** The correct fit and proper positioning of the brace must be checked during handover to the patient. The brace is seated correctly when the lower edge of the brace extends just below the groin.
- 7) **Instruct the patient:** Explain to the patient how to don the brace and ensure that it is seated correctly. Practice donning and doffing the brace with the patient.

4.3 Cleaning

NOTICE

Use of improper cleaning agents

Damage to the product due to use of improper cleaning agents

- ▶ Only clean the product with the approved cleaning agents.

Clean the brace regularly:

- 1) Close all hook-and-loop closures.
- 2) Recommendation: Use a laundry bag or net.

- 3) Wash the brace in warm water at **40 °C** with standard mild detergent. Do not use fabric softener. Rinse well.
- 4) Allow to air dry. Do not expose to direct heat sources (e.g. sunlight, stove or radiator).

5 Disposal

Dispose of the product in accordance with national regulations.

6 Legal Information

All legal conditions are subject to the respective national law of the country of use and may vary accordingly.

6.1 Local Legal Information

Legal information that applies **exclusively** to specific countries is written in the official language of the respective country of use after this chapter.

6.2 Liability

The manufacturer shall be liable in the event that the product is used in accordance with the descriptions and instructions in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorized modification of the product.

6.3 CE Conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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