


8E38=8*, 8E39=8*, 8E43=8*

DE Gebrauchsanweisung (Fachpersonal)	3
EN Instructions for use (qualified personnel)	23
FR Instructions d'utilisation (Personnel spécialisé)	43
IT Istruzioni per l'uso (Personale tecnico specializzato)	65
ES Instrucciones de uso (Personal técnico especializado)	85
PT Manual de utilização (Pessoal técnico)	106
NL Gebruiksaanwijzing (Vakmensen)	127
SV Bruksanvisning (Fackpersonal)	148
DA Brugsanvisning (Faguddannet personale)	167
NO Bruksanvisning (Fagpersonell)	187
HU Használati utasítás (szakszemélyzet)	207
TR Kullanma talimatı (Uzman personel)	226
EL Οδηγίες χρήσης (Τεχνικό προσωπικό)	245
RU Руководство по применению (Квалифицированный персонал)	267
JA 取扱説明書 (有資格担当者)	289

INFORMATION

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INFORMATION

Date of last update: 2022-03-23

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

The products "8E38=8*", 8E39=8* and 8E43=8*" are referred to as the product/prosthesis/terminal device in the following.

These instructions for use provide you with important information on the use, adaptation and handling of the product.

Only put the product into use in accordance with the information contained in the accompanying documents supplied.

According to the manufacturer (Otto Bock Healthcare Products GmbH), the patient is the operator of the product according to the IEC 60601-1:2005/A1:2012 standard.

2 Product description

2.1 Function

The SensorHand Speed is a prosthetic hand with myoelectric control, featuring a particularly high gripping speed combined with an innovative, sensitive control concept.

The SensorHand Speed is equipped with the SUVA sensors grip stabilisation system, the FlexiGrip function and a programming function. Proportional control with DMC (dynamic mode control) allows the patient to regulate the gripping speed and gripping force in proportion to the level of the muscle signal. When the strength of the muscle signal changes, the gripping speed and gripping force immediately adjust to the changed muscle signal.

The coaxial bushing of the 8E38=8* or 8E43=8 SensorHand Speed is marked with an orange ring, and the connecting cable of the 8E39=8* SensorHand Speed is marked with an orange sleeve.

Various control programmes make optimal adaptation to the individual needs and abilities possible.

SUVA sensors

The SUVA sensors integrated in the thumb (see fig. 1) recognise when a gripped object changes its position and is therefore at risk of slipping out of the user's grasp. The system then increases the gripping force automatically and continuously, starting from the originally specified initial gripping force, until the object is in a stable position again.

FlexiGrip function

The FlexiGrip function allows a grasped object to be turned or shifted in the hand without having to loosen the grip via electrode signals and subsequently re-grip the object.

The SensorHand Speed follows the changes in the position of the gripped objects like a natural hand would. This makes the grip flexible.

The FlexiGrip function can be stopped at any time with a brief "OPEN" impulse.

2.1.1 Wrist joint versions

The products "8E38=8*", 8E39=8* and 8E43=8*" are distinguished by different variants of the wrist joint:

8E38=* (quick-disconnect wrist)

Makes it easy to separate the terminal device from the prosthetic socket. The terminal device can be removed quickly when needed via a 360° rotating movement and replaced by a different terminal device with the same quick-disconnect unit.

8E39=* (lamination ring)

Low-profile connection for users with a long forearm or transcarpal amputation. The hand can be rotated against a constant friction, which can be adjusted during the fitting process.

The 9S110=* lamination ring has to be ordered separately.

The 13E190 or 13E190=150 switch block is required when this terminal device is used.


8E43=* (flexion joint)

The flexion joint with quick-disconnect wrist makes positioning possible for the user at 20° or 40° flexion, in the neutral position and at 20° or 40° extension. The terminal device can be removed quickly when needed via a 360° rotating movement and replaced by a different terminal device with the same quick-disconnect unit.

2.1.2 Unlocking/locking the wrist joint (8E43=*)

The individual flexion and extension of the wrist joint can be locked in five different positions (in 20° increments each).

- 1) Press the release button in the direction indicated by the arrow.
- 2) Move the terminal device to the desired position while holding the release button. From the relaxed hand position, it engages at 20° and 40° in each direction.
- 3) Let go of the release button to lock the terminal device in the respective position.

2.2 Combination possibilities

This product can be combined with the following Ottobock components:

Power supply (rechargeable battery)

- 757B35=1 / 757B35=3 / 757B35=5 MyoEnergy Integral
- 757B20/757B21 EnergyPack

Battery chargers

The following battery chargers can be used depending on the rechargeable battery:

- 757L20 charger (including 757L16* power supply)
- 757L35 battery charger

Lamination rings

- 10S1=* lamination ring (for 8E38=*, 8E43=*)
- 9S110=* lamination ring (for 8E39=*)

Active rotation

- 10S17 electric wrist rotator
- 13E205 MyoRotronic

Passive rotation

- 9E169 coaxial plug
- 10S4 coupling piece

Elbow

- ErgoArm Hybrid plus: 12K44=*
- ErgoArm Electronic plus: 12K50=*
- DynamicArm: 12K100N=*/12K110N=*

3 Intended use

3.1 Indications for use

The product is intended **exclusively** for upper limb exoprosthetic fittings.

3.2 Conditions of use

The product is intended **exclusively** for use on **one** user. Use of the product by another person is not approved by the manufacturer.

The product was developed for everyday use and must not be used for unusual activities. These unusual activities include, for example, sports which involve excessive strain and/or shocks to the wrist joint (pushups, downhill racing, mountain biking, ...) or extreme sports (free climbing, paragliding, etc.). Furthermore, the product should not be used to operate motor vehicles, heavy equipment (e.g. construction machines), industrial machines or motor-driven equipment.

Permissible ambient conditions are described in the technical data (see page 42).

3.3 Indications

- Amputation level below-elbow, above-elbow and shoulder disarticulation
- For unilateral or bilateral amputation
- Dysmelia of the forearm or upper arm
- The patient must be able to understand usage and safety messages and put them into practice.
- The patient must fulfil the physical and mental requirements for perceiving visual/acoustic signals and/or mechanical vibrations.

3.4 Contraindications

- All conditions which contradict or go beyond the specifications listed in the section on "Safety" and "Indications for use".

3.5 Qualification

The fitting of a patient with the product may only be carried out by O&P professionals who have been authorised with the corresponding Ottobock training.

4 Safety

4.1 Explanation of warning symbols

 WARNING	Warning regarding possible serious risks of accident or injury.
 CAUTION	Warning regarding possible risks of accident or injury.
 NOTICE	Warning regarding possible technical damage.

4.2 Structure of the safety instructions

 WARNING
The heading describes the source and/or the type of hazard
The introduction describes the consequences in case of failure to observe the safety instructions. Consequences are presented as follows if more than one consequence is possible:
> E.g.: Consequence 1 in the event of failure to observe the hazard
> E.g.: Consequence 2 in the event of failure to observe the hazard
▶ This symbol identifies activities/actions that must be observed/carried out in order to avert the hazard.

4.3 General safety instructions

 WARNING
Non-observance of safety notices
Personal injury/damage to the product due to using the product in certain situations.
▶ Observe the safety notices and the stated precautions in this accompanying document.

 WARNING
Operating a vehicle with the prosthesis
Accidents due to unexpected actions of the prosthesis.
▶ The prosthesis should not be used for the operation of motor vehicles and heavy equipment (e.g. construction machines).

 WARNING
Using the prosthesis while operating machines
Injury due to unexpected prosthesis behaviour.
▶ The prosthesis should not be used to operate industrial machines or motor-driven equipment.

 WARNING
Operating the prosthesis near active implanted systems
Interference with active implantable systems (e.g., pacemaker, defibrillator, etc.) due to electromagnetic emissions of the prosthesis.
▶ When operating the prosthesis in the immediate vicinity of active implantable systems, ensure that the minimum distances stipulated by the manufacturer of the implant are observed.

- ▶ Make sure to observe any operating conditions and safety instructions stipulated by the manufacturer of the implant.

⚠ WARNING

Use of damaged power supply unit, adapter plug or battery charger

Risk of electric shock due to contact with exposed, live components.

- ▶ Do not open the power supply unit, adapter plug or battery charger.
- ▶ Do not expose the power supply unit, adapter plug or battery charger to extreme loading conditions.
- ▶ Immediately replace damaged power supply units, adapter plugs or battery chargers.

⚠ WARNING

Skin contact with leaking lubricants due to mechanical defects

Injury due to skin irritation.

- ▶ Avoid contact between leaking lubricants and the mouth, nose and eyes.
- ▶ The product must be inspected by an authorised Ottobock Service Center.

⚠ CAUTION

Signs of wear on the product

Injury due to faulty control or malfunction of the product

- ▶ In the interest of patient safety and in order to maintain operating reliability, the terminal device has to be inspected by an authorised Ottobock Service Centre in case of noticeable restrictions of functionality.
- ▶ Note that the functionality of the terminal device may be limited if the battery charge level gets too low.

⚠ CAUTION

Independent manipulation of the product

Injury due to malfunction and resulting unexpected prosthesis actions.

- ▶ Manipulations to the product other than the tasks described in these instructions for use are not permitted.
- ▶ The product and any damaged components may only be opened and repaired by authorised, qualified Ottobock personnel.

⚠ CAUTION

Use of a damaged product

Injury due to loss of product functionality.

- ▶ Prior to use, conduct an external visual inspection to verify that all parts of the product are undamaged.
- ▶ In case of damage, have the product repaired promptly.

⚠ CAUTION

Penetration of product with dirt and humidity

Injury due to unexpected product behaviour or malfunction.

- ▶ Ensure that neither solid particles nor liquids can penetrate into the product.

4.4 Information on Alignment/Adjustment

CAUTION

Operator errors during the adjustment process with the adjustment software

Injury due to unexpected product behaviour.

- ▶ Participation in an Ottobock product training course is required prior to initial use. Additional product training courses may become necessary to qualify for software updates.
- ▶ Transfer the setting changes to the terminal device before you check the settings on the patient.
- ▶ Use the online help which is integrated into the software.

CAUTION

Incorrect electrode settings/electrode assignment

Injury due to unexpected product behaviour.

- ▶ The electrodes are to be placed on intact skin only and with as much electrode-skin contact as possible. In the case of strong interference from electronic devices, the position of the electrodes should be checked and changed if necessary. If the interference cannot be eliminated or if you do not achieve the expected results by adjustment or selection of the appropriate control programme, please contact the Ottobock branch responsible for your country.
- ▶ Set the electrode sensitivity as low as possible in order to reduce interference from powerful electromagnetic radiation (e.g. visible or concealed theft prevention systems at the entrance/exit of stores), metal detectors/body scanners for persons (e.g. in airports) or other sources of strong electromagnetic interference (e.g. high-voltage lines, transmitters, transformer stations, computer tomographs, magnetic resonance tomographs, etc.).
- ▶ Make sure that the electrode connection positions correspond to physiological opening and closing for the corresponding muscle groups.

CAUTION

Incorrect electrode settings due to muscle fatigue

Injury due to faulty control or malfunction of the product.

- ▶ Allow the patient to rest during the adjustment of the electrodes.

CAUTION

Insufficient skin contact of the electrodes

Injury due to unexpected product behaviour as a result of uncontrolled activation of the prosthesis component.

- ▶ The electrodes are to be placed on intact skin only and with as much electrode-skin contact as possible.
- ▶ Make sure that skin contact of the electrodes continues to be given when carrying heavy loads.
- ▶ If the product cannot be controlled correctly using muscle signals, turn off the entire prosthesis and contact your O&P professional.

CAUTION

Failure to observe the instructions for use of all prosthesis components used

Injury due to unexpected product behaviour.

- ▶ Observe the instructions for use of the prosthesis components used.

⚠ CAUTION

Use of unsuitable prosthetic components

Injury due to unexpected product behaviour.

- ▶ Use the product only in combination with components listed in the section "Combination possibilities" (see page 26).

⚠ CAUTION

Battery not connected correctly

Injury due to faulty control or malfunction of the product.

- ▶ Make sure the battery connections are not switched.
- ▶ Check the battery connections by pulling slightly on the cable.

⚠ CAUTION

Use of silicone spray when putting on the prosthetic glove

Injury due to a loss of grip resulting from insufficient adhesion of the prosthetic glove on the hand.

- ▶ Only use 640F18 donning spray to put on the prosthetic glove. Note the instructions for use for the prosthetic glove.

4.5 Information on Proximity to Certain Areas

⚠ CAUTION

Insufficient distance to HF communication devices (e.g. mobile phones, Bluetooth devices, WiFi devices)

Injury due to unexpected behaviour of the product caused by interference with internal data communication.

- ▶ Therefore, keeping a minimum distance of 30 cm from HF communication devices is recommended.

⚠ CAUTION

Operating the product in very close proximity to other electronic devices

Injury due to unexpected behaviour of the product caused by interference with internal data communication.

- ▶ Do not operate the product in the immediate vicinity of other electronic devices.
- ▶ Do not stack the product with other electronic devices during operation.
- ▶ If simultaneous operation cannot be avoided, monitor the product and verify proper use in the existing setup.

⚠ CAUTION

Proximity to sources of strong magnetic or electrical interference (e.g. theft prevention systems, metal detectors)

Injury due to unexpected behaviour of the product caused by interference with internal data communication.

- ▶ Avoid remaining in the vicinity of visible or concealed theft prevention systems at the entrance/exit of stores, metal detectors/body scanners for persons (e.g. in airports) or other sources of strong magnetic and electrical interference (e.g. high-voltage lines, transmitters, transformer stations, computer tomographs, magnetic resonance tomographs, etc.).
- ▶ When walking through theft prevention systems, body scanners or metal detectors, watch for unexpected behaviour of the product.

4.6 Information on Use

CAUTION

Improper use

Injury due to faulty operation or malfunction of the product.

- ▶ Instruct the patient in the proper use of the product.

CAUTION

Improper product care

- > Injuries due to faulty control/malfunction of the product or damage to the mechanical components
- > Damage or breakage due to brittleness of plastics caused by the use of acetone, petrol or similar solvents.
- ▶ Clean the product only as described in the "Cleaning" section (see page 41).
- ▶ Do not clean the product under running water.
- ▶ When using a prosthetic glove, also note the instructions for use of the glove.

CAUTION

Grasping objects with incorrect gripping forces

Injury due to unexpected product behaviour.

- ▶ Note that the gripping has to be controlled manually depending on the consistency (soft/hard) of the object being grasped.

CAUTION

Risk of pinching between the fingertips

Injury due to pinching of body parts.

- ▶ Ensure that no body parts are between the fingertips when using the product.
- ▶ Ensure that no body parts are between the fingertips when closing the hand.
- ▶ When closing the hand, ensure that fingers and other body parts are not in the area of the finger joints.
- ▶ Make sure the product is switched off for cleaning.

CAUTION

Unintentional unlocking of the terminal device

Risk of injury due to the terminal device disconnecting from the forearm (e.g. while carrying objects).

- ▶ Ensure that the connection is carried out correctly when connecting the hand to the prosthetic socket or components.
- ▶ When using a quick-disconnect wrist unit, position the terminal device so that slight twisting cannot disconnect the terminal device from the forearm.

5 Scope of delivery

- 1 pc. 8E38=8* SensorHand Speed
- or
- 1 pc. 8E39=8* SensorHand Speed
- or
- 1 pc. 8E43=8* SensorHand Speed
 - 1 pc. Instructions for use (qualified personnel)
 - 1 pc. Instructions for use (user)

6 Preparing the product for use

6.1 Charging the battery

The following information is found in the instructions for use of the respective batteries or elbow components:

- Handling the battery
- Querying the charge level
- Feedback (beeps and vibration signals)

6.2 Adjusting the electrodes

The ideal position of the electrodes is determined with the 757M11=X-CHANGE MyoBoy, see 647G265=ALL_INT instructions for use.

Program 1, 5 and 6:

Select the desired program variant in the PAULA program. Adjust each electrode so the patient can hold the respective muscle signal above the HIGH value for about two seconds (see fig. 7, (1)).

Program 2:

Select the desired program variant in the PAULA program. Adjust each electrode so the patient can hold the respective muscle signal above the LOW value for about two seconds (see fig. 7, (2)).

Program 3:

In the PAULA program, select the SensorHand Speed in the AutoControl mode. Adjust the electrode so the patient can hold the muscle signal above the ON value for about two seconds, see fig. 8.

Program 4:

In the PAULA program, select the SensorHand Speed in the VarioControl mode. Adjust the electrode so the patient can hold the muscle signal above the HIGH value for about two seconds, see fig. 8.

INFORMATION

Only one muscle signal can be viewed with programs 3 and 4, see fig. 7 and see fig. 8.

6.3 Control programs

Six programs are provided for optimal adaptation to the patient. These programs are selected using coding plugs of different colours (see fig. 4).

The SensorHand Speed is equipped with a black coding plug and preconfigured with program 1 at the factory. When the black coding plug is used, the programs are selected with the 757T13 MyoSelect.

Adjusting the speed of the SensorHand Speed with the 757T13 MyoSelect can be carried out with any coding plug; see the 647G131 instructions for use.

Changing the coding plug, see fig. 4

- 1) Remove the lock ring.
- 2) Push the inner hand back so the control electronics are exposed.
- 3) Change the coding plug.
- 4) Push the inner hand forward.
- 5) Install the lock ring.
- 6) Remove and reinstall the battery.

INFORMATION: This causes the control device to recognise the new coding plug.

6.3.1 Program overview

Program 1: DMC plus Sensorik

Coding plug, white: 13E184=1

Activation	Open	Close	Indication
Two electrodes	Muscle signal via the electrode Speed: Proportional	Muscle signal via the electrode Speed: Proportional	For patients with two strong muscle signals

Program 2: AutoControl LowInput

Coding plug, red: 13E184=2

Activation	Open	Close	Indication
Two electrodes	Muscle signal via the electrode Speed: Proportional	Muscle signal via the electrode, digital (short signal at any level) Speed: Constant	For patients with two weak muscle signals
One electrode and any MyoBock switch	Muscle signal via the electrode Speed: Proportional	Signal via the switch Speed: Constant	For patients with only one muscle with a weak muscle signal
Any MyoBock switch	The hand opens as long as the OPEN side of the switch is activated. Speed: Constant	The hand closes via the CLOSE side of the switch. Speed: Constant	For patients with muscle signals that are too weak or no muscle signals

Program 3: AutoControl

Coding plug, green: 13E184=3

Activation	Open	Close	Indication
One electrode	Fast, sustained muscle signal via the electrode Speed: Constant	Very slow muscle relaxation via the electrode: hand stops in the open position Quick muscle relaxation via the electrode: hand closes Speed: Constant	For patients with only one muscle with a very weak muscle signal
Any MyoBock switch	The hand opens as long as the switch is activated. Speed: Constant	The hand closes automatically as soon as the switch is released. Speed: Constant	For patients with muscle signals that are too weak or no muscle signals

Program 4: VarioControl

Coding plug, blue: 13E184=4

Activation	Open	Close	Indication
One electrode	Speed and level of muscle contraction on the electrode Speed: Proportional	Speed and level of muscle relaxation on the electrode Speed: Proportional	For patients with one muscle and a strong muscle signal or tendency to co-contraction
One linear control element	Speed and level of pull on the linear control element Speed: Proportional	Speed of relaxing the pull on the linear control element Speed: Proportional	For patients with muscle signal that is too weak or no muscle signal

Program 5: VarioDual

Coding plug, yellow: 13E184=5

Activation	Open	Close	Indication
Two electrodes	Speed and level of muscle contraction on the first electrode Speed: Proportional	Speed and level of muscle relaxation on the first electrode Speed: Proportional Gripping force proportional to strength of the muscle signal on the second electrode	For patients with two strong muscle signals

Program 6: DMC plus Sensorik

Coding plug, purple: 13E184=6

Activation	Open	Close	Indication
Two electrodes	Muscle signal via the electrode Speed: Proportional	Muscle signal via the electrode Speed: Proportional	For patients with two strong muscle signals SUVA sensors and FlexiGrip can be deactivated

6.3.2 Program description

INFORMATION

Both the automatic gripping force readjustment and the FlexiGrip function can be stopped at any time with a brief "OPEN" impulse.

INFORMATION

All gripping force specifications (e.g. 100 N ~ 10 Kp) are non-binding guiding values to illustrate the functional principle of the various operating modes.

6.3.2.1 Program 1: DMC plus Sensorik

Control with two electrodes

This control variant corresponds to DMC plus control with the integrated "virtual hand switch" (dynamic mode control) but also features the "SUVA sensors" grip stabilisation system. The level of the gripping speed or gripping force is determined by the strength of the electrode signal (resulting from the muscle contraction). After gripping with maximum force, the activation threshold in the OPEN direction is raised to a higher value ("virtual hand switch"). This increase reduces the

risk of opening the hand with unintended muscle signals. The grip is more reliable, for example, when holding cutlery.

Opening: Proportional via the OPEN electrode
 Closing: Proportional via the CLOSE electrode

Example 1: With a low muscle signal, the lowest gripping force (10 N) is built up to grasp an object. When the sensors recognise that the position of the object is changing, a readjustment is automatically performed at up to 1.5 times the initial gripping force (15 N) as needed. FlexiGrip takes effect starting at 20 N. When there is no longer a load, the SensorHand Speed grips with the previous gripping force again.

Example 2: A higher muscle signal produces a higher gripping force and, when the position of the gripped object changes, a readjustment is performed up to the maximum gripping force (130 N) as needed. FlexiGrip takes effect when the load on the hand exceeds 130 N (gripping force of the hand and externally applied force). When there is no longer a load, the SensorHand Speed grips with the previous gripping force again.

Open	Close
Muscle signal via the electrode	Muscle signal via the electrode
Speed: Proportional, 15 mm/s to 300 mm/s	Speed: Proportional, 15 mm/s to 300 mm/s

Initial gripping force	Automatic gripping force readjustment	FlexiGrip function
Proportional: 0 N bis 100 N	Proportional: Up to max. 1.5 times the initial gripping force, e.g. if the initial gripping force is 10 N, gripping force readjustment up to max. 15 N	Depending on the initial gripping force, takes effect slightly above the maximum gripping force readjustment in each case Min. from 20 N Max. from 130 N

The gripping force can be increased up to the maximum gripping force (100 N) at any time as needed via a stronger muscle signal, independently of the automatic gripping force readjustment.

6.3.2.2 Program 2: AutoControl - LowInput

Control with: Two electrodes, one electrode and one switch or one switch

The hand closes at the fastest speed and grips an object with the lowest gripping force (10 N). When the sensors recognise that the position of the object is changing, a readjustment is automatically performed up to the maximum gripping force (130 N) as needed. FlexiGrip takes effect starting at the maximum gripping force. When there is no longer a load, the SensorHand Speed grips with the previous gripping force again.

Control with two electrodes

Opening: Proportional via the OPEN electrode.
 Closing: At maximum speed with short muscle signal at any level over the ON threshold on the CLOSE electrode.

Open	Close
------	-------

Muscle signal via the electrode Speed: Proportional, 15 mm/s to 300 mm/s	Muscle signal via the electrode Speed: Constant, 300 mm/s
--	---

Initial gripping force	Automatic gripping force readjustment	FlexiGrip function
10 N	Up to max. 130 N	Takes effect starting at the maximum gripping force

Control with one electrode and one switch

Opening: Proportional via the OPEN electrode.
Closing: At maximum speed by briefly activating the switch.

Open	Close
Muscle signal via the electrode Speed: Proportional, 15 mm/s to 300 mm/s	Signal via the CLOSE side of the switch Speed: Constant, 300 mm/s

Initial gripping force	Automatic gripping force readjustment	FlexiGrip function
10 N	Up to max. 130 N	Takes effect starting at the maximum gripping force

Control with one switch

This program can be used in combination with any MyoBock switch.

Opening: At maximum speed as long as the OPEN side of the switch is activated. The hand then remains open.
Closing: At maximum speed by activating the CLOSE contact of the switch.

Open	Close
The hand opens as long as the OPEN side of the switch is activated Speed: Proportional, 15 mm/s to 300 mm/s	Signal via the CLOSE side of the switch Speed: Constant, 300 mm/s

Initial gripping force	Automatic gripping force readjustment	FlexiGrip function
10 N	Up to max. 130 N	Takes effect starting at the maximum gripping force

Connect the following cables to the 9E169 coaxial plug.

- see fig. 6 (1): battery connection cable
- see fig. 6 (2): electrode
- see fig. 6 (3): red/white cable of the 13E99 connecting cable (closing contacts)

6.3.2.3 Program 3: AutoControl

Control with: One electrode or one linear control element or one switch

The hand closes at the fastest speed and grips an object with the lowest gripping force (10 N). When the sensors recognise that the position of the object is changing, a readjustment is automatically and continuously performed up to the respective gripping force required (max. 130 N).

FlexiGrip takes effect when the load on the hand exceeds 130 N (gripping force of the hand and externally applied force). When there is no longer a load, the SensorHand Speed grips with the previous gripping force again.

Control with one electrode

- Opening: At maximum speed with fast, sustained muscle signal via the electrode.
- Closing: At maximum speed with quick relaxation of the muscle.
- Stopping: Relaxing the muscles very slowly via the electrode causes the hand to stop in the open position.
- Example 1: After opening, relax the muscle at very slow speed. The opening position remains unchanged.
- Example 2: After opening, relax the muscle at the fastest speed. The hand closes automatically at the fastest speed and starts grasping the object with a gripping force of 10 N. When the sensors recognise that the position of the object is changing, a readjustment is automatically performed up to the maximum gripping force (130 N) as needed.

Open	Close
Fast, sustained muscle signal via the electrode. Speed: Constant, 300 mm/s	Relaxing the muscles very slowly via the electrode: Hand stops in the open position. Relaxing the muscles quickly via the electrode: Hand closes Speed: Constant, 300 mm/s

Initial gripping force	Automatic gripping force readjustment	FlexiGrip function
10 N	Up to max. 130 N	Takes effect starting at the maximum gripping force

Control with one linear control element

- Opening: At maximum speed with a quick pull on the linear control element.
- Closing: At maximum speed by quickly relaxing the pull on the linear control element.
- Stopping: Relaxing the pull on the linear control element very slowly causes the hand to stop in the open position.
- Example 1: After opening, relax the pull on the linear control element very slowly. The opening position remains unchanged.
- Example 2: After opening, relax the pull on the linear control element at the fastest speed. The hand closes automatically at the fastest speed and starts grasping the object with a gripping force of 10 N. When the sensors recognise that the position of the object is changing, a readjustment is automatically performed up to the maximum gripping force (130 N) as needed.

Open	Close
Pull on the linear control element at high speed. Speed: Constant, 300 mm/s	Relaxing the pull on the linear control element very slowly: Hand stops in the open position Relaxing the pull on the linear control element quickly:

	Hand closes Speed: Constant, 300 mm/s
--	---

Initial gripping force	Automatic gripping force readjustment	FlexiGrip function
10 N	Up to max. 130 N	Takes effect starting at the maximum gripping force

Control with one switch

Opening: At maximum speed as long as the switch is activated.

Closing: After releasing the switch, the hand closes automatically at the fastest speed and starts grasping the object with a gripping force of 10 N.

Open	Close
Opens as long as the switch is activated. Speed: Constant, 300 mm/s	Closes automatically as soon as the switch is released. Speed: Constant, 300 mm/s

Connect the following cables to the 9E169 coaxial plug.

- see fig. 5 (1): battery connection cable
- see fig. 5 (2): red/white cable of the 13E99 connecting cable (opening contacts)

6.3.2.4 Program 4: VarioControl

Control with: One electrode or one linear control element

With this program, the opening speed is determined by the level and speed of the muscle contraction. The closing speed depends on the reduction of muscle tension. FlexiGrip takes effect starting at the maximum gripping force. When there is no longer a load, the SensorHand Speed grips with the previous gripping force again.

Control with one electrode

Opening: Proportional. The opening speed is determined by the speed and strength of the muscle contraction.

Closing: Proportional. The closing speed is determined by the speed and strength with which the muscles are relaxed. This also specifies the level of the maximum gripping force readjustment.

Stopping: Relaxing the muscles very slowly via the electrode causes the hand to stop in the open position.

Example 1: After opening, relax the muscle at slow speed. The hand is closed analogue to the duration of muscle relaxation, at slow speed. The object is gripped with low force (10 N). There is no automatic gripping force readjustment.

Example 2: After opening, relax the muscle at the fastest speed. The hand closes automatically at the fastest speed and starts grasping the object with a gripping force of 10 N. When the sensors recognise that the position of the object is changing, a readjustment is automatically performed up to the maximum gripping force (130 N) as needed.

Open	Close
------	-------

Via speed and strength of the muscle contraction on the electrode. Speed: Proportional, 15 mm/s to 300 mm/s	Via speed and strength of muscle relaxation on the electrode. Speed: Proportional, 15 mm/s to 300 mm/s
---	--

Initial gripping force	Automatic gripping force readjustment	FlexiGrip function
10 N	At low to moderate closing speed: None	Effective from 20 N
10 N	At moderate to high closing speeds: Up to max. 130 N	Takes effect starting at the maximum gripping force

Control with one linear control element

- Opening: Proportional. The opening speed is determined by the speed and strength of the pull on the linear control element.
- Closing: Proportional. The closing speed is determined by the speed at which the pull on the linear control element is relaxed. This also specifies the level of the maximum gripping force readjustment.
- Stopping: Relaxing the pull on the linear control element very slowly causes the hand to stop in the open position.
- Example 1: After opening, relax the pull on the linear control element slowly. The hand is closed analogue to the duration of muscle relaxation, at slow speed. The object is gripped with low force (10 N). There is no automatic gripping force readjustment.
- Example 2: After opening, relax the pull on the linear control element at the fastest speed. The hand closes at the fastest speed and starts grasping the object with a gripping force of 10 N. When the sensors recognise that the position of the object is changing, a readjustment is automatically performed up to the maximum gripping force (130 N) as needed.

Open	Close
Via speed and strength of the pull on the linear control element Speed: Proportional, 15 mm/s to 300 mm/s	Via speed at which the pull on the linear control element is relaxed Speed: Proportional, 15 mm/s to 300 mm/s

Initial gripping force	Automatic gripping force readjustment	FlexiGrip function
10 N	At low to moderate closing speed: None	At low to moderate closing speed: Effective from 15 N
10 N	At moderate to high closing speeds: Up to max. 130 N	At moderate to high closing speeds: Effective starting at the maximum gripping force

6.3.2.5 Program 5: VarioDual

Control with two electrodes

With this program, the opening speed is determined by the level and speed of the muscle contraction. The closing speed up to reaching the minimum gripping force of approx. 10 N depends on the speed of muscle relaxation. The gripping force is determined by the subsequent or simultaneous muscle signal on the second electrode. Depending on the initial gripping force, FlexiGrip takes effect slightly above the maximum gripping force readjustment in each case. When there is no longer a load, the SensorHand Speed grips with the previous gripping force again.

Electrode 1

- Opening: Proportional. The opening speed is determined by the speed and strength of the muscle contraction.
- Closing: Proportional. The closing speed is determined by the speed and strength with which the muscles are relaxed. The gripping force is approx. 10 N.
- Stopping: Relaxing the muscles very slowly via the electrode causes the hand to stop in the open position.

Electrode 2

- Gripping: The build-up of the gripping force is determined by the level of the muscle signal on the second electrode. The maximum gripping force is approx. 100 N.
- Example 1: After opening, relax the muscle at any speed. Closing is proportional to the speed of muscle relaxation. The object is grasped with the minimum gripping force (10 N). FlexiGrip takes effect starting at 20 N. When there is no longer a load, the SensorHand Speed grips again with the minimum gripping force.
- Example 2: After grasping as in example 1, the object is to be gripped with a higher gripping force. To do so, produce a muscle signal on the second electrode. A gripping force between 10 N and 100 N can be built up proportionally. When the position of the gripped object changes, the gripping force is increased by up to approx. 1.5 times the specified gripping force. FlexiGrip takes effect at about twice the value of the specified gripping force up to max. 130 N; when there is no longer a load, the SensorHand Speed grips again with the original gripping force.

Open	Close
Via speed and strength of the muscle contraction on electrode 1 Speed: Proportional, 15 mm/s to 300 mm/s	Via speed and strength of muscle relaxation on electrode 1 Speed: Proportional, 15 mm/s to 300 mm/s Build-up of gripping force: Gripping force depending on the level of the muscle signal on electrode 2. Gripping force: Proportional, 10 N to 100 N

Initial gripping force	Automatic gripping force readjustment	FlexiGrip function
10 N Proportional: 10 N to 100 N	During closing: No gripping force readjustment During build-up of gripping force: Proportional, up to max. 1.5 times the initial gripping force	From 20 N depending on the initial gripping force, takes effect slightly above the maximum gripping force readjustment in each case Min. from 20 N Max. from 130 N

6.3.2.6 Program 6: DMC plus Sensorik, can be deactivated

Control with two electrodes

This control variant corresponds to program 1, but the “SUVA sensors” and FlexiGrip function can be temporarily deactivated.

Activating and deactivating the “SUVA sensors” and FlexiGrip function

The “SUVA sensors” can be deactivated for gripping very soft and pliable objects, such as very soft foams or tweezers. To do so, open the SensorHand Speed to the stop and keep it open with a muscle signal at any level. At the same time, apply some pressure to the “SUVA sensors” (Fig. 1), e.g. by pushing on the edge of a table. A short vibration signal confirms deactivation. Repeat the same process to activate the “SUVA sensors”. Two short vibration signals confirm activation

of the “SUVA sensors”.

INFORMATION: Note that when the “SUVA sensors” are deactivated, the gripping force is not readjusted automatically and gripped objects may slip from the grasp. After the battery is inserted, vibration signals inform you of the current mode.

One vibration signal: Sensors are deactivated

Two vibration signals: Sensors are activated

7 Use

7.1 Terminal device on/off

The terminal device can be switched off when a grip is to be maintained for a longer period of time, while holding cutlery, writing or during longer passive pauses, for example, while travelling by rail or air. This prevents unintentional opening of the hand caused by inadvertent tensing of the muscles or extreme electrical interference. The operating time of the prosthesis battery is also extended.

The ON/OFF button can be operated by pressing on the prosthetic glove.

Area	Function
Back of the hand	ON (see fig. 2)
Thumb	OFF (see fig. 3)

7.2 Opening the terminal device in an emergency

The integrated slip clutch allows the terminal device to be opened regardless of the active control signals.

8 Storage

If the System Electric Hand will not be used, store the System Electric Hand in the open position to protect the sensors and mechanism.

9 Cleaning

- 1) Switch off the product before cleaning.
- 2) Clean the product with a damp cloth and mild soap when needed.
Make sure that no liquids get into the product and product components.
- 3) Dry the product with a lint-free cloth and allow it to air dry fully.

INFORMATION

When using a prosthetic glove, note the cleaning information in the instructions for use for the prosthetic glove.

10 Maintenance

Performing regular maintenance (service inspections) every 24 months is recommended to prevent injuries and maintain the quality of the product.

In general, all products are subject to compliance with the maintenance intervals during the warranty period. This is the only way to maintain full warranty cover.

Additional services such as repairs may be provided in the course of maintenance. These additional services may be provided free of charge or can be billable according to an advance cost estimate, depending on the extent and validity of the warranty.

11 Legal information

11.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

11.2 Trademarks

All product names mentioned in this document are subject without restriction to the respective applicable trademark laws and are the property of the respective owners.

All brands, trade names or company names may be registered trademarks and are the property of the respective owners.

Should trademarks used in this document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

11.3 CE conformity

Otto Bock Healthcare Products GmbH hereby declares that the product is in compliance with applicable European requirements for medical devices.

The product meets the requirements of the RoHS Directive 2011/65/EU on the restriction of the use of certain hazardous substances in electrical and electronic devices.

The full text of the regulations and requirements is available at the following Internet address: <http://www.ottobock.com/conformity>

12 Technical data

Ambient conditions	
Storage (with and without packaging)	+5 °C/+41 °F to +40 °C/+104 °F Max. 85% relative humidity, non-condensing
Transport (with and without packaging)	-20 °C/-4 °F to +60 °C/+140 °F Max. 90% relative humidity, non-condensing
Operation	-5 °C/+23 °F to +45 °C/+113 °F Max. 95% relative humidity, non-condensing

Reference number	8E38=8*, 8E39=8*	8E43=8*
Opening width	100 mm	
Proportional speed	15–300 mm/s	
Proportional gripping force	0–100 N	
Weight (including system inner hand)	462 g	517 g
Operating voltage	Approx. 7.4 V	
Static current	2 mA	
Lifetime	5 years	

Rechargeable battery reference number	757B20	757B21	757B35 =1	757B35 =3	757B35 =5
Temperature range for charging [°C]	+5 to +40				
Capacity [mAh]	900 ¹ / 950 ¹	680 ¹ / 800 ¹	600	1150	3450
Rechargeable battery lifetime [years]	2				
Product behaviour during the charging process	The product is non-functional				

Rechargeable battery reference number	757B20	757B21	757B35 =1	757B35 =3	757B35 =5
Operating time of the product with fully charged battery [grip cycles]	Approx. 2500– 3000	Approx. 2000– 2500 (depend- ing on the bat- tery capacity)	Approx. 1000– 2000	Approx. 4000	Approx. 10000
Charging times (rechargeable battery fully dis- charged) [hours]	Approx. 3.5	Approx. 3	Approx. 2.5		Approx. 3
Nominal voltage [V]	Approx. 7.4				
Battery technology	Li-Ion		Li-Po	Li-Ion	
Approved battery chargers	757L20		757L35		

¹ See labeling on battery

13 Appendices

13.1 Symbols Used

Manufacturer

In some jurisdictions it is not permissible to dispose of these products with unsorted household waste. Disposal that is not in accordance with the regulations of your country may have a detrimental impact on health and the environment. Please observe the instructions of your national authority pertaining to return and collection.

Declaration of conformity according to the applicable European directives

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