





50P20 Manu Arexa

50P21 Manu Arexa Pollex

DE Gebrauchsanweisung	5	FI Käyttöohje.....	35
EN Instructions for use	8	PL Instrukcja użytkowania.....	38
FR Instructions d'utilisation.....	11	HU Használati utasítás	42
IT Istruzioni per l'uso	14	CS Návod k použití.....	45
ES Instrucciones de uso	17	SK Návod na používanie	48
PT Manual de utilização.....	20	ET Kasutusjuhis	51
NL Gebruiksaanwijzing.....	23	RU Руководство по применению	54
SV Bruksanvisning	26	JA 取扱説明書.....	57
DA Brugsanvisning	29	ZH 使用说明书.....	60
NO Bruksanvisning.....	32		

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INFORMATION

Date of last update: 2022-03-31

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

These instructions for use provide important information on the fitting and application of the 50P20 Manu Arexa and 50P21 Manu Arexa Pollex wrist braces.

2 Intended use

2.1 Indications for use

The brace is intended **exclusively** for orthotic fittings of the upper limbs and **exclusively** for contact with intact skin.

The brace must be used in accordance with the indications.

2.2 Indications

- Sprains of the wrist joint or the thumb saddle joint

- Rheumatoid arthritis
- Wrist drop due to paresis or paralysis of the radial nerve

50P20 Manu Arexa additional indications

- Tendinopathy of the wrist joint
- Ligament injury, ligament strain in the wrist
- Wrist instability
- Impingement syndrome
- Irritation in the wrist area
- Carpal tunnel syndrome

50P21 Manu Arexa Pollex additional indications

- Fractures of the distal forearm or carpal bone
- Ligament rupture (partial rupture)
- Partial arthrodesis
- Preoperative/postoperative/posttraumatic conditions that require immobilization of the wrist
- Irritation in the wrist and thumb areas
- Night positioning

The indication must be determined by the physician.

2.3 Contraindications

2.3.1 Absolute Contraindications

None known.

2.3.2 Relative Contraindications

The following indications require consultation with a physician: skin diseases/injuries, inflammation, prominent, swollen scars, reddening and hyperthermia in the relevant area; pronounced varicose veins, especially with impaired blood return, lymph drainage disorders and unclear soft tissue swelling distal to the device; sensory and circulatory disorders in the area of the forearm and hand, e.g. associated with diabetic neuropathy.

2.4 Lifetime

The product is designed for a maximum lifetime of **1 year**.

2.5 Mechanism of Action

The brace stabilizes the wrist.

3 Safety

3.1 Explanation of warning symbols

⚠ CAUTION

Warning regarding possible risks of accident or injury.

NOTICE

Warning regarding possible technical damage.

3.2 General safety instructions

⚠ CAUTION

Reuse on other persons and improper cleaning

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ The product may be used by one person only.
- ▶ Clean the product regularly.

⚠ CAUTION

Contact with heat, embers or fire

Risk of injury (such as burns) and risk of product damage

- ▶ Keep the product away from open flames, embers and other sources of heat.

NOTICE

Use of a worn or damaged product

Limited effectiveness

- ▶ Before each use, check the product for functional reliability and for possible wear or damage.
- ▶ Do not continue using a product that is no longer functional, or that is worn or damaged.

NOTICE

Improper molding or application

Damage to the product due to overloading of the material and improper fit of the product due to breakage of load-bearing components

- ▶ The product may only be molded by qualified personnel.
- ▶ Do not make any improper changes to the product.
- ▶ Always apply the product according to the information in the instructions.

NOTICE

Contact with oils, salves, lotions or other products that contain oils or acids

Insufficient stabilization due to loss of material functionality

- ▶ Do not expose the product to oils, salves, lotions or other products that contain oils or acids.

4 Handling

INFORMATION

- ▶ The daily duration of use and period of application are generally determined by the physician.
- ▶ The initial fitting and application of the product must be carried out by qualified personnel.
- ▶ Instruct the patient in the handling and care of the product.
- ▶ Instruct the patient to see a physician immediately if any exceptional changes are noted (e.g. worsening of the complaint).

4.1 Selecting the Size

- 1) Measure the wrist circumference.
- 2) Determine the brace size (see size chart).

4.2 Fitting and Application

⚠ CAUTION

Incorrect or excessively tight application

Skin irritation due to overheating, local pressure points due to tight application to the body

- ▶ Ensure that the product is applied properly and fits correctly.
- ▶ When in doubt, do not continue to use the product if skin irritation occurs.

The 50P20 and 50P21 braces are pre-configured for the right hand upon delivery (see fig. 1, see fig. 2). The metal splints must be modified for use on the left hand (2772886). Two metal splints – a palmar spoon-shaped splint and a dorsal splint – are located in the outer pockets of the brace. Shape the splints by hand to customize them for the user.

- 1) Open the brace.
- 2) Remove the palmar splint (A) from the outer pocket of the brace (see fig. 6, see fig. 4).
- 3) Mold the splint so it fits the anatomical shape of the patient's palm.
- 4) Insert the splint into the pocket with the outer spoon-shaped curve facing the palm.
- 5) **50P20: The dorsal splint (B1) can optionally be used on the ulnar side (B2).** Remove the splint from the outer pocket, shape it by hand, and place it in the dorsal or ulnar pocket (see fig. 4).
- 6) **50P21:** Remove the dorsal splint (B1) from the outer pocket, shape it by hand, and place it in the pocket.
- 7) Pull the brace over the wrist, guiding the thumb through the thumb opening (see fig. 5). Tighten the lacing and fasten the hook-and-loop closure.
- 8) Check the fit of the brace.
- 9) **Optional:** Trim the hook-and-loop closure with scissors to better customize the fit of the brace (see fig. 7).

→ The brace fits snugly without pinching.

4.3 Removal

- 1) Open the hook-and-loop closure of the brace.
- 2) Take off the brace.

5 Cleaning

NOTICE

Use of improper cleaning agents

Damage to the product due to use of improper cleaning agents

- ▶ Only clean the product with the approved cleaning agents.

- 1) Remove all metal splints.
- 2) Fasten all hook-and-loop closures.
- 3) Recommendation: Use a laundry bag or net.
- 4) Wash the brace in warm water at **40 °C** (104 °F) with standard mild detergent. Do not use fabric softener. Rinse thoroughly.
- 5) Allow to air dry. Do not expose to direct heat sources (e.g. sunlight, stove or radiator).
- 6) Reinsert the metal splints.

6 Disposal

Dispose of the product in accordance with national regulations.

7 Legal Information

All legal conditions are subject to the respective national law of the country of use and may vary accordingly.

7.1 Liability

The manufacturer shall be liable in the event that the product is used in accordance with the descriptions and instructions in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorized modification of the product.

7.2 CE Conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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