## 6Y44, 6Y45 Compression Liner

<b>DE</b> Gebrauchsanweisung .....	2	<b>CS</b> Návod k použití.....	49
<b>EN</b> Instructions for use .....	5	<b>RO</b> Instrucțiuni de utilizare.....	52
<b>FR</b> Instructions d'utilisation.....	9	<b>HR</b> Upute za uporabu .....	56
<b>IT</b> Istruzioni per l'uso .....	13	<b>SL</b> Navodila za uporabo .....	59
<b>ES</b> Instrucciones de uso .....	16	<b>SK</b> Návod na používanie .....	62
<b>PT</b> Manual de utilização.....	20	<b>BG</b> Инструкция за употреба .....	66
<b>NL</b> Gebruiksaanwijzing.....	24	<b>TR</b> Kullanma talimatı .....	70
<b>SV</b> Bruksanvisning .....	27	<b>EL</b> Οδηγίες χρήσης .....	73
<b>DA</b> Brugsanvisning .....	31	<b>RU</b> Руководство по применению .....	77
<b>NO</b> Bruksanvisning .....	34	<b>JA</b> 取扱説明書.....	81
<b>FI</b> Käyttöohje.....	38	<b>ZH</b> 使用说明书.....	85
<b>PL</b> Instrukcja użytkownika.....	41	<b>KO</b> 사용 설명서.....	88
<b>HU</b> Használati utasítás .....	45		

## 1 Produktbeschreibung

Deutsch

### INFORMATION

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#### **INFORMATION**

Date of last update: 2022-03-03

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

#### **INFORMATION**

- ▶ The patient must be fitted with the liner no later than the date printed on the inside of the packaging or on the label (MHD/BBDD).

### **1.1 Function**

The liners for postoperative compression therapy after an amputation are intended to reduce oedema, alleviate pain, shape the residual limb and speed up the rehabilitation process. This is intended to enable the patient to wear a prosthesis as soon as possible. Compression decreases from distal to proximal due to the liner shape and the decreasing wall thickness.

These instructions for use apply to the following liners:

- 6Y44, 6Y45

## 1.2 Combination possibilities

**6Y44:** The product is not suitable for use in a prosthesis.

**6Y45:** The product may be used in an interim prosthesis.

## 2 Intended use

### 2.1 Indications for use

The product is intended for compression therapy after an amputation, to prepare for wearing a prosthesis.

The product must **not** be used for compression therapy for partial hand or partial foot amputations.

### 2.2 Contraindications

- Peripheral arterial occlusive disease (Doppler index less than 0.6 in non-amputated area)
- Decompensated heart failure
- Septic phlebitis
- Phlegmasia cerulea dolens
- Weeping dermatosis
- Severe sensory disturbances
- Intolerance and allergic reactions to the materials used

### 2.3 Lifetime

The product is a wear part, which means it is subject to normal wear and tear.

The maximum permissible number of reprocessing cycles limits the product lifetime (see page 8).

## 3 Safety

### 3.1 Explanation of warning symbols

**WARNING** Warning regarding possible serious risks of accident or injury.

**CAUTION** Warning regarding possible risks of accident or injury.

**NOTICE** Warning regarding possible technical damage.

## 3.2 General safety instructions

**WARNING!**

### Risk of severe injuries

- ▶ Store the product out of reach of children. There is a risk of suffocation if the product is pulled over the mouth and nose.

**CAUTION!**

### Risk of injury and risk of product damage

- ▶ Keep the product away from open flames, embers and other sources of heat.
- ▶ Make sure the product is applied properly and fits correctly. Improper or excessively tight application can cause pressure points and constriction of blood vessels and nerves.
- ▶ Do not continue to use the product if it does not fit properly.

### Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ Reprocess the product before using it on another person (see page 8).
- ▶ Observe the cleaning instructions.

## 4 Scope of delivery

Quantity	Designation	Reference number
1	Instructions for use	–
1	Liner	–

## 5 Information for use

Compression therapy usually begins 5 to 7 days after amputation. The timing is determined by the attending doctor. Compression in the early phase is important because the greatest change in the residual limb circumference occurs in the initial weeks after amputation.

### Course of compression therapy

The course of compression therapy is determined by the attending doctor. Initially, the bandaged residual limb is measured and the appropriate liner is selected. Ottobock recommends wearing the liner twice daily. Start with a short time and extend this daily.

#### Proposed wearing time (respectively morning and afternoon)

Day 1: 1 hour

Day 3: 3 hours

Day 2: 2 hours

From day 4 on: 4 hours

### Changing the liner size

As the oedema decreases, the liner size has to be adjusted to ensure continuous compression. The residual limb should therefore be measured regularly and the values noted. This allows the progress of therapy to be tracked and the timing for the next liner size to be determined.

## 6 Use

#### NOTICE

### Mechanical overload

Impaired functionality due to mechanical damage

- ▶ Check the product for damage prior to each use.
- ▶ Do not use the product if its functionality has been impaired.
- ▶ Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

### 6.1 Wound dressing

The wound on the residual limb must not be left exposed. The wound area has to be fully covered.

- ▶ Cover the entire wound area with an unmovable dressing (occlusive dressing, for example).

### 6.2 Size selection

The compression effect of the liner depends on the chosen size. Therefore, the liner size should correspond as closely as possible to the measured residual limb circumference. A shorter wearing time can make it easier for the patient to get used to the liner if compression is uncomfortable. A larger liner should be chosen only if it is not possible to wear the appropriate liner.

- 1) Measure the residual limb circumference **4 cm** proximal of the residual limb end.
- 2) Select the liner size that corresponds to the resulting measurement.  
**If the measurement is between two sizes:** Select the smaller size.
- 3) Check the fit of the liner. With a larger liner, ensure that the liner stretches and thereby generates compression.

### 6.3 Optional: Adapting the liner

The liner may be shortened to adapt to the patient's individual needs. This may limit the ability to reuse the liner on other patients.

- > **Recommended materials:** Conventional straight edge scissors (no serrated scissors)
- ▶ **INFORMATION: For a transfemoral interim prosthesis, the end of the liner must always be outside the prosthetic socket.**  
Trim the liner to the required length and contour. Make sure the cut edges are smooth so the material does not tear.

### 6.4 Putting on and taking off the liner

#### Donning

- 1) Check the liner for damage.
- 2) Turn the liner inside-out and place it on the end of the residual limb.

- 3) Unroll the liner onto the residual limb without displacing the soft tissue and without leaving wrinkles or air pockets.
- 4) Check the fit and the alignment of the liner.

### Doffing

- 1) Roll up the liner and remove it from the residual limb.
- 2) **Recommendation:** Clean the liner (see page 8).

## 6.5 Cleaning and Care

#### NOTICE

#### Use of improper cleaning agents

Damage to the product due to use of improper cleaning agents

- ▶ Only clean the product with the approved cleaning agents.

#### INFORMATION

Wearing a liner requires daily intensive skin care. Ottobock recommends the Derma skin care products.

- 1) Turn the liner inside out.
- 2) Hand wash the liner in warm water (**30 °C**) using neutral detergent (e.g. 453H10=1 Ottobock Derma Clean).
- 3) Thoroughly rinse the liner with clear water until all detergent residues have been removed.
- 4) Turn the liner again so the outer side is on the outside.
- 5) Put a towel into the liner and allow the liner to dry lying flat on a drying rack. Alternatively, use a drying stand (TT: 5Y100=1, TF: 5Y100=2).

## 7 Reprocessing

#### INFORMATION

**Provide the instructions for sterilisation of the product to the person performing this task.**

The liner must be reprocessed before using it on another patient. The liner is sterilised to eliminate all pathogens from the previous patient.

Subsequent non-sterile storage of the liner is permissible. The process is defined in the DIN EN ISO 17665-1:2006-11 standard.

### The liner may be reprocessed six times.

> **Prerequisite:** The liner has been cleaned.

- 1) Fill the liner with paper so it stays open during the process.
- 2) Put the liner into a bag and seal the bag.
- 3) Prepare and label the sterilisation test set.
- 4) **INFORMATION: Ottobock does not recommend sterilisation with ethylene oxide.**  
Perform sterilisation in the autoclave. The following procedure has been tested:
  - **121 °C, 15 psi, 30 minutes**
  - **134 °C, fractionated vacuum, moist heat, min. 5 minutes**
- 5) Store the liner according to hospital standards.

## 8 Disposal

In some jurisdictions it is not permissible to dispose of the product with unsorted household waste. Improper disposal can be harmful to health and the environment. Observe the information provided by the responsible authorities in your country regarding return, collection and disposal procedures.

## 9 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

### 9.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

## 9.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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