

28U90

DE Gebrauchsanweisung	4
EN Instructions for use	8
FR Instructions d'utilisation	11
IT Istruzioni per l'uso	16
ES Instrucciones de uso	20
PT Manual de utilização	24
NL Gebruiksaanwijzing	28
SV Bruksanvisning	32
DA Brugsanvisning	36
FI Käyttöohje	39
PL Instrukcja użytkowania	43
HU Használati utasítás	47
CS Návod k použití	51
RO Instrucțiuni de utilizare	55
BG Инструкция за употреба	59
EL Οδηγίες χρήσης	64
RU Руководство по применению	68
ZH 使用说明书	72

	Artikelnummer / Article Number	Seite / Side	Schuhgröße/ Shoe Size	Farbe/ Colour	Farbe Polster/ Colour Pad
	28U90=L35-37-0	links/left	L35-37	Transparent/transparent	Beige/beige
	28U90=L35-37-7	links/left	L35-37	Schwarz/black	Schwarz/black
	28U90=R35-37-0	rechts/right	R35-37	Transparent/transparent	Beige/beige
	28U90=R35-37-7	rechts/right	R35-37	Schwarz/black	Schwarz/black
	28U90=L37-39-0	links/left	L37-39	Transparent/transparent	Beige/beige
	28U90=L37-39-7	links/left	L37-39	Schwarz/black	Schwarz/black
	28U90=R37-39-0	rechts/right	R37-39	Transparent/transparent	Beige/beige
	28U90=R37-39-7	rechts/right	R37-39	Schwarz/black	Schwarz/black
	28U90=L39-41-0	links/left	L39-41	Transparent/transparent	Beige/beige
	28U90=L39-41-7	links/left	L39-41	Schwarz/black	Schwarz/black
	28U90=R39-41-0	rechts/right	R39-41	Transparent/transparent	Beige/beige
	28U90=R39-41-7	rechts/right	R39-41	Schwarz/black	Schwarz/black
	28U90=L41-44-0	links/left	L41-44	Transparent/transparent	Beige/beige
	28U90=L41-44-7	links/left	L41-44	Schwarz/black	Schwarz/black
	28U90=R41-44-0	rechts/right	R41-44	Transparent/transparent	Beige/beige
	28U90=R41-44-7	rechts/right	R41-44	Schwarz/black	Schwarz/black
	28U90=L44-47-0	links/left	L44-47	Transparent/transparent	Beige/beige
	28U90=L44-47-7	links/left	L44-47	Schwarz/black	Schwarz/black
	28U90=R44-47-0	rechts/right	R44-47	Transparent/transparent	Beige/beige
28U90=R44-47-7	rechts/right	R44-47	Schwarz/black	Schwarz/black	

Material	PA, PP
-----------------	--------

INFORMATION

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INFORMATION

Date of last update: 2022-05-18

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

These instructions for use contain important information for fitting and setting up your 28U90 ankle foot orthosis.

2 Intended use

2.1 Indications for use

The orthosis is intended **exclusively** for orthotic fittings of the lower limbs and **exclusively** for contact with intact skin.

The orthosis must be used in accordance with the indications.

2.2 Indications

- Dorsal flexor weakness, (e.g. after stroke, traumatic brain injury, in multiple sclerosis, neuromuscular atrophy, or peroneal paralysis)

Indications must be determined by the physician.

2.3 Contraindications

2.3.1 Absolute Contraindications

- Moderate to severe spasticity in the lower leg
- Moderate to severe oedema
- Leg ulcers
- Moderate to severe foot deformities

2.3.2 Relative Contraindications

The following indications require consultation with a physician: skin diseases/injuries, inflammation, prominent scars that are swollen, reddening and hyperthermia of the fitted limb/body area; lymphatic flow disorders, including unclear soft tissue swelling distal to the body area to which the medical device will be applied; sensory and circulatory disorders in the lower leg and the foot area.

2.4 Restrictions for use

The orthosis must not be used for the prevention of contractures (shortening) etc.

The orthosis must not be used for sports activities that involve jumping, sudden movements or rapid step sequences (e.g. basketball, badminton, riding sports). Sports activities should generally be discussed with the patient.

2.5 Lifetime

The product is designed for a maximum lifetime of **6 months**.

2.6 Effects

The ankle foot orthosis provides lifting of the foot and limits plantar flexion.

3 Safety

3.1 Explanation of warning symbols

⚠ WARNING Warning regarding possible serious risks of accident or injury.

⚠ CAUTION Warning regarding possible risks of accident or injury.

NOTICE Warning regarding possible technical damage.

3.2 General safety instructions

⚠ CAUTION

Reuse on other persons and improper cleaning

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ The product may be used by one person only.
- ▶ Clean the product regularly.

⚠ CAUTION

Contact with heat, embers or fire

Risk of injury (such as burns) and risk of product damage.

- ▶ Keep the product away from open flames, embers and other sources of heat.

⚠ CAUTION

Mechanical damage to the product

Injuries due to changes in or loss of functionality

- ▶ Use caution when working with the product.
- ▶ Check the product for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, discontinue use of the product and have it checked by authorised, qualified personnel.

NOTICE

Contact with oils, salves, lotions or other products that contain oils or acids

Insufficient stabilisation due to loss of material functionality

- ▶ Do not expose the product to oils, salves, lotions or other products that contain oils or acids.

4 Use

⚠ WARNING

Operation of motor vehicles

Risk of accident due to limited body function

- ▶ The product has not been approved by the manufacturer for the operation of motor vehicles.
- ▶ Observe the applicable legal regulations for the operation of motor vehicles and have your driving ability examined by an authorised test centre.

INFORMATION

- ▶ The daily duration of use and period of application are generally determined by the physician.
- ▶ The initial fitting and application of the product must be carried out by qualified personnel.
- ▶ Instruct the patient in the handling and care of the product.

- ▶ Instruct the patient to see a physician immediately if any exceptional changes are noted (e.g. worsening of the complaint).

4.1 Size Selection

- ▶ Select the orthosis size based on the shoe size (see sizing table).

4.2 Adaptation

⚠ CAUTION

Improper moulding or application

Injuries or damage to the orthosis due to overloading of the material and improper fit of the orthosis due to breakage of load-bearing components

- ▶ Do not make any improper changes to the orthosis.
- ▶ Always apply the orthosis according to the information in the instructions.

⚠ CAUTION

Incorrect or excessively tight application

Risk of local pressure and constriction of blood vessels and nerves due to improper or excessively tight application

- ▶ Ensure that the product is applied properly and fits correctly.

NOTICE

Use without suitable footwear

Restriction of functionality due to insufficient stabilisation

- ▶ Only use the orthosis with closed footwear, taking the allowable effective heel height into account.

Adapting the orthosis size:

- 1) **If necessary:** Cut or sand down the orthosis to size according to shoe size.
- 2) Then, deburr the worked edges and polish out any notches.

Adapting the orthosis:

- 1) **If necessary:** Thermoplastically re-form the orthosis by applying gentle heat to a limited extent.
NOTICE! Re-forming of the orthosis in cold conditions or at temperatures above 80°C is not permitted.
- 2) Attach the calf pad and the hook-and-loop strap (see fig. 2).
INFORMATION: The hook-and-loop strap can be fastened medially or laterally.
- 3) Once fitting is complete, reattach the hook-and-loop strap to the calf pad (see fig. 3).

Final inspection:

- ▶ Conduct trial walking with the patient. Testing on inclines, ramps and stairs is mandatory. Fine tune as required.

4.3 Application

NOTICE

Mechanical overload

Impaired functionality due to mechanical damage

- ▶ Check the product for damage prior to each use.
- ▶ Do not use the product if its functionality has been impaired.
- ▶ Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

1. Open the hook-and-loop strap.
2. Insert the orthosis into the shoe.
3. Put the orthosis on together with the shoe.
4. Tighten and close the hook-and-loop strap.

5 Cleaning

NOTICE

Use of improper cleaning agents

Damage to the product due to use of improper cleaning agents

- ▶ Only clean the product with the approved cleaning agents.

Clean the orthosis regularly:

- ▶ **Allowable cleaning agent:** pH neutral soap (e.g. 453H10 Derma Clean)

Orthosis frame

- 1) Wipe with a damp cloth as needed.
- 2) Allow to air dry. Do not expose to direct heat sources (e.g. sunlight, stove or radiator).

Textile pads

- 1) Remove the textile pads from the orthosis.
- 2) Fasten all hook-and-loop closures.
- 3) Hand wash the textile pads in warm water at **30 °C** with standard mild detergent. Rinse thoroughly.
- 4) Allow to air dry. Do not expose to direct heat sources (e.g. sunlight, stove or radiator).

6 Disposal

Dispose of the product in accordance with national regulations.

7 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

7.1 Local Legal Information

Legal information that applies **exclusively** to specific countries is written in the official language of the respective country of use after this chapter.

7.2 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

7.3 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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