
## 5R1=1, 5R1=2, 5R1=6, 5R1=6-H

<b>DE</b> Gebrauchsanweisung .....	3	<b>CS</b> Návod k použití.....	24
<b>EN</b> Instructions for use .....	4	<b>RO</b> Instrucțiuni de utilizare.....	26
<b>FR</b> Instructions d'utilisation.....	6	<b>HR</b> Upute za uporabu .....	28
<b>IT</b> Istruzioni per l'uso .....	8	<b>SL</b> Navodila za uporabo .....	29
<b>ES</b> Instrucciones de uso.....	9	<b>SK</b> Návod na používanie .....	31
<b>PT</b> Manual de utilização.....	11	<b>BG</b> Инструкция за употреба .....	32
<b>NL</b> Gebruiksaanwijzing .....	13	<b>TR</b> Kullanma talimatı .....	34
<b>SV</b> Bruksanvisning .....	15	<b>EL</b> Οδηγίες χρήσης .....	36
<b>DA</b> Brugsanvisning .....	16	<b>RU</b> Руководство по применению .....	38
<b>NO</b> Bruksanvisning .....	18	<b>JA</b> 取扱説明書.....	39
<b>FI</b> Käyttöohje.....	19	<b>ZH</b> 使用说明书.....	41
<b>PL</b> Instrukcja użytkowania.....	21	<b>KO</b> 사용 설명서.....	42
<b>HU</b> Használati utasítás .....	23		

# 1 Produktbeschreibung

Deutsch

## INFORMATION

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#### **INFORMATION**

Date of last update: 2022-05-30

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

#### 1.1 Construction and Function

The 5R1=1, 5R1=2, 5R1=6 and 5R1=6-H socket attachment blocks are bonded to the distal end of the prosthetic socket and then laminated together with it. They serve to connect the prosthetic socket to a socket adapter.

#### 1.2 Combination possibilities

This prosthetic component is compatible with Ottobock's system of modular connectors. Functionality with components of other manufacturers that have compatible modular connectors has not been tested.

### 2 Intended use

#### 2.1 Indications for use

The product is intended exclusively for lower limb exoprosthetic fittings.

#### 2.2 Area of application

- **5R1=1, 5R1=2, 5R1=6:** Approved for a body weight of up to **125 kg**.
- **5R1=6-H:** Approved for a body weight of up to **150 kg**.

#### 2.3 Environmental conditions

##### **Storage and transport**

Temperature range -20 °C to +60 °C (-4 °F to +140 °F), relative humidity 20 % to 90 %, no mechanical vibrations or impacts

##### **Allowable environmental conditions**

**Temperature range:** -10 °C to +45 °C (14 °F to 113 °F)

**Moisture:** relative humidity: 20% to 90%, non-condensing

##### **Unacceptable environmental conditions**

**Chemicals/liquids:** fresh water, salt water, perspiration, urine, acids, soapsuds, chlorine water

**Solids:** dust, sand, highly hygroscopic particles (e. g. talcum)

## 2.4 Lifetime

This product was tested by the manufacturer with 3 million load cycles. Depending on the user's activity level, this corresponds to a maximum lifetime of 5 years.

## 3 General safety instructions

### CAUTION!

### Risk of injury and risk of product damage

- ▶ Comply with the product's field of application and do not expose it to excessive strain (see page 4).
- ▶ Note the combination possibilities/combinations exclusions in the instructions for use of the products.
- ▶ Do not expose the product to prohibited environmental conditions.
- ▶ Check the product for damage if it has been exposed to prohibited environmental conditions.
- ▶ Do not use the product if it is damaged or in a questionable condition. Take suitable measures (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop).
- ▶ Observe the maximum lifetime of the product.
- ▶ To prevent mechanical damage, use caution when working with the product.
- ▶ If you suspect the product is damaged, check it for proper function and readiness for use.
- ▶ Do not use the product if its functionality is restricted. Take suitable measures (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop).

### Signs of changes in or loss of functionality during use

Among other factors, changes in functionality can be indicated by an altered gait pattern, a change in the positioning of the prosthetic components relative to each other and by the development of noises.

## 4 Scope of delivery

see fig. 1, item	Quantity	Designation	Reference number
–	1	Instructions for use	–
1	1	Socket attachment block	5R1=1
2			5R1=2
3			5R1=6-H
6			5R1=6
4	1	Lamination dummy	4X6
5	4	Countersunk head screw	501S41=M6x-25

## 5 Preparing the product for use

### CAUTION!

#### Incorrect alignment or assembly

Risk of injury due to damaged prosthetic components

- ▶ Observe the alignment and assembly instructions.

### INFORMATION

The layout described in this document was approved for the maximum product user body weight. The prosthetist assumes full responsibility for any change to the layout.

### 5.1 Fabricating the prosthetic socket

#### First lamination process

- > **Required materials:** 99B81 PVA bag, 616G6 Dacron® felt, 623T3 Perlon stockinette, 616G12 carbon fibre cloth, 617H119 Orthocryl lamination resin 80:20 PRO
- 1) Pull a soaked PVA bag over the model.
- 2) Place a layer of Dacron® felt over the entire model.
- 3) Place a layer of Dacron® felt around the proximal edge of the model.
- 4) Pull 2 layers of Perlon stockinette over the model.
- 5) Place 3 layers of carbon fibre cloth (e.g. 15 cm x 15 cm) on the distal end of the model, offsetting the fibre alignment.
- 6) Pull 2 layers of Perlon stockinette over the model.
- 7) Pull a soaked PVA bag over the model.
- 8) Complete the lamination process with lamination resin.

- 9) Once the lamination resin has cured, remove the PVA bag.

#### Foaming and fastening the socket attachment block

- > **Required materials:** 617H32 Pedilen rigid foam 300, 617H21 Orthocryl sealing resin
- 1) Sand down the surplus lamination resin on the distal end of the model.
- 2) Foam the distal end of the model with rigid foam.
- 3) Trim the rigid foam to the maximum possible extent and sand it smooth in the desired position.
- 4) Fasten the socket attachment block to the rigid foam with sealing resin.
- 5) Mill a transition between the socket attachment block and the rigid foam (see fig. 2). While doing this, do not sand the shape or size of the adapter connection surface, which is determined by the socket attachment block.
- 6) Insert the lamination dummy in the threaded holes of the socket attachment block (see fig. 3).
- 7) Apply sealing resin to the outside of the prosthetic socket and to the lamination dummy (see fig. 4).

#### Second lamination process

- > **Required materials:** 99B81 PVA bag, 623T3 Perlon stockinette, 616B1=50 carbon fibre cloth strap, 616G12 carbon fibre cloth, 616G15 woven carbon fibre stockinette, 617H119 Orthocryl lamination resin 80:20 PRO
- 1) Pull a Perlon stockinette over the model.
- 2) Apply 2 strips of carbon fibre cloth strap in a cross over the socket attachment block up to the middle of the model.
- 3) Place 2 layers of carbon fibre cloth (e.g. 15 cm x 15 cm) on the distal end of the model, offsetting the fibre alignment.
- 4) Pull a Perlon stockinette over the model.
- 5) Pull 2 layers of woven carbon fibre stockinette over the model. Do not allow the stockinette to protrude beyond the socket edge.
- 6) Pull 2 layers of Perlon stockinette over the model.
- 7) Pull a soaked PVA bag over the model.
- 8) Complete the lamination process with lamination resin.
- 9) Once the lamination resin has cured, remove the PVA bag.

## 5.2 Final assembly

### CAUTION!

#### Improper assembly of the screw connections

Risk of injury due to breakage or loosening of the screw connections

- ▶ Clean the threads before every installation.
- ▶ Apply the specified torque values.
- ▶ Follow the instructions regarding the length of the screws and about how to secure the screws.

- > **Required materials:** sanding tool, 710D20 torque wrench

- 1) **NOTICE! Do not sand down the laminate around the edge of the product.**  
Sand down the laminate on the distal end of the prosthetic socket to the surface of the lamination dummy and remove the lamination dummy.
- 2) **CAUTION! Do not use thread lock (e.g. Loctite®) to mount the socket adapter. Chemical thread locks damage the laminate.**  
Fasten the socket adapter to the socket attachment block using the countersunk head screws (see fig. 5):  
Screw in 2 countersunk head screws on the posterior side (12 Nm).  
Screw in 2 countersunk head screws on the anterior side (12 Nm).

## 6 Cleaning

- 1) Clean the product with a damp, soft cloth.
- 2) Dry the product with a soft cloth.
- 3) Allow to air dry in order to remove residual moisture.

## 7 Maintenance

- ▶ A visual inspection and functional test of the prosthetic components should be performed after the first 30 days of use.
- ▶ Inspect the entire prosthesis for wear during normal consultations.
- ▶ Conduct annual safety inspections.

## 8 Disposal

In some jurisdictions it is not permissible to dispose of the product with unsorted household waste. Improper disposal can be harmful to health and the environment. Observe the information provided by the responsible authorities in your country regarding return, collection and disposal procedures.

## 9 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

### 9.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

### 9.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

### 9.3 Warranty

The manufacturer warrants this device from the date of purchase. The warranty covers defects that can be proven to be a direct result of flaws in the material, production or construction and that are reported to the manufacturer within the warranty period.

Further information on the warranty terms and conditions can be obtained from the competent manufacturer distribution company.

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