4142 Manu 3D 4143 Manu 3D Basic 4146 Manu 3D Pollex

DE Gebrauchsanweisung	5
EN Instructions for use	9
FR Instructions d'utilisation	14
IT Istruzioni per l'uso	19
ES Instrucciones de uso	24
NL Gebruiksaanwijzing	29
SV Bruksanvisning	33
DA Brugsanvisning	38
NO Bruksanvisning	43
FI Käyttöohje	47
PL Instrukcja użytkowania	52
RU Руководство по применению	57
JA 取扱説明書	62
ZH 使用说明书	66

4142

4143

4146

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INFORMATION

Date of last update: 2022-07-12

- Please read this document carefully before using the product and observe the safety notices.
- Instruct the user in the safe use of the product.
- Please contact the manufacturer if you have questions about the product or in case of problems.

- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

These instructions for use provide important information on the fitting and application of the 4142 Manu 3D, 4143 Manu 3D Basic and 4146 Manu 3D Pollex wrist supports.

2 Intended use

2.1 Indications for use

The support is intended **exclusively** for treating the upper limbs and **exclusively** for contact with intact skin.

The support must be used according to the indication.

2.2 Indications

4142 Manu 3D

- Irritation in the wrist area
- Carpal tunnel syndrome
- Osteoarthritis of the wrist
- Subluxation/dislocation
- Tendovaginitis in the hand
- Sprains

4143 Manu 3D Basic

- Irritation in the wrist area
- Osteoarthritis of the wrist
- Subluxation/dislocation
- Tendovaginitis in the hand
- Sprains

4146 Manu 3D Pollex

- De Quervain's tenosynovitis
- Irritation in the hand region
- Irritation of the carpometacarpal and the metacarpophalangeal joint
- Collateral ligament injuries of the metacarpophalangeal joint (e.g. skier's thumb)
- Osteoarthritis of the wrist
- Subluxation/dislocation
- Tendovaginitis in the hand and thumb areas
- Sprains

The indication must be determined by the physician.

2.3 Contraindications

2.3.1 Absolute Contraindications

None known.

2.3.2 Relative Contraindications

The following indications require consultation with a physician: skin diseases/injuries, inflammation, prominent, swollen scars, reddening and hyperthermia in the relevant area; pronounced varicose veins, especially with impaired blood return, lymph drainage disorders and unclear soft tissue swelling distal to the device; sensory and circulatory disorders in the area of the forearm and hand, e.g. associated with diabetic neuropathy.

2.4 Lifetime

The product is designed for a maximum lifetime of **6 months**.

2.5 Mechanism of Action

The support stabilizes the wrist.

In addition for 4146 Manu 3D Pollex:

The support braces the CMC and MCP joint of the thumb.

3 Safety

3.1 Explanation of warning symbols

 **CAUTION** Warning regarding possible risks of accident or injury.

 **NOTICE** Warning regarding possible technical damage.

3.2 General safety instructions

 **CAUTION**

Contact with heat, embers or fire

Risk of injury (such as burns) and risk of product damage

- ▶ Keep the product away from open flames, embers and other sources of heat.

 **CAUTION**

Reuse on other persons and improper cleaning

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ The product may be used by one person only.
- ▶ Clean the product regularly.

NOTICE

Improper use and changes

Change in or loss of functionality as well as damage to the product

- ▶ Use this product with care and only for its intended purpose.
- ▶ Do not make any improper changes to the product.

NOTICE

Contact with oils, salves, lotions or other products that contain oils or acids

Insufficient stabilization due to loss of material functionality

- ▶ Do not expose the product to oils, salves, lotions or other products that contain oils or acids.

4 Handling

INFORMATION

- ▶ The daily duration of use and period of application are generally determined by the physician.
- ▶ The initial fitting and application of the product must be carried out by qualified personnel.
- ▶ Instruct the patient in the handling and care of the product.
- ▶ Instruct the patient to see a physician immediately if any exceptional changes are noted (e.g. worsening of the complaint).

4.1 Selecting the Size

- 1) Measure the wrist circumference.
- 2) Determine the size of the support (see size chart).

4.2 Fitting and Application

CAUTION

Incorrect or excessively tight application

Skin irritation due to overheating, local pressure points due to tight application to the body

- ▶ Ensure that the product is applied properly and fits correctly.

- ▶ When in doubt, do not continue to use the product if skin irritation occurs.

NOTICE

Use of a worn or damaged product

Limited effectiveness

- ▶ Before each use, check the product for functional reliability and for possible wear or damage.
- ▶ Do not continue using a product that is no longer functional, or that is worn or damaged.

- > Open the hook-and-loop closures in the area of the wrist joint.
- 1) Pull the support over the thumb (see fig. 1).
- 2) Fasten all hook-and-loop closures in the area of the wrist joint (see fig. 2).
- 3) **Only 4143 and 4146:** Readjust the hook-and-loop closure in the area of the thumb (see fig. 3).

5 Cleaning

NOTICE

Use of improper cleaning agents

Damage to the product due to use of improper cleaning agents

- ▶ Only clean the product with the approved cleaning agents.

Clean the support regularly:

- 1) Take out all removable splints.
- 2) Fasten all hook-and-loop closures.
- 3) Recommendation: Use a laundry bag or net.
- 4) Wash the support in warm water at **40 °C** with standard mild detergent. Do not use fabric softener. Rinse thoroughly.
- 5) Allow to air dry. Do not expose to direct heat sources (e.g. sunlight, stove or radiator).
- 6) Reinsert all splints.

6 Disposal

Dispose of the product in accordance with national regulations.

7 Legal Information

All legal conditions are subject to the respective national law of the country of use and may vary accordingly.

7.1 Liability

The manufacturer shall be liable in the event that the product is used in accordance with the descriptions and instructions in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorized modification of the product.

7.2 CE Conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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