

## 28F10

<b>DE</b> Gebrauchsanweisung .....	7
<b>EN</b> Instructions for use .....	12
<b>FR</b> Instructions d'utilisation .....	18
<b>IT</b> Istruzioni per l'uso .....	24
<b>ES</b> Instrucciones de uso .....	30
<b>PT</b> Manual de utilização .....	35
<b>NL</b> Gebruiksaanwijzing .....	41
<b>SV</b> Bruksanvisning .....	47
<b>DA</b> Brugsanvisning .....	52
<b>PL</b> Instrukcja użytkowania .....	58
<b>HU</b> Használati utasítás .....	64
<b>CS</b> Návod k použití .....	69
<b>SK</b> Návod na používanie .....	75
<b>TR</b> Kullanma talimatı .....	80
<b>EL</b> Οδηγίες χρήσης .....	86
<b>RU</b> Руководство по применению .....	92
<b>ZH</b> 使用说明书 .....	98


24

**Material**

PP, PA, PUR, PU, NBR/SBR, EVA, PES, PVC

**INFORMATION**

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### **INFORMATION**

Date of last update: 2022-10-27

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

These instructions for use contain important information on the fitting and application of the 2BF10 heel relief orthosis designed by Dr Settner and master O&P professional Münch.

Ottobock recommends that you attend a practical seminar which provides training on the fitting of the orthosis. Detailed steps are shown in the heel relief orthosis video and can be found using the following link:  
<http://videoguides.ottobock.com/orthotics>

## 2 Product description

### 2.1 Components

Item	Component	Item	Component
1	Calf pad	9	Compensation sole
2	Adhesive strips (calf pad)	10	Self-adhesive hook-and-loop strap
3	Forefoot cap	11	Heel wedge
4	Hook-and-loop closure (instep)	12	Longitudinal arch piece (left/right)
5	Setting nuts and fastening screws (instep hook-and-loop closure)	13	Heel container
6	Fastening screws (forefoot cap)	14	Hook-and-loop closure (lower shin)
7	Compression pads	15	Hook-and-loop closure (upper shin)
8	Plantar pad (foam)		

### 2.2 Design

The heel relief orthosis is made of glass fibre-reinforced polypropylene.

The heel container (13) is designed to enclose the dorsal calf and support the longitudinal arch (12) of the foot, thus enabling suspended positioning of the calcaneus in the orthosis. Plantar foot support is provided by the moulded foot piece for the longitudinal arch and the compensation sole (9), which can be affixed to the heel container. The heel wedge (11) can be used to compensate for the height of the shoe on the unaffected side.

To facilitate check-ups and for improved air circulation, the orthosis has an opening at the back of the heel container. The weight-bearing pad set for the heel, which includes a plantar pad (8) and compression pads (7), is applied according to the therapy plan and allows the weight borne by the calcaneus to be built up gradually.

The hook-and-loop closures (4, 14 and 15) and the forefoot cap (3) hold the foot in place. The zip fastener on the forefoot cap makes it easy to apply and remove the orthosis and also protects the toes. The enclosed fastening screws and setting nuts (5) are used to affix the hook-and-loop closure (4) to the orthosis. The hook-and-loop closures at shin level (14 and 15) can also be riveted or sewn together to increase stability. The calf pad (1) can be cut to size and serves as padding within the heel container.

The underside of the orthosis features a non-skid profile for increased safety when walking.

## 3 Intended use

### 3.1 Indications for use

The orthosis is intended **exclusively** for orthotic fittings of the lower limbs and **exclusively** for contact with intact skin.

The orthosis must be used in accordance with the indications.

### 3.2 Indications

- Functional follow-up treatment for calcaneal fractures (unilateral or bilateral), regardless of the type of fracture and course of primary treatment.
- Weight-bearing relief following arthrodesis of the lower ankle

Indications must be determined by the physician.

### 3.3 Contraindications

#### 3.3.1 Absolute Contraindications

- External fixator (wire osteosynthesis)

#### 3.3.2 Relative Contraindications

The following indications require consultation with a physician: skin diseases/injuries, inflammation, prominent scars that are swollen, reddening and hyperthermia of the fitted limb/body area; lymphatic flow disorders, including unclear soft tissue swelling distal to the body area to which the medical device will be applied; sensory and circulatory disorders in the lower leg and the foot area.

#### 3.4 Lifetime

The product is designed for a maximum lifetime of **6 months**.

#### 3.5 Effects

The orthosis relieves the calcaneus by providing targeted support to the longitudinal arch beneath the os naviculare as well as a form-fitting enclosure for the metatarsals and the calf. Accelerated bone healing is achieved by applying the basic principle of the Ilizarov technique (compression and distraction of the fractured areas). An almost natural rollover of the injured foot with the orthosis exerts a positive influence on the neuromuscular structures and calcareous salt content of the bone. The functioning muscle pump may also contribute to thrombosis prophylaxis.

#### 3.6 Treatment and therapy plan

The individual therapy plan is determined by the physician and is independent of the type of fracture and treatment approach.

After the swelling subsides or the stitches are removed: fitting of the orthosis with insertion of the plantar pad (4 hrs.)

Days 8–12	Full weight-bearing without crutches
Week 4	1st x-ray evaluation (lateral view of upper ankle, special/lateral view of calcaneus, upper ankle under full weight bearing/Hallgrímsson imaging)
Week 6	Insert 1st compression pad
Week 8	2nd x-ray evaluation Insert 2nd compression pad
Week 10	Insert 3rd compression pad
Week 11	If necessary, take measurements for orthopaedic footwear or insoles (within 5–6 days)
Week 12	Conclusion of medical treatment / test for fitness for work for worker's compensation cases

(The time indications are from the BG Unfallklinik (hospital for accident cases) Duisburg, Germany, and are calculated from the date of the injury)

## 4 Safety

### 4.1 Explanation of warning symbols

 <b>CAUTION</b>	Warning regarding possible risks of accident or injury.
 <b>NOTICE</b>	Warning regarding possible technical damage.

## 4.2 General safety instructions

### ⚠ CAUTION

#### Reuse on other persons and improper cleaning

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ The product may be used by one person only.
- ▶ Clean the product regularly.

### ⚠ CAUTION

#### Contact with heat, embers or fire

Risk of injury (such as burns) and risk of product damage.

- ▶ Keep the product away from open flames, embers and other sources of heat.

### NOTICE

#### Contact with oils, salves, lotions or other products that contain oils or acids

Insufficient stabilisation due to loss of material functionality

- ▶ Do not expose the product to oils, salves, lotions or other products that contain oils or acids.

### NOTICE

#### Exposure of the product to unsuitable environmental conditions

Damage, brittleness or destruction due to improper handling

- ▶ Avoid storage in condensing ambient humidity.
- ▶ Avoid contact with abrasive substances (e.g. sand, dust).
- ▶ Do not expose the product to temperatures below  $-10\text{ °C}$  or above  $+50\text{ °C}$  (e.g. sauna, excessive sunlight, drying on a radiator).

## 5 Handling

### INFORMATION

- ▶ The daily duration of use and period of application are generally determined by the physician.
- ▶ The initial fitting and application of the product must be carried out by qualified personnel.
- ▶ Instruct the patient in the handling and care of the product.
- ▶ Instruct the patient to see a physician immediately if any exceptional changes are noted (e.g. worsening of the complaint).

### 5.1 Size Selection

- ▶ Select the orthosis size based on the shoe size (see sizing table).

Size	Shoe size	Size Heel container	Size Longitudinal arch piece	Size Forefoot cap
S	36–41	Container size S (36–41)	Longitudinal arch piece size S (36–41)	Forefoot cap size S (36–41)
M	42–46	Container size M–L (42–50)	Longitudinal arch piece size M (42–46)	Forefoot cap size M (42–46)
L	47–50	Container size M–L (42–50)	Longitudinal arch piece size L (47–50)	Forefoot cap size L (36–41)

## 5.2 Adaptation

- 1) Place the longitudinal arch piece (12) in the heel container (13).

**INFORMATION: To determine the placement, affix the longitudinal arch piece to a piece of paper. Position the patient on the longitudinal arch piece, align it with the navicular bone (os naviculare) and mark the contour of the foot using a pencil (see fig. 2). Then cut out the contour (see fig. 3) and place it in the heel container (see fig. 4) so the end of the heel meets the mark in the heel container.**

- 2) Mark the position of the longitudinal arch piece in the heel container (see fig. 5).
- 3) Fix the longitudinal arch piece along the mark inside the heel container.
- 4) Insert the plantar pad (8) in the heel container (see fig. 6).
- 5) Carry out an initial weight bearing test with the longitudinal arch piece (see fig. 7, see fig. 8) and the forefoot cap (3) in place (see fig. 9).

**INFORMATION: When the heel is under maximum weight, use the opening cut in the heel container to check the degree of relief provided by the orthosis.**

- 6) **Optional:** Adapt the heel container to the individual patient using milling (see fig. 11) or heat (see fig. 12), then widening it as necessary.
- 7) **Optional:** Fix the the compensation sole (9) with hook-and-loop strap (10) to further increase the height of the longitudinal arch.
- 8) Position the plantar pad (8) in the heel area so that it maintains continual contact with the heel.

**INFORMATION: During the course of therapy, compression pads (7) are inserted under the plantar pad according to the therapy plan and serve the gradual build-up of weight borne by the calcaneus.**

- 9) Mount the hook-and-loop closure and the guide loop (4) at the instep using the supplied setting nuts and screws (5) (see fig. 13). Do this by attaching the setting nut from the inside and applying the strap from the outside (see fig. 14).

**INFORMATION: The longitudinal arch piece must be removed temporarily before mounting the hook-and-loop closure.**

- 10) **Optional:** Cut the calf pad (1) to the required size (see fig. 15).
- 11) Attach the adhesive strips (2) (see fig. 16) and fix the calf pad in the heel container.
- 12) **Optional:** Cut the sides of the forefoot cap (see fig. 19). The two pieces can then be riveted or sewn together.
- 13) Push the forefoot cap over the heel container until the desired shoe size is set (scaling displayed on the side of the heel container).
- 14) Fix the forefoot cap to the heel container using screws (6) (see fig. 20).
- 15) Position the foot in the orthosis (see fig. 21).
- 16) Close all hook-and-loop closures one after the other, ensuring they exert even pressure (see fig. 22).
- 17) Close the zip on the forefoot cap (see fig. 23).
- 18) Carry out a weight bearing test with the fully assembled orthosis.

**INFORMATION: Use a strip of paper in the opening cut in the heel container to check the degree of relief provided by the orthosis (see fig. 24). To do so, place the strip of paper between the heel and the plantar pad. When the calcaneus is bearing the user's full weight, pull the paper strip out of the orthosis in one even stroke. There should be slight resistance as you pull the paper strip out of the orthosis, but it should not tear.**

## 5.3 Application

### CAUTION

#### **Direct skin contact with the product**

Skin irritation due to chafing or sweat formation

- ▶ Do not wear this product directly on the skin.

### **⚠ CAUTION**

#### **Improper moulding or application**

Injuries or damage to the orthosis due to overloading of the material and improper fit of the orthosis due to breakage of load-bearing components

- ▶ Do not make any improper changes to the orthosis.
- ▶ Always apply the orthosis according to the information in the instructions.

### **⚠ CAUTION**

#### **Incorrect or excessively tight application**

Risk of local pressure and constriction of blood vessels and nerves due to improper or excessively tight application

- ▶ Ensure that the product is applied properly and fits correctly.

### **NOTICE**

#### **Use of a worn or damaged product**

Limited effectiveness

- ▶ Before each use, check the product for functional reliability and for possible wear or damage.
- ▶ Do not continue using a product that is no longer functional, or that is worn or damaged.

- 1) Open the hook-and-loop closures and the zip fastener.
- 2) Place the foot in the orthosis (see fig. 21). The toes should be placed in the orthosis first.
- 3) Close the hook-and-loop closures on the instep and lower leg, ensuring they exert even pressure (see fig. 22).
- 4) Close the zip fastener (see fig. 23).

## **6 Cleaning**

### **NOTICE**

#### **Use of improper cleaning agents**

Damage to the product due to use of improper cleaning agents

- ▶ Only clean the product with the approved cleaning agents.

- ▶ **Allowable cleaning agent:** pH neutral soap (e.g. 453H10 Derma Clean)

#### **Orthosis frame**

- 1) Wipe with a damp cloth as needed.
- 2) Allow to air dry. Do not expose to direct heat sources (e.g. sunlight, stove or radiator).

#### **Straps and pads**

- 1) Remove the straps and pads from the orthosis.
- 2) Close all hook-and-loop closures.
- 3) Hand-wash the straps and pads in warm water at 30 °C (86 °F) using neutral detergent. Rinse thoroughly.
- 4) Allow to air dry. Do not expose to direct heat sources (e.g. sunlight, stove or radiator).

## **7 Disposal**

Dispose of the product in accordance with national regulations.

## **8 Legal information**

All conditions are subject to the respective national laws of the country of use and may vary accordingly.

## 8.1 Local Legal Information

Legal information that applies **exclusively** to specific countries is written in the official language of the respective country of use after this chapter.

## 8.2 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

## 8.3 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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