

170K1=80 Free Walk, 170K1=120 Free Walk










| | |
|--|-----|
| DE Gebrauchsanweisung (Fachpersonal) | 5 |
| EN Instructions for use (qualified personnel) | 14 |
| FR Instructions d'utilisation (Personnel spécialisé) | 23 |
| IT Istruzioni per l'uso (Personale tecnico specializzato) | 33 |
| ES Instrucciones de uso (Personal técnico especializado) | 43 |
| PT Manual de utilização (Pessoal técnico) | 53 |
| NL Gebruiksaanwijzing (Vakmensen) | 63 |
| SV Bruksanvisning (Fackpersonal) | 72 |
| DA Brugsanvisning (Faguddannet personale) | 81 |
| NO Bruksanvisning (Fagpersonell) | 90 |
| FI Käyttöohje (Ammattihenkilöstö) | 99 |
| HU Használati utasítás (szakszemélyzet) | 108 |
| TR Kullanma talimatı (Uzman personel) | 117 |


1

2

INFORMATION

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INFORMATION

Date of last update: 2022-10-27

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

These instructions for use provide you with important information on the processing of the 170K1=80 and 170K1=120 Free Walk orthotic system.

2 Product description

2.1 Available sizes

The orthosis is available in two sizes. The choice of size is based on the patient's body weight.

| Article number | Max. body weight [kg] |
|----------------|-----------------------|
| 170K1=80 | 80 |
| 170K1=120 | 120 |

2.2 Components

| Scope of delivery (see fig. 1) | | | | |
|--|----------------------|---|---------------------------------|----------------------------------|
| Item | Quantity [Piece (s)] | Component | 170K1=L/R-80* Article number | 170K1=L/R-120* Article number |
| ¹⁾ Colour: beige (0), black (7) | | | | |
| 1 | 4 | Y-hook-and-loop ¹⁾ | 170D29=38.* | |
| 2 | 2 | Hook-and-loop closure ¹⁾ | 170Z4=600.* | |
| 3 | 2 | Hook-and-loop closure ¹⁾ | 170Z4=400.* | |
| 4 | 4 | Strap guide loop with tube connection ¹⁾ | 170D18=120.* | |
| 5 | 1 | Support tube with pad | - | |
| 6 | 1 | Base tube | - | |
| 7 | 1 | Support tube with pad | - | |
| 8 | 1 | Joint protector ¹⁾ | 170Z120=* | |
| 9 | 1 | 3-phase switch ¹⁾ | 170D71=* | |
| 10 | 1 | Support tube with pad | - | |
| 11 | 1 | Base tube | - | |
| 12 | 1 | Support tube with pad | - | |
| 13 | 1 | Ankle joint, upper part | - | |
| 14 | 1 | Ankle joint, lower part | 170X17N=* | |
| 15 | 1 | Foot stirrup set | | |
| 16 | 4 | Pad button | 170D69 | |
| 17 | 4 | Strap guide loop ¹⁾ | 170D20=* | |
| 18 | 3 | Foam pads ¹⁾ | 170D100=1.* | |
| 19 | 1 | Foam pads ¹⁾ | 170D100=2.* | |
| 20 | 1 | Shrinkable tubing | 616Z9=19x61 | |

| Maintenance sets (see fig. 2) | | | | |
|-------------------------------|----------------------|--|---------------------------------|----------------------------------|
| Item | Quantity [Piece (s)] | Component | 170K1=L/R-80* Article number | 170K1=L/R-120* Article number |
| 1 | 1 | Maintenance set for lock | 170Z99=2 | |
| | 1 | Steel pull cable from maintenance set for lock | 170Z96 | |
| 2 | 1 | Maintenance set for Free Walk | 170X18 | |
| 3 | 1 | Foot stirrup set | 170X17N=* | |

| 170X17N=* foot stirrup set (see fig. 2, maintenance set 3) | | | | |
|---|-----------------------------|--------------------------|-------------------------------------|--------------------------------------|
| Item | Quantity [Piece (s)] | Component | 170K1=L/R-80* Article number | 170K1=L/R-120* Article number |
| 1 | 1 | Threaded plate | | 170D28 |
| 2 | 1 | Foot stirrup, upper part | | 170F14 |
| 3 | 1 | Foot stirrup, lower part | | 170F15 |
| 4 | 2 | Setting nuts | | 170D33 |
| 5 | 2 | Screws | | 501F9=M4x4 |
| 6 | 1 | Ankle joint, lower part | | 170C11=L/R |
| 7 | 1 | Screw | | 501F9=M6x10-1 |

2.3 Design

The Ottobock FreeWalk orthoses include an orthotic knee joint with lock. Three modes can be selected on the orthotic knee joint:

- **FreeWalk mode:** The orthotic knee joint locks in the stance phase and unlocks in the swing phase.
- **Permanent lock:** The orthotic knee joint is permanently locked.
- **Free motion mode:** The orthotic knee joint is permanently unlocked.

3 Intended use

3.1 Indications for use

The orthosis is intended **exclusively** for orthotic fittings of the lower limbs and **exclusively** for contact with intact skin.

The orthosis must be used in accordance with the indications.

3.2 Indications

- Flaccid paresis or paralysis of the lower limb
 - Adequate muscle strength of the hip or knee extensors or hyperextension in the knee joint
- Indications must be determined by the physician.

3.3 Contraindications

3.3.1 Absolute Contraindications

- Ankle movement < 10°
- Knee flexion contractures > 10°
- Unstable knee joint with compensating movement in varus position
- Valgus deformity > 15°

3.3.2 Relative Contraindications

- Valgus deformity 0° to 15° (fitting only with 170D30 medial support)
- Spasticity resulting in involuntary movements in the ankle joint

3.4 Lifetime

The product is designed for a lifetime of **3 years** when used as intended and assembled professionally.

3.5 Qualification

Patients may be fitted with the product only by trained qualified personnel. The qualified personnel must be familiar with the handling of the various techniques, materials, machines and tools.

4 Safety

4.1 Explanation of warning symbols

⚠ WARNING Warning regarding possible serious risks of accident or injury.

⚠ CAUTION Warning regarding possible risks of accident or injury.

NOTICE Warning regarding possible technical damage.

4.2 General safety instructions

⚠ CAUTION

Mechanical damage to the product

Injuries due to changes in or loss of functionality

- ▶ Use caution when working with the product.
- ▶ Check the product for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, discontinue use of the product and have it checked by authorised, qualified personnel.

⚠ CAUTION

Contact with heat, embers or fire

Risk of injury (such as burns) and risk of product damage.

- ▶ Keep the product away from open flames, embers and other sources of heat.

NOTICE

Contact with oils, salves, lotions or other products that contain oils or acids

Insufficient stabilisation due to loss of material functionality

- ▶ Do not expose the product to oils, salves, lotions or other products that contain oils or acids.

NOTICE

Contact with saltwater or water containing chlorine/soap, fresh water or abrasive substances (e.g. sand)

Damage and premature product wear

- ▶ Following any contact with the substances identified above, promptly clean the product in accordance with the section "Cleaning".

NOTICE

Exposure of the product to unsuitable environmental conditions

Damage, brittleness or destruction due to improper handling

- ▶ Avoid storage in condensing ambient humidity.
- ▶ Avoid contact with abrasive substances (e.g. sand, dust).
- ▶ Do not expose the product to temperatures below -10 °C (14 °F) or above +60 °C (140 °F) (e.g. sauna, excessive sunlight, drying on a radiator).

5 Preparing the product for use

CAUTION

Improper moulding or application

Injuries or damage to the orthosis due to overloading of the material and improper fit of the orthosis due to breakage of load-bearing components

- ▶ Do not make any improper changes to the orthosis.
- ▶ Always apply the orthosis according to the information in the instructions.

CAUTION

Incorrect alignment, assembly or adjustment

Risk of injury due to change in or loss of functionality

- ▶ Assembly, adjustment and maintenance operations may only be completed by qualified personnel.
- ▶ Observe the alignment, assembly and adjustment instructions.

NOTICE

Thermal overloading of product components

Loss of function due to improper thermal treatment

- ▶ Prior to the thermal treatment of the product, remove all temperature-critical components (e.g. plastic parts).

5.1 Assembly

Adapting the upper part of the foot stirrup

> **Required tools:** Felt tip pen, hacksaw, file, sandpaper

- 1) Shape the insole or inner sandal.
- 2) Position the upper part of the foot stirrup together with the lower part of the foot stirrup and shaped insole on the patient's foot. In doing so, the Torx screw of the orthotic ankle joint is at the height of the lateral malleolus (see fig. 3).
- 3) Mark the height of the lower part of the foot stirrup on the upper part of the foot stirrup using a felt tip pen (see fig. 4).
- 4) Shorten the upper part of the foot stirrup between two bores so that **at least two bores and no more than three bores remain below the marking** (see fig. 5).
- 5) Round and deburr the cut edge on the upper part of the foot stirrup.
- 6) **Optional:** To avoid damaging the shoe, shorten the protruding section of the upper part of the foot stirrup to two bores if needed.

Mounting the foot stirrup

> **Required tools:** Torx T20 screwdriver

- 1) Insert the ankle joint setting nuts into the slotted hole from the inside of the lower part of the foot stirrup, noting the rotation lock (see fig. 6).
- 2) **INFORMATION: Position the scale of the upper part of the foot stirrup laterally.** Position the upper part of the foot stirrup on the outside of the lower part of the foot stirrup and secure it using the two ankle joint setting nuts (see fig. 7).
- 3) Insert the screws and tighten them (see fig. 8).

Assembling the orthosis

> **Required tools:** Torx T20, T30 screwdrivers.

- 1) Loosen the joint screw of the lower part of the ankle joint and remove the threaded plate.
- 2) Place the foot stirrup onto the lower part of the ankle joint (see fig. 9).
- 3) Reinsert the threaded plate (see fig. 10).
- 4) Install the foot stirrup on the lower part of the ankle joint with the joint screw (see fig. 11).

- 5) Put the orthosis on together with the shoe and insert.
- 6) Check that the length of the foot stirrup is correct. If necessary, adjust the foot stirrup using the bores in the upper part of the foot stirrup and the slotted hole in the lower part of the foot stirrup so that the joint screw of the orthotic ankle joint is at the height of the lateral malleolus.
- 7) Tighten the screws to secure the foot stirrup (see fig. 8).
- 8) Mark the position of the foot stirrup in the a-p direction and the exterior rotation on the insert.
- 9) Drill the holes for riveting in the insert.
- 10) Rivet the insert to the foot stirrup.

5.2 Trial fitting

After assembly, the orthosis can be fitted and the fit and function can be checked on the patient.

5.2.1 Functional test

To verify technically correct functionality, take the foot section in one hand and the femoral section in the other:

Free Walk mode

- > The 3-phase switch is in the upper position.
 - ▶ Push the orthotic ankle joint in the dorsiflexion direction and extend the orthotic knee joint.
 - The lock disengages automatically and the orthotic knee joint can be flexed.

Permanent locking

- > The orthotic knee joint is set to Free Walk mode.
 - 1) Slide the 3-phase switch to the lower position with the knee extended.
 - 2) Push the orthotic ankle joint in the dorsiflexion direction and extend the orthotic knee joint.
 - The orthotic knee joint is permanently locked.

Free motion mode

- > The orthotic knee joint is set to permanently locked mode.
 - 1) Slide the 3-phase switch to the upper position.
 - 2) Push the orthotic ankle joint in the dorsiflexion direction and extend the orthotic knee joint.
 - 3) Flex the orthotic knee joint.
 - 4) Slide the 3-phase switch to the lower position with the knee flexed.
 - The orthotic knee joint can be flexed freely at all times.

Please contact Ottobock Customer Service in case of questions.

5.2.2 Verification of fit

Checking the fit in the frontal plane

- > **Required tools:** 170W2 tool set, 170W13 aluminium clamping jaws, 170W11 frontal contour bending tool
 - 1) Check the fit of the base tube on the leg contour.
 - 2) Mount the aluminium clamping jaws on a vice.
 - 3) Adapt the base tube in the vice. Use the frontal contour bending tool for this purpose.

NOTICE! Avoid crushing the base tube by excessive shaping.

Checking the fit of the support tubes

Support tubes that have too small a diameter may result in medial pressure areas. Support tubes with too large a diameter could obstruct sitting, especially in the thigh area.

In either case, the patient's individual wishes should be taken into account while checking the fit of the support tubes.

- > **Required tools:** 170W13 aluminium clamping jaws, 170W12* pad holder bending tool, 170W18 bending iron
 - 1) Loosen the screw on the strap guide loop with tube connection by two turns.
 - 2) Push the screw in until a distinct clicking is heard. Use an Allen key to push it in.
 - 3) Unscrew the strap guide loop with tube connection.
 - 4) Remove the support tube. Loosen the two screws in order to do so.

- 5) Mount the aluminium clamping jaws on a vice.
- 6) Secure the support tube. Use the pad holder bending tool in order to do so (see fig. 12).
- 7) **NOTICE! Avoid crushing the support tube by excessive clamping.**
- 8) Correct the shape of the support tubes. Use the bending iron to do so.
- 9) Install the support tubes. In order to do so, secure the screws with Loctite® 241 and tighten them alternately until the tensioning wedge is flush with no gaps (tightening torque **4 Nm**).

Checking the fit in the sagittal plane

The fit of the orthosis in the sagittal plane is of crucial importance for the locking and unlocking function of the orthotic knee joint.

- > To verify technical functionality, have the patient sit on a chair. If necessary, help the patient unlock the orthotic knee joint so they can sit with the knee flexed.
 - ▶ Take the calf in your hand and bring it slowly into the extended position.
 - If the anatomical limit in the knee joint is reached before the orthosis extension stop, then the orthosis must be aligned with slightly more flexion.

The corrections for the standing position can be performed by bending the femoral support tubes.

- 1) Disassemble the axis of the orthotic knee joint.
 - The upper and lower parts are separated from each other.
- 2) Correct the fit in the sagittal plane in **2°** increments.

INFORMATION: Corrections of 2° or 4° may result in significant differences in the fit.

5.2.3 Checking freedom of movement in the ankle

NOTICE

Improper settings

Loss of functionality due to damaged steel cable

- ▶ Do not use the set screw of the steel cable on the ankle joint bottom part to adjust the lock.

- 1) Check the freedom of movement of the orthotic ankle joint.

INFORMATION: The orthotic ankle joint must have at least 5° mobility up to the dorsiflexion stop in order to ensure automatic unlocking of the orthotic knee joint.
- 2) Use the joint screw of the orthotic ankle joint to adjust the neutral position according to the scale on the upper part of the foot stirrup (see fig. 11).
 - **If the lock of the orthotic knee joint does not open at the end of the stance phase:** Adjust the foot stirrup in increments in plantar flexion.
 - **If the lock of the orthotic knee joint opens in the mid-stance phase:** Adjust the foot stirrup in increments in dorsal extension.
- 3) **INFORMATION: In a standing position, ensure that the orthotic ankle joint is not too far in dorsal extension or plantar flexion, as this can result in unintentional unlocking of the orthotic knee joint.**
- 4) **Optional:** Bend the foot stirrup with insert with the help of the 711S4=3x3.5 bending iron.

5.3 Finishing the orthosis

CAUTION

Unallowable use of lubricants

Risk of injury and loss of functionality as well as damage to the product

- ▶ The product is designed to be free of lubricants. Avoid contact with lubricants (e.g. grease, oil).
- ▶ After contact with lubricants, promptly clean the product or parts of the product with a degreasing cleaning agent (e.g. acetone or isopropyl alcohol).

> **Required tools:** Torque wrench

- 1) Secure the joint screw of the orthotic ankle joint with Loctite 241 and tighten it (tightening torque **14 Nm**).

- 2) Secure the screws of the foot stirrup with Loctite 241 and tighten them (tightening torque **2.5 Nm**).
- 3) **Optional:** Trim all hook-and-loop straps of the closure straps so sufficient material remains for attaching the Y-hook-and-loop closure.
- 4) Attach the Y-hook-and-loop closure.
- 5) **Optional:** Sew the Y-hook-and-loop closure to the loop strap.
- 6) To protect the orthotic ankle joint, apply the shrinkable tubing (see fig. 1, item 20) over the upper part of the ankle joint (item 13), lower part of the ankle joint (item 14) and foot stirrup (item 15).

6 Handling

WARNING

Operation of motor vehicles

Risk of accident due to limited body function

- ▶ The product has not been approved by the manufacturer for the operation of motor vehicles.
- ▶ Observe the applicable legal regulations for the operation of motor vehicles and have your driving ability examined by an authorised test centre.

CAUTION

Excessive strain on load-bearing components

Injuries due to changes in or loss of functionality

- ▶ Only use the product for the defined area of application.
- ▶ If the product has been exposed to extreme strain (e.g. due to falling), take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

CAUTION

Use on more than one patient

Skin irritation, formation of eczema or infections

- ▶ Only use the orthosis on one patient.

CAUTION

Improper handling

Risk of falling due to accidentally opening or engaging the lock

- ▶ Avoid accidentally opening or engaging the lock.
- ▶ Instruct the patient in the correct handling of the orthosis joints.

CAUTION

Incorrect or excessively tight application

Risk of local pressure and constriction of blood vessels and nerves due to improper or excessively tight application

- ▶ Ensure that the product is applied properly and fits correctly.

CAUTION

Using orthosis without padding

Risk of local pressure points, superficial injuries or skin irritation due to skin contact with rigid or sharp orthosis components

- ▶ Only use the orthosis with undamaged padding.

6.1 Trial walking

The Ottobock Free Walk orthosis deviates from conventional orthosis designs. This means that patients who have been using a conventional orthosis for a long time in particular will have to learn the principle.

For further information, please consult document 646A214 on therapeutic application and gait training.

7 Cleaning

CAUTION

Improper cleaning

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ Clean the product regularly.
- ▶ **Inform the patient.**

Clean the orthosis regularly:

- 1) Remove the straps from the orthosis.
- 2) Fasten all hook-and-loop closures.
- 3) Hand-wash the straps and pads in warm water at **30 °C (86 °F)** using neutral detergent. Rinse thoroughly.
- 4) Allow to air dry. Do not expose to direct heat sources (e.g. sunlight, stove or radiator).

8 Maintenance

The manufacturer requires at least a semi-annual inspection of the product to verify functionality and check for wear.

The wear and tear parts of the orthotic knee joints are to be replaced by the 170Z99=2, 170X18 and 170X17N=L/R maintenance sets (see fig. 2).

Pads and closures for replacement are shown in the product overview and scope of delivery table (see fig. 1).

Replacing pads

- 1) Loosen the screw on the strap guide loop with tube connection by two turns.
- 2) Push the screw in until a distinct clicking is heard. Use an Allen key to push it in.
- 3) Unscrew the strap guide loop with tube connection.
- 4) Remove and replace the pads.
- 5) Remove the strap guide loops from the end of the support tube.

Replacing the lock

> **Required materials and tools:** 636K13 Loctite, 2 mm Allen key, torque wrench adjustable to 0.5 Nm

- 1) Remove the retaining ring from the knee axis in the orthotic knee joint.
- 2) Push the knee axis out to the inside.
→ The orthosis is disassembled in the orthotic knee joint.
- 3) Loosen the lock screw of the steel cable with the Allen key.
- 4) Remove the lock screw and pressure ball from the orthotic ankle joint.
- 5) Remove the retaining ring of the lock axis.
- 6) Push the lock axis out to the inside.
- 7) Move the lock with steel cable up and out of the orthotic knee joint and replace the lock.
- 8) Insert the flange bushings into the lock by hand and press them in lightly using the vice.
- 9) After pressing them in, finish the flange bushings with a 6H7 reamer.
- 10) Knot the pull cable at the lock.

INFORMATION: Note that the knot has to be at least 5 cm away from the lock.

- 11) Put the lower part of the orthosis onto the workbench with the strap guide loops down.

- 12) Insert the steel cable into the hole in the orthotic knee joint and push it through the base tube.
INFORMATION: Note that the steel cable comes out the back of the base tube but before the orthotic ankle joint.
- 13) With one hand, pull the steel cable further out of the orthotic ankle joint and insert the lock in the orthotic knee joint with the other hand.
INFORMATION: Verify that it fits properly in the bearing bushings.
- 14) Position the lower spring of the lock so that it does not obstruct the steel cable.
- 15) Insert the lock axis from the medial side into the orthotic knee joint and slide on the retaining ring.
- 16) Bring the orthotic ankle joint forward to the stop and thread the steel cable from above through the orthotic ankle joint.
- 17) Insert the upper part of the orthosis with new bearing washers into the orthotic knee joint.
- 18) Slide the knee axis from the medial side into the orthotic knee joint and slide on the retaining ring.
→ Now the orthotic knee joint must be locked.
- 19) Insert the pressure ball into the orthotic ankle joint.
- 20) Secure the lock screw with Loctite 241 and install it. Do not tighten the lock screw.
- 21) Hold the steel cable with one hand while tensioning it slightly and lift the lock by **about 1 mm**. Tighten the lock screw at the same time (tightening torque **0.5 Nm**).
→ The lock has to unlock due to dorsiflexion in the orthotic ankle joint and extension in the orthotic knee joint.
- 22) Shorten the pull cable. The pull cable has to project **2 cm** out of the ankle joint.

9 Disposal

Dispose of the product in accordance with national regulations.

10 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

10.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

10.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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