
451F10, 451F11, 451F12, 451F13


DE Gebrauchsanweisung	3	CS Návod k použití.....	37
EN Instructions for use	5	RO Instrucțiuni de utilizare.....	39
FR Instructions d'utilisation.....	8	HR Upute za uporabu	42
IT Istruzioni per l'uso	11	SL Navodila za uporabo	44
ES Instrucciones de uso	13	SK Návod na používanie	47
PT Manual de utilização.....	16	BG Инструкция за употреба	49
NL Gebruiksaanwijzing.....	19	TR Kullanma talimatı	52
SV Bruksanvisning	21	EL Οδηγίες χρήσης	55
DA Brugsanvisning	24	RU Руководство по применению	57
NO Bruksanvisning	26	JA 取扱説明書.....	60
FI Käyttöohje.....	29	ZH 使用说明书.....	62
PL Instrukcja użytkowania.....	31	KO 사용 설명서.....	65
HU Használati utasítás	34		

1

2

1 Produktbeschreibung

Deutsch

INFORMATION

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INFORMATION

Date of last update: 2022-11-08

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

1.1 Function

The 451F10, 451F11, 451F12, 451F13 residual limb compression socks are intended to reduce or prevent oedemas on the residual limb. They are worn directly on the skin.

2 Intended use

2.1 Indications for use

The product is intended for compression therapy after an amputation.

- **451F11, 451F12:** The product is used after a transfemoral amputation or knee disarticulation.
- **451F10, 451F13:** The product is used after a transtibial or Syme amputation.

2.2 Contraindications

- Peripheral arterial occlusive disease (Doppler index less than 0.6 in non-amputated area)
- Decompensated heart failure
- Septic phlebitis
- Phlegmasia cerulea dolens
- Weeping dermatosis
- Severe sensory disturbances
- Intolerance and allergic reactions to the materials used

2.3 Lifetime

The product is a wear part, which means it is subject to normal wear and tear.

3 General safety instructions

Risk of injury and risk of product damage

- ▶ Keep the product away from open flames, embers and other sources of heat.
- ▶ The product may only be used by one person; reuse on another person is prohibited.
- ▶ Observe the cleaning instructions.

Risk of product damage and limited functionality

- ▶ Check the product for damage and readiness for use prior to each use.
- ▶ Do not use the product if its functionality is restricted. Take suitable measures (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop).
- ▶ Ensure the product is functioning reliably and check for wear and tear as well as damage before each use.
- ▶ Do not continue using a product that is no longer functional or that is worn or damaged.
- ▶ Do not expose the product to prohibited environmental conditions.
- ▶ Check the product for damage if it has been exposed to prohibited environmental conditions.
- ▶ Do not use the product if it is damaged or in a questionable condition. Take suitable measures (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop).

4 Scope of delivery

Quantity	Designation	Reference number
1	Instructions for use	–
1	Residual limb sock	–

5 Use

5.1 Information for use

The course of compression therapy (for example, the timing, compression class, duration of use) is determined by the attending doctor.

- As the oedema decreases, the size of the residual limb sock has to be adjusted to ensure continuous compression
- Regular examination by a doctor is essential to ensure treatment success.

5.2 Donning and doffing the residual limb sock

Donning

- > **Prerequisite:** The residual limb sock is dry.
- > **Prerequisite:** Remove all sharp objects from the hands and wrists (such as rings, bracelets or watches).
- 1) **Recommendation:** Put on thin gloves to prevent damage due to rough or sharp fingernails.
- 2) Turn the residual limb sock inside out, almost to the end.
- 3) Slowly slip the sock over the residual limb without twisting it.
- 4) Use a massaging motion to evenly distribute the tension of the material. This optimises the effect.
- 5) **If the sock has a fastening strap:** Fasten the strap around the hips and make sure it is not too tight.

Doffing

- > **Prerequisite:** Remove all sharp objects from the hands and wrists (such as rings, bracelets or watches).
- 1) Roll the residual limb sock up from the open end to the closed end and take it off the residual limb.
- 2) **Recommendation:** Clean the residual limb sock (see page 7).

5.3 Cleaning and Care

- 1) **If present:** Fasten all closures
- 2) **NOTICE! To avoid damaging the product, do not wring out or use the spin cycle.**
Wash the product by hand or in the washing machine (gentle cycle) (temperature: **40 °C (104<nbsp>/>°F)**, mild detergent, no fabric softener).
- 3) Rinse the product with clear water.
- 4) Roll the product in a towel and squeeze out water.
- 5) Leave the product to dry lying flat on a hand towel or drying rack.

6 Disposal

Do not dispose of the product with unsorted household waste. Improper disposal can be harmful to health and the environment.

Observe the information provided by the responsible authorities in your country regarding return, collection and disposal procedures.

7 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

7.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

7.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

8 Technical data

451F11, 451F12						
Size	XS	S	M	L	XL	XXL
Compression class	451F11: CCL2, 451F12: CCL1					
Length [cm]	20/25/30/35					
Circumference, proximal [cm]	41 to 44	44 to 48	48 to 52	52 to 56	56 to 60	60 to 64
Circumference, distal [cm]	29 to 31	31 to 34	34 to 37	37 to 40	40 to 43	43 to 46

451F10, 451F13						
Size	XS	S	M	L	XL	XXL
Compression class	451F10: CCL2, 451F13: CCL1					
Length [cm]	30/38/46					
Circumference, proximal [cm]	39 to 41	41 to 44	44 to 47	47 to 50	50 to 53	60 to 64
Circumference, knee [cm]	29 to 31	31 to 34	34 to 37	37 to 40	40 to 43	43 to 46
Circumference, distal [cm]	27 to 29	29 to 32	32 to 35	35 to 38	38 to 41	41 to 44

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