12K41, 12K42, 12K44, 12K50

DE	Gebrauchsanweisung (Fachpersonal)	2
EN	Instructions for use (qualified personnel)	23
FR	Instructions d'utilisation (Personnel spécialisé)	44
IT	Istruzioni per l'uso (Personale tecnico specializzato)	66
ES	Instrucciones de uso (Personal técnico especializado)	88
PT	Manual de utilização (Pessoal técnico)	110
NL	Gebruiksaanwijzing (Vakmensen)	132
SV	Bruksanvisning (Fackpersonal)	153
HU	Használati utasítás (szakszemélyzet)	174
RU	Руководство по применению (Квалифицированный персонал)	195
JA	取扱説明書 (有資格担当者)	217

1 Vorwort

INFORMATION

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INFORMATION

Date of last update: 2022-10-19

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- ▶ Please keep this document for your records.

The "12K41=* ErgoArm, 12K42=* ErgoArm plus, 12K44=* ErgoArm Hybrid plus, 12K50=* ErgoArm Electronic plus" product is referred to simply as the product/elbow joint below.

These instructions for use provide you with important information on the use, adaptation and handling of the product.


Only put the product into use in accordance with the information contained in the accompanying documents supplied.

The product may not be transferred to the patient without prior instruction.

2 Product description

2.1 Construction and Function

2.1.1 Design

1. Easy Plug
2. Strap clamp
3. Bowden cable (mechanical lock)
4. Dial
5. Humeral rotation feature adjusting screw

The product is a mechanically/electronically driven elbow joint. It helps the wearer perform everyday tasks in combination with other prosthetic components (see page 25).

2.1.2 Function

Various control programmes make optimal adaptation to the individual needs and abilities possible.

The ErgoArm models differ in the following characteristics:

Models	AFB	Easy Plug	Lock	Slip-stop
12K41=* ErgoArm	-	-	Mechanical	Mechanical
12K42=* ErgoArm plus	x	-	Mechanical	Mechanical
12K44=* ErgoArm Hybrid plus	x	x	Mechanical	Mechanical
12K50=* ErgoArm Electronic plus	x	x		


Models	AFB	Easy Plug	Lock	Slip-stop
			Mechanical or Electronic	Mechanical or Electronic

2.1.2.1 Definition of terms

AFB (Automatic Forearm Balance)

The AFB flexion aid stores the energy that is released when the prosthetic arm is extended and uses it to assist with flexion. The complex mechanism of the flexion aid controls the dynamics of forearm movement and equalises the leverage forces that develop via a compensation gear. As a result, the patient only has to exert minimal force on the flexion cable or through body momentum. The compensation coefficient can be easily adjusted to the individual weight of the prosthetic forearm and different clothing.

Easy Plug

1. Battery connection
2. CLOSE electrode connection
3. OPEN electrode connection
4. Switch connection

The Easy Plug is the control electronics for the 12K44=* and 12K50=* elbow joints.

The connecting cables of the electrodes, switches and batteries are connected to the Easy Plug. Its interior cable routing minimises the risk of cable breaks.

Lock

The product is either locked mechanically or electronically (see page 23). The internal ratchetless lock can be unlocked or locked in any position, even under a load. The product has a load capacity of up to 230 N in the locked state (this value may decline over the course of several years due to wear and tear). The lock will slip under higher loads.

Mechanical lock

A strong pull on the cable (audible clicking sound) locks or unlocks the joint.

Thanks to the integrated compensation mechanism, the lock can be released with little effort, even under load. The integrated strain relief reduces the risk of damage when the pull cable is pulled too hard.

Electronic lock

The lock can be released and engaged using electrode signals or switch signals (depending on the chosen switching version).

Slip-stop function

The slip-stop function makes it possible to lower the forearm in a controlled manner without having to fully release the lock.

Slip-stop function, mechanical

A gentle pull on the cable (no audible click) unlocks the product so the forearm can be moved continuously. When the pull on the cable stops, the product is locked again.

Slip-stop function, electronic

The slip-stop function is activated by the electrode signals (depending on the chosen control version). When the electrode signal stops, the product is locked again.

2.2 Combination possibilities

Wrist joints

- MovoWrist Flex: 10V39=*
- Ball ratchet wrist joint: 10V8
- Ottobock wrist joint: 10V18=*
- Robo-Wrist, mechanical wrist joint: 10V41=*
- Electric wrist rotator: 10S17 (only for 12K44=* and 12K50=*)
- MyoRotronic: 12E205 (only for 12K44=* and 12K50=*)
- VC hand: 8K27=*
- Ottobock system hand -passive-: 8K18=*
- Ottobock system hand -passive-: 8K19=*
- All-purpose hook: 10A3
- All-purpose ring: 10A4
- All-purpose hook: 10A12=*
- MovoHook 2Grip: 10A71=*
- MovoHook 2Grip: 10A81=*
- System Electric Hand: 8E38=* (only for 12K44=* and 12K50=*)
- System Electric Greifer: 8E33=* (only for 12K44=* and 12K50=*)
- bebionic hand EQD: 8E70=* (only for 12K44=* and 12K50=*)

Terminal devices

- VO hand: 8K22=*
- VO hand: 8K23=*
- VC hand: 8K26=*

3 Intended use

3.1 Indications for use

The product is intended **exclusively** for upper limb exoprosthetic fittings.

3.2 Conditions of use

The product is designed for use by unilateral or bilateral amputees.

The product is intended **exclusively** for use on adults.

The product was developed for everyday use and must not be used for unusual activities. These unusual activities include, for example, extreme sports (free climbing, parachuting, paragliding, etc.).

3.3 Contraindications

- All conditions which contradict or go beyond the specifications listed in the section on "Safety" and "Indications for use".

3.4 Qualification

The fitting of a patient with the product may only be carried out by O&P professionals who have been authorised with the corresponding Ottobock training.

4 Safety

4.1 Explanation of warning symbols

 WARNING	Warning regarding possible serious risks of accident or injury.
 CAUTION	Warning regarding possible risks of accident or injury.
 NOTICE	Warning regarding possible technical damage.

4.2 Structure of the safety instructions

CAUTION

The heading describes the source and/or the type of hazard

The introduction describes the consequences in case of failure to observe the safety instructions. Consequences are presented as follows if more than one consequence is possible:

- > E.g.: Consequence 1 in case of failure to observe the hazard
- > E.g.: Consequence 2 in case of failure to observe the hazard
- ▶ This symbol identifies activities/actions that must be observed/carried out in order to avert the hazard.

4.3 General safety instructions

WARNING

Non-observance of safety notices

Personal injury/damage to the product due to using the product in certain situations.

- ▶ Observe the safety notices and the stated precautions in this accompanying document.

WARNING

Use of the product while operating a vehicle or machinery

- > Risk of accident due to unexpected behaviour of the product.
- > Risk of injury due to faulty control or malfunction of the product.
- ▶ Inform the patient that observing the applicable national legal regulations is mandatory when operating motor vehicles and machinery.
- ▶ Inform the patient that their driving ability must be examined and confirmed by an authorised test centre.
- ▶ Inform the patient that the product must be switched off before operating a vehicle or machinery.

WARNING

Charging the prosthesis without taking it off

Risk of electric shock due to defects in the power supply unit or in the battery charger.

- ▶ For safety reasons, remove the prosthesis prior to charging it.

CAUTION

Changing prosthetic components when switched on

Injury due to faulty control or malfunction of the prosthesis system.

- ▶ Before changing prosthetic components (such as the terminal device), remove the rechargeable battery from the battery receptacle or switch off the prosthesis system.

CAUTION

Overloading due to unusual activities

Injury due to unexpected product behaviour as a result of a malfunction.

- ▶ The product was developed for everyday use and must not be used for unusual activities. These unusual activities include, for example, sports with excessive strain and/or shocks to the wrist joint (pushups, downhill mountain biking, ...) or extreme sports (free climbing, paragliding, etc.).
- ▶ Careful handling of the product and its components not only increases their service life but, above all, ensures the patient's personal safety!
- ▶ If the product and its components have been subjected to extreme loads (e.g. due to a fall, etc.), then the product must be inspected for damage immediately. If necessary, forward the product to an authorised Ottobock Service Centre.

⚠ CAUTION**Proximity to sources of strong magnetic or electrical interference (e.g. theft prevention systems, metal detectors)**

Injury due to unexpected behaviour of the product caused by interference with internal data communication.

- ▶ Avoid remaining in the vicinity of visible or concealed theft prevention systems at the entrance/exit of stores, metal detectors/body scanners for persons (e.g. in airports) or other sources of strong magnetic and electrical interference (e.g. high-voltage lines, transmitters, transformer stations, computer tomographs, magnetic resonance tomographs, etc.).
- ▶ When walking through theft prevention systems, body scanners or metal detectors, watch for unexpected behaviour of the product.

⚠ CAUTION**Signs of wear and tear on the product components**

Injury due to faulty control or malfunction of the product

- ▶ In the interest of preventing patient injuries and to maintain product quality, we recommend performing regular service.
- ▶ Please see the service pass for more information about the service intervals. This is included in the purchased extended warranty packages.

⚠ CAUTION**Incorrect electrode settings due to muscle fatigue**

Injury due to faulty control or malfunction of the product.

- ▶ Allow the patient to rest during the adjustment of the electrodes.

⚠ CAUTION**Penetration of dirt and humidity into the product**

Injury due to faulty control or malfunction of the product.

- ▶ Ensure that neither solid particles nor liquids can penetrate into the product.

⚠ CAUTION**Use of inappropriate components**

Risk of injury due to malfunction of the product.

- ▶ Operate the product only with components recommended by the manufacturer. You can find the list of components in the "Combination possibilities" section.
- ▶ Operate the product only with options that have been recommended by the manufacturer. You can find the list of components in the "Scope of delivery and accessories" section.

⚠ CAUTION**Independent user manipulation of system components**

Injury due to faulty control or malfunction of the product.

- ▶ Manipulations to the product other than the tasks described in these instructions for use are not permitted.
- ▶ The product and any damaged components may only be opened and repaired by authorised, qualified Ottobock personnel.

⚠ CAUTION**Manual unlocking of elbow lock under load**

Injury by release of elbow lock under load.

- ▶ Particular caution should be exercised when unlocking the elbow lock while lifting heavy loads.
- ▶ Be careful when unlocking the lock under such conditions due to the possibility of injury.

⚠ CAUTION

Damage to connecting cable due to kinking or small radii

Risk of injury due to unexpected product behaviour as a result of a malfunction.

- ▶ Make sure there are no tight angles or kinks when installing the cables.
- ▶ Replace the connecting cable immediately if it is damaged.

⚠ CAUTION

Charging the product with damaged power supply unit/charger/charger cable

Injury due to unexpected behaviour of the product caused by insufficient charging.

- ▶ Check the power supply unit, charger and charger cable for damage before use.
- ▶ Replace any damaged power supply unit, charger or charger cable.

⚠ CAUTION

Sudden flexion of the product due to high level of compensation

Risk of injury due to sudden flexion of the prosthesis.

- ▶ Note the level of compensation when taking off heavy clothing.
- ▶ Note the setting for the level of compensation when the prosthesis is not being worn.
- ▶ The prosthesis must always be flexed for donning and doffing.

⚠ CAUTION

Incorrect donning and doffing of the prosthesis

Risk of injury due to sudden flexion of the prosthesis.

- ▶ Switch off the prosthesis before donning or doffing.
- ▶ The prosthesis must always be flexed for donning and doffing.

NOTICE

Coating, gluing or painting the prosthesis

Damage or fracture due to chemical processes.

- ▶ The prosthesis must not be coated, glued or painted.

NOTICE

Damage to the battery

Damage to the battery due to using both connection options.

- ▶ Only use one of the two connections for the battery (connection on the forearm cable or Easy Plug connection).

NOTICE

Improper product care

Damage to the product due to the use of incorrect cleaning agents.

- ▶ Only clean the product with a damp cloth and mild soap (e.g. 453H10=1-N Ottobock DermaClean).

5 Scope of Delivery and Accessories

12K41 ErgoArm plus scope of delivery

- 1 pc. 12K41=* ErgoArm
- 1 pc. Lamination protection cover (disc with arrow)
- 1 pc. Cable lock
- 2 pc. bracket cover
- 1 pc. 60x2.5 O-ring
- 1 pc. 647G1475 Instructions for use (user)
- 1 pc. 647H437 Instructions for use (qualified personnel)

12K42 ErgoArm plus scope of delivery

- 1 pc. 12K42=* ErgoArm plus
- 1 pc. Lamination protection cover (disc with arrow)
- 1 pc. Cable lock
- 2 pc. bracket cover
- 1 pc. 60x2.5 O-ring
- 1 pc. 647G1475 Instructions for use (user)

- 1 pc. 647H437 Instructions for use (qualified personnel)

12K44 ErgoArm Hybrid plus scope of delivery

- 1 pc. 12K44=* ErgoArm Hybrid plus
- 1 pc. Lamination protection cover (disc with arrow)
- 1 pc. Cable lock
- 2 pc. bracket cover
- 1 pc. 60x2.5 O-ring
- 1 pc. Lamination protection cover (hollow ball)
- 1 pc. 647G1475 Instructions for use (user)
- 1 pc. 647H437 Instructions for use (qualified personnel)

12K50 ErgoArm Electronic plus scope of delivery

- 1 pc. 12K50=* ErgoArm Electronic plus
- 1 pc. Lamination protection cover (disc with arrow)
- 1 pc. Cable lock
- 2 pc. bracket cover
- 1 pc. 60x2.5 O-ring
- 1 pc. Lamination protection cover (hollow ball)
- 1 pc. spherical cord lock
- 1 pc. 647G1475 Instructions for use (user)
- 1 pc. 647H437 Instructions for use (qualified personnel)

Accessories

- Alignment aid for ErgoArm: 743A23
- Clamp stopple set: 21A207
- Adapter: 13Z68
- Battery connecting cable: 13E51=2
- Battery receptacle set, screwable: 757Z184=2, 757Z190=2
- Battery receptacle set: 757Z184=1, 757Z190=1

Options for programme settings of the 12K50

- 1 pc. 757T13 MyoSelect
 - Coding plug set: 13E182
- or

6 Preparation for use

6.1 Required machines, tools and materials

- 633F11 silicone grease
- 633F30 Special grease
- 701P250=1 Belt sander (1x 230 V, 50/60 Hz)
- 756B2 Electric jigsaw (1x 230 V N/PE, 50 Hz)
- Set of screwdrivers
- 636W60 Adhesive (Loctite 243)
- Pattex Repair
- 89Z1 Tweezers
- Knives

6.2 Laminating instructions

INFORMATION

Laminate the socket according to the state of the art (see 646T333 technical information).

6.3 Removing the lamination ring

Complete the following steps to remove the lamination ring:

1) Set the compensation to the minimum level using the dial.

2) Secure the lifter strap.

3) Turn the friction screw until the lamination ring can be unscrewed from the elbow ball.


- 4) Completely unscrew the lamination ring from the elbow ball.
- 5) Carefully remove the electronics housing from the lamination ring.

NOTICE
Do not remove the lifter strap
The lifter strap must not be removed from the elbow ball. Fixing the lifter strap in place (e.g. with a haemostat) is no longer necessary.

6.4 Sealing the product

Complete the following steps to seal the product:

- 1) Clean the lamination ring so it is free of grease and dirt.


- 2) Lubricate the groove, inner edges and thread of the lamination ring with special grease.

Additional steps for the products 12K44 and 12K50

NOTICE
Misalignment during installation of the electronics housing
Misalignment during installation results in the risk of loss of contact or short-circuiting of the spring contacts. The mechanical strength can also be reduced due to thread misalignment. The circuit board and the edge of the lamination ring must be parallel to each other.

Complete the following additional steps to seal the product:

- 1) Thoroughly lubricate the plugs of the connecting cables (electrodes, switch, battery) with silicone grease.
- 2) Remove the blanking covers from the electronics housing.
- 3) Connect the plugs to the electronics housing.
- 4) Thoroughly seal the plugs and remaining blanking covers from above with silicone grease.
- 5) Push the electronics housing, including the connected cables, into the lubricated lamination ring and press firmly.

INFORMATION: A loud popping sound must be audible.

- 6) Check the fit of the electronics housing.
INFORMATION: The circuit board of the electronics housing and the edge of the lamination ring are flush/parallel to each other.


- 7) Connect the desired programming plug to the electronics housing if necessary.

6.5 Removing the elbow joint

Removing the elbow joint for the 12K41 product

Complete the following additional steps to remove the elbow joint:

- 1) Extend the arm.

- 2) Loosen both joint caps.
- 3) Loosen the mounting screws from the elbow joint.
- 4) Carefully pull the elbow joint off the forearm.


Removing the elbow joint for the 12K42, 12K44, 12K50 products

Complete the following steps to remove the elbow joint:

- 1) Flex the arm.
- 2) Set the compensation to the minimum level using the dial.
- 3) Extend the arm again.
- 4) Lock the elbow joint in this position to fix the removal and installation position.

- 5) Loosen both joint caps.
- 6) Pry out the cable plugs.

- 7) Disconnect the cables.

- 8) Loosen the mounting screws from the elbow joint.
- 9) Unscrew the strap clamp.
- 10) Pull out the strap.
- 11) Carefully pull the elbow joint off the forearm.


- 12) Screw the elbow joint into the upper arm.

6.6 Removing the AFB flexion aid

Complete the following steps to remove the AFB flexion aid:

- 1) Unscrew the two mounting screws on the forearm.

- 2) Pull out the AFB flexion aid.

6.7 Shortening the forearm

INFORMATION

Wrap the elbow area with plastic wrap to prevent swarf from entering it during cutting.

Complete the following steps to shorten the forearm:

- 1) Add 10 mm to the length that was determined and mark this on the forearm.

INFORMATION: The forearm length is measured from the distal end of the lamination ring to the middle of the elbow axle.

- 2) Shorten the forearm at the distal end to the patient's individual dimensions (trim or grind).

INFORMATION: The minimum length of the forearm must be 220 mm after shortening (measured from the elbow axle).

- 3) Sand all sides of the cut edges at a right angle.
- 4) Deburr and slightly round the inside and outside edges.

6.8 Optional: Positioning the battery receptacle

Complete the following steps to position the battery receptacle:


- 1) Position the battery receptacle.

INFORMATION: Make sure the battery receptacle does not collide with the AFB flexion aid and electric wrist rotator.

- 2) Connect the plug of the battery connecting cable to the circuit board of the AFB flexion aid.

6.9 Installing the AFB flexion aid

Complete the following steps to install the AFB flexion aid:

- 1) Insert the AFB flexion aid in the forearm.

- 2) Tighten the mounting screws on the forearm.

6.10 Converting the cable lock

Complete the following steps to convert the cable lock:

17

- 1) Pull the ball cap off the elbow ball.
- 2) Disengage the telescopic spring.
- 3) Remove the ladder lock buckle.

18

- 4) Pull out the telescopic spring.

19

- 5) Remove the dummy and connect it again on the other side.


20

- 6) Thread the telescopic spring through on the other side.
- 7) Engage the telescopic spring on the protruding metal axle.
- 8) Put the ball cap back on.
- 9) Reinstall the ladder lock buckle.

6.11 Finishing the switch cable

Complete the following additional steps to finish the switch cable:

- 1) Thread the switch cable through the cable lock.
- 2) Secure the end of the switch cable with a button.
- 3) Tighten the switch cable and make a loop.

6.11.1 Shortening the pull cable (only for 12K50)

Complete the following steps to shorten the pull cable:

- 1) Trim the pull cable to the desired length.
- 2) Tie off the end of the pull cable with the cord lock.

⚠ CAUTION

Removing the pull cable

Injuries due to loss of product functionality.

- ▶ For safety reasons, the pull cable must not be removed entirely.

6.12 Installing the elbow joint

NOTICE

Use of the product with a damaged flat cable

Malfunction due to incorrect or interrupted data transfer.

- ▶ Avoid bending or pressing the flat cable multiple times in the same position.

Complete the following steps to install the elbow joint:

- 1) Push the elbow joint into the flat toothed disc on the forearm.

INFORMATION: Do not pinch the existing cables (12K44, 12K50). The cables are exposed and may be damaged if they do not have cover caps.

- 2) Apply Loctite 243 to the two elbow joint mounting screws and tighten them evenly.
- 3) **12K44/12K50:** Lubricate the wiring connectors with silicone grease.
- 4) **12K44/12K50:** Reconnect both cables.
- 5) Extend the arm.
- 6) Apply Pattex Repair to the brackets.

- 7) Put the joint cap on the brackets.
INFORMATION: Do not pinch the existing cables (12K44, 12K50).

Additional steps: installing the elbow joint for the 12K42, 12K44, 12K50 products

- 1) Thread the strap under axle.
INFORMATION: Make sure the lamination ring is correctly positioned.
INFORMATION: There needs to be space for the strap clamp between the strap clamp guide and the lamination ring.
INFORMATION: Ensure the strap is straight and not twisted.
- 2) Slide the strap clamp under the lamination ring.

- 3) Secure the strap clamp with the screw.

6.13 Setting the preflexion

Complete the following steps to set the preflexion:

- 1) Flex the arm.
- 2) Set the compensation to the minimum level using the dial.
- 3) Extend the forearm.

INFORMATION: The forearm has to hang down freely.

INFORMATION: The preflexion angle is initially set at about 10°.

- 4) Secure the strap cam on the extension stop.
- 5) Flex the forearm until the strap forms a loop.
- 6) Lock the joint with the cable lock.
- 7) Loosen the screw on the strap clamp.
- 8) Shorten or lengthen the strap.

INFORMATION: Lengthening the strap decreases the preflexion. Shortening the strap increases the preflexion.

INFORMATION: A 3 mm change in the length of the strap changes the flexion angle by about 5°. The free swing range is changed by the same amount as the flexion angle.

- 9) Tighten the screw on the strap clamp.
- 10) Free the strap cam.

INFORMATION: For the flexion aid to function, the strap cam has to be up against the extension stop when the forearm is extended.

6.14 Fastening the flexion cable

Complete the following steps to fasten the flexion cable:

- Fasten the flexion cable inside the marked area.

INFORMATION: The end of the flexion pull cable must not be drawn into the socket when the forearm is released.

6.15 Installing the elbow joint in the upper arm socket

Complete the following steps for installation in the upper arm socket:

- 1) Apply silicone grease to the exterior thread of the elbow joint.
- 2) Screw the elbow joint into the upper arm.

7 Handling

7.1 Adjusting the level of compensation

CAUTION

Incorrect donning and doffing of the prosthesis

Risk of injury due to sudden flexion of the prosthesis.

- Switch off the prosthesis before donning or doffing.
- The prosthesis must always be flexed for donning and doffing.

CAUTION**Sudden flexion of the product due to high level of compensation**

Risk of injury due to sudden flexion of the prosthesis.


- ▶ Note the level of compensation when taking off heavy clothing.
- ▶ Note the setting for the level of compensation when the prosthesis is not being worn.
- ▶ The prosthesis must always be flexed for donning and doffing.

Complete the following steps to set the level of compensation:

- 1) Flex the arm.
 - 2) Use the dial to adjust the compensation.
- INFORMATION: The adjustment range is limited by stops.**
INFORMATION: The adjustment mechanism is protected by a slip clutch.

There are three different compensation ranges:

**Range 1 = low compensation level**

Permits the arm to swing freely while walking.

Range 2 = moderate compensation level

Balance is positively supported when flexing the arm. Compensation is reduced again during extension.

Range 3 = constant compensation level

If set correctly, the weight of the forearm is balanced by the compensation level.

Range 4 = low compensation level

Low compensation until flexion stop.

7.2 Adjusting the upper arm rotation joint

Complete the following steps to rotate and adjust the upper arm rotation joint:

- 1) Loosen the adjusting screw.
- 2) Move the upper arm rotation joint to the desired position.
- 3) Re-tighten the adjusting screw.

7.3 Adjusting the electrodes (only for 12K50)

The electrode setting depends on the chosen electric hand and the components that are used (e.g. 12E205 MyoRotronic). Observe the supplied instructions for use for the electrodes, electric hand and other components when making settings.

7.4 Program settings (only for 12K50)

Program settings with the MyoSelect

Complete the following steps to make program settings with the MyoSelect:

- 1) Remove the battery or turn off the MyoEnergy Integral.
- 2) Remove the electric wrist rotator and MyoRotronic if present.
- 3) Connect the 9E169 coaxial plug to the electronics housing.
- 4) Connect the MyoSelect with the 9E169 coaxial plug and 12K50.
- 5) Select the appropriate program in the MyoSelect.
- 6) Disconnect the 9E169 coaxial plug from the electronics housing.
- 7) Reinstall the electric wrist rotator and MyoRotronic if applicable.
- 8) Install the battery or switch on the MyoEnergy Integral.

Program settings using the programming plugs

Complete the following steps to make program settings with the programming switch:

- ▶ Connect the appropriate program switch (see page 40) to the electronics housing.

7.5 Program switching (only for 12K50)

The following switching methods can be selected:

- Switch
- Electrodes (generating a co-contraction)

7.5.1 Switching modes

The following control programs are available for switching:

Control program	Electronic elbow lock	Slip-stop function	Switching between hand and rotation
1 (white)	Press and release switch: - Lock - Unlock	No	All MyoRotronic programs
2 (red)	Co-contraction	No	Program 1 and program 4 of the MyoRotronic
3 (green)	Press and hold switch (elbow mode) - CLOSE electrode = lock - OPEN electrode = unlock Release switch (hand mode)	Yes	All MyoRotronic programs
4 (blue)	Press and release switch (elbow mode) - CLOSE electrode = lock - OPEN electrode = unlock Press and release switch (hand mode)		
5 (yellow)	Press and release switch (elbow mode) - CLOSE electrode = lock - OPEN electrode = unlock No electrode signal for 10 seconds (hand mode) Press and release switch (hand mode)		
6 (purple)	Co-contraction (elbow mode)	Yes	Program 1 and program 4 of the MyoRotronic

Control program	Electronic elbow lock	Slip-stop function	Switching between hand and rotation
	- CLOSE electrode = lock - OPEN electrode = unlock Co-contraction (hand mode)	Yes	Program 1 and program 4 of the MyoRotronic
7 (orange)	Co-contraction (elbow mode) - CLOSE electrode = lock - OPEN electrode = unlock No electrode signal for 10 seconds (hand mode) Co-contraction (hand mode)		

8 Charging the battery

The following points must be observed when charging the battery:

- The capacity of a fully charged battery is sufficient for one full day.
- We recommend charging the product once a day when used by the patient on a daily basis.
- The battery should be charged for at least 3 hours prior to initial use.

8.1 757L20 battery charger and 757B2* EnergyPack

NOTICE

Deep discharge of the battery

Irreparable damage to the product and resulting loss of functionality in case of irregular and incomplete charging of the battery.

- ▶ In case of long-term storage, fully charge the battery every 4 to 6 months.

- 1) Plug the charging station with power supply unit into the wall socket.
→ The green LED on the power supply unit and the green LED in the middle of the LED row light up.
- 2) Insert one or two rechargeable batteries into the charging station.
- 3) The green LED for the respective charging socket lights up ● and the charging process starts.
- 4) When the rechargeable battery is fully charged, the green LED for the respective charging socket flashes ●.
- 5) Remove the rechargeable battery after the charging process is complete.

8.2 757L35 battery charger and 757B35=* rechargeable battery

- 1) Connect the charging plug to the charging receptacle on the prosthetic socket.
→ A correct connection between the battery charger and the prosthesis is indicated by feedback (Status signals).
- 2) The charging process starts.
→ The prosthesis is switched off automatically.
- 3) Disconnect the prosthesis when the charging process is complete.

9 Cleaning

- 1) Clean the product with a damp cloth and mild soap (e.g. Ottobock 453H10=1-N Derma Clean) when needed. Ensure that no liquid penetrates into the system component(s).
- 2) Dry the product with a lint-free cloth and allow it to air dry fully.

10 Maintenance

Performing regular maintenance (service inspections) every 24 months is recommended to prevent injuries and maintain the quality of the product.

The grace period is no more than one month before or three months after maintenance is due.

In general, all products are subject to compliance with the maintenance intervals during the warranty period. This is the only way to maintain full warranty cover.

Additional services such as repairs may be provided in the course of maintenance. These additional services may be provided free of charge or can be billable according to an advance cost estimate, depending on the extent and validity of the warranty.

The following components must always be sent in for maintenance and repairs:

The product, battery charger and power supply. The shipping container for the loaner unit you receive must be reused for sending back the components requiring inspection.

11 Legal information

11.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

11.2 Trademarks

All product names mentioned in this document are subject without restriction to the respective applicable trademark laws and are the property of the respective owners.

All brands, trade names or company names may be registered trademarks and are the property of the respective owners.

Should trademarks used in this document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

11.3 CE conformity

Otto Bock Healthcare Products GmbH hereby declares that the product is in compliance with European requirements for medical devices.

The product meets the requirements of the RoHS Directive 2011/65/EU on the restriction of the use of certain hazardous substances in electrical and electronic devices.

The full text of the regulations and requirements is available at the following Internet address: <http://www.ottobock.com/conformity>

12 Technical data

Ambient conditions	
Storage (with and without packaging)	+5 °C/+41 °F to +40 °C/+104 °F Max. 85% relative humidity, non-condensing
Transport (with and without packaging)	-20 °C/-4 °F to +60 °C/+140 °F Max. 90% relative humidity, non-condensing
Operation	-5 °C/+23 °F to +45 °C/+113 °F Max. 95% relative humidity, non-condensing

Reference number	12K41=*	12K42=*	12K44=*	12K50=*
Current draw (in mA)	-	-	-	3.5
Power supply	-	-	-	757B20 757B21 757B25=* 757B35=*

Reference number	12K41=* 12K42=* 12K44=* 12K50=*
Max. load in locked state (N)	230
Flexion angle (in °)	15°–145°
Service life	5 years

13 Appendices

13.1 Symbols Used

Declaration of conformity according to the applicable European directives

Serial number (YYYY WW NNN)
 YYYY – year of manufacture
 WW – week of manufacture
 NNN – sequential number

In some jurisdictions it is not permissible to dispose of these products with unsorted household waste. Disposal that is not in accordance with the regulations of your country may have a detrimental impact on health and the environment. Please observe the instructions of your national authority pertaining to return and collection.

Medical device

Manufacturer

13.2 Display of the current charge level

13.2.1 757B2* EnergyPack

LED indicator	Event
	Charge level above 50% (illuminated in green)
	Charge level below 50% (flashing alternating green and orange)
	Charge level below 5% (illuminated in orange)

When inserting/engaging the rechargeable battery in the prosthesis, a charge level indicator is activated on the battery for a few seconds.

13.2.2 757B35=* Battery

The battery charge level can be queried at any time.

- 1) With the prosthesis switched on, press the charging receptacle button and hold for less than one second.
- 2) The LED display on the charging receptacle provides information on the current battery charge level (Status signals).

Charging receptacle	Event
	Battery fully charged (illuminated in green)
	Battery 50% charged (illuminated in yellow)
	Battery drained (illuminated in orange)

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